





Please read these instructions carefully before beginning the installation. Failure to understand and follow installation instructions may result in injury to installer and/or end user and may void the warranty. If you have any questions call Sunrise Medical Technical support at 800-333-4000.

Included Parts:

- 1. 1 ea Legrest Harness
- 2. 4 ea Zip-tie

Tools needed:


1. Wire cutter

A. Preparing for installation:

- 1. Remove the rear shroud by detaching the Velcro® at the sides and lifting off from the blocks at the top of the frame.
- 2. Cut the zip-tie (A) holding the wires to the plate, Detail 1a.
- 3. Unplug all the wires from the top and bottom rows of the cxsm. Ensure that the labels on the wires are not disturbed. Fig 2.


B. Installing the Harness:


- 1. Route connectors from the Power Center Mount towards the wire clips using a smooth radius (do not kink wires), Fig 3.
- 2. Zip-tie the harness to the seat behind the legrest connectors, Fig 4.
- 3. Route the harness through the harness clips under the seat. Fig 5a and 5b.
- 4. CAUTION! Continue routing harness between the seat pan and the tilting links, Fig 6.
- 5. Plug the bottom row of harnesses into the cxsm first in this order: Leg Lift (yellow label) first, and Leg Extend (white label) last. Then plug in the feedback harnesses (top row) in this order: Recline (label 1) first, Tilt or Tilt/Lift (label 2) second, and Leg (label 3) last. Make sure the graphic or number label on the wires matches the corresponding graphic or number labels on the cxsm. Note- Not all chairs will have the same amount of wires. Yours may differ from the one shown. Figures 7 & 8.
- 6. Ensure that the wires are tucked in out of the way, Fig 9.
- 7. Use supplied zip-tie (A) to secure wires to the plate. Do not zip-tie Recline Position 1 top row. Detail 1a & 9.
- 8. Reinstall the rear shroud by lining up the shroud to the blocks on the top of the frame and then pushing into the Velcro® at the sides.












©2017 Sunrise Medical (US) LLC




LEG EXTEND (White)

LEG LIFT (Yellow)

LIFT (Green)

TILT (Purple)


RECLINE (Red)



LEGS (3)

TILT or TILT/LIFT (2)

RECLINE (1)



Sunrise Medical 2842 Business Park Ave Fresno, CA 93727 USA In Canada (800) 263-3390

