J2 Deep Contour Cushion

User Instruction Manual & Warranty

J2 Dyb-kontur pude

Brugsanvisning og Garanti

J2 Djupkonturdyna

Bruksanvisning & garanti

J2 Deep Contour Cushion

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The J2 Deep Contour Cushion is the newest member of the J2 cushion family. The J2 Deep Contour Cushion was developed for clinicians to use with clients who are at extremely high risk for skin breakdown. Designed for superior skin protection and stability, the J2 Deep Contour Cushion is the maintenance free alternative you've been waiting for.

The J2 Deep Contour Cushion Components

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



J2 Deep Contour underlag Velcro-striper #3

Caution

Prior to prolonged sitting, any cushion should be tried for a few hours at a time while a clinician inspects your skin to ensure that red pressure spots are not developing. You should regularly check for skin redness. The clinical indicator for tissue breakdown is skin redness. If your skin develops redness, discontinue the use of the cushion immediately and see your doctor or therapist.

INSTRUCTIONS FOR J2 DEEP CONTOUR CUSHION

Caution

Installing a cushion on a wheelchair may affect the center of gravity of the wheelchair and may cause the wheelchair to tip backwards, potentially resulting in injury. Always assess for the potential need for anti-tippers or amputee axle adapter brackets to be added to the wheelchair to help increase stability.

The J2 Deep Contour Cushion is designed to help reduce the incidence of pressure sores. However, no cushion can completely eliminate sitting pressure or prevent pressure sores. The J2 Deep Contour Cushion is not a substitute for good skin care including; proper diet, cleanliness, and regular pressure reliefs.

In addition, prior to prolonged sitting, any cushion should be tried for a few hours at a time while a clinician inspects your skin to ensure that red pressure spots are not developing. You should regularly check for skin redness. The clinical indicator for tissue breakdown is skin redness. If your skin develops redness, discontinue the use of the cushion immediately and see your doctor or therapist.

Do not leave the cushion outside overnight at temperatures below 40° F (5° C). Allow the cushion to warm to room temperature before using. Sitting on extremely cold or hot surfaces can cause skin damage.

Avoid sharp objects or exposure to excessive heat or open flame.

Fitting The J2 Deep Contour Cushion

Select the cushion to match your hip width when in the seated position. The cushion length should extend to within 1" - 2" (2.5cm - 5cm) from the back of the knees. Place the cushion to the back of the wheelchair with the material identification tags and fluid pad toward the rear. When fit correctly, your hips should be to the back of the chair. Your ischials (seat bones) should be centered on the

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Jay Flow™ Fluid pad in the seat well of the cushion. Proper footrest adjustment may enhance sitting comfort and help lower peak sitting pressures. Adjust the footrests so your legs rest firmly but comfortably on the cushion.

Note: Do not use the wheelchair's seat width as a guide in determining the appropriate cushion width. Improper cushion fit may be detrimental to the performance of a pressure-distributing cushion and its support of the body.

Avoid using the cushion on overly stretched seat upholstery. Stretched upholstery does not provide adequate support for the cushion and improper positioning may result. Replace stretched upholstery or purchase the J2 Deep Contour Solid Seat Insert or order the J2 Deep Contour Cushion with the Solid Drop Seat. The Solid Drop Seat **Option** is attached to the cushion for a one piece design that replaces the wheelchair's seat upholstery.

Avoid using overly stretched back rest upholstery. Stretched back upholstery does not provide adequate support for the body and improper positioning may result. Replace stretched back rest upholstery or purchase a back support product from Sunrise that is appropriate for you. Contact an authorized supplier for a free demonstration of these products.

Checking for bottoming out on the J2 Cushion

Bottoming out occurs on a J2 Cushion when you displace the fluid out from underneath your seat bones, which leaves you sitting on the foam base. It sometimes occurs on very thin individuals, people using recliner wheelchairs, people who slouch when sitting, or if using a cushion that is **too wide**.

When bottoming out occurs, increased pressure is placed onto the ischials and coccyx increasing the risk for skin breakdown.

To check for bottoming out, sit on the cushion **without the cover** for a minimum of two minutes. Transfer up and off the cushion (or have someone help you transfer), trying not to disturb the fluid underneath you. Push down in the depressions on the pad where your ischials (seat bones) and coccyx (tailbone) were. You should have to push through at least 1/2" (1cm) of fluid before you feel the firm cushion base below.

If the cushion is properly positioned and the footrests are properly adjusted, and there is not at

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

least the minimum 1/2" (1cm) of fluid, the cushion is bottoming out and should not be used. If you are bottoming out, discontinue use of the cushion and see your clinician. Usually bottoming out is easily solved by using fluid supplement pads (part# F119). Call your local authorized supplier to see if this is appropriate for you.


Note: We recommend using a J2 Deep Contour Solid Seat Insert or order the J2 Deep Contour Cushion with the Solid Drop Seat. The Solid Drop Seat Option is attached to the cushion for a one piece design that replaces the wheelchair's seat upholstery. It will optimize the performance of your J2 Deep Contour Cushion and help provide better positioning, and increase the life of the cushion.

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

ASSEMBLING AND ADJUSTING THE SOLID DROP SEAT

Located on the bottom of the Solid Drop Seat are 10 pairs of holes, five pairs each side of the platform. Each hook has two slots that line up with one pair of holes. Locate where the hooks should best be placed on the platform so that the installation doesn't interfere with hardware on the wheelchair seat frame. There are four hooks. Each side requires installing two. Align the slots in the hooks with one pair of holes. Slide a washer onto a bolt and fasten the bolt through the slot in the hook into the Solid Drop Seat. Each hook requires two bolts. Do the same for all four hooks. Leave all hooks slightly loose. Place the solid seat on the seat rail frame and adjust the width. Tighten as many hook bolts as possible prior to removing Solid Drop Seat from chair (this will preserve their alignment). Remove the Solid Drop Seat and finish tightening the bolts to secure the hooks.

Strong platform promotes level sitting surface

Five hook attachment locations accommodate most wheelchair frames

Each hook adjusts 1/2" in width for slightly irregular wheelchair widths

Cushion is notched to drop up to 1 1/2"

1/2", 1" or 1 1/2" drop hooks may be used to change seat to floor height

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Secure hook retainers and stopper

You should now be ready to secure the restraining hardware onto your wheelchair rails.

Snap the two hook retainers toward the back of the wheelchair seat rail.

The rear hooks will slide under these as shown below. These will prevent your cushion from accidentally popping out of the wheelchair. Snap the hook stopper towards the front of the wheelchair seat rails. The front hooks rest behind the stopper to prevent your cushion from sliding forward.

Retainers and stoppers should be secured tightly to wheelchair seat rails with the nuts and bolts provided. Point the heads of the bolts to the outside of the wheelchair.



Note: Enclosed are three lengths of 8-32 bolts. Using bolt longer than required can interfere with installation and removal of seat. Use:

- the 8-32 x 3/4" when the stopper and/or retainers are mounted on 7/8" tubing.
- the 8-32 x 1" when the stopper and/or retainers are mounted on 1" diameter tubing.
- the 8-32 x 1-1/4" when the stopper and/or retainers are mounted over "rivnut" threaded insert on 1" diameter tubing.

To install cushion, simply slide rear hooks under the rear hooks retainers and drop front hook behind hook stopper.

Note: If hook retainers or stopper cannot be installed due to interference with hardware on your wheelchair, you might need to select an alternate location for your brackets and hooks.

Use hook stopper 1/4" forward of one front hook

Use hook retainers behind rear hooks

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Length adjustment

For a longer seat length, position Solid Drop Seat 1" (2 cm) forward of the rear vertical posts of your wheelchair. If front hooks are less than 1" (2 cm) from front end of seat rails, then remove front hooks and move them to the alternate location directly behind (see picture).

Fitting Hints:

- Once you have determined the position of your Solid Drop Seat on your wheelchair, you may find that the hooks interfere with the hardware that is particular to your wheelchair. If this is the case, reposition those particular hardware hooks to an alternate location.

Easy Maintenance and Cleaning

Remove the cover and turn it inside out. Machine wash on warm cycle and tumble dry on low heat. Do not use industrial dryers. Do not steam autoclave. Wipe to clean the fluid pad with warm water and towel dry. Do not submerge in water. Water may be absorbed into the foam through the air vents on the top of the pad. Wipe to clean the foam base and towel dry.

Note: When the cover is removed from the cushion for cleaning, inspect the fluid pad and foam base for unusual wear. Feel the fluid for any inconsistencies. If inconsistencies are found, knead the fluid pad to original consistency and/or call Jay customer service at Sunrise Medical at 1-800-648-8282.

Reassembling the J2 Deep Contour Cushion after cleaning

Follow cleaning and maintenance instructions as described in "Easy maintenance and cleaning", then reassemble as follows after cleaning. Attach the fluid pad to the foam base. Attach the rear three Velcro® strips first, then the front two Velcro strips. Reinstall the cover with the material I.D. tags to the back. If necessary, reinstall the four hooks on the solid drop seat.

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Märk: When the cover is removed from the cushion for cleaning, inspect the fluid pad and foam base for unusual wear. Feel the fluid for any inconsistencies. If inconsistencies are found, knead the fluid pad to original consistency and/or call Jay customer service at Sunrise Medical at 1-800-648-8282.

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J2 Deep Contour Cushion Warranty

Each Jay fluid cushion is carefully inspected and tested to provide peak performance. Every Jay fluid cushion is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to the cushion's removable cover. Claims and repairs should be processed through the nearest authorized supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded. There are not warranties which extend beyond the description on the face hereof. Remedies for breach of express warranties herein are limited to repair or replacement of the goods. In no event shall damages for breach of any warranty include any consequential damages or exceed the cost of non-conforming goods sold.

This product is manufactured to comply with the 'Medical Device Directive' 93/42/EEC.

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This product is manufactured to comply with the 'Medical Device Directive' 93/42/EEC.

ENGLAND

Sunrise Medical Ltd.
High Street
Wollaston
West Midlands D48 4PS
England
44-1384-44-66-84

DENMARK

G2 Wheelchairs and Gear
Jegstrupvej 56
DK 8361 Hasselager Denmark
DK 8361

SWEDEN

Sunrise Medical AB
Britta Sahlgrens gata 8A
SE-421 31 VÄSTRA FRÖLUNDA

NORWAY

Norsk Rehab
1450 Nesoddtangen
Oslo, Norway 1450

Sunrise Home Healthcare Group**Mobility Products Division**

7477 East Dry Creek Parkway • Longmont, Colorado 80503 USA

(303) 218-4500 or (800) 648-8282

In Canada (800) 263-3390

Manufactured under one or more of the following patents plus other patents pending: U.S. Patent # 4,588,229; 4,660,238; 4,726,624; 4,761,843; 4,842,330; 5,018,790; 5,255,404; 5,352,023; 5,362,543; 5,369,829; 5,390,384; 5,395,162; 5,457,833; 5,490,299; 5,397,517, 5,687,436; and other foreign patents.

© 1998, Sunrise Medical 11.98 Rev. A
XT2405S