


**!** Please read these instructions carefully before beginning the installation. Failure to understand and follow installation instructions may result in injury to installer and/or end user and may void the warranty. If you have any questions call Sunrise Medical Technical support at 800.333.4000

Available Pad sizes:




**NOTE:** Before mounting the laterals be sure that the correct mounting slots (A) are available. Do not tighten any hardware on this installation until you have completed all of the adjustments needed. Always consult a professional when making adjustments, and adding components to any seating system. Additionally, note that the lateral can be installed on either the inside (Figure 2) or outside (Figure 2A) of the backshell.

**I. Installation should be:** **2 2A**

- Remove the lateral pad and cover (C).
- Determine whether the swing-away lateral should be mounted internally or externally.
- Insert the nut plates (B) into the lateral slots (A). Position the JAY lateral (F) and hand-thread the smaller button head screws (D) and washers into the nutplate on the backshell and lateral pad. (Figure 2 and 2A)
- Once lateral pads and hardware have been mounted loosely on both sides, re-attach the covers to the back and lateral pad for the adjustment step.

**Tools needed for installation**

- 1 x 3mm hex key (supplied)
- 1 x 5mm hex key (supplied)



**2. Adjusting the lateral pads should be:** **2 2A 3**

- With user seated, and in a recommended posture, use the 3mm hex key to loosen the screws (D). Slide the nutplate and adjustable bracket up or down to adjust the height of the laterals. Snug the screws to hold the adjustment.
- Adjust the length of the laterals by sliding the lateral pad (C) and nutplate horizontally within the bracket slot (E) snug screws (D).
- The JAY swing-away lateral can further be adjusted by using the 5mm hex key and angle adjusting at the 6mm screws (G, Figure 3)
- Once height, length and angle adjustments are completed, tighten the screws to 20 in-lbs.

**3. Checkout**

Sunrise recommends that all components, and fasteners be checked monthly for wear or breakage. **Tighten loose fasteners to 20 in-lbs.** Broken components must be replaced immediately. If you discover a broken component, discontinue use until a suitable part can be found by an authorized Sunrise Medical Supplier

