



Evolution™ Wave™

Adjustable Air-Foam Floatation skin protection cushions with a solid positioning base —available in three styles.

SEATING SYSTEMS



At VARILITE, innovation means superior solutions like our new Wave Series, a family of contoured bases for positioning support. Therapists can now combine Evolution and Evolution PSV™ adjustable Air-Foam Floatation skin protection cushions with a solid positioning base. Three different Wave designs are available: CPW (Contoured Positioning Wedge), LPB (Lateral Positioning Base), or CPB (Contoured Positioning Base).



Evolution™ Wave™

- CPW (Contoured Positioning Wedge)**
 - lateral thigh support to prevent thigh abduction
 - medial thigh support to prevent thigh adduction
- LPB (Lateral Positioning Base)**
 - lateral pelvic support to maintain pelvis in center of seat
 - lateral thigh support to prevent thigh abduction
 - solid base of support for postural stability
- CPB (Contoured Positioning Base)**
 - lateral pelvic support to maintain pelvis in center of seat
 - lateral thigh support to prevent thigh abduction
 - medial thigh support to prevent thigh adduction
 - solid base of support for postural stability

VARILITE Wave Positioning Support Comparison Chart

	Lateral Pelvic Support	Posterior Pelvic Support	Lateral Thigh Support	Medial Thigh Support	Firm Base of Support
CPW (Contoured Positioning Wedge)			•	•	
LPB (Lateral Positioning Base)	•	•	•		•
CPB (Contoured Positioning Base)	•	•	•	•	•

Each Wave positioning support fits securely beneath the Evolution cushion. Both cushion and Wave positioning support fit inside the tailored Evolution cover.

Evolution Wave Series are available with the following options:

- Cover: Mesh or Incontinence
- Valve: Standard or PSV (Pressure Setting Valve)

See price list for size availability and HCPCS coding.



4000 1st Avenue South
Seattle, WA 98134
Phone: (800) 827-4548
Toll Free Fax: (800) 363-7739
Fax: (206) 343-5795
www.varilite.com

Authorized EU Representative:
Cascade Designs, Limited
Dwyer Road, Midleton, Co. Cork
Republic of Ireland

© 2005 Cascade Designs, Inc.



Evolution™ Wave™ Series

Each Wave Design Enhances the Multi-Stiffness Foam of the Evolution



Evolution Wave CPW (Contoured Positioning Wedge)

Evolution Wave CPW provides exceptional positioning of the lower extremities for enhanced sitting posture.




1. Wave Contoured Positioning Wedge

Lateral thigh supports position the lower extremities and prevent unwanted abduction of the thighs. Medial thigh support provides the correct amount of abduction to prevent unwanted adduction of the thighs. The Wave CPW increases the weight-bearing load on the back of the thighs and decreases the load on the buttocks, without changing the seat-to-back angle. The closed-cell foam CPW will not absorb fluid and is resistant to bacteria.

2. Air-Foam Floatation Cushion

3. Cover

Evolution Wave CPW



Evolution Wave LPB (Lateral Positioning Base)

Evolution Wave LPB provides a solid base of support and a high degree of pelvic and lower extremity lateral positioning.




1. Wave Lateral Positioning Base

Lateral pelvic supports, located posterior to the trochanters, provide lateral stability that maintains the pelvis in the center of the seat. Lateral thigh supports position the lower extremities and prevent unwanted abduction of the thighs. The Wave LPB provides a firm foundation for postural stability and discourages chair upholstery from slinging. The closed-cell foam LPB will not absorb fluid and is resistant to bacteria.

2. Air-Foam Floatation Cushion

3. Cover

Evolution Wave LPB



Evolution Wave CPB (Contoured Positioning Base)

Evolution Wave CPB provides a solid base of support and a high degree of pelvic and lower extremity lateral and medial positioning.




1. Wave Contoured Positioning Base

Lateral pelvic supports, located posterior to the trochanters, provide lateral stability that maintains the pelvis in the center of the seat. Lateral thigh supports position the lower extremities and prevent unwanted abduction of the thighs. Medial thigh support provides the correct amount of abduction to prevent unwanted adduction of the thighs. The Wave CPB provides a firm foundation for postural stability and discourages chair upholstery from slinging. The closed-cell foam CPB will not absorb fluid and is resistant to bacteria.

2. Air-Foam Floatation Cushion

3. Cover

Evolution Wave CPB



*Cover is shown with cut-away for illustration purposes only.