



# Breezy 600

## EIB60, July 2014

ORDERING INFORMATION PHONE: (800) 333-4000, FAX: (800) 300-7502, WEBSITE: www.sunrisemedical.com

### ACCOUNT INFORMATION

Account No: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

Drop Ship/Ship To: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_  
 Zip Code: \_\_\_\_\_ Phone Number: \_\_\_\_\_

#### BREEZY 600

Hangers not included in base price

086WM20 K0004 Breezy 600 \$1,095.00

#### Frame Color

086C6 Forest Green Std  
 086C16 Silver Texture NC  
 086C1 Blue NC  
 086C2 Black NC

#### SEATING

##### Seat to Floor Height Range

Standard (15" - 20 1/2", see matrix) Std  
 Ultra Hemi (13" - 18", see matrix) NC

##### Seat Width

086W5 14" Std  
 086W7 16" Std  
 086W9 18" Std  
 086W42 E2201 20" \$140.00

##### Seat Depth

086FD7 16" Std  
 086FD9 18" NC

##### Seat Height

086HT1 13" \$95.00  
 086HT2 14" \$95.00  
 086HT3 14 1/2" \$95.00  
 086HT4 15" \$95.00  
 086HT5 15 1/2" \$95.00  
 086HT6 16" \$95.00  
 086HT7 16 1/2" \$95.00  
 086HT8 17" NC  
 086HT9 17.5" NC  
 086HT10 18" NC  
 086HT11 18.5" NC  
 086HT12 19" <sup>12</sup> NC  
 086HT13 19.5" <sup>12</sup> NC  
 086HT14 20" <sup>12</sup> NC  
 086HT16 20.5" <sup>12</sup> NC

##### Seat Sling

086SL3 Padded Seat Sling \$15.00  
 086SL4 Non-padded Seat Sling <sup>2</sup> Std  
 086SL6 Omit Seat Sling -\$15.00

##### Cushion

To order any JAY cushions on-chair, please see JAY price list.

086JC01 Jay Basic Price List  
 086JC02 Jay Easy Price List  
 086JC03 Jay Soft Combi P Price List  
 086JC04 Jay J2 Price List  
 086JC05 Jay J2 Deep Contour Price List  
 086JC06 Jay Lite Price List  
 086JC07 Jay J3 Price List  
 086CU1 2" Cushion \$120.00

#### BACKREST (16"-19", 8° w/ Push Handles)

Jay J3 Back <sup>15</sup>  
 086JB01 Please see the J3 Back order form for part #'s.  
 086B12 Adjustable w/tall back canes <sup>1</sup> \$15.00  
 086BU14 Omit Back Upholstery -\$15.00

#### ARMRESTS

##### Flip Back Armrests

Fixed Height, Desk Length Std  
 086AA5 Fixed Height, Full Length NC  
 086AA6 E0973 Flip Back Height-Adjustable, Desk \$295.00  
 E0973 Flip Back Height-Adjustable, Full \$295.00

#### FOOTREST

##### Hanger Type

086H4 60° Swing-Away Hemi <sup>1,3</sup> \$140.00  
 086H131 60° Swing-Away <sup>4</sup> \$140.00  
 086H2 70° Swing-Away <sup>4,5</sup> \$140.00  
 086H24 K0053 Articulating Legrests <sup>4,6,7,8,9</sup> \$350.00  
 086H14 E0990 Elevating Legrests <sup>4,10</sup> \$295.00

##### Footrest Extension Tubes

086E1 Short (13"-15 1/2") <sup>4,11</sup> NC  
 086E2 Medium (15 1/2" - 18") <sup>4,6</sup> Std  
 086E44 Long (18" - 20 1/2") <sup>4,6,12</sup> \$35.00  
 086E21 Omit Ext. Tubes & Footplates NC

##### Footplates

086F1 Composite w/Heel Loops <sup>13</sup> Std  
 086F5 Aluminum <sup>14</sup> NC  
 086F3 E0951 Adult Angle-Adj. w/Heel Loops <sup>15</sup> \$185.00  
 K0040

#### CASTERS

5" Polyurethane NC  
 8" Polyurethane <sup>16</sup> Std  
 8" Pneumatic <sup>6</sup> \$85.00

##### Footnotes:

- <sup>1</sup> Not available when using Standard frame.
- <sup>2</sup> For use with cushion.
- <sup>3</sup> Only available with standard hemi extension tube.
- <sup>4</sup> Not available with Ultra Hemi frame.
- <sup>5</sup> Not available with 8" pneumatic casters.
- <sup>6</sup> Not available with 20" rear wheels when using Standard frame.
- <sup>7</sup> Only available with 8" casters.
- <sup>8</sup> Not available with 14" seat width.
- <sup>9</sup> Pull-to-lock wheel locks required.
- <sup>10</sup> Not available with composite footplates.
- <sup>11</sup> Not available with articulating legrest.
- <sup>12</sup> Must use 1-1/2" stem.
- <sup>13</sup> Not available with ELR.
- <sup>14</sup> Only available with ELR.
- <sup>15</sup> To order a J3 Back please select Omit Back Upholstery.  
Please note: The J3 Back requires back canes or push handles.
- <sup>16</sup> Not available with 70° hangers due to interference

PROPER USE OF HCPCS CODES DOES NOT ENSURE COVERAGE OR PAYMENT. FOR COVERAGE INFORMATION, VERIFY THE APPROPRIATE PAYER'S COVERAGE POLICY.

**Caster Options**

<input type="checkbox"/>	1-1/2" Fork Stem (Double Bearing)	\$35.00
<input type="checkbox"/>	086CP13 Double Bearing Caster Plate	\$35.00

**AXLES AND WHEELS**

**Axle**

<input type="checkbox"/>	086AL4 Threaded Offset Axles	Std
<input type="checkbox"/>	086AX43 Quick-Release Rear Axles	\$95.00
<input type="checkbox"/>	86OA2 E0958 One-Arm Drive <sup>7,8</sup>	\$850.00
<input type="checkbox"/>	086HM2 Left Mount	
<input type="checkbox"/>	086HM1 Right Mount	

**Anti-Tip Tubes**

<input type="checkbox"/>	086AT1 E0971 Rear Anti-Tip Tubes	\$140.00
--------------------------	----------------------------------	----------

**Rear Wheel Size**

<input type="checkbox"/>	086WS4 20"	NC
<input type="checkbox"/>	086WS5 22"	NC
<input type="checkbox"/>	086WS6 24"	Std

**Rear Wheel Tires**

<input type="checkbox"/>	086RT4 Low-profile Polyurethane <sup>1,5</sup>	Std
<input type="checkbox"/>	086RT29 Full-profile Polyurethane <sup>2,5</sup>	\$50.00
<input type="checkbox"/>	086RT35 E2211 E2212 Pneumatic <sup>5</sup>	\$165.00
<input type="checkbox"/>	086RT33 E2213 Pneumatic w/Airless Insert <sup>5,6</sup>	\$215.00

**Hand Rims**

<input type="checkbox"/>	086HR5 Composite <sup>5</sup>	Std
<input type="checkbox"/>	086HR2 Aluminum - painted <sup>5</sup>	NC
<input type="checkbox"/>	086HR1 Aluminum - anodized <sup>5</sup>	\$0.00
<input type="checkbox"/>	086HR4 Plastic Coated <sup>5</sup>	\$95.00
<input type="checkbox"/>	086HR7 Omit Hand Rims <sup>5</sup>	NC

**Projections <sup>3,5</sup>**

<input type="checkbox"/>	086HR8 Oblique - 8	\$325.00
--------------------------	--------------------	----------

**Wheel Locks**

<input type="checkbox"/>	086WL1 High-Push Breezy-Style	Std
<input type="checkbox"/>	086WL2 High-Pull Breezy-Style	NC
<input type="checkbox"/>	086WL22 High-Push Quickie-Style	\$40.00
<input type="checkbox"/>	086WL21 High-Pull Quickie-Style	\$40.00

**Wheel Lock Options**

<input type="checkbox"/>	086EH1 E0961 6" Ext. Handles with Bungee	\$75.00
<input type="checkbox"/>	086GA1 E0974 Grade Aids <sup>4,9</sup>	\$150.00

**ACCESSORIES**

**Backpack & Seat Pouch**

<input type="checkbox"/>	086A1 Utility Back Pack (Black/Charcoal)	\$55.00
<input type="checkbox"/>	086A8 Seat Pouch (Black/Charcoal)	\$45.00

**Pelvic Positioning Belt**

<input type="checkbox"/>	086A12 E0978 Aircraft Buckle (2" wide)	\$110.00
<input type="checkbox"/>	086A142 E0978 1.5" Auto Buckle	\$70.00

**Wheelchair Tray Table - Tuffak**

<input type="checkbox"/>	086A31 E0950 Medium (15"-17")	\$250.00
<input type="checkbox"/>	086A32 E0950 Large (18"-20")	\$250.00

**Miscellaneous**

<input type="checkbox"/>	086A82 Caddy	\$45.00
<input type="checkbox"/>	086A51 E2207 Crutch Holder	\$195.00
<input type="checkbox"/>	086A83 Tool Kit	\$45.00

See Sunrise Seating Solutions on reverse side.

TOTAL CHAIR PRICE: \$ \_\_\_\_\_

**FOOTNOTES:**

- <sup>1</sup> Not available with Projection Hand Rims.
- <sup>2</sup> "No Charge" option when used with Projection Hand Rims.
- <sup>3</sup> Not available with Low-profile polyurethane tires. Full-profile Polyurethane tires available as "No Charge" option. Requires 24" rear wheels.
- <sup>4</sup> Not available with Polyurethane tires.
- <sup>5</sup> Not available with 24" Mag SLT
- <sup>6</sup> Must order Pneumatic tire with Flat Free Insert
- <sup>7</sup> Not available with Ultra Hemi frame.
- <sup>8</sup> Not available with Quick-Release Rear Axles.
- <sup>9</sup> Grade Aids require Quickie style wheel locks



The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Standard** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
15.0"	5" lp	N/A	N/A	N/A	N/A
16.0"	5" lp	5" lp	N/A	N/A	N/A
16.5"	5" lp	5" lp	N/A	N/A	N/A
17.0"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
17.5"	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp
18.0"	N/A	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem
18.5"	N/A	8" lp *20" full poly not recommended	8" lp	8" lp	8" lp
19.0"	N/A	N/A	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
19.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem
20.0"	N/A	N/A	N/A	N/A	8" lp & 1-1/2" stem
20.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem

The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Ultra Hemi** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
13.0"	5" lp	N/A	N/A	N/A	N/A
14.0"	5" lp	5" lp	N/A	N/A	N/A
14.5"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
15.0"	8" lp	8" lp	8" lp	8" lp	N/A
15.5"	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem
16.0"	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu
16.5"	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem
17.0"	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
17.5"	N/A	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu
18.0"	N/A	N/A	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem