



Breezy 600

EIB60, July 2014

ORDERING INFORMATION PHONE: (800) 333-4000, FAX: (800) 300-7502, WEBSITE: www.sunrisemedical.com

ACCOUNT INFORMATION

Account No: _____
 Date: _____
 P.O. #: _____
 Buyer: _____
 Marked For: _____

Drop Ship/Ship To: _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____ Phone Number: _____

BREEZY 600

Hangers not included in base price

086WM20 K0004 Breezy 600 \$1,095.00

Frame Color

086C6 Forest Green Std
 086C16 Silver Texture NC
 086C1 Blue NC
 086C2 Black NC

SEATING

Seat to Floor Height Range

Standard (15" - 20 1/2", see matrix) Std
 Ultra Hemi (13" - 18", see matrix) NC

Seat Width

086W5 14" Std
 086W7 16" Std
 086W9 18" Std
 086W42 E2201 20" \$140.00

Seat Depth

086FD7 16" Std
 086FD9 18" NC

Seat Height

086HT1 13" \$95.00
 086HT2 14" \$95.00
 086HT3 14 1/2" \$95.00
 086HT4 15" \$95.00
 086HT5 15 1/2" \$95.00
 086HT6 16" \$95.00
 086HT7 16 1/2" \$95.00
 086HT8 17" NC
 086HT9 17.5" NC
 086HT10 18" NC
 086HT11 18.5" NC
 086HT12 19" ¹² NC
 086HT13 19.5" ¹² NC
 086HT14 20" ¹² NC
 086HT16 20.5" ¹² NC

Seat Sling

086SL3 Padded Seat Sling \$15.00
 086SL4 Non-padded Seat Sling ² Std
 086SL6 Omit Seat Sling -\$15.00

Cushion

To order any JAY cushions on-chair, please see JAY price list.

086JC01 Jay Basic Price List
 086JC02 Jay Easy Price List
 086JC03 Jay Soft Combi P Price List
 086JC04 Jay J2 Price List
 086JC05 Jay J2 Deep Contour Price List
 086JC06 Jay Lite Price List
 086JC07 Jay J3 Price List
 086CU1 2" Cushion \$120.00

BACKREST (16"-19", 8° w/ Push Handles)

Jay J3 Back ¹⁵
 086JB01 Please see the J3 Back order form for part #'s.
 086B12 Adjustable w/tall back canes ¹ \$15.00
 086BU14 Omit Back Upholstery -\$15.00

ARMRESTS

Flip Back Armrests

Fixed Height, Desk Length Std
 086AA5 Fixed Height, Full Length NC
 086AA6 E0973 Flip Back Height-Adjustable, Desk \$295.00
 E0973 Flip Back Height-Adjustable, Full \$295.00

FOOTREST

Hanger Type

086H4 60° Swing-Away Hemi ^{1,3} \$140.00
 086H131 60° Swing-Away ⁴ \$140.00
 086H2 70° Swing-Away ^{4,5} \$140.00
 086H24 K0053 Articulating Legrests ^{4,6,7,8,9} \$350.00
 086H14 E0990 Elevating Legrests ^{4,10} \$295.00

Footrest Extension Tubes

086E1 Short (13"-15 1/2") ^{4,11} NC
 086E2 Medium (15 1/2" - 18") ^{4,6} Std
 086E44 Long (18" - 20 1/2") ^{4,6,12} \$35.00
 086E21 Omit Ext. Tubes & Footplates NC

Footplates

086F1 Composite w/Heel Loops ¹³ Std
 086F5 Aluminum ¹⁴ NC
 086F3 E0951 Adult Angle-Adj. w/Heel Loops ¹⁵ \$185.00
 K0040

CASTERS

5" Polyurethane NC
 8" Polyurethane ¹⁶ Std
 8" Pneumatic ⁶ \$85.00

Footnotes:

- ¹ Not available when using Standard frame.
- ² For use with cushion.
- ³ Only available with standard hemi extension tube.
- ⁴ Not available with Ultra Hemi frame.
- ⁵ Not available with 8" pneumatic casters.
- ⁶ Not available with 20" rear wheels when using Standard frame.
- ⁷ Only available with 8" casters.
- ⁸ Not available with 14" seat width.
- ⁹ Pull-to-lock wheel locks required.
- ¹⁰ Not available with composite footplates.
- ¹¹ Not available with articulating legrest.
- ¹² Must use 1-1/2" stem.
- ¹³ Not available with ELR.
- ¹⁴ Only available with ELR.
- ¹⁵ To order a J3 Back please select Omit Back Upholstery.
Please note: The J3 Back requires back canes or push handles.
- ¹⁶ Not available with 70° hangers due to interference

PROPER USE OF HCPCS CODES DOES NOT ENSURE COVERAGE OR PAYMENT. FOR COVERAGE INFORMATION, VERIFY THE APPROPRIATE PAYER'S COVERAGE POLICY.

Caster Options

<input type="checkbox"/>	1-1/2" Fork Stem (Double Bearing)	\$35.00
<input type="checkbox"/>	086CP13 Double Bearing Caster Plate	\$35.00

AXLES AND WHEELS

Axle

<input type="checkbox"/>	086AL4 Threaded Offset Axles	Std
<input type="checkbox"/>	086AX43 Quick-Release Rear Axles	\$95.00
<input type="checkbox"/>	86OA2 E0958 One-Arm Drive ^{7,8}	\$850.00
<input type="checkbox"/>	086HM2 Left Mount	
<input type="checkbox"/>	086HM1 Right Mount	

Anti-Tip Tubes

<input type="checkbox"/>	086AT1 E0971 Rear Anti-Tip Tubes	\$140.00
--------------------------	----------------------------------	----------

Rear Wheel Size

<input type="checkbox"/>	086WS4 20"	NC
<input type="checkbox"/>	086WS5 22"	NC
<input type="checkbox"/>	086WS6 24"	Std

Rear Wheel Tires

<input type="checkbox"/>	086RT4 Low-profile Polyurethane ^{1,5}	Std
<input type="checkbox"/>	086RT29 Full-profile Polyurethane ^{2,5}	\$50.00
<input type="checkbox"/>	086RT35 E2211 E2212 Pneumatic ⁵	\$165.00
<input type="checkbox"/>	086RT33 E2213 Pneumatic w/Airless Insert ^{5,6}	\$215.00

Hand Rims

<input type="checkbox"/>	086HR5 Composite ⁵	Std
<input type="checkbox"/>	086HR2 Aluminum - painted ⁵	NC
<input type="checkbox"/>	086HR1 Aluminum - anodized ⁵	\$0.00
<input type="checkbox"/>	086HR4 Plastic Coated ⁵	\$95.00
<input type="checkbox"/>	086HR7 Omit Hand Rims ⁵	NC

Projections ^{3,5}

<input type="checkbox"/>	086HR8 Oblique - 8	\$325.00
--------------------------	--------------------	----------

Wheel Locks

<input type="checkbox"/>	086WL1 High-Push Breezy-Style	Std
<input type="checkbox"/>	086WL2 High-Pull Breezy-Style	NC
<input type="checkbox"/>	086WL22 High-Push Quickie-Style	\$40.00
<input type="checkbox"/>	086WL21 High-Pull Quickie-Style	\$40.00

Wheel Lock Options

<input type="checkbox"/>	086EH1 E0961 6" Ext. Handles with Bungee	\$75.00
<input type="checkbox"/>	086GA1 E0974 Grade Aids ^{4,9}	\$150.00

ACCESSORIES

Backpack & Seat Pouch

<input type="checkbox"/>	086A1 Utility Back Pack (Black/Charcoal)	\$55.00
<input type="checkbox"/>	086A8 Seat Pouch (Black/Charcoal)	\$45.00

Pelvic Positioning Belt

<input type="checkbox"/>	086A12 E0978 Aircraft Buckle (2" wide)	\$110.00
<input type="checkbox"/>	086A142 E0978 1.5" Auto Buckle	\$70.00

Wheelchair Tray Table - Tuffak

<input type="checkbox"/>	086A31 E0950 Medium (15"-17")	\$250.00
<input type="checkbox"/>	086A32 E0950 Large (18"-20")	\$250.00

Miscellaneous

<input type="checkbox"/>	086A82 Caddy	\$45.00
<input type="checkbox"/>	086A51 E2207 Crutch Holder	\$195.00
<input type="checkbox"/>	086A83 Tool Kit	\$45.00

See Sunrise Seating Solutions on reverse side.

TOTAL CHAIR PRICE: \$ _____

FOOTNOTES:

- ¹ Not available with Projection Hand Rims.
- ² "No Charge" option when used with Projection Hand Rims.
- ³ Not available with Low-profile polyurethane tires. Full-profile Polyurethane tires available as "No Charge" option. Requires 24" rear wheels.
- ⁴ Not available with Polyurethane tires.
- ⁵ Not available with 24" Mag SLT
- ⁶ Must order Pneumatic tire with Flat Free Insert
- ⁷ Not available with Ultra Hemi frame.
- ⁸ Not available with Quick-Release Rear Axles.
- ⁹ Grade Aids require Quickie style wheel locks



The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Standard** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
15.0"	5" lp	N/A	N/A	N/A	N/A
16.0"	5" lp	5" lp	N/A	N/A	N/A
16.5"	5" lp	5" lp	N/A	N/A	N/A
17.0"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
17.5"	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp
18.0"	N/A	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem
18.5"	N/A	8" lp *20" full poly not recommended	8" lp	8" lp	8" lp
19.0"	N/A	N/A	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
19.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem
20.0"	N/A	N/A	N/A	N/A	8" lp & 1-1/2" stem
20.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem

The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Ultra Hemi** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
13.0"	5" lp	N/A	N/A	N/A	N/A
14.0"	5" lp	5" lp	N/A	N/A	N/A
14.5"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
15.0"	8" lp	8" lp	8" lp	8" lp	N/A
15.5"	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem
16.0"	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu
16.5"	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem
17.0"	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
17.5"	N/A	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu
18.0"	N/A	N/A	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem