


MULTI-FUNCTION, LIGHT WEIGHT SYSTEM:

- BATH SLIDER
- ROLL-IN SHOWER CHAIR
- TRANSFER CHAIR
- BEDSIDE COMMODE
- OVER-THE-TOILET COMMODE




Shown with optional adjustable bridge.

FEATURES

- Compact 18" wide tub base transfer frame fits inside bathtub
- Locking system for added safety
- Height adjustments to fit over most toilets and tubs
- Removable 13" bridge
- Lightweight aluminum construction
- Flip-up armrests and swing-away footrests
- Latex and silicone free
- Includes commode pail, locking casters, and footrests

BENEFITS

- Reduces the need for multiple pieces of equipment
- Designed for independent use, or with assistance
- Shower door or curtain can be drawn when bridge is removed



*All dimensions are based on footprint unless noted.

Optional Bridges

- Fixed length: 7", 9" or 11"
- Adjustable length: 13½" - 17¾" or 17½" - 21¼"