



J3 Cushion

ACCOUNT #: _____ SHIP TO: _____
 DATE: _____ NAME: _____
 P.O. #: _____ ADDRESS: _____
 BUYER: _____ CITY/ STATE: _____
 MARKED FOR: _____ ZIP CODE: _____ PHONE #: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

To order a J3 Cushion, a J3 SMART Part Number must be built. Follow the steps below to build your own unique Smart Part #.

Example

J3C S 1816 B FF S P
 J3 CUSHION STEP 1 STEP 2 STEP 3 STEP 4 STEP 5 STEP 6



For additional information, please refer to the J3 Cushion Quick Reference Guide. Part # 111903.

Selected Component

J3 CUSHION BASE	Select	Part #	Charge
J3 Cushion (K0734/K0735)	<input type="checkbox"/>	J3C	\$308



STEP 1

CUSHION HEIGHT	Select	Part #	Upcharge
Standard (3.75" unloaded state)	<input type="checkbox"/>	S	-
Deep (4.25" unloaded state)	<input type="checkbox"/>	D	ADD \$30

Deep height designed for maximum immersion for a higher risk user with lower mass & tissue integrity



STEP 2

CUSHION WIDTH	Select	Part #	Upcharge
Standard (12"-21")	<input type="checkbox"/>	Select box to the right	-
Bariatric Cushion (22"-24") [K0735 for 22" wide and greater]	<input type="checkbox"/>	D	Add \$85



Width in Inches

CUSHION DEPTH	Select	Part #	Upcharge
Standard (12"-21")	<input type="checkbox"/>	Select box to the right	-
Bariatric Cushion (22"-24")	<input type="checkbox"/>	D	-



Depth in Inches

WIDTH
 (shaded boxes indicated custom measurement)

	12	14	15	16	17	18	19	20	21	22	24
12	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
14	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
15	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
17	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
18	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
19	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
21	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

STEP 3

PELVIC LOADING AREA (PLA) SIZE	Select	Part #	Upcharge
A - Small (Standard with 12" - 15" widths*)	<input type="checkbox"/>	A	-
B - Medium (Standard with 16" - 19" widths*)	<input type="checkbox"/>	B	-
C - Large (Standard with 20" - 24" widths*)	<input type="checkbox"/>	C	-

Optimal PLA sizing can help provide lower Ischial pressures, improved stability and comfort for the user.



A, B, or C

ASIS WIDTH	PLA SIZE
<9.75"	A
9.75" -11.75"	B
>11.75"	C



*Unless another size is requested

To determine the "Pelvic Loading Area" (PLA) size for Trochanter Loading, locate and measure(using a Jcaliper or ART Caliper) the distance between the two ASIS on the Pelvis. Once the measurement is taken, use the chart above to select the PLA size.

Please note the following 12" or 13" wide cushions not available in deep height.

B and C PLA sizes are not available on 12"-14" wide cushions.

C PLA sizes are not available in 15"-20" wide cushions



J3 Cushion

STEP 4

PELVIC LOADING AREA (PLA) INSERT

Select	Part #	Upcharge	
<input type="checkbox"/>	FF	-	<input type="text"/>
<input type="checkbox"/>	FV	ADD \$10	

*Field Variable Insert comes standard with two 3D fluid inserts to be used for volume adjustments or obliquity changes.

FF or FV

STEP 5

COVER OPTION

Select	Part #	Upcharge	
<input type="checkbox"/>	S	-	<input type="text"/>
<input type="checkbox"/>	I	-	
<input type="checkbox"/>	M	-	

S, I, or M

STEP 6

(Optional)

POSITIONING COMPONENTS

(K0736/K0737)

Select	Part #	Up charge	
<input type="checkbox"/>	P	ADD \$20	<input type="text"/>

P if applicable

OPTIONAL CUSTOMIZATIONS

JAY supports our customers by providing a wide range of customizations product to benefit the end-user, many of which that are no-charge.

To order a modification, simply contact our customer service team and discuss the users needs and what options may be available.

Modification	Description	Up charge
Increased Height	Adds total height to cushion.	NC
Angle Cut Posterior/Anterior Slope	Removed material from cushion to create forward or rearward angle.	NC
Angle Wedge	Adds material to create forward or rearward angle	NC
Underfills	Remove fluid to create ideal immersion	NC
Overfills	Add fluid to prevent user from bottoming out in cushion	NC
Obliquity Underfills or Overfills	Adds or removes fluid in one of the 3 fluid pad segments to handle obliquities	NC
Notch for Rails	Remove material from cushion to lower cushion between the wheelchair seat rails	NC
Notch for back cane	Add or remove material to adjust the location the cushion lines up against the	NC
Leg Length Cut	Cuts different lengths into the cushion with a custom cover.	NC
Contracture Cut	Standard angle to remove on the J3 Cushion is 20 degrees. Can increase this up to 45 degrees	NC
Leg Trough Cut (1" std)	For improved positioning, use this modification to increase the depth of the cushion to support the femurs	NC
Custom Drop Seat	Adds a solid hinge system to be able to place the seat cushion below the seat rails and effectively lower seat to floor height	ADD \$75
Embroidery	Add personalization for identification purposes. (Hint: Great for large facilities or product identification during cleaning)	ADD \$30
Change Length or Width	Customizes the size of the cushion by removing or adding material with a custom cover.	NC
Change PLA Location	For Instance: Add 2" of foam to rear, remove from front to adjust PLA forward for sacral sitters.	NC

* If there is a modification that is not listed above, contact customer service and we will do our best to accommodate your needs.

J3 SMART PART

J3C						
-----	--	--	--	--	--	--

POSITIONING COMPONENTS
(If applicable)

PRICE

NOTES