

matrx

CONTOUR BACK

The MaTRx Contour back provides comfortable back support with a limited amount of adjustability for users with basic positioning needs.

Product Highlights:

- ▶ Moisture resistant, breathable cover
- ▶ Onyx cover upgrade available
- ▶ Mildly contoured for gentle midline support
- ▶ Ergonomically contoured aluminum back pan
- ▶ Fixed hardware for secure non-removable attachment
- ▶ Height adjustable on back canes small range of depth recline adjustment
- ▶ Fits 7/8" or 1" diameter back of canes (inserts available at no extra charge for 3/4" back canes)
- ▶ Widths: 15" to 22"
- ▶ Heights: 12", 15", 17"
- ▶ 300 lb weight capacity
- ▶ 2 year Warranty



ONYX COVER UPGRADE

PRODUCT OPTIONS

PELVIC/THORACIC LATERAL SUPPORTS



- ▶ Designed especially for MaTRx Contour or High Backs providing infinite/easy adjustability up and down
- ▶ 3 pad sizes:
Small 3 1/2" x 5"
Medium 3 3/4" x 5 3/4"
Large 4 1/2" x 6 1/4"

SWINGAWAY LATERALS SUPPORTS



- ▶ Designed especially for MaTRx Contour or High Backs providing infinite/easy adjustability and obstruction-free transfers
- ▶ 3 pad sizes:
Small 3 1/2" x 5"
Medium 3 3/4" x 5 3/4"
Large 4 1/2" x 6 1/4"

EXTENDED SHEAR REDUCTION MOUNTING HARDWARE (ESR)



- ▶ Allows the Contour or High Back to be mounted directly to Motion Concepts TRx power seating

PRODUCT REFERENCE



MaTRx Contour Back 12" height MaTRx Contour Back 15" height MaTRx Contour Back 17" height

Weight capacity	300 lb	300 lb	300 lb
Back width	15" to 18"	15" to 22"	15" to 22"
Fits w/c widths	15" to 18"	15" to 22"	15" to 22"
Lateral Option	✓	✓	✓
Contour depth	3"	3"	3"

Please see Price Lists/Order Forms for ADP/Medicare Codes

CONTOUR BACK WEIGHTS - *Weights include mounting hardware*

MaTRx-PB Contour Back 12" H

SIZE	KG	LB
1512	1.0	2.25
1612	1.1	2.5
1712	1.2	2.75
1812	1.4	3.0
1912	1.5	3.25
2012	1.6	3.5
2112	1.7	3.75
2212	1.8	4.0

MaTRx-PB Contour Back 15" H

SIZE	KG	LB
1515	1.2	2.75
1615	1.4	3.0
1715	1.5	3.25
1815	1.6	3.5
1915	1.7	3.75
2015	1.8	4.0
2115	1.9	4.25
2215	2.0	4.5

MaTRx-PB Contour Back 17" H

SIZE	KG	LB
1517	1.5	3.25
1617	1.6	3.5
1717	1.7	3.75
1817	1.8	4.0
1917	1.9	4.25
2017	2.0	4.5
2117	2.2	4.75
2217	2.3	5.0

COVER WEIGHTS

MaTRx-PB Contour Back 12" H

SIZE	KG	LB
1512	0.15	0.3
1612	0.15	0.3
1712	0.15	0.3
1812	0.15	0.3
1912	0.15	0.3
2012	0.15	0.3
2112	0.15	0.3
2212	0.15	0.3

MaTRx-PB Contour Back 15" H

SIZE	KG	LB
1515	0.17	0.4
1615	0.17	0.4
1715	0.17	0.4
1815	0.17	0.4
1915	0.17	0.4
2015	0.17	0.4
2115	0.17	0.4
2215	0.17	0.4

MaTRx-PB Contour Back 17" H

SIZE	KG	LB
1517	0.19	0.42
1617	0.19	0.42
1717	0.19	0.42
1817	0.19	0.42
1917	0.19	0.42
2017	0.19	0.42
2117	0.19	0.42
2217	0.19	0.42