



# KidWalk / KidWalk II

## Dynamic Mobility System



# Order Form

## Price List

Effective 10-1-08

Patent Protected 6,832,770 B1

Client: \_\_\_\_\_ Age: \_\_\_\_\_  
 Height: \_\_\_\_\_ Weight: \_\_\_\_\_ Dx: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Facility: \_\_\_\_\_  
 Therapist: \_\_\_\_\_ Phone: \_\_\_\_\_

PO: \_\_\_\_\_ Date: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_

### KidWalk / KidWalk II Dynamic Mobility System

**Standard Features Include:** Auto centering pelvic guides; Lateral weight shift; Vertical weight relief with adjustable load compensating spring assembly; Angle adjustable mast; Anti-tippers; Wheel locks; Quick release wheels; Padded seat; Adjustable axle position; Light up rear casters; Two piece frame for easy transport.

**KidWalk**  
**\$2,795.00**  
Dynamic Mobility System  
 Seat height 9" - 22"  
 Client Weight up to 65 pounds

**Specifications:**  
 Base: 24" wide by 27" long  
 Overall length: 33" - 35"  
 Mast Height: 33.5"  
 Upholstery: Black nylon  
 Wheels: 20" Mag  
 Casters: 4" Lighted



**KidWalk II**  
**\$3,395.00**  
Dynamic Mobility System  
 Seat height 15" - 28"  
 Client Weight up to 90 pounds

**Specifications:**  
 Base: 27.5" wide by 32" long  
 Overall length: 36" - 38"  
 Mast Height: 37.25"  
 Upholstery: Black nylon  
 Wheels: 24" Mag  
 Casters: 5" Lighted

**Frame Color:**  Razzberry  Shocker Green  Candy Red  Candy Blue  Black

**KidWalk Hip Support Configuration Selection: \***

Standard Configuration:  Width 5.5" - 9"  Depth 3.5" x 6"  
 Optional Configurations:  Width 7" - 12"  Depth 5" - 7"  
 Flip Up Pad Brackets:  Straight (std)  10 Degree Angle  
 \*If no selections are marked the standard configuration will be supplied.

**KidWalk II Hip Support Configuration Selection: \***

Standard Configuration:  Width 7" - 12"  Depth 5" x 7"  
 Optional Width Configurations:  5.5" - 9"  8.5" - 15"  
 Optional Depth Configurations:  3.5" - 6" (with non-extended hip pads)  
 Flip Up Pad Brackets:  Straight (std)  10 Degree Angle

**Upper Body Support Options:**

**Lateral Chest System:** Includes lateral side supports, swingaway anterior support pads and lumbar support pad .....\$ 445.00

<b>KidWalk</b>	<b>KidWalk II</b>
<input type="checkbox"/> Std 6.5" - 10.5" W x 4" - 6.5" D	<input type="checkbox"/> Std 8.5" - 14.5" W x 5.5" - 8.5" D
<input type="checkbox"/> Optional Width Adj. 8" - 14.5"	<input type="checkbox"/> Optional Width Adj. 6.5" - 10.5"
<input type="checkbox"/> Optional Depth Adj. 5.5" - 8.5"	<input type="checkbox"/> Optional Depth Adj. 4" - 6.5"

Flip Up Chest Pad Brackets:  10 Degree Angle (std)  Straight

**Upper Body Soft Support Wrap:** For the client that needs minimal upper body support .....\$ 239.00

**Headrest Options:**

**Headpad:** 3" x 4" flat headpad .....\$ 195.00  
 **Contoured Small Headrest:** 4"H x 6"W (includes mount).....\$ 299.00  
 **Contoured Medium Headrest:** 5"H x 8"W (includes mount).\$ 345.00  
 **1/2 Square Tube Mount:** (For use with OEM headrests)....\$ 95.00  
 **i2i Flip Down Headrest Mount:** .....\$ 295.00  
 **Specific Headrest Adapter:** (Specify headrest type) ..... per quote

**Lower Extremity Positioning Options:**

**Extended Hip Pads:** (Standard on 5" - 7" depth) .....\$ 95.00  
 **Dynamic Leg Wraps:** Assist in supplying feedback to the client while encouraging movement in the lower extremities.....\$ 125.00

**Seat Options:**

**Wide Seat: 3" x 8":** Large padded seat .....\$ 175.00  
 **Extended Seats:** Longer sides assist in controlling adduction and scissoring: Specify length of sides below ..... \$ 275.00

Short 4" length  
 Medium 6" length  
 Long 8" length

**Leg Straps:** Restricts movement of the lower extremities (Requires 6" or 8" extended seat).....\$ 175.00  
 **Additional Standard Seat:** (When ordered with unit).....\$ 100.00

**Frame Options:**

**Cruiser Wheels:** Allows use as an omni directional walker for clients up to 40" tall (**Not for use on KidWalkII**) .....\$ 299.00  
 **Non-Lighted Caster Wheels:** .....\$ no charge  
 **Sliding Ankle Prompts:** Straps hold ankles to help keep the client's feet from crossing.....\$ 259.00  
 **Directional Locks:** Locks out caster swivel.....\$ 249.00

**Wheel Options:**

**Anti-Rollback Device:**.....\$ 295.00  
 **20" Spoke Wheels (KidWalk only):**.....\$ 395.00  
 **24" Spoke Wheels (KidWalk II only):**.....\$ 395.00  
 **Spoke Protectors 20":**.....\$ 175.00  
 **Spoke Protectors 24":**.....\$ 195.00

