

Innovative mono-frame chair offering maximum strength and performance at the lowest possible weight.



Nothing Beats a Quickie™

You drive our design. In development, we consider how you interact with the chair and how both can work as one to maintain maximum efficiency, performance, style and fit. Before you receive the chair, the product undergoes rigorous testing above and beyond the required testing necessary for federal guidelines. Through our focus on these core elements: design, development and testing; these chairs are designed by you and for you.



Quickie GT
with XTR suspension



Back-release lanyard



Natural Fit™ Handrim

Lightweight Aluminum Chair

- Weighs in at less than 21.5 bs.

LIFESTYLE OPTIONS

Give your chair some attitude

- Choose from 16 standard color options and four Kolorfusion™ colors
- Anodized parts
- Folding backrest with back-release lanyard

Without sacrificing comfort

- Natural Fit™ Handrims
- Folding, lockdown, angle-adjustable backrest
- GT with XTR suspension

Kolorfusion™ colors



Anodized parts



SPECIFICATIONS

HCPCS Code	K0005
Frame Style	Rigid - open
Transport Weight	12.4 lbs.
Overall Weight	21.3 lbs.
User Weight Capacity	265 lbs.
Seat Width	12" to 20"
Seat Depth	12" to 22"
Back Angle Adjustment	85° to 102°
Front Seat To Floor Height	16" to 21"
Rear Seat to Floor Height	13" to 20"
Frame Angle	75°, 80°, 85°
WARRANTY	
Warranty	Limited Lifetime