



Quickie GP/GPV September 2016

EIGP70

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base Model				
<input type="checkbox"/>	024WM23	K0005	Quickie GP-70°	\$1,895.00
<input type="checkbox"/>	024WM23	K0005	Quickie GPV-70° (NA 12"-13")	\$1,895.00

Frame Angle				
<input type="checkbox"/>	024FA6		60° (NA 18 1/4" and 19 1/4" frame length)	NC
<input type="checkbox"/>	024FA7		70°	NC
<input type="checkbox"/>	024FA1		75° (NA with rigid footrest)	NC
<input type="checkbox"/>	024FA8		80° (NA with rigid footrest)	NC

Attitude Bike (Optional)				
<input type="checkbox"/>	024IN04		Attitude Bike bundle NA w/ GPV	See Order Form

2. SIZE / SEATING

Seat Width				
<i>Seat Width 1/2" Narrower</i>				
<input type="checkbox"/>	024W3		12"	NC
<input type="checkbox"/>	024W4		13"	NC
<input type="checkbox"/>	024W5		14"	NC
<input type="checkbox"/>	024W6		15"	NC
<input type="checkbox"/>	024W7		16"	NC
<input type="checkbox"/>	024W8		17"	NC
<input type="checkbox"/>	024W9		18"	NC
<input type="checkbox"/>	024W10		19"	\$125.00
<input type="checkbox"/>	024W11	E2201	20"	\$500.00
<input type="checkbox"/>	024W12	E2201	21"	\$500.00
<input type="checkbox"/>	024W13	E2201	22"	\$500.00

Frame Length				
<input type="checkbox"/>	024FL30		18-1/4" (Avail. on 12" Sling Depth)	NC
<input type="checkbox"/>	024FL31		19-1/4" (Avail. on 13" Sling Depth)	NC
<input type="checkbox"/>	024FL32		20-1/4" (Avail. on 14", 15" Sling Depth)	NC
<input type="checkbox"/>	024FL33		21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
<input type="checkbox"/>	024FL34		22-1/4" (Avail. on 16", 17" Sling Depth)	NC
<input type="checkbox"/>	024FL35		23-1/4" (Avail. on 17", 18" Sling Depth)	NC
<input type="checkbox"/>	024FL36		24-1/4" (Avail. on 18" Sling Depth)	NC
<input type="checkbox"/>	024FL37		25-1/4" (Avail. on 19" Sling Depth)	NC
<input type="checkbox"/>	024FL38		26-1/4" (Avail. on 20" Sling Depth)	NC
<input type="checkbox"/>	024FL39		27-1/4" (Avail. on 21" Sling Depth)	NC
<input type="checkbox"/>	024FL40		28-1/4" (Avail. on 22" Sling Depth)	NC

Seat Depth				
<i>21 & 22" seat depths only available with 16" or wider chairs</i>				
<input type="checkbox"/>	024SD5		12"	NC
<input type="checkbox"/>	024SD6		13"	NC
<input type="checkbox"/>	024SD7		14"	NC
<input type="checkbox"/>	024SD8		15"	NC
<input type="checkbox"/>	024SD9		16"	NC
<input type="checkbox"/>	024SD10		17"	NC
<input type="checkbox"/>	024SD11		18"	NC
<input type="checkbox"/>	024SD72		19"	\$75.00
<input type="checkbox"/>	024SD13	E2203	20"	\$550.00
<input type="checkbox"/>	024SD14	E2203	21"	\$550.00
<input type="checkbox"/>	024SD33	E2204	22"	\$900.00

Frame Height				
<i>See Appendix 11 for height chart. Short frame NA w/Xtender</i>				
<input type="checkbox"/>	024FH2		Short (front 9.25"-rear 8.25") NA with Xtender Power Assist	\$110.00
<input type="checkbox"/>	024FH5		Standard (front 10" - rear 9")	NC
<input type="checkbox"/>	024FH7		Tall (front 11" - rear 10")	\$110.00

Cushion				
<input type="checkbox"/>	024CU9		Reinforced Seat Sling	\$30.00
<input type="checkbox"/>	024CU1	K0669	2" Cushion	\$120.00
<input type="checkbox"/>	024CU2		Omit Cushion	NC
<input type="checkbox"/>	024CU4		Omit Cushion and Seat Sling	-\$15.00

JAY Bundled Cushions - Optional				
<i>To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.</i>				
<input type="checkbox"/>	024JC20		JAY Bundled Cushion Cushion Name or Part #: Width: _____ Depth: _____	Price List

3. CASTERS

Caster Wheels			
<input type="checkbox"/>	024CW1		3" Micro NC
<input type="checkbox"/>	024CW3		4" Low Profile Polyurethane NC
<input type="checkbox"/>	024CW31		4" Micro Lighted \$70.00
<input type="checkbox"/>	024CW6		5" Soft Roll Semi-Pneumatic \$85.00
<input type="checkbox"/>	024CW5		5" Low Profile Polyurethane NC
<input type="checkbox"/>	024CW7		5" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i> \$85.00
<input type="checkbox"/>	024CW33	E2214 E2215	6" Pneumatic \$75.00
<input type="checkbox"/>	024CW34		6" Polyurethane NC
<input type="checkbox"/>	024CW35		6" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i> \$85.00
<input type="checkbox"/>	024CW21	E2214 E2215	8" Pneumatic <i>Only available on 60° frame option</i> \$75.00

Caster Wheels - Continued			
<input type="checkbox"/>	024CW32		8" Polyurethane <i>Only available on 60° frame option</i> NC
<input type="checkbox"/>	024CW36	E2214 E2215	8" x 2" Pneumatic <i>Only available on 60° frame option</i> \$145.00
<input type="checkbox"/>	024CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert <i>Only available on 60° frame option</i> \$225.00

Caster Fork			
<input type="checkbox"/>	024CS5		3/4" Longer Fork Stem Bolt NC
<input type="checkbox"/>	024CS9		1-1/2" Longer Fork Stem Bolt NC
<input type="checkbox"/>	024CS3	K0073	Caster Pin Locks \$195.00

4. REAR WHEELS

Xtender Power Assist NEW!			
<i>RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME</i>			
<input type="checkbox"/>	024XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL \$7,895.00
<input type="checkbox"/>	024XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM \$8,695.00
<input type="checkbox"/>	024XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL \$7,895.00
<input type="checkbox"/>	024XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM \$8,695.00
<i>Please select one of the following handrim options</i>			
<input type="checkbox"/>	024XH5		Plastic Coated NC
<input type="checkbox"/>	024XH6		Stainless Steel NC

Xtender - Other Options			
<input type="checkbox"/>	024XA10		Spare NiMH Power Pack - NICKEL \$1,395.00
<input type="checkbox"/>	024XA11		Spare Li-ion Power Pack - LITHIUM \$1,795.00

Rear Wheel Size			
<i>Spinergy NA w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options</i>			
<input type="checkbox"/>	024RW8		22" Mag NC
<input type="checkbox"/>	024RW9		22" Spoke \$200.00
<input type="checkbox"/>	024RW1		24" Mag NC
<input type="checkbox"/>	024RW2		24" Spoke NC
<input type="checkbox"/>	024RW11		26" Spoke <i>(NA with Xtender Power Assist)</i> \$65.00
<input type="checkbox"/>	024RW42		24" Spinergy (540) SPOX Black \$795.00
<input type="checkbox"/>	024RR25		Omit Rear Wheel and Axle -\$55.00
Intended Size for Set Up : _____			

Rear Wheel Tire			
<input type="checkbox"/>	024RT1	E2211 E2212	Pneumatic NC
<input type="checkbox"/>	024RT2	E2211 E2213	Pneumatic Airless Insert <i>NA with 26" wheels & Xtender</i> \$80.00
<input type="checkbox"/>	024RT3	E2218	Full Poly <i>N/A with 26" wheels</i> NC
<input type="checkbox"/>	024RT6	E2211 E2212	High-Pressure Clincher \$95.00

Hand Rim			
<input type="checkbox"/>	024HR1		Aluminum Anodized NC
<input type="checkbox"/>	024HR4		Plastic Coated <i>-NA on all rims</i> \$95.00
<input type="checkbox"/>	024HR7		Omit Hand Rim NC

Projection Hand Rim			
<i>Projections NA w/ plastic coated handrims. NA w/ Spoke guards</i>			
<input type="checkbox"/>	024HR8		8 Oblique (24"/26") \$325.00

Axle			
<input type="checkbox"/>	024AX1		Quick Release Steel NC
<input type="checkbox"/>	024AX11		Quad Quick Release \$125.00

Axle Plate			
<input type="checkbox"/>	024AL2		Standard Adjustable NC
<input type="checkbox"/>	024AL13		Double Slotted Adjustable \$60.00
<input type="checkbox"/>	024AL6	E0959	Amputee Adjustable <i>NA with padded swing-away armrests or stroller handles</i> \$175.00

5. WHEEL LOCKS

Wheel Lock			
<input type="checkbox"/>	024WL1		High Mount Push NC
<input type="checkbox"/>	024WL2		High Mount Pull NC
<input type="checkbox"/>	024WL5		High Mount Scissor <i>NA w/ extension handles, grade aides</i> NC
<input type="checkbox"/>	024WL12		Do Not Mount Wheel Locks NC
<input type="checkbox"/>	024WL10		Omit NC

Extension Handle			
<input type="checkbox"/>	024EH1	E0961	6" Pair \$75.00

Grade Aid			
<input type="checkbox"/>	024GA1	E0974	Grade Aid <i>NA w/ poly, high-pressure clincher tires, or Xtender</i> \$150.00

6. FOOTREST

Footplates

Heel loop or leg strap not included, See Footplate Options Section

1. GP FOOTREST OPTIONS

<input type="checkbox"/>	024F56	K0040	Rigid NA w/Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform NA w/Ext. Tubes	\$95.00
<input type="checkbox"/>	024F1	K0040	Composite Footplate	NC
<input type="checkbox"/>	024F3	K0040	Adult Angle Adj.	\$185.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$190.00

2. GPV FOOTREST OPTIONS

<input type="checkbox"/>	024F56	K0040	Rigid NA w/Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform NA w/Ext. Tubes	\$95.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$190.00

Footrest Extension Tubes

<input type="checkbox"/>	024E1		Short (14" to 16-1/2")	NC
<input type="checkbox"/>	024E2		Medium (16-1/2" to 19")	NC
<input type="checkbox"/>	024E5		Long (19" to 21-1/2")	NC
<input type="checkbox"/>	024E16	K0037	High Mount Only avail. with adult angle adjustable footplates	\$55.00

Footplate Options (Ext Tube & Direct Mount)

<input type="checkbox"/>	024HL1	E0951	Heel Loop	\$30.00
<input type="checkbox"/>	024HL5	K0038	Leg Strap	\$35.00
<input type="checkbox"/>	024HL4		Omit Heel Loop	NC
<input type="checkbox"/>	024HL6		Omit Leg Strap	NC
<input type="checkbox"/>	024F7		Omit Footplate & Ext Tube	-\$15.00
<input type="checkbox"/>	024IG2		Impact Guard Neoprene	\$85.00

7. BACKREST

Back Type

<input type="checkbox"/>	024BT3		Non-Folding Backrest (Fixed) <i>NA with push handles or adj. upholstery</i>	
<input type="checkbox"/>	024BB1		8" Back Height	NC
<input type="checkbox"/>	024BB2		9" Back Height	NC
<input type="checkbox"/>	024BB3		10" Back Height	NC
<input type="checkbox"/>	024BB4		11" Back Height	NC
<input type="checkbox"/>	024BB5		12" Back Height	NC
<input type="checkbox"/>	024BB6		13" Back Height	NC
<input type="checkbox"/>	024BB7		14" Back Height	NC
<input type="checkbox"/>	024BB8		15" Back Height	NC
<input type="checkbox"/>	024BB9		16" Back Height	NC
<input type="checkbox"/>	024BT1		Folding Backrest (Adjustable)	
<input type="checkbox"/>	024BA1		11" - 14" Back Height	\$65.00
<input type="checkbox"/>	024BA2		12" - 16" Back Height	\$65.00
<input type="checkbox"/>	024BA3		14" - 18" Back Height	\$65.00
<input type="checkbox"/>	024BA4		16" - 20" Back Height	\$65.00

Back Upholstery

<input type="checkbox"/>	024BU1		Standard	NC
<input type="checkbox"/>	024BU8		Adjustable	\$175.00
<input type="checkbox"/>	024BU16		Reinforced (NA w/ adj. upholstery)	\$30.00
<input type="checkbox"/>	024BU14		Omit Back Upholstery	-\$15.00

Backrest Options

<input type="checkbox"/>	024BT33		Folding w/ Integral Push Handles	NC
<input type="checkbox"/>	024A85		Bolt-On Push Handle (Folding only)	\$115.00
<input type="checkbox"/>	024SR2		Regular Stroller Handle <i>NA with amputee axle plate</i>	\$125.00
<input type="checkbox"/>	024SR3		Tall Stroller Handle <i>NA with amputee axle plate</i>	\$125.00

JAY Back (Optional)

<input type="checkbox"/>	024JB01		J3 Back (See J3 Order Form)	Price List
--------------------------	---------	--	-----------------------------	------------

8. ARMREST / SIDEGUARDS

Single Post

<input type="checkbox"/>	024AA9	E0973	Single Post Ht. Adjust. Low <i>Must be removed to fold back down</i>	\$285.00
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Contour (14")	NC
<input type="checkbox"/>	024AA8	E0973	Single Post Ht. Adjust. Std <i>Must be removed to fold back down</i>	\$285.00
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Contour (14")	NC

Padded Swing Away Armrest

All Padded Swing-Away NA with amputee axle plate or Xtender Power Assist

<input type="checkbox"/>	024AA19	E0973	Padded S-A Medium	\$210.00
<input type="checkbox"/>	024AA16	E0973	Padded S-A Adult	\$210.00

Side Guard

<input type="checkbox"/>	024SG3	K0108	Plastic Regular	\$125.00
<input type="checkbox"/>	024SG4	K0108	Plastic Kids	\$125.00

9. KEY OPTIONS AND ACCESSORIES

Rear Anti-Tip

<input type="checkbox"/>	024AT1	E0971	Rear Anti-Tip	\$140.00
--------------------------	--------	-------	---------------	----------

Backpack

<input type="checkbox"/>	024A1		Quickie Backpack Black	\$55.00
<input type="checkbox"/>	024A2		Quickie Kids Backpack Black	\$50.00

Seat Pouch

<input type="checkbox"/>	024A8		Quickie Seat Pouch Black	\$45.00
--------------------------	-------	--	--------------------------	---------

Positioning Belts

<input type="checkbox"/>	024A142	E0978	Belt 1.5" Auto Buckle w/Clamps	\$70.00
<input type="checkbox"/>	024A12	E0978	Belt 2" Aircraft Buckle w/Clamps	\$110.00
<input type="checkbox"/>	024A13	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$120.00

Spoke Guards

NA w/ Xtender Power Assist; Only available w/ Mag and Spoke wheels

<input type="checkbox"/>	024A23	K0065	Black Spoke Guards	\$115.00
<input type="checkbox"/>	024A24	K0065	Clear Spoke Guards	\$115.00

Crutch Holder

<input type="checkbox"/>	024A53	E2207	Left Side	\$195.00
<input type="checkbox"/>	024A52	E2207	Right Side	\$195.00

Extra Leg Strap

<input type="checkbox"/>	024A55	K0038	Extra Leg Strap Black	\$30.00
--------------------------	--------	-------	-----------------------	---------

Miscellaneous

<input type="checkbox"/>	024A83		Tool Kit	\$45.00
--------------------------	--------	--	----------	---------

Tray Table

<input type="checkbox"/>	024A31	E0950	Tray Table	\$250.00
--------------------------	--------	-------	------------	----------

10. COLOR

Standard Frame Colors

<input type="checkbox"/>	024C140	Aztec Gold	NC
<input type="checkbox"/>	024C2	Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	024C141	Black Cherry	NC
<input type="checkbox"/>	024C10	Black Opal	NC
<input type="checkbox"/>	024C38	Blue Opal	NC
<input type="checkbox"/>	024C36	Candy Blue	NC
<input type="checkbox"/>	024C14	Candy Purple	NC
<input type="checkbox"/>	024C35	Candy Red	NC
<input type="checkbox"/>	024C142	Electric Blue	NC
<input type="checkbox"/>	024C143	Evergreen	NC
<input type="checkbox"/>	024C126	Glow	NC
<input type="checkbox"/>	024C144	Green Apple	NC
<input type="checkbox"/>	024C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	024C148	Mauve Pink	NC
<input type="checkbox"/>	024C103	Rootbeer	NC
<input type="checkbox"/>	024C114	Sparkle Silver	NC
<input type="checkbox"/>	024C145	Sunrise Orange	NC
<input type="checkbox"/>	024C113	Titanium Color Paint	NC
<input type="checkbox"/>	024C8	Yellow	NC

Matte Frame Colors

<input type="checkbox"/>	024C116	Matte Black	NC
<input type="checkbox"/>	024C151	Matte Black Cherry	NC
<input type="checkbox"/>	024C150	Matte Electric Blue	NC
<input type="checkbox"/>	024C149	Matte Evergreen	NC
<input type="checkbox"/>	024C152	Matte Purple	NC

Custom Frame Color

<input type="checkbox"/>	024C112	Desert Camouflage	\$295.00
<input type="checkbox"/>	024C43	Mossy Oak Camouflage	\$295.00
<input type="checkbox"/>	024C42	Stars & Stripes	\$295.00
<input type="checkbox"/>	024C108	Zebra (Choose above base color) Std with White unless above base paint is chosen NA w/ Black, Matte Black, Black Opal or Titanium.	\$295.00

11. APPENDIX

GP Series Rear STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

*** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			3"	4" Poly	5" Poly	6" Poly	6" Pneu	8" Poly	8" Pneu	8x2
Fork	Spacer	Fork Hole	4" Lighted	5" Soft Roll	6" x 1.5					
			4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00	x
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50	x
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50	x
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00	x
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50	x
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50	x
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top	x	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	x	21.00

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

Frame Dimensions

