

Account Information

Company Name: _____
Address: _____
City, State Zip: _____
Phone: _____
Fax: _____
Email: _____
(please provide email for shipping notification)

Mark For: _____
Ship To: _____
Address: _____
City, State Zip: _____
Phone: _____
Purchase Order: _____

Order Information

Buddy Roamer Size 1

Lower Leg 12-27"
Under Arm 23-40"
Hip Width 7-10"
Trunk Width 7-10"
Headrest Height 27-51"
Max User's Weight 110 lb
Unit Width 25"
Unit Length 31"

BASIC SYSTEM INCLUDES rolling base with 8" rear wheels and casters, hip support, telescoping headrest, adjustable dynamic spring, rear wheel locks and anti tippers. **(Must select one chest support and head support)**

Buddy Roamer Size 2

Lower Leg 16-32"
Under Arm 27-43"
Hip Width 8-14"
Trunk Width 8-14"
Headrest Height 31-55"
Max User's Weight 110 lb
Unit Width 25"
Unit Length 31"

BASIC SYSTEM INCLUDES rolling base with 8" rear wheels and casters, hip support, telescoping headrest, adjustable dynamic spring, rear wheel locks and anti tippers. **(Must select one chest support and head support)**

- ☐ **ROA 01** Buddy Roamer Size 1 4,600
- choose one**
- ☐ **ROA 01A** Lateral Chest Support \$402
- ☐ **BRSA-008** Lateral Chest Support, with shoulder straps \$894
- choose one**
- ☐ **ROA 007S** Standard Headrest 4" W x 6" H \$180
- ☐ **ROA 007L** Wider Headrest 9.5" W x 6" H \$234



- ☐ **ROA 02** Buddy Roamer Size 2 4,650.00
- choose one**
- ☐ **ROA 02B** Lateral Chest Support \$402
- ☐ **BRSA-008** Lateral Chest Support, with shoulder straps \$894
- choose one**
- ☐ **ROA 007S** Standard Headrest 4" W x 6" H \$180
- ☐ **ROA 007L** Wider Headrest 9.5" W x 6" H \$234

Options

Buddy Roamer Options, size 1 and size 2

- ☐ **ROA 009** Hand Holds- individual handle type \$1,400
- ☐ **ROA 009B** Hand Holds-single horizontal front bar \$1,400
(always must be removed for transfers)
- ☐ **ROA 012** Abduction Pommel (replaces sit bar) \$595
(11" L x 12" H x 2.25" W)
- ☐ **ROA 011** No Roll Back Option \$480
- ☐ **ROA 002A** Front Castor Directional Swivel Lock \$120
- ☐ **ROA 005** Custom Hip/Chest Lateral Width as per quote

How to Measure



Gait is dynamic! Your walker must be too!