

## Companion Model Comparison



\*Shown with optional wrap



	Companion	Companion Plus
<b>Motor</b>	250W hub motor	350W hub motor with NEW more powerful controller
<b>Brakes</b>	Drum brakes	Drum brakes with NEW parking brake function
<b>Tire</b>	Solid rubber tire	NEW non-marking tire
<b>Battery</b>	Lithium-ion battery 36V 7.5Ah	NEW lithium-ion battery 36V 7.5Ah with improved range
<b>Range</b>	Up to 25km/16mi on a full charge	Up to 30km/18.6mi on a full charge
<b>Speed</b>	Three-speed control - max. 24kmh/15mph	NEW Five-speed control - max. 26kmh/16mph
<b>Power Switch</b>	ON/OFF keyed power switch	ON/OFF keyed power switch
<b>Turning Radius</b>	180 degree handle angle range Fully adjustable handlebars with rubber grip handles	180 degree handle angle range Fully adjustable handlebars with NEW soft multi-color grip handles
<b>Display</b>	2-inch LCD display - battery charge status, trip distance, and speed	NEW programmable display with 5 speed settings, speedometer, odometer, trip distance, and battery status
<b>Weight</b>		
Companion Unit	9.4kgs/20lbs	9.4kgs/20lbs
Adaptor	0.8kg/1.7lbs	0.8kgs/1.7lbs
Rigid Mount	1.6kg/3.6lbs	1kg/2.2lbs
Folding Mount	17-19in 1.55kg/3.5lbs	17-19in 1.55kg/3.5lbs