



Quickie QX

February 2023

EIF6

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base Model

☐	196WM94	K0004	Quickie QX (Fixed Axle) 300 lbs For 14" sizes use E1236	\$2,006.00
--------------------------	---------	-------	--	------------

Frame Style

See Seat Height Matrix on page 2

☐	196FS9	Standard Frame	NC
☐	196FS30	Hemi Frame	NC

2. SIZE / SEATING

Frame Width

☐	196W5		14"	NC
☐	196W6		15"	NC
☐	196W7		16"	NC
☐	196W8		17"	NC
☐	196W9		18"	NC
☐	196W10		19"	NC
☐	196W11	E2201	20"	\$559.00

Frame Width - Built 4 Me Options

Not available with EXO standard or Tension adjustable upholstery, Not available with Ballistic Nylon or 3DX Tension adjustable upholstery.

☐	196W59	E2201	B4ME - 21" Width	\$614.00
☐	196W60	E2201	B4ME - 22" Width	\$614.00

Seat Sling

☐	196SL1		Standard Seat Sling	NC
☐	196SL6		Omit Seat Sling	-\$15.00

Seat Depth

Frame Depths equals Seat Depths.

☐	196SD7		14"	NC
☐	196SD8		15"	NC
☐	196SD9		16"	NC
☐	196SD10		17"	NC
☐	196SD11		18"	NC
☐	196SD71		19"	NC
☐	196SD13	E2203	20"	\$614.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

JAY Bundled Cushion		Price List
Cushion Name or Part #		
☐ 196JC20	Width: _____ Depth: _____	

Personalized Seating

☐	196QC01	Surefit	
☐	196QC02	Configurfit	
☐	196QC03	JayFit	



Quickie QX, Ultra Lightweight Folding Wheelchair



Quickie QX featuring the Jay X2™ cushion

3. CASTERS / SEAT HEIGHT

Front Seat Height						
See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger						
12.5" & 13" only available w/ Alum Angle Adj Footplate						
☐	196SH34	12.5"	NC	☐	196SH14	16.0"
☐	196SH35	13.0"	NC	☐	196SH15	16.5"
☐	196SH40	13.5"	NC	☐	196SH16	17.0"
☐	196SH41	14.0"	NC	☐	196SH17	17.5"
☐	196SH42	14.5"	NC	☐	196SH18	18.0"
☐	196SH12	15.0"	NC	☐	196SH20	18.5"
☐	196SH13	15.5"	NC	☐	196SH21	19.0"
				☐	196SH23	19.5"

Caster Wheels				
☐	196CW3		4" Low Profile Polyurethane	NC
☐	196CW44	E2219	4" x 1.25" Semi Pneumatic	\$61.00
☐	196CW5		5" Low Profile Polyurethane	NC
☐	196CW7	E2219	5" x 1.5" Semi Pneumatic	\$96.00
☐	196CW34		6" Polyurethane	NC
☐	196CW35	E2219	6" x 1.5" Semi Pneumatic	\$96.00
☐	196CW32		8" Polyurethane	NC
☐	196CW39	E2219	8" x 1.5" Semi Pneumatic	\$96.00

Caster Fork				
☐	196CF71		3.5" (standard stem bolt only)	NC
☐	196CF12		5.25"	NC
☐	196CF17		6"	NC
☐	196CF23		7"	NC

4. REAR WHEELS / REAR SEAT HEIGHTS

Rear Seat Height				
See Rear Seat Height Matrix. 12.5" & 18.5" RSH not available with fixed axle.				
☐	196RH9		12"	NC
☐	196RH11		13"	NC
☐	196RH12		13.5"	NC
☐	196RH13		14"	NC
☐	196RH14		14.5"	NC
☐	196RH15		15"	NC
☐	196RH16		15.5"	NC
☐	196RH17		16"	NC
☐	196RH18		16.5"	NC
☐	196RH19		17"	NC
☐	196RH26		17.5"	NC
☐	196RH20		18"	NC
☐	196RH21		19"	NC

Rear Wheel Size				
See Rear Tire & Wheel Matrix for availability				
☐	196RW3		20" Mag Wheel	NC
☐	196RW8		22" Mag Wheel	NC
☐	196RW1		24" Mag Wheel	NC
☐	196RW103		20" Mag Wheel 5 Spoke	NC
☐	196RW104		22" Mag Wheel 5 Spoke	NC
☐	196RW102		24" Mag Wheel 5 Spoke	NC
☐	196RW90		20" Lite Spoke Wheel	\$224.00
☐	196RW91		22" Lite Spoke Wheel	\$224.00
☐	196RW92		24" Lite Spoke Wheel	NC
☐	196RW94		26" Lite Spoke Wheel	\$224.00
☐	196RW25		Omit	NC

Caster Stem Bolt				
☐	196CS1		Standard	NC
☐	196CS5		3/4"	NC
☐	196CS9		1 1/2"	NC

Front Seat Height Matrix								
Front Seat Height Matrix - Standard Stem Bolt								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	19.5	17.0	18.5
Seat squeeze and seat depth can affect FSH availability								
For 3/4" and 1 1/2" Stem bolt matrix see page 5								

Rear Wheel Tire				
See Rear Tire & Wheel Matrix for availability				
☐	196RT1		Pneumatic	NC
☐	196RT2	E2211 E2213	Pneumatic w/ Airless Insert	\$96.00
☐	196RT3		Full Poly	NC
☐	196RT39	E2211 E2212	High-Pressure Clincher Black	\$106.00

Handrim				
Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic N/A w/ 25, 26"				
☐	196HR1		Aluminum Anodized	NC
☐	196HR4		Plastic	\$123.00
☐	196HR22		Natural Fit - Standard Grip	\$526.00
☐	196HR23		Natural Fit - Super Grip	\$585.00

Seat Squeeze Matrix		
Frame Depth	Max Seat Squeeze	
14"	2.0"	
15" - 16"	2.5"	
17" - 18"	2.5"	
19" - 20"	3.0"	

Rear Seat Height Matrix			
Rear Wheel Size	Available Rear Seat Heights		
	Pneumatic Tire	Full Poly Tire	HP Clincher Tire
20"	12" - 16"	12" - 16"	N/A
22"	13.5" - 17.5"	13" - 17"	13" - 17"
24"	14" - 18"	13.5" - 17.5"	13.5" - 17.5"
26"	15" - 19"	N/A	15" - 19"
12.5" and 18.5" RSH not available with QX. Values above assume flat seat hights. RSH availability may decrease in value with increased seat angle.			

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price																																										
4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)																																																			
Rear Tire & Wheel Matrix					One Arm Drive																																														
<table><tr><th rowspan="2">Tire</th><th colspan="3">Rear Wheel</th></tr><tr><th>5 Spoke Mag</th><th>Lite Spoke</th><th>Mag Wheel</th></tr><tr><td>Pneumatic</td><td>20,22,24</td><td>20,22,24,26</td><td>20,22,24</td></tr><tr><td>Pneumatic w/ Airless Insert</td><td>20,22,24</td><td>20,22,24,26</td><td>20,22,24</td></tr><tr><td>Full Poly</td><td>20,22,24</td><td>20,22,24</td><td>20,22,24</td></tr><tr><td>High-Pressure Clincher Blk</td><td>N/A</td><td>22,24,26</td><td>N/A</td></tr><tr><td colspan="4">20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm</td></tr></table>					Tire	Rear Wheel			5 Spoke Mag	Lite Spoke	Mag Wheel	Pneumatic	20,22,24	20,22,24,26	20,22,24	Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24	Full Poly	20,22,24	20,22,24	20,22,24	High-Pressure Clincher Blk	N/A	22,24,26	N/A	20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm				Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim. 1AD option is only available on 14" - 19" frame widths. <table><tr><td>☐</td><td>196AL16</td><td>E0958</td><td>One Arm Drive Adj. Axle Plate</td><td>\$950.00</td></tr><tr><td></td><td>196HM1</td><td></td><td>☐ Right Hand Mount</td><td>NC</td></tr><tr><td></td><td>196HM2</td><td></td><td>☐ Left Hand Mount</td><td>NC</td></tr></table>					☐	196AL16	E0958	One Arm Drive Adj. Axle Plate	\$950.00		196HM1		☐ Right Hand Mount	NC		196HM2		☐ Left Hand Mount	NC
Tire	Rear Wheel																																																		
	5 Spoke Mag	Lite Spoke	Mag Wheel																																																
Pneumatic	20,22,24	20,22,24,26	20,22,24																																																
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24																																																
Full Poly	20,22,24	20,22,24	20,22,24																																																
High-Pressure Clincher Blk	N/A	22,24,26	N/A																																																
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm																																																			
☐	196AL16	E0958	One Arm Drive Adj. Axle Plate	\$950.00																																															
	196HM1		☐ Right Hand Mount	NC																																															
	196HM2		☐ Left Hand Mount	NC																																															
					<table><tr><th colspan="4">Axle</th></tr><tr><td>☐</td><td>196AX1</td><td></td><td>Quick Release Steel</td><td>NC</td></tr><tr><td>☐</td><td>196AX11</td><td></td><td>Quad Quick Release Steel N/A on Spoke</td><td>\$141.00</td></tr><tr><td>☐</td><td>196AX45</td><td></td><td>Threaded Axle</td><td>\$111.00</td></tr></table>					Axle				☐	196AX1		Quick Release Steel	NC	☐	196AX11		Quad Quick Release Steel N/A on Spoke	\$141.00	☐	196AX45		Threaded Axle	\$111.00																							
Axle																																																			
☐	196AX1		Quick Release Steel	NC																																															
☐	196AX11		Quad Quick Release Steel N/A on Spoke	\$141.00																																															
☐	196AX45		Threaded Axle	\$111.00																																															

5. WHEEL LOCKS					Extension Handle				
Wheel Lock					Ext. handles only available with high mount push or pull wheel locks				
☐	196WL1		High Mount Push	NC	☐	196EH7	E0961	Extension Handle (Pair)	\$84.00
☐	196WL2		High Mount Pull	NC					
☐	196WL30		Compact	\$157.00					
☐	196WL24		Unilateral Wheel Lock Right	\$335.00					
☐	196WL23		Unilateral Wheel Lock Left	\$335.00					
☐	196WL10		Omit Wheel Lock	NC					
☐	196WL12		Do Not Mount	NC					

6. HANGERS AND FOOTRESTS					Extension Tubes				
Front Hangers					See Lower Leg Length Matrix for range				
☐	196H145		60° Swing In/Out Hanger	NC	☐	196E1		Short	NC
☐	196H139		70° Swing In/Out Hanger	NC	☐	196E2		Medium	NC
☐	196H141		80° Swing In/Out Hanger	NC	☐	196E5		Long	NC
☐	196H159	K0053	Articulating Hanger N/A w/ Alum Angle Adj Front Mount footplates.	\$413.00	Footplate				
☐	196H160	E0990	Elevating Hanger - Ext Mount N/A w/ Alum Angle Adj Front Mount footplates.	\$352.00	Heel loop not included, see Footplate Option section below				
☐	196H161	E0990	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$352.00	☐	196F1		Composite Footplate	NC
☐	196H144		Omit Hanger & Footplate	NC	☐	196F271	K0040 E0951	Composite Angle Adj-(16"-22"w)	\$202.00
Articulating / Elevating Hanger Cover					☐	196F3	K0040	Aluminum Angle Adj. Footplate	\$212.00
☐	196DI47		Neoprene Cover	\$380.00	☐	196F5		Aluminum Only available w/Elevating Legrest	NC
☐	196DI48		Gel Pad Cover	\$447.00	Hanger & Footplate Options				
					☐	196HL12	E0951	Adjustable Heel Loop	\$40.00
					☐	1968A02		Gel Cover For Hangers (Pair)	\$123.00

7. BACKREST / SEATING					Push Handles				
Back Type					Integral Push Handles, Folding Height Adjustable & Omit 8° for Std Backrest, Straight for Angle Adj. Backrest				
☐	196BT10		Standard	NC	☐	196PH6		Integral Push Handles	NC
☐	196BT100		Angle Adjustable Back	\$184.00	☐	196PH5		Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height	\$111.00
Backrest Height					☐	196PH17		Folding Push Handle with Ht. Adj. Back	\$252.00
☐	196B5		Med (13"-17")	NC	☐	196PH2		Omit Push Handles	NC
☐	196B9		Tall (16"-20")	NC					

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING (CONT)									
Back Upholstery					Military Service Patches				
☐	196BU1		Standard Upholstery	NC	☐	196VU1		US Marines Patch	\$111.00
☐	196BU88		Standard Upholstery EXO	\$56.00	☐	196VU2		US Army Patch	\$111.00
☐	196BU70	E2611	Backrest Tension Adj Ballistic Nylon <i>*Reference 186BU70 for billing</i>	\$352.00	☐	196VU3		US Navy Patch	\$111.00
☐	196BU71	E2611	Backrest Tension Adj 3DX Vented <i>*Reference 186BU71 for billing</i>	\$370.00	☐	196VU4		US Air Force Patch	\$111.00
☐	196BU89	E2611	Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i>	\$370.00	☐	196VU12		US Coast Guard	\$111.00
☐	196BU14		Omit Back Upholstery	-\$15.00	☐	196VU5		Purple Heart	\$111.00
JAY Backs					☐	196VU8		Purple Heart with Marines Patch	\$173.00
☐	196JB01		JAY Back (see Order forms)	Price List	☐	196VU9		Purple Heart with Army Patch	\$173.00
					☐	196VU10		Purple Heart with Navy Patch	\$173.00
					☐	196VU11		Purple Heart with Air Force Patch	\$173.00
					☐	196VU13		Purple Heart with Coast Guard	\$173.00
8. ARMREST / SIDEGUARDS									
Armrest - Padded Swing Away					Armpad - For Height Adjustable				
Some limitations when used with the Angle Adjustable Back in the rear most setting					☐	196AM23		Desk Length (10")	NC
Please circle a setting, will default to 11" if no selection is made					☐	196AM24		Full Length - Full (14")	NC
☐	196AA18	E0973	Padded Swing-Away STD (9" to 13")	\$252.00					
		Setting	☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"						
Armrest - Height Adjustable									
Flip Back HA Classic N/A w / 14"-15" depth									
☐	196AA104	E0973	Single Post Ht. Adj - Std (9.75"-14.75")	\$319.00					
☐	196AA105	E0973	Single Post Ht Adj - Low (8.25"-12.25")	\$319.00					
☐	196AA65	E0973	Flip Back Ht Adj LITE (10"-15")	\$352.00					
☐	196AA103	E0973	Dual Post Flip Back Ht Adj	\$330.00					
☐	196AA159	E0973	Flip Back Ht Adj Classic (11"-15")	\$352.00					
☐	196AA14		Omit Armrest	NC					
9. KEY OPTIONS									
Anti Tips					Pack Accessories				
☐	196AT1	E0971	Anti Tips (Pair)	\$168.00	☐	196A1		Backpack - Black/Charcoal Adult	\$84.00
☐	196AT86	E0971	Anti Tips (Pair) Low Profile	\$168.00	☐	1961A01		Seat Pouch Black	\$56.00
Positioning Belt									
☐	196A141	E0978	1.5" Auto Buckle w/ Clamps	\$78.00					
☐	196A804		Positioning Belt Slide Clamps	\$56.00					
JAY Positioning									
☐	196K1	E1020	JAY Lower Extremity S/A support	\$441.00					
			☐ Left ☐ Right Position						
☐	196K2		Omit Hanger	-\$10.00					
			☐ Left ☐ Right						
☐	196K7		Pad Size	NC					
			☐ 6" x 8" ☐ 8" x 10"						
10. COLOR									
Standard Frame Colors					Matte Frame Colors				
☐	196C140		Aztec Gold	NC	☐	196C116		Matte Black	NC
☐	196C2		Black (Not avail as base w/Zebra)	NC	☐	196C151		Matte Black Cherry	NC
☐	196C141		Black Cherry	NC	☐	196C150		Matte Electric Blue	NC
☐	196C10		Black Opal	NC	☐	196C149		Matte Evergreen	NC
☐	196C38		Blue Opal	NC	☐	196C152		Matte Purple	NC
☐	196C36		Candy Blue	NC					
☐	196C14		Candy Purple	NC					
☐	196C35		Candy Red	NC					
☐	196C142		Electric Blue	NC					
☐	196C143		Evergreen	NC					
☐	196C44		Glossy White	NC					
☐	196C126		Glow	NC					
☐	196C144		Green Apple	NC					
☐	196C147		Hot Sparkle Pink	NC					
☐	196C148		Mauve Pink	NC					
☐	196C103		Rootbeer	NC					
☐	196C114		Sparkle Silver	NC					
☐	196C145		Sunrise Orange	NC					
☐	196C113		Titanium Color Paint	NC					
☐	196C8		Yellow	NC					
					Personalized Frame Color				
					☐	196C112		Desert Camouflage	\$441.00
					☐	196C43		Mossy Oak Camouflage	\$441.00
					☐	196C42		Stars & Stripes	\$441.00
					☐	196C108		Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Blk, Matte Black, Blk Opal, Titanium)	\$441.00
Additional Options Through Built for Me									
To inquire about any unique personalized Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011									
									



APPENDIX

Front Seat Height Matrix Cont.

Front Seat Height Matrix - 3/4" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5	
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5	
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5	
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5	
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5	

Seat squeeze and seat depth can affect FSH availability


Front Seat Height Matrix - 1 1/2" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
6 x 1.5 Semi Pneu	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5	
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5	

Seat squeeze and seat depth can affect FSH availability

Lower Leg Length Range Matrix

Extension Mount Hanger Leg Lengths					
Hanger	Ext Tube	Seat-to-Footrest Range			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A
	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A

All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).



Min Allowable Front Seat Heights (for Extension Mount Hangers)

Hanger	Ext Tube	Min Allowable FSH			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	14.0	13.5	12.5	N/A
	Medium	14.5	14.5	12.5	N/A
	Long	16.0	16.0	13.5	N/A
70° Swing In/Out	Short	14.0	14.0	13.0	N/A
	Medium	15.0	15.0	13.5	N/A
	Long	16.5	16.5	15.0	N/A
80° Swing In/Out	Short	14.0	14.0	13.5	N/A
	Medium	14.5	14.5	14.5	N/A
	Long	16.5	16.5	15.5	N/A

ALR & ELR Hanger Leg Lengths

Hanger	Extension	Seat Sling to Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	Short	8	11	10
	Med	10	13	12
	Long	12	15	14

Rear Wheel & Handrim Matrix

Handrim	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Aluminum Anodized	20,22,24	20,22,24,26	20,22,24
Plastic Coat	20,22,24	22,24	20,22,24
Natural Fit - Standard Grip	24	22,24,26	22,24
Natural Fit - Super Grip	24	22,24,26	22,24



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA