



Quickie QXi

February 2023

EIF6

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ ST: _____ ZIP: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ ST: _____ ZIP: _____
Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base Model

Weight capacity 300 lbs

☐	196WM78	K0005	Quickie QXi (Adjustable Axle) For 14" sizes use E1236	\$2,084.00
--------------------------	---------	-------	--	------------

Frame Style

See Seat Height Matrix on page 2

☐	196FS9		Standard Frame	NC
☐	196FS30		Hemi Frame	NC

Transit

☐	196TO38	K0108	Transit System WC19	\$359.00
--------------------------	---------	-------	---------------------	----------

2. SIZE / SEATING

Frame Width

☐	196W5		14"	NC
☐	196W6		15"	NC
☐	196W7		16"	NC
☐	196W8		17"	NC
☐	196W9		18"	NC
☐	196W10		19"	NC
☐	196W11	E2201	20"	\$559.00



Frame Width - Built 4 Me Options

Not available with EXO standard or Tension adjustable upholstery, Not available with Ballistic Nylon or 3DX Tension adjustable upholstery.

☐	196W59	E2201	B4ME - 21" Width	\$614.00
☐	196W60	E2201	B4ME - 22" Width	\$614.00

Seat Sling

☐	196SL1		Standard Seat Sling	NC
☐	196SL6		Omit Seat Sling	-\$15.00

Seat Depth

Frame Depths equals Seat Depths.

☐	196SD7		14"	NC
☐	196SD8		15"	NC
☐	196SD9		16"	NC
☐	196SD10		17"	NC
☐	196SD11		18"	NC
☐	196SD71		19"	NC
☐	196SD13	E2203	20"	\$614.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	196JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width: _____	Depth: _____	

Personalized Seating

☐	196QC01	Surefit	
☐	196QC02	Configurfit	
☐	196QC03	JayFit	



Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair



Quickie QXi featuring the Jay ION™ cushion

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
3. CASTERS / SEAT HEIGHT									

Front Seat Height							
See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger							
12.5" & 13" only available w/ Alum Angle Adj Footplate							
☐	196SH34	12.5"	NC	☐	196SH14	16.0"	NC
☐	196SH35	13.0"	NC	☐	196SH15	16.5"	NC
☐	196SH40	13.5"	NC	☐	196SH16	17.0"	NC
☐	196SH41	14.0"	NC	☐	196SH17	17.5"	NC
☐	196SH42	14.5"	NC	☐	196SH18	18.0"	NC
☐	196SH12	15.0"	NC	☐	196SH20	18.5"	NC
☐	196SH13	15.5"	NC	☐	196SH21	19.0"	NC
				☐	196SH23	19.5"	NC

Caster Wheels			
☐	196CW3	4" Low Profile Polyurethane	NC
☐	196CW44	E2219 4" x 1.25" Semi Pneumatic	\$61.00
☐	196CW5	5" Low Profile Polyurethane	NC
☐	196CW7	E2219 5" x 1.5" Semi Pneumatic	\$96.00
☐	196CW34	6" Polyurethane	NC
☐	196CW35	E2219 6" x 1.5" Semi Pneumatic	\$96.00
☐	196CW32	8" Polyurethane	NC
☐	196CW39	E2219 8" x 1.5" Semi Pneumatic	\$96.00

Caster Fork			
☐	196CF71	3.5" (standard stem bolt only)	NC
☐	196CF12	5.25"	NC
☐	196CF17	6"	NC
☐	196CF23	7"	NC

Caster Stem Bolt			
☐	196CS1	Standard	NC
☐	196CS5	3/4"	NC
☐	196CS9	1 1/2"	NC

Front Seat Height Matrix

Front Seat Height Matrix - Standard Stem Bolt								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	19.5	17.0	18.5
Seat squeeze and seat depth can affect FSH availability								
For 3/4" and 1 1/2" Stem bolt matrix see page 5								

4. REAR WHEELS / REAR SEAT HEIGHTS

Rear Seat Height			
See Rear Seat Height Matrix			
☐	196RH9	12"	NC
☐	196RH10	12.5"	NC
☐	196RH11	13"	NC
☐	196RH12	13.5"	NC
☐	196RH13	14"	NC
☐	196RH14	14.5"	NC
☐	196RH15	15"	NC
☐	196RH16	15.5"	NC
☐	196RH17	16"	NC
☐	196RH18	16.5"	NC
☐	196RH19	17"	NC
☐	196RH26	17.5"	NC
☐	196RH20	18"	NC
☐	196RH27	18.5"	NC
☐	196RH21	19"	NC

Rear Wheel Size			
See Rear Tire & Wheel Matrix for availability			
☐	196RW3	20" Mag Wheel	NC
☐	196RW8	22" Mag Wheel	NC
☐	196RW1	24" Mag Wheel	NC
☐	196RW103	20" Mag Wheel 5 Spoke	NC
☐	196RW104	22" Mag Wheel 5 Spoke	NC
☐	196RW102	24" Mag Wheel 5 Spoke	NC
☐	196RW90	20" Lite Spoke Wheel	\$224.00
☐	196RW91	22" Lite Spoke Wheel	\$224.00
☐	196RW92	24" Lite Spoke Wheel	NC
☐	196RW93	25" Lite Spoke Wheel	\$224.00
☐	196RW94	26" Lite Spoke Wheel	\$224.00
☐	196RW25	Omit	NC

Rear Seat Height Matrix		
Rear Wheel Size	Std Frame	Hemi Frame
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"
24" Rear Wheel	15.0" - 18.5"	13.5" - 16.5"
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.0"
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"

Rear Wheel Tire			
See Rear Tire & Wheel Matrix for availability			
☐	196RT1	Pneumatic	NC
☐	196RT2	E2211 E2213 Pneumatic w/ Airless Insert	\$96.00
☐	196RT3	Full Poly	NC
☐	196RT39	E2211 E2212 High-Pressure Clincher Black	\$106.00
☐	196RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$157.00

Handrim			
Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic N/A w/ 25, 26"			
See Appendix for Rear wheel / Handrim matrix			
☐	196HR1	Aluminum Anodized	NC
☐	196HR4	Plastic	\$123.00
☐	196HR22	Natural Fit - Standard Grip	\$526.00
☐	196HR23	Natural Fit - Super Grip	\$585.00
☐	196HR32	Natural Fit - No Thumb Grip	\$455.00
☐	196HR26	The Surge - Oval w/Gription Strip	\$472.00
☐	196HR27	The Surge LT - Oval w/Gription Strip	\$508.00
☐	196HR28	Q-Grip Neoprene Coated	\$414.00

Seat Squeeze Matrix	
Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)

Rear Tire & Wheel Matrix

Tire	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Pneumatic	20,22,24	20,22,24,26	20,22,24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24
Full Poly	20,22,24	20,22,24	20,22,24
High-Pressure Clincher Blk	N/A	22,24,26	N/A
Schwalbe Marathon	N/A	24,25,26	N/A
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm			

One Arm Drive

Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim. 1AD option is only available on 14" - 19" frame widths.

☐	196AL16	E0958	One Arm Drive Adj. Axle Plate	\$950.00
	196HM1		☐ Right Hand Mount	NC
	196HM2		☐ Left Hand Mount	NC

Axle

☐	196AX1		Quick Release Steel	NC
☐	196AX11		Quad Quick Release Steel N/A on Spoke	\$141.00
☐	196AX45		Threaded Axle	\$111.00

Adjustable Axle Options

☐	196AS4		Adjustable sleeve	NC
☐	196AS5		Fixed Sleeve 0°	NC
☐	196AS6		Fixed Sleeve 3°	NC

Center of Gravity

Measure from front of back post to center of rear axle

☐	196CC2		1"	NC
☐	196CC3		1.5"	NC
☐	196CC4		2"	NC
☐	196CC5		2.5"	NC
☐	196CC6		3"	NC

5. WHEEL LOCKS

Wheel Lock

☐	196WL1		High Mount Push	NC
☐	196WL2		High Mount Pull	NC
☐	196WL30		Compact	\$157.00
☐	196WL24		Unilateral Wheel Lock Right	\$335.00
☐	196WL23		Unilateral Wheel Lock Left	\$335.00
☐	196WL10		Omit Wheel Lock	NC
☐	196WL12		Do Not Mount	NC

Extension Handle

Ext. handles only available with high mount push or pull wheel locks

☐	196EH7	E0961	Extension Handle (Pair)	\$84.00
--------------------------	--------	-------	-------------------------	---------

6. HANGERS AND FOOTRESTS

Front Hangers

Footplate mounts with Extension Tube using clamping hardware

☐	196H145		60° Swing In/Out Hanger	NC
☐	196H139		70° Swing In/Out Hanger	NC
☐	196H141		80° Swing In/Out Hanger	NC
☐	196H159	K0053	Articulating Hanger N/A w/ Alum Angle Adj Front Mount footplates.	\$413.00
☐	196H160	E0990	Elevating Hanger - Ext Mount N/A w/ Alum Angle Adj Front Mount footplates.	\$352.00
☐	196H161	E0990	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$352.00
☐	196H144		Omit Hanger & Footplate	NC

Articulating / Elevating Hanger Cover

☐	196DI47		Neoprene Cover	\$380.00
☐	196DI48		Gel Pad Cover	\$447.00

Extension Tubes

See Lower Leg Length Matrix for range

☐	196E1		Short	NC
☐	196E2		Medium	NC
☐	196E5		Long	NC

Footplate

Heel loop not included, see Footplate Option section below

☐	196F1		Composite Footplate	NC
☐	196F271	K0040 E0951	Composite Angle Adj-(16"-22"w)	\$202.00
☐	196F3	K0040	Aluminum Angle Adj. Footplate	\$212.00
☐	196F5		Aluminum Only available w/Elevating Legrest	NC

Hanger & Footplate Options

☐	196HL12	E0951	Adjustable Heel Loop	\$40.00
☐	1968A02		Gel Cover For Hangers (Pair)	\$123.00

7. BACKREST / SEATING

Back Type

☐	196BT10		Standard	NC
☐	196BT100		Angle Adjustable Back	\$184.00

Backrest Height

☐	196B5		Med (13"-17")	NC
☐	196B9		Tall (16"-20")	NC

Push Handles

Integral Push Handles, Folding Height Adjustable & Omit is 8° for Standard Backrest, Straight for Angle Adj. Backrest

☐	196PH6		Integral Push Handles	NC
☐	196PH5		Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height	\$111.00
☐	196PH17		Folding Push Handle with Ht. Adj. Back	\$252.00
☐	196PH2		Omit Push Handles	NC

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING (CONT)									

Back Upholstery			
☐	196BU1	Standard Upholstery	NC
☐	196BU88	Standard Upholstery EXO	\$56.00
☐	196BU70	E2611 Backrest Tension Adj Ballistic Nylon *Reference 186BU70 for billing	\$352.00
☐	196BU71	E2611 Backrest Tension Adj 3DX Vented *Reference 186BU71 for billing	\$370.00
☐	196BU89	E2611 Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$370.00
☐	196BU14	Omit Back Upholstery	-\$15.00

JAY Backs			
☐	196JB01	JAY Back (see Order forms)	Price List

Military Service Patches			
☐	196VU1	US Marines Patch	\$111.00
☐	196VU2	US Army Patch	\$111.00
☐	196VU3	US Navy Patch	\$111.00
☐	196VU4	US Air Force Patch	\$111.00
☐	196VU12	US Coast Guard	\$111.00
☐	196VU5	Purple Heart	\$111.00
☐	196VU8	Purple Heart with Marines Patch	\$173.00
☐	196VU9	Purple Heart with Army Patch	\$173.00
☐	196VU10	Purple Heart with Navy Patch	\$173.00
☐	196VU11	Purple Heart with Air Force Patch	\$173.00
☐	196VU13	Purple Heart with Coast Guard	\$173.00

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away			
Some limitations when used with the Angle Adjustable Back in the rear most setting Please circle a setting, will default to 11" if no selection is made			
☐	196AA18	E0973 Padded Swing-Away STD (9" to 13") Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"	\$252.00

Armrest - Height Adjustable			
Flip Back HA Classic N/A w / 14"-15" depth			
☐	196AA104	E0973 Single Post Ht. Adj - Std (9.75"-14.75")	\$319.00
☐	196AA105	E0973 Single Post Ht Adj - Low (8.25"-12.25")	\$319.00
☐	196AA65	E0973 Flip Back Ht Adj LITE (10"-15")	\$352.00
☐	196AA103	E0973 Dual Post Flip Back Ht Adj	\$330.00
☐	196AA159	E0973 Flip Back Ht Adj Classic (11"-15")	\$352.00
☐	196AA14	Omit Armrest	NC

Armpad - For Height Adjustable			
☐	196AM23	Desk Length (10")	NC
☐	196AM24	Full Length - Full (14")	NC

Armrest Options			
☐	1968A01	Gel Cover For Swing Away Arms (pair)	\$168.00

Side Guards			
☐	196SG33	K0108 Plastic - Standard (8.5" High)	\$162.00
☐	196SG29	K0108 Plastic - Low (6.5" High)	\$162.00

9. KEY OPTIONS

Anti Tips			
☐	196AT1	E0971 Anti Tips (Pair)	\$168.00
☐	196AT86	E0971 Anti Tips (Pair) Low Profile	\$168.00

Positioning Belt			
☐	196A141	E0978 1.5" Auto Buckle w/ Clamps	\$78.00
☐	196A804	Positioning Belt Slide Clamps	\$56.00

JAY Positioning			
☐	196K1	E1020 JAY Lower Extremity S/A support ☐ Left ☐ Right Position	\$441.00
☐	196K2	Omit Hanger ☐ Left ☐ Right	-\$10.00
☐	196K7	Pad Size ☐ 6" x 8" ☐ 8" x 10"	NC

Pack Accessories			
☐	196A1	Backpack - Black/Charcoal Adult	\$84.00
☐	1961A01	Seat Pouch Black	\$56.00

Spoke Guards			
☐	1964A01	K0065 Quickie Black Spoke Guard	\$128.00
☐	1964A02	K0065 Quickie Clear Spoke Guard	\$128.00

Miscellaneous			
☐	1968A03	Caddy	\$56.00
☐	196A2100	E2207 Crutch Holder	\$218.00
☐	196A999	Quickie Water Bottle & Holder	\$61.00

10. COLOR

Standard Frame Colors			
☐	196C140	Aztec Gold	NC
☐	196C2	Black (Not avail as base w/Zebra)	NC
☐	196C141	Black Cherry	NC
☐	196C10	Black Opal	NC
☐	196C38	Blue Opal	NC
☐	196C36	Candy Blue	NC
☐	196C14	Candy Purple	NC
☐	196C35	Candy Red	NC
☐	196C142	Electric Blue	NC
☐	196C143	Evergreen	NC
☐	196C44	Glossy White	NC
☐	196C126	Glow	NC
☐	196C144	Green Apple	NC
☐	196C147	Hot Sparkle Pink	NC
☐	196C148	Mauve Pink	NC
☐	196C103	Rootbeer	NC
☐	196C114	Sparkle Silver	NC
☐	196C145	Sunrise Orange	NC
☐	196C113	Titanium Color Paint	NC
☐	196C8	Yellow	NC

Matte Frame Colors			
☐	196C116	Matte Black	NC
☐	196C151	Matte Black Cherry	NC
☐	196C150	Matte Electric Blue	NC
☐	196C149	Matte Evergreen	NC
☐	196C152	Matte Purple	NC

Personalized Frame Color			
☐	196C112	Desert Camouflage	\$441.00
☐	196C43	Mossy Oak Camouflage	\$441.00
☐	196C42	Stars & Stripes	\$441.00
☐	196C108	Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Blk, Matte Black, Blk Opal, Titanium)	\$441.00

Additional Options Through Built for Me			
To inquire about any unique personalized Quickie products, contact: Phone: (800) 333-4000, Fax: (800) 333 9011			



APPENDIX


Front Seat Height Matrix Cont.

Front Seat Height Matrix - 3/4" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5	
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5	
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5	
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5	
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5	
Seat squeeze and seat depth can affect FSH availability									

Front Seat Height Matrix - 1 1/2" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
6 x 1.5 Semi Pneu	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5	
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5	
Seat squeeze and seat depth can affect FSH availability									

Lower Leg Length Range Matrix

Extension Mount Hanger Leg Lengths					
Hanger	Ext Tube	Seat-to-Footrest Range			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A
	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A
All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).					



Min Allowable Front Seat Heights (for Extension Mount Hangers)					
Hanger	Ext Tube	Min Allowable FSH			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	14.0	13.5	12.5	N/A
	Medium	14.5	14.5	12.5	N/A
	Long	16.0	16.0	13.5	N/A
70° Swing In/Out	Short	14.0	14.0	13.0	N/A
	Medium	15.0	15.0	13.5	N/A
	Long	16.5	16.5	15.0	N/A
80° Swing In/Out	Short	14.0	14.0	13.5	N/A
	Medium	14.5	14.5	14.5	N/A
	Long	16.5	16.5	15.5	N/A

ALR & ELR Hanger Leg Lengths				
Hanger	Extension	Seat Sling to Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	Short	8	11	10
	Med	10	13	12
	Long	12	15	14

Handrim	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Aluminum Anodized	20,22,24	20,22,24,25,26	20,22,24
Plastic Coat	20,22,24	22,24	20,22,24
Natural Fit - Standard Grip	24	22,24,25,26	22,24
Natural Fit - No Thumb Grip	24	22,24,25,26	22,24
Natural Fit - Super Grip	24	22,24,25,26	22,24
Surge - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Surge LT - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Q Grip Neoprene Coated	N/A	24,25	N/A



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA