



Quickie 2

February 2023

EIQ2X

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓ Option # HCPCS Definition Price ✓ Option # HCPCS Definition Price

1. FRAME

Base model

E1236 coding still applies for EIQ2X with sizes 14" or less.

☐	002WM32	K0005	Quickie 2	\$2,363.00
☐	002WM109	K0005	Quickie 2 - Transit included	\$2,722.00

Transit Option

N/A on Stroller Handles (Removable Push Handles) or Xtender

☐	002TO46	K0108	Transit System WC19	\$359.00
--------------------------	---------	-------	---------------------	----------

Frame Style

☐	002FS31		Standard Swing Away	NC
☐	002FS30		Hemi Swing Away	NC

Rear Frame Type (300 lbs. Weight Capacity)

☐	002RS1		Rear Side Frame Straight	NC
☐	002RS2		Rear Side Frame Rounded	NC
☐	002RS3		Depth Adjustable Back (DAB) 3.0" of rearward adjustment. N/A w/ swing-away or flip-back armrest.	\$218.00

Heavy Duty Package

☐	002HV1		Heavy Duty (18"-24" Wide only) (350 lbs. weight capacity). Includes Reinforced Back Upholstery & Stabilizer Bar.	\$559.00
--------------------------	--------	--	--	----------

2. SIZE / SEATING

Frame Width

☐	002W3		12"	NC
☐	002W4		13"	NC
☐	002W5		14"	NC
☐	002W6		15"	NC
☐	002W7		16"	NC
☐	002W8		17"	NC
☐	002W9		18"	NC
☐	002W10		19"	NC
☐	002W11	E2201	20"	\$559.00
☐	002W12	E2201	21" (Only Avail. w/ Heavy Duty)	\$559.00
☐	002W13	E2201	22" (Only Avail. w/ Heavy Duty)	\$559.00
☐	002W57	E2201	B4ME - 23" Width (Heavy Duty)	\$782.00
☐	002W58	E2202	B4ME - 24" Width (Heavy Duty)	\$782.00

Frame Depth

Seat Depth can be ordered up to 3" smaller for growth. Std frame available with 13"-22" frame Depth & HD frame available with 16"-22" frame Depth. Hemi available with 16"-18" & 20"-22" frame depth. Please refer to matrix in page 8 of Appendix for frame Depth Vs seat Depth matrix. Growth is achieved to the front of the depth.

☐	002FD4		13"	NC
☐	002FD5		14"	NC
☐	002FD6		15"	NC
☐	002FD7		16"	NC
☐	002FD8		17"	NC
☐	002FD9		18"	NC
☐	002FD10		19" Steel reinforced	NC
☐	002FD36		20" Steel reinforced	NC
☐	002FD51		B4ME - 21" Frame Depth	\$670.00
☐	002FD52		B4ME - 22" Frame Depth	\$670.00

Seat Depth

12" seat depth not available with Hemi frame

☐	002SD5		12"	NC
☐	002SD6		13"	NC
☐	002SD7		14"	NC
☐	002SD8		15"	NC
☐	002SD9		16"	NC
☐	002SD10		17"	NC
☐	002SD11		18"	NC
☐	002SD72		19"	\$84.00
☐	002SD13	E2203	20"	\$614.00
☐	002SD105	E2203	B4ME - 21" Depth	\$670.00
☐	002SD106	E2204	B4ME - 22" Depth	\$670.00
☐	002SD107	E2203	B4ME - 21" Depth (Heavy Duty)	\$782.00
☐	002SD108	E2204	B4ME - 22" Depth (Heavy Duty)	\$782.00

Seat Sling				
☐	002SL1		Standard Seat Sling	NC
☐	002SL6		Omit Seat Sling	-\$15.00

JAY Cushions (Optional)				
To order any JAY cushions on-chair, please see the JAY price list.				
Cushion name or part number must be filled out with desired width and depth				
☐	002JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part # _____		
		Width: _____	Depth: _____	
Cushion				
☐	002CU1	K0669	2" Cushion	\$168.00

3. CASTERS / SEAT HEIGHT


Front Seat Height							
Please see matrix on Page 9 for additional information. 13.5"-14" Hemi only							
21.5"-22" not available in Hemi Frame							
☐	002SH40	13.5"	NC	☐	002SH18	18.0"	NC
☐	002SH41	14.0"	NC	☐	002SH20	18.5"	NC
☐	002SH42	14.5"	NC	☐	002SH21	19.0"	NC
☐	002SH12	15.0"	NC	☐	002SH23	19.5"	NC
☐	002SH13	15.5"	NC	☐	002SH24	20.0"	NC
☐	002SH14	16.0"	NC	☐	002SH26	20.5"	NC
☐	002SH15	16.5"	NC	☐	002SH27	21.0"	NC
☐	002SH16	17.0"	NC	☐	002SH29	21.5"	NC
☐	002SH17	17.5"	NC	☐	002SH36	22.0"	NC

Caster Wheels				
See matrix on page 9.				
☐	002CW1		3" Micro	NC
☐	002CW3		4" Low Profile Polyurethane N/A w/ HD package	NC
☐	002CW31		4" Micro Lighted	\$78.00
☐	002CW72		4" x 1" Aluminum Polyurethane	\$123.00
☐	002CW44	E2219	4" x 1.25" Semi Pneumatic	\$61.00
☐	002CW64		4" x 1.5" Aluminum Soft Roll	\$207.00
☐	002CW5		5" Low Profile Polyurethane	NC
☐	002CW51		5" Micro Lighted	\$96.00
☐	002CW6	E2219	5" Soft Roll (Grey) N/A w/ Frog Legs	\$101.00
☐	002CW65		5" x 1" Aluminum Polyurethane	\$123.00
☐	002CW7	E2219	5" x 1.5" Semi-Pneumatic	\$96.00
☐	002CW66		5" x 1.5" Aluminum Soft Roll	\$218.00
☐	002CW34		6" Polyurethane	NC
☐	002CW33	E2214 E2215	6" Pneumatic	\$84.00
☐	002CW35	E2219	6" x 1.5" Semi Pneumatic	\$96.00
☐	002CW68		6" x 1.5" Aluminum Soft Roll	\$229.00
☐	002CW16		7" Polyurethane	NC
☐	002CW32		8" Polyurethane	NC
☐	002CW21	E2214 E2215	8" Pneumatic	\$96.00
☐	002CW39	E2219	8" x 1.5" Semi Pneumatic	\$96.00
☐	002CW36	E2214 E2215	8" x 2" Pneumatic	\$162.00
☐	002CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$252.00

Caster Housing				
Integral bearing forks only available with Integral caster housing				
☐	002CP2		Standard	NC
☐	002CP14		Integral Caster Housing N/A with caster Pin locks	\$84.00

Caster Fork				
Please see matrix on Page 9 for additional information. Frog Legs not available with Integral caster housing				
☐	002CF46		Multi Position Fork-Short	NC
☐	002CF47		Multi Position Fork- Tall	NC
☐	002CF62		3" Fork	NC
☐	002CF12		5 1/4" Fork	NC
☐	002CF17		6 " Fork	NC
☐	002CF23		7" Fork	NC
☐	002CF29		8" x 2" Fork	NC
☐	002CF65		Multiple Position Integral caster housing only	NC
☐	002CF69		Multiple Position Single Sided Fork Integral caster housing only	\$189.00
☐	002CF53	E1015	Frog Legs Suspension (0-100 lbs)	\$486.00
☐	002CF54	E1015	Frog Legs Suspension (101-190 lbs)	\$486.00
☐	002CF55	E1015	Frog Legs Suspension (191-265 lbs)	\$486.00

Caster Options				
☐	002OC1	K0073	Caster Pin Locks (N/A w/8"x2") NA w/ Xtender, Frog Legs or Transit	\$218.00



✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR WHEELS / SEAT HEIGHT									

Rear Seat Height						
☐	002RH11	13.0"	NC	☐	002RH26	17.5" NC
☐	002RH12	13.5"	NC	☐	002RH20	18.0" NC
☐	002RH13	14.0"	NC	☐	002RH27	18.5" NC
☐	002RH14	14.5"	NC	☐	002RH21	19.0" NC
☐	002RH15	15.0"	NC	☐	002RH37	19.5" NC
☐	002RH16	15.5"	NC	☐	002RH22	20.0" NC
☐	002RH17	16.0"	NC	☐	002RH38	20.5" NC
☐	002RH18	16.5"	NC	☐	002RH23	21.0" NC
☐	002RH19	17.0"	NC	☐	002RH39	21.5" NC
				☐	002RH40	22.0" NC

Axle Plate			
☐	002AL102	Standard Adjustable (All in-one)	NC
☐	002AL6	E0959 Amputee Adjustable N/A w/ Padded Armrest	\$196.00

Axle			
☐	002AX1	Quick Release Steel	NC
☐	002AX11	Quad Quick Release 5 Spoke Mag & Mountain Wheel only	\$141.00

One Arm Drive			
N/A w/ Xtender, 12", 23" & 24" wide, amputee axle plates, and small tray table. N/A w/ combination of Padded SA arms and Anti tips. Only available w/ 20", 22", 24" mag wheels. Max weight capacity of 300 lbs. Only available in 6 spoke Mags			

☐	002AL16	E0958	One Arm Drive Adj. Axle Plate				\$950.00	
		☐	002OA1	Small Handrim (N/A w/ 20" wheels) ☐ RH position ☐ LH position				NC
		☐	002OA2	Standard Handrim ☐ RH position ☐ LH position				NC
Select Wheel Size:		☐	20"	☐	22"	☐	24"	

Xtender Power Assist	
RSH of 14.5"-20" for 22" rear wheels and RSH 15.5"-21" for 24" wheels. Front seat height must be greater than rear. N/A w/ Transit Option Max weight capacity of 265 lbs. Available 12"-22" widths.	

Xtender Power Assist				
☐	002XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$8,991.00
☐	002XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$9,885.00
☐	002XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$8,991.00
☐	002XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$9,885.00

Please select one of the following handrim options			
☐	002XH5	Plastic Coated	NC
☐	002XH6	Stainless Steel	NC

Xtender - Other Options			
☐	002XA10	Spare NiMH Power Pack - NICKEL	\$1,559.00
☐	002XA11	Spare Li-ion Power Pack - LITHIUM	\$2,006.00

Xtender Wheel Options			
☐	002XA8	SPARE 22" Lite Spoke Pneu Alum	NC
☐	002XA9	E2211 E2213 SPARE 22" Lite Spoke Pneu Plastic	\$224.00
☐	002XA2	SPARE 24" Lite Spoke Pneu Alum	NC
☐	002XA4	E2211 E2213 SPARE 24" Lite Spoke Pneu Plastic	\$224.00

Rear Wheel Type			
Not applicable w/ Xtender Power Assist Option. Ultra Lightweight Wheels only available w/ Alum Anodized Handrim. Handrims cannot be omitted on Lite Spoke Wheels.* Omit Rear wheel will default chair configuration set up to a 24" wheel			
☐	002RW103	20" Mag Wheel 5 Spoke	NC
☐	002RW104	22" Mag Wheel 5 Spoke	NC
☐	002RW102	24" Mag Wheel 5 Spoke	NC
☐	002RW90	20" Lite Spoke Wheel	\$224.00
☐	002RW91	22" Lite Spoke Wheel	\$224.00
☐	002RW92	24" Lite Spoke Wheel	NC
☐	002RW93	25" Lite Spoke Wheel	NC
☐	002RW94	26" Lite Spoke Wheel	\$224.00
☐	002RW77	22" Ultra Lightweight Spoke Wheel	\$279.00
☐	002RW78	24" Ultra Lightweight Spoke Wheel	\$279.00
☐	002RW28	22" Spinergy SPOX Black	\$944.00
☐	002RW95	24" Spinergy SPOX Black	\$944.00
☐	002RW96	25" Spinergy SPOX Black	\$944.00
☐	002RW97	26" Spinergy SPOX Black	\$944.00
☐	002RW32	22" Spinergy LX Black	\$1,185.00
☐	002RW98	24" Spinergy LX Black	\$1,185.00
☐	002RW99	25" Spinergy LX Black	\$1,185.00
☐	002RW40	26" Spinergy LX Black	\$1,185.00
☐	002RW50	24" Mountain	\$279.00
☐	002RW25	Omit Rear Wheel and Axle *	-\$55.00
☐	002RW49	Omit Rear Wheel (Use w/24" Xtender)	NC
☐	002RW111	Omit Rear Wheel (Use w/22" Xtender)	NC

Rear Wheel Tire			
☐	002RT1	Pneumatic	NC
☐	002RT2	E2211 E2213 Pneumatic Airless Insert N/A with Xtender.	\$96.00
☐	002RT24	E2211 E2213 Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$109.00
☐	002RT3	Full Poly	NC
☐	002RT38	E2218 SHOX Performance Poly	\$106.00
☐	002RT25	E2218 Primo Orion Poly	\$121.00
☐	002RT39	E2211 E2212 High-Pressure Clincher Black	\$106.00
☐	002RT27	E2211 E2212 Knobby V-TRAK	\$106.00
☐	002RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$157.00
☐	002RT12	E2211 E2212 Mountain Tire <i>Only available on 24" Mountain Rear Wheel</i>	\$106.00


✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR WHEELS / SEAT HEIGHT (CONT)									

Rear Seat Height Availability Matrices

Rear tire, side frame, and axle plate selections dictate seat height availability. *Not all combinations are possible within acceptable seat height ranges

Rear Tire	Min RSH	Max RSH
20" Pneumatic, Poly	13"	19.0"
22" Pneumatic, Poly	14.0"	20.0"
22" HP Clincher Blk	13.5"	19.5"
22" Schwalbe	13.5"	19.5"
24" Pneumatic, Poly	15"	21.0"
24" HP Clincher Blk	14.5"	20.5"
24" Knobby	15"	21.0"
24" Schwalbe	14.5"	20.5"
24" SHOX Poly	14.5"	20.5"
24" Mountain	15.0"	21.0"
25" HP Clincher Blk	15.0"	21.0"
25" Schwalbe	15.0"	21.0"
25" SHOX Poly	14.5"	20.5"
26" Pneumatic	16.0"	22.0"
26" HP Clincher Blk	15.5"	21.5"
26" Schwalbe	15.5"	21.5"

Rear Side Frame	Min RSH	Max RSH
Curved	13"	21.0"
Straight	13	22.0"
DAB	13	22.0"



Axle Plate	Min RSH	Max RSH
Std Adj (All in one)	13"	22.0"
Amputee	14"	22.0"

Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Tire					
Pneumatic	20,22,24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		24,26
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Pressure Clinch Blk			22,24,26	22,24	22,24,26
Knobby V-TRAK			24		24
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Perform. Poly	20,22,24		24,25		24,25
20"= 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590					

Handrim

Natural Fits available in 22", 24", and 25". Handrim selection N/A w/ Xtender. Ultra lightweight wheels only available w/ alum anodized handrim. Plastic coated limited to certain rims and sizes. Natural Fit w/ 5 spoke MAG only available in 24"

☐	002HR1	Aluminum Anodized	NC
☐	002HR4	Plastic Coated - N/A on every rim	\$123.00
☐	002HR22	Natural Fit - Standard Grip	\$526.00
☐	002HR25	Natural Fit - No Thumb Grip	\$455.00
☐	002HR23	Natural Fit - Super Grip	\$585.00
☐	002HR30	Natural Fit LT- Standard Grip	\$526.00
☐	002HR31	Natural Fit LT- Super Grip	\$585.00
☐	002HR32	Natural Fit LT - No Thumb Grip	\$455.00
☐	002HR26	The Surge - Oval w/Gription Strip	\$472.00
☐	002HR27	The Surge LT - Oval w/Gription Strip	\$508.00
☐	002HR28	Q-Grip Neoprene Coated	\$414.00
☐	002HR24	Spinergy Flex Rim Available on 24" & 25" Spinergy LX wheel	\$782.00
☐	002HR7	Omit Handrim Option may not be available on all wheels	NC
☐	002HR8	8 Oblique N/A w/ Xtender. Only available with 24" or 26" 6 spoke Mag or Lite Spoke rims.	\$363.00

Camber

☐	002CI1	0°	NC
☐	00217CI	2° (default)	NC
☐	002CI3	4°	NC

Center of Gravity Preset

Amputee Axle Plate will come preset at 1.0" Center of Gravity (CG), range of adjustment is -3.0" to -1.75" and 1.0" to 3.5" CG.

☐	002CC2	1.0" (default)	NC
☐	002CC3	1.5"	NC
☐	002CC4	2.0"	NC
☐	002CC5	2.5"	NC
☐	002CC6	3.0"	NC
☐	002CC7	3.5"	NC

One Arm Drive CG

☐	002CC30	1AD Forward position	NC
☐	002CC31	1AD Reward position	NC
☐	002CC32	1AD Center position	NC

5. WHEEL LOCKS

Wheel Lock

Ergo & compact N/A w/ 13"-14" frame depth or amputee axle plate. Additional No Gos may exist depending on wheel size and position

☐	002WL46	High Mount Push	NC
☐	002WL50	High Mount Pull	NC
☐	002WL75	Ergo Scissor	\$123.00
☐	002WL64	Compact	\$157.00
☐	002WL24	Unilateral Wheel Lock Right	\$335.00
☐	002WL23	Unilateral Wheel Lock Left	\$335.00
☐	002WL10	Omit	NC

Wheel Lock Option

☐	002WL12	Do Not Mount	NC
--------------------------	---------	--------------	----

Extension Handle

N/A w/ ergo scissor or compact wheel locks

☐	002EH7	E0961	Extension Handle (Pair)	\$84.00
--------------------------	--------	-------	-------------------------	---------

Grade Aid

Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options

☐	002GA1	E0974	Grade Aid	\$168.00
--------------------------	--------	-------	-----------	----------

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
6. HANGER & FOOTREST									

Hangers

Front mount hangers only available w/ aluminum angle adj footplates.
ALR/ELR Hangers N/A w/ Push To Lock wheel locks.

☐	002H145		60° Swing In/Out Hanger	NC
☐	002H139		70° Swing In/Out Hanger	NC
☐	002H141		80° Swing In/Out Hanger	NC
☐	002H154		60° Swing In/Out Hanger Hemi	NC
☐	002H155		70° Swing In/Out Hanger Hemi	NC
☐	002H140		70° Swing In/Out Front Mount	NC
☐	002H142		80° Swing In/Out Front Mount	NC
☐	002H158		90° Swing In/Out Front Mount	NC
☐	002H159	K0053	Articulating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$413.00
☐	002H160	E0990	Elevating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$352.00
☐	002H161	E0990	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$352.00
☐	002H32		Omit Hanger & Footplate	-\$15.00

Articulating / Elevating Hanger Cover

☐	002DI47		Neoprene Cover	\$380.00
☐	002DI48		Gel Pad Cover	\$447.00

Footplates

Heel loop or leg strap not included, See Footplate Options Section

☐	002F1		Composite (13"-22"w)	NC
☐	002F271	K0040	Composite Angle Adj-(16"-22"w)	\$202.00
☐	002F3	K0040	Alum Angle Adj - Ext Mt (12"-22"w)	\$212.00
☐	002F8	K0040	Alum Angle Adj - Front Mt(12"-22"w)	\$212.00
☐	002F5		Alum Footplate - Fixed (14"-22"w) Available w/ Ext Mount ELR or ALR	NC
☐	002F48	K0040	Platform Flip Up Angle Adj (12"-22"w) N/A w/ Hemi Hangers,ELR/ALR,Front Mnt	\$212.00
☐	002F7		Omit Footplate & Ext Tube	-\$15.00

Footplate Options (Ext Tube & Front Mount)

☐	002HL1	E0951	Heel Loop	\$45.00
☐	002HL5	K0038	Leg Strap	\$45.00

Footrest Extension Tubes

Required w/ extension tube mounted hangers. See leg length matrix

☐	002E61		Short Hemi	NC
☐	002E1		Short	NC
☐	002E2		Medium	NC
☐	002E5		Long	NC
☐	002E16		High Mount N/A w/ one-arm drive option or ELR/ALR or 90° front mount. Only available w/angle adj footplates	\$84.00
☐	002E23		2" - Front Mount Only N/A w/ 6", 7" or 8" casters with 90°hangers, or ELR's	\$61.00
☐	002E24		4" - Front Mount Only N/A w/6", 7" or 8" casters w/ 90°hangers or ELR/ALR	\$61.00

Lower Leg Length Setting (Extension Mount Only)

Recommend 2" less than FSTFH. (*)denotes frame, 9" & 9.5" Hemi Only

☐	002LO01	9.0"(*)	NC	☐	002LO13	15.0"	NC
☐	002LO02	9.5"(*)	NC	☐	002LO14	15.5"	NC
☐	002LO03	10.0"	NC	☐	002LO15	16.0"	NC
☐	002LO04	10.5"	NC	☐	002LO16	16.5"	NC
☐	002LO05	11.0"	NC	☐	002LO17	17.0"	NC
☐	002LO06	11.5"	NC	☐	002LO18	17.5"	NC
☐	002LO07	12.0"	NC	☐	002LO19	18.0"	NC
☐	002LO08	12.5"	NC	☐	002LO20	18.5"	NC
☐	002LO09	13.0"	NC	☐	002LO21	19.0"	NC
☐	002LO10	13.5"	NC	☐	002LO22	19.5"	NC
☐	002LO11	14.0"	NC	☐	002LO23	20.0"	NC
☐	002LO12	14.5"	NC	☐	002LO24	20.5"Std. Only	NC

Z-finity Hangers/Footrest System

(Maximum weight capacity for Z-Finity footrests is 165 lbs)

Some interference may exist in certain configurations. N/A with Heavy Duty

☐	002H162		Z-finity Hanger	\$391.00
--------------------------	---------	--	-----------------	----------

Z-finity Footrest Mount

☐	002DI22		Direct Hanger Mt/Fixed Frame - High	NC
☐	002DI23		Direct Hanger Mt/Fixed Frame - Low	NC
☐	002DI24		Extension Tube Mount	NC

Z-finity Footrest Pivot Arm

☐	002DI20		Short 2"	NC
☐	002DI21		Long 3 1/4"	NC

Z-finity Footplate

Some interference may exist in certain configurations.

☐	002F320	K0040	Multi Angle Adjustable (Pair) (not avail. w/8"-9" widths)	\$207.00
☐	002F321	K0040	Angle Adjustable Platform Footboard	\$239.00
☐	002F6		Omit Footplate	NC

Z-finity ExtensionTube Length (When Choosing Ext. Tube Mount)

☐	002E62		Hemi Short	NC
☐	002E1		Short	NC
☐	002E2		Medium	NC
☐	002E5		Long	NC

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING									

Back Type			
☐	002BT10	Standard	NC
☐	002BT20	Angle Adjustable Back (3-5" Pivot)	\$207.00
☐	002BT22	Angle Adjustable Back (5"-7" Pivot)	\$207.00
☐	002BT6	Depth Adjustable (Required with DAB)	NC

Back Type - Build 4 Me Options			
Must choose DAB Back type with B4ME options			
☐	002BT127	E2398 B4ME - Dynamic Back	\$1,671.00
		☐ Fixed Back post	NC
		☐ Std (19") ☐ Med (22") ☐ Tall (25")	
		☐ Height Adjustable Stroller handles	NC
☐	002BT128	E1226 B4ME - Reclining Backrest	\$1,224.00
		☐ Std (18"-23") ☐ Tall (19"-26")	
		☐ Med (19") ☐ Tall (22")	

Back Post Options			
Omit push handles N/A w/ stabilizer bar. Folding HA push handle N/A w/8 deg or Low back height.			

Standard Back type options:			
☐	002BP1	Push Handle Straight	NC
☐	002BP3	Push Handle 8° Med or Tall only. N/A w/ Angle Adjustable Backrest	\$40.00
☐	002BP2	Omit Push Handle Straight (N/A w/HD)	NC
☐	002BP4	Omit Push Handle 8° (N/A w/HD) Med or Tall only	\$40.00
☐	002PH17	Folding Push Handle with Ht. Adj. Back N/A w/Angle Adj Back & Heavy D.	\$252.00

Angle Adj Back type options:			
☐	002BP5	Angle Adj. Straight Push Handles	NC
☐	002BP6	Angle Adj. Omit Push Handles, N/A w/HD	NC

Back Height (Std or Angle Adjustable*)			
*For Angle Adj Back see the angle Adj Backrest chart in the right column of the page.			
☐	002B1	Low (9"-12" Std Back height range)	NC
☐	002B5	Medium (12"-15" Std Back height range)	NC
☐	002B9	Tall (16"-19" Std Back height range)	NC

Back Height Setting (Std or Angle Adjustable)					
EXO upholstery not available with 21" back height					
☐	002UX4	9"	NC	☐	002UX18 16" NC
☐	002UX6	10"	NC	☐	002UX20 17" NC
☐	002UX8	11"	NC	☐	002UX22 18" NC
☐	002UX10	12"	NC	☐	002UX24 19" NC
☐	002UX12	13"	NC	☐	002UX26 20" (*) NC
☐	002UX14	14"	NC	☐	002UX28 21" (*) NC
☐	002UX16	15"	NC	(*)20" & 21" only w/ Ang Adj Back	

Stroller Handles (Removable Push Handles)			
Only available w/ omit push handles. N/A w/ DAB, amputee axle plate, push handles, or when angle adjustable back is selected, flip-back or flip-back height adjustable armrests, and JAY backs.			
☐	002SR2	Regular Stroller Handles	\$151.00
☐	002SR3	Tall Stroller Handles	\$151.00

Stabilizer Bar			
N/A with Folding Push Handle, or with Ht. Adj. Back or DAB Stroller handles			
☐	002SZ1	Stabilizer Bar Standard with HeavyDuty Package and Recline	\$118.00


Backrest Upholstery (Standard or Angle Adjustable)			
Std EXO upholstery and Tension adj binding is always black. EXO Tension Adj Uph N/A with Heavy Duty. Max back height 20"for EXO.			
☐	002BU1	Standard Upholstery	NC
☐	002BU88	Standard Upholstery EXO	\$56.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue \$15.00			
☐	002BU70	E2611 Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$352.00
☐	002BU89	E2611 Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$370.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue \$15.00			
☐	002BU71	E2611 Backrest Tension Adjust 3DX Vented *Reference 186BU71 for billing	\$370.00
☐	002BU82	E2611 Jay GO Backrest N/A on Heavy Duty	\$455.00
☐	002BU14	Omit Back Upholstery	-\$15.00

Standard Backrest Option			
☐	002OD4	Reinforced Upholstery Only avail w/ Std Back Upholstery. Std w/HD package	\$56.00

Angle Adjustable Backrest Pivot Height			
Standard push handles only. Some interference w/ Padded SA arms. N/A w/ flip-back or flip-back height adjustable armrests (classic).			
☐	002AH1	Angle Adjustable Pivot 3"	NC
☐	002AH2	Angle Adjustable Pivot 4"	NC
☐	002AH3	Angle Adjustable Pivot 5"	NC
☐	002AH4	Angle Adjustable Pivot 6"	NC
☐	002AH5	Angle Adjustable Pivot 7"	NC

Preset Angle Adjustable Backrest (-15° to +3°) (Optional)			
☐	002BD5	Back Angle - 15°	NC
☐	002BD4	Back Angle - 12°	NC
☐	002BD3	Back Angle - 6°	NC
☐	002BD2	Back Angle - 0°	NC
☐	002BD1	Back Angle +3°	NC

Angle Adjustable Backrest Chart			
Back Tube Height	Low	Medium	Tall
Pivot Regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot Regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot Regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot Tall set at 7"	16" 17" 18"	19" 20" 21"	N/A
Adjustable Pivot Point Measurement from Seat Sling			
If back height option is not chosen, middle option will be present			



✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING (cont)									

DAB Back Height				
☐	002B66		Depth Adjustable 12"	NC
☐	002B67		Depth Adjustable 14"	NC
☐	002B68		Depth Adjustable 16"	NC
☐	002B69		Depth Adjustable 18"	NC
☐	002B70		Depth Adjustable 20"	NC

DAB Back Post Options				
<i>Diagram of DAB in Appendix on page 9. SA Stroller Handle height range is 21.25"-28.5"</i>				
☐	002BP33		DAB - Swing Away Stroller Handles	\$111.00
☐	002BP27		DAB - Post Standard (19")	NC
☐	002BP28		DAB - Post Medium (22")	NC
☐	002BP29		DAB - Post Tall (25")	NC

DAB Backrest options				
<i>DAB Solid Back NA w/23" & 24" frame widths</i>				
☐	002BU28		Omit DAB Solid Back and Hardware	NC
☐	002BU29		Omit DAB Solid Back Include Hardware	\$90.00
☐	002BU23		Depth Adjustable Solid Back	\$218.00

8. ARMREST / SIDE GUARDS

Padded Swing Away				
<i>N/A w/ DAB, Amputee axle. Some limitations may exist with integral push handles, some RSTFHs and wheel sizes. **Now available w/Transit.</i>				
☐	002AA18	E0973	Padded Swing-Away STD (9" to 13") Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"	\$252.00

Flip Back				
<i>N/A w/ DAB, Xtender, 12" - 14" frame depth, Classic Style armpad in highest/lowest axle RSH on rounded rear frame, or angle adjustable back.4 Point (FS) J3 mounting bracket will cause interference with Flip Back Height Adj armrest. Flip back N/A with amputee axle plate for seat heights requiring the highest axle plate position</i>				
☐	002AA6	E0973	Flip Back Height Adj(Classic 8.5"-13.5")	\$352.00
☐	002AA5		Flip Back (incl. Desk Length Pad)	NC
☐	002AA65	E0973	Flip Back Height Adj (LITE 8.5"-13.5")	\$352.00

Flip Up (For use with Depth Adjustable Back)				
<i>Requires an Armpad or Foam from the Arm Pad section below</i>				
☐	002AA7	E0973	Length Adjustable Locking Flip-up	\$308.00
☐	002AA1	E0973	Adj. Locking Flip-up (Std. w/ foam)	\$296.00

9. KEY OPTIONS & ACCESSORIES

Anti-Tip				
☐	002AT1	E0971	Rear Anti-Tip (NC w/ Xtender)	\$168.00

Backpack				
☐	002A1		Quickie Backpack Black	\$84.00
☐	002A2		Quickie Kids Backpack Black	\$67.00

Seat Pouch				
☐	0021A01		Seat Pouch Black /Charcoal	\$56.00

JAY Positioning				
☐	002K1	E1020	JAY Lower Extremity S/A support ☐ Left ☐ Right Position	\$441.00
☐	002K2		Omit Hanger ☐ Left ☐ Right	-\$10.00
☐	002K7		Pad Size ☐ 6" x 8" ☐ 8" x 10"	NC

Military Service Patches (Optional)				
☐	002VU2		Army Patch	\$111.00
☐	002VU4		Air Force Patch	\$111.00
☐	002VU12		Coast Guard Patch	\$111.00
☐	002VU1		Marines Patch	\$111.00
☐	002VU3		Navy Patch	\$111.00
☐	002VU5		Purple Heart Patch	\$111.00
☐	002VU9		Purple Heart with Army Patch	\$173.00
☐	002VU11		Purple Heart with Air Force Patch	\$173.00
☐	002VU13		Purple Heart with Coast Guard Patch	\$173.00
☐	002VU8		Purple Heart with Marines Patch	\$173.00
☐	002VU10		Purple Heart with Navy Patch	\$173.00

JAY Back (Optional)				
☐	002JB01		J3 Back (See J3 Order Form)	Price List

Height Adjustable				
☐	002AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$319.00
☐	002AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$319.00

Omit Armrest				
☐	002AA14		Omit Armrest	NC

Armpad for Height Adjustable				
<i>Ergo armpad only available with Single post Height adj & Flip back armrest</i>				
<i>Ergo not available with Classic Flip Back armrest</i>				
☐	002AM23		Desk Length - Classic (10")	NC
☐	002AM24		Full Length - Classic (14")	NC
☐	002AM9		Armpad Foam - Short (Flip Up armrests)	NC
☐	002AM10		Armpad Foam - Long (Flip Up armrests)	NC
☐	002AM27		Ergo Pad (9.5"-13.5")	NC

Side Guard				
<i>Plastic side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist.</i>				
☐	002SG33	K0108	Plastic Standard (8.5" High)	\$162.00
☐	002SG29	K0108	Plastic Low (6.5" High)	\$162.00

Additional Cover & Pack Accessories				
☐	0028A02		Gel Cover Hanger/Frame	\$123.00
☐	0028A01		Gel Cover Swing Away Arm	\$168.00

Positioning Belts				
☐	002A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$78.00
☐	002A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$123.00
☐	002A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$133.00
☐	002A120		Positioning Belt Slide Clamps	\$56.00

Spoke Guards				
<i>Only available w/Mountain, and Lite Spoke wheels,</i>				
<i>Not available with 5 spoke MAG wheels or Projection handrims</i>				
☐	0024A02	K0065	Quickie Clear Spoke Guards	\$128.00
☐	0024A01	K0065	Quickie Black Spoke Guards	\$128.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
9. KEY OPTIONS & ACCESSORIES (cont)									
Tray Tables					Miscellaneous				
☐	002A30	E0950	Small (13"-14"W)	\$285.00	<i>O2 holder N/A with Xtender or 12" - 13" frame width. Crutch holder N/A with EXO Tension Adj upholstery and Omit push handle</i>				
☐	002A31	E0950	Medium (15"-17"W)	\$285.00	☐	0028A03		Caddy (N/A w/ front mount hangers)	\$56.00
☐	002A32	E0950	Large (18"-20"W)	\$285.00	☐	002A2100	E2207	Crutch Holder	\$218.00
☐	002A33	E0950	X-Large (21"-22"W)	\$285.00	☐	002LS14	K0038	Extra Leg Strap Black	\$40.00
					☐	002IG2		Impact Guard Neoprene	\$101.00
					☐	002A73	K0105	IV Holder	\$308.00
					☐	002A145	E2208	O2 Holder	\$347.00
					☐	0028A21		Tool Kit	\$61.00
					☐	002A999		Quickie Water Bottle & Holder	\$61.00

10. COLOR

Standard Frame Colors				
☐	002C140		Aztec Gold	NC
☐	002C2		Black <i>(Not avail as base w/Zebra)</i>	NC
☐	002C141		Black Cherry	NC
☐	002C10		Black Opal	NC
☐	002C38		Blue Opal	NC
☐	002C36		Candy Blue	NC
☐	002C14		Candy Purple	NC
☐	002C35		Candy Red	NC
☐	002C142		Electric Blue	NC
☐	002C143		Evergreen	NC
☐	002C44		Glossy White	NC
☐	002C126		Glow	NC
☐	002C144		Green Apple	NC
☐	002C147		Hot Sparkle Pink	NC
☐	002C148		Mauve Pink	NC
☐	002C103		Rootbeer	NC
☐	002C114		Sparkle Silver	NC
☐	002C145		Sunrise Orange	NC
☐	002C113		Titanium Color Paint	NC
☐	002C8		Yellow	NC

Personalized Frame Color				
☐	002C112		Desert Camouflage	\$441.00
☐	002C43		Mossy Oak Camouflage	\$441.00
☐	002C42		Stars & Stripes	\$441.00
☐	002C108		Zebra <i>(Choose above base color) Std with White unless above base paint is chosen N/A w/ Black</i>	\$441.00

Anodized Color Package				
Standard Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.				
☐	002EB07		Black - Caster only	\$56.00
☐	002EB08		Blue	\$111.00
☐	002EB13		Gold	\$111.00
☐	002EB09		Lime Green	\$111.00
☐	002EB10		Purple	\$111.00
☐	002EB11		Red	\$111.00
☐	002EB12		Ti Grey	\$111.00

Additional Options Through Built for Me				
To inquire about any unique personalized Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011				

Matte Frame Colors				
☐	002C116		Matte Black	NC
☐	002C151		Matte Black Cherry	NC
☐	002C150		Matte Electric Blue	NC
☐	002C149		Matte Evergreen	NC
☐	002C152		Matte Purple	NC



Commonly Requested Options Include:

- ✓ Personalized Anodizing
- ✓ Xtender Battery on right side
- ✓ Dual Mount Wheel Locks - Attendant Foot Lock
- ✓ Hub Wheel Locks on Mag or Spoke Wheels
- ✓ Colored handrims




11. APPENDIX

Frame Depth matrix		Seat Depth											
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
Frame Depth	13"	S	S										
	14"	S	S	S									
	15"	S	S	S	S								
	16"		S	S	S	S							
	17"			Hemi*	S	S	S	S					
	18"				S	S	S	S					
	19"					S	S	S	S				
	20"						S	S	S	S			
	21"							S	S	S	S		
	22"								S	S	S	S	

"S" = Standard Frame Style

"H" = Hemi Frame Style

*Only (6) Frame Depths avail. on Hemi Style Frame



Lower Leg Length Matrix

2" & 4" extensions for front mount hangers shipped separate, not mounted. Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers. All values shown apply to Standard frame style, subtract 1.5" to get Hemi equivalent. More than 2" will require on field adjustment.

Hanger	Extension	Seat Sling to Top of Footrest in inches		Min Allowable FSH
		Min	Max	
60° Hemi Swing In/Out	Short Hemi	10	13	14
	Short	-	-	-
	Med	12.5	15.5	15.5
	Long	14.5	17.5	17.5
70° Hemi Swing In/Out	Short Hemi	11	14	14
	Short	-	-	-
	Med	13.5	16.5	16
	Long	15.5	18.5	18
60° Swing In/Out	Short	11.5	15	15
	Med	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	12	16.5	15
	Med	12.5	18.5	15.5
	Long	14.5	20.5	17.5
80° Swing In/Out	Short	13.5	17	15.5
	Med	13.5	19	15.5
	Long	15.5	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
80° Swing In/Out w/Platform Flip Up	Short	13.5	14.5	15.5
	Long	15	19	17
70° Front Mount	No Ext	6	11.5	14
	2"	6	14	16
	4"	6	16	18
80° Front Mount	No Ext	6	12	14
	2"	6	14.5	16.5
	4"	6	16.5	18.5
90° Front Mount	No Ext	5.5	11.5	14
	Short	15	16.5	17
	Med	15	18.5	17
ALR	Long	16.5	20.5	18.5
	Short	12	16.5	15.5
	Med	12	18.5	15.5
ELR Ext. Mount	Long	13.5	20.5	17
	No Ext	5.5	9.5	13

Bckrest Height Matrix

Back Type	Ang-Adj Pivot Ht Setting	Push Handles	Backrest Height	Avail. Ht. Range
Standard - (Fixed)	N/A	Integral	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"
		Omit	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"
Angle Adjustable (3"-5" Pivot)	3"	Integral	Low	12"-14"
			Medium	15"-17"
			Tall	17"-19"
	4"	Integral	Low	13"-15"
			Medium	16"-18"
			Tall	18"-20"
	5"	Integral	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	3"	Omit	Low	12"-14"
			Medium	15"-17"
			Tall	17"-19"
	4"	Omit	Low	13"-15"
			Medium	16"-18"
			Tall	18"-20"
	5"	Omit	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
Angle Adjustable (5"-7" Pivot)	5"	Integral	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	6"	Integral	Low	15"-17"
			Medium	18"-20"
			Tall	N/A
	7"	Integral	Low	16"-18"
			Medium	18"-20"
			Tall	N/A
	5"	Omit	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	6"	Omit	Low	15"-17"
			Medium	18"-20"
			Tall	N/A
	7"	Omit	Low	16"-18"
			Medium	18"-20"
			Tall	N/A
Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	N/A	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"

* Values Given are Estimates. Actual seat to floor height may vary based on seat depth and front seat to floor height selection.

All 3" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0
	1 1/2"	-	-	-	18.0	18.5
Multi Position Tall	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
3"	Std	-	-	-	-	15.0
	1/2"	-	-	-	-	15.5
	1"	-	-	-	-	16.0
	1 1/2"	-	-	-	-	16.5
5 1/4"	Std	-	-	-	-	17.5
	1/2"	-	-	-	-	18.0
	1"	-	-	-	-	18.5
	1 1/2"	-	-	-	-	19.0
6"	Std	-	-	-	-	18.0
	1/2"	-	-	-	-	18.5
	1"	-	-	-	-	19.0
	1 1/2"	-	-	-	-	19.5
7"	Std	-	-	-	-	19.0
	1/2"	-	-	-	-	19.5
	1"	-	-	-	-	20.0
	1 1/2"	-	-	-	-	20.5
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	16.5	17.0
	1 1/2"	-	-	-	18.0	18.5
Integral Multi Pos Forks (Low)	Std	-	-	-	15.5	16.0
	1/2"	-	-	-	16.0	16.5
	1"	-	-	-	16.5	17.0
Integral Multi Pos Forks (Std)	Std	-	-	-	16.5	17.0
	1"	-	-	-	17.0	17.5

All 6" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	19.5	20.0
Multi Position Tall	Std	18.5	19.0	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
	1 1/2"	-	-	-	21.5	22.0
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	18.0	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
6"	Std	-	18.0	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	20.5	21.0
7"	Std	-	19.0	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
	1 1/2"	-	-	-	21.5	22.0
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	18.5
	1 1/2"	-	-	-	-	20.0
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5

All 4" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	16.5	17.0	17.5
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	18.5	19.0
Multi Position Tall	Std	-	-	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	20.5	21.0
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	17.0	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.0	19.5
6"	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	19.5	20.0
7"	Std	-	-	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	20.5	21.0
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	16.5	17.0	17.5
	1 1/2"	-	-	18.0	18.5	19.0
Integral Multi Pos Forks (Low)	Std	-	-	-	16.0	16.5
	1/2"	-	-	-	16.5	17.0
	1"	-	-	-	17.0	17.5
Integral Multi Pos Forks (Std)	Std	-	-	16.5	17.0	17.5
	1"	-	-	-	17.5	18.0

All 5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	16.5	17.0	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
	1 1/2"	-	-	-	19.0	19.5
Multi Position Tall	Std	-	18.5	19.0	19.5	20.0
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
	1 1/2"	-	-	-	19.5	20.0
6"	Std	-	-	18.0	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
7"	Std	-	-	19.0	19.5	20.0
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	17.5	18.0
	1 1/2"	-	-	-	19.0	19.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	17.0	17.5	18.0
	1"	-	-	-	18.0	18.5

All 7" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	19.0
	1/2"	-	-	-	-	19.5
	1"	-	-	-	-	20.0
	1 1/2"	-	-	-	-	20.5
Multi Position Tall	Std	19.0	19.5	20.0	20.5	21.0
	1/2"	-	-	-	21.0	21.5
	1"	-	-	-	21.5	22.0
	1 1/2"	-	-	-	22.0	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
	1 1/2"	-	-	-	20.5	21.0
6"	Std	-	-	19.0	19.5	20.0
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
7"	Std	-	19.5	20.0	20.5	21.0
	1/2"	-	-	-	21.0	21.5
	1"	-	-	-	21.5	22.0
	1 1/2"	-	-	-	22.0	-
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	19.0
	1"	-	-	-	-	19.5

All 8" x 1.5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	-	20.0
	1"	-	-	-	-	20.5
	1 1/2"	-	-	-	-	21.0
6"	Std	-	-	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
	1 1/2"	-	-	-	21.5	22.0
7"	Std	-	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
	1 1/2"	-	-	-	-	-
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

11. APPENDIX

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"

* Values Given are Estimates. Actual seat to floor height may vary based on seat depth and front seat to floor height selection.

All 8" x 2" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
6"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
7"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
8" x 2"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
Frog Legs	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

All Other 8" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Multi Position Tall	Std	19.5	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	20.0
	1/2"	-	-	-	-	20.5
	1"	-	-	-	-	21.0
6"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
7"	Std	-	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA