



Quickie 5^R

February 2023

EIR10

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ ST: _____ ZIP: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ ST: _____ ZIP: _____
Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base Model

☐	129WM86	K0005	Quickie 5 ^R (weight capacity 265 lbs) Use E1235 for sizes 14" or less	\$2,648.00
--------------------------	---------	-------	---	------------

Frame Angle (see Figure 1)

☐	129FA1	75°	NC
☐	129FA8	80°	NC
☐	129FA2	85°	NC

Frame Inset Per Side (See Figure 2)

☐	129FF1	0"	NC
☐	129FF2	1" (N/A w/ 12" wide)	NC
☐	129FF3	2" (N/A w/ 12", 13", or 14" wide)	NC

Frame Inset

Footrest Width

0

Subtract 0 from Frame Width

1"

Subtract 2" from Frame Width

2"

Subtract 4" from Frame Width

Transit

Transit not available w/ XTR suspension & Xtender.

☐	129TO45	K0108	Transit System	\$359.00
--------------------------	---------	-------	----------------	----------

Frame Style (See Figure 3)

XTR Suspension N/A with Xtender

☐	129FS34	Classic Frame	NC
☐	129FS35	Performance Frame	NC
☐	129FS36	XTR Suspension Frame Only available in Performance Frame Set Up	\$1,000.00
☐	129SK1	Light Shock (up to 180 lbs.)	NC
☐	129SK3	Heavy Shock (up to 265 lbs.)	NC

XTR available squeeze and RSH

Seat Depth	Squeeze (FSH - RSH)
14"	2.0" or 3.0"
15"	2.0" or 3.5"
16"	2.5" or 3.5"
17"	2.5" or 4.0"
18"	2.5" or 4.0"
19"	2.5" or 4.5"
20"	3.0" or 4.5"

XTR Tire - Available RSH	
Rear Tire	RSH
20" SHOX	13.5" thru 15.0"
20" all others	14.0" thru 15.5"
22" SHOX	14.5" thru 16.0"
22" Pneumatic	15.5" thru 16.5"
22" all others	15.0" thru 16.0"
24" SHOX	15.5" thru 16.5"
24" all others	16.0" thru 17.0"
25" SHOX	16.0" thru 17.0"
25" Marathon	16.0" thru 17.5"
26" all tires	17.0" thru 18.0"

Figure 1




Figure 2





Figure 3



Classic Frame - Traditional seat sling to front frame bend angle dimensions.

Performance Frame - Seat Sling to front frame bend angle is 1.5" shorter than the classic position.


Front of back cane post to front edge of seat sling will remain the same for both Classic and Performance.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
2. SIZE / SEATING									

Depth Adjustable Backrest (Optional)

See diagram Depth Adjust Backrest extends the frame to the rear of the seat sling, select a factory setting. N/A on 19" and 20" Seat Depth or XTR frame.

☐	129DP02	0 Setting	Backrest at standard position, field adjustable to -1" or -2"	\$184.00
☐	129DP03	-1 Setting	Backrest -1" offset, field adjustable to flush or -2"	\$184.00
☐	129DP04	-2 Setting	Backrest -2" offset, field adjustable to flush or -1"	\$184.00



Seat Width

See Figure 2. Measure from outside of seat tube to outside of opposite seat tube. XTR and J3 Back N/A w/ 12" and 13" widths

☐	129S3	12"	NC
☐	129S4	13"	NC
☐	129S5	14"	NC
☐	129S6	15"	NC
☐	129S7	16"	NC
☐	129S8	17"	NC
☐	129S9	18"	NC
☐	129S10	19"	NC
☐	129S48	E2201	20" \$559.00

Seat Depth

See Figure .1 Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered. XTR frame N/A with 12" and 13" depth

☐	129SD5	12"	NC
☐	129SD6	13"	NC
☐	129SD7	14"	NC
☐	129SD8	15"	NC
☐	129SD9	16"	NC
☐	129SD10	17"	NC
☐	129SD11	18"	NC
☐	129SD12	19"	NC
☐	129SD13	E2203	20" \$614.00

Seat Sling

☐	129SL1	Standard Seat Sling	NC
☐	129SL6	Omit Seat Sling	-\$15.00
☐			

Cushion

☐	129CU1	K0669	2" Cushion	\$168.00
--------------------------	--------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	129JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width:	Depth:

3. CASTERS / SEAT HEIGHT

Front Seat Height

See Figures 4 and 5

☐	129SH14	16.0"	NC	☐	129SH20	18.5"	NC
☐	129SH15	16.5"	NC	☐	129SH21	19.0"	NC
☐	129SH16	17.0"	NC	☐	129SH23	19.5"	NC
☐	129SH17	17.5"	NC	☐	129SH24	20.0"	NC
☐	129SH18	18.0"	NC	☐	129SH26	20.5"	NC
☐				☐	129SH27	21.0"	NC

Front to Rear Seat Height Differences

Frame Depth	Max Seat Sqz
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Caster Wheels

☐	129CW1	3" Micro	NC
☐	129CW3	4" Low Profile Polyurethane	NC
☐	129CW31	4" Micro Lighted	\$78.00
☐	129CW72	4" x 1" Aluminum Polyurethane	\$123.00
☐	129CW44	E2219 4" x 1.25" Semi Pneumatic	\$61.00
☐	129CW64	4" x 1.5" Aluminum Soft Roll	\$207.00
☐	129CW5	5" Low Profile Polyurethane	NC
☐	129CW51	5" Micro Lighted	\$96.00
☐	129CW6	E2219 5" Soft Roll (Grey) N/A w/ Frog Legs	\$101.00
☐	129CW65	5"x 1" Aluminum Polyurethane	\$123.00
☐	129CW7	E2219 5"x 1.5" Semi Pneumatic	\$96.00
☐	129CW66	5"x 1.5" Aluminum Soft Roll	\$218.00
☐	129CW34	6" Polyurethane	NC
☐	129CW35	E2219 6"x 1.5" Semi Pneumatic	\$96.00
☐	129CW68	6"x 1.5" Aluminum Soft Roll	\$229.00

Caster Fork

See Figures 5 and 6. Figure 6 in appendix pg 6.

☐	129CF65	Multi Position Fork	NC
☐	129CF69	Single Side Fork	\$189.00

Caster Fork Suspension

☐	129CF53	E1015	Frog Legs 0 weight to 100 lbs	\$486.00
☐	129CF54	E1015	Frog Legs 101 lbs weight to 190 lbs	\$486.00
☐	129CF55	E1015	Frog Legs 191 lbs weight to 265 lbs	\$486.00




Figure 5

Available FSH per Caster Fork

Caster Size	Short (Multi, Single Sided)	Standard (Multi, Single Sided)	Frog Legs
3"	16"-20"	16.5"-21"	17"-21"
4"	16"-20.5"	16.5"-21"	17"-21"
5"	17,18,19,20,21"	17"-21"	18"-21"
6"x1"	N/A	18"-21"	19,20,21"
6"x1.5"	N/A	18"-21"	19,20,21"

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR WHEELS / SEAT HEIGHT									

Rear Seat Height

Please reference RSH-Axle plate-wheel diagram on bottom of the page.

☐	129RH11	13.0"	NC
☐	129RH12	13.5"	NC
☐	129RH13	14.0"	NC
☐	129RH14	14.5"	NC
☐	129RH15	15.0"	NC
☐	129RH16	15.5"	NC
☐	129RH17	16.0"	NC
☐	129RH18	16.5"	NC
☐	129RH19	17.0"	NC
☐	129RH26	17.5"	NC
☐	129RH20	18.0"	NC
☐	129RH27	18.5"	NC
☐	129RH21	19.0"	NC
☐	129RH37	19.5"	NC
☐	129RH22	20.0"	NC

Camber

0,3,6 available with std frame. 0,4,8 available with XTR frame

☐	129CI1	0°	NC
☐	129CI2	3°	NC
☐	129CI3	4°	NC
☐	129CI4	6°	NC
☐	129CI5	8°	NC

Axle Plate - Wheel - Available RSH					
Axle Plate	20" Wheel	22" Wheel	24" Wheel	25" Wheel	26" Wheel
Std Short	13"-17"	13.5"-18"	14.5"-19"	15"-19.5"	15.5"-20"
Std Long	13"-18.5"	13.5"-19.5"	14.5"-20"	15"-20"	15.5"-20"
Amp Short	13.5"-17"	14.5"-18"	15.5"-19"	16"-19.5"	16.5"-20"
Amp Long	13.5"-19"	14.5"-20"	15.5"-20"	16"-20"	16.5"-20"
Susp	14"-17.5"	15"-18.5"	16"-19.5"	16.5"-20"	17"-20"
Xtender	N/A	17.5"-19.5"	18.5"-20"	N/A	N/A

*Note: other factors such as CG or Transit may limit RSH

Xtender Power Assist

RSH 17.5"-19.5" with 22" wheels, RSH 18.5"-20" with 24" wheels 0 camber, N/A w/ Amputee Axle. 22" & 24" wheels only. Not available with XTR or standard rear suspension. Anti-Tips no charge with Xtender. Dual Drive mode available w/B4ME

☐	129XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$8,991.00
☐	129XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$9,885.00
☐	129XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$8,991.00
☐	129XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$9,885.00

Please select one of the following handrim options

☐	129XH5	Plastic Coated	NC
☐	129XH6	Stainless Steel	NC

Xtender - Other Options

☐	129XA10	Spare NiMH Power Pack - NICKEL	\$1,559.00
☐	129XA11	Spare Li-ion Power Pack - LITHIUM	\$2,006.00

Xtender Wheel Options

☐	129XA8	SPARE 22" Lite Spoke Pneu Alum	NC
☐	129XA9	E2211 E2213 SPARE 22" Lite Spoke Pneu Plastic	\$224.00
☐	129XA2	SPARE 24" Lite Spoke Pneu Alum	NC
☐	129XA4	E2211 E2213 SPARE 24" Lite Spoke Pneu Plastic	\$224.00

Rear Wheel Size

Ultra Lightweight Wheels only available with aluminum anodized handrim

Contact Built 4 Me for Spinergy color options

☐	129RW103	20" Mag Wheel 5 Spoke	NC
☐	129RW104	22" Mag Wheel 5 Spoke	NC
☐	129RW102	24" Mag Wheel 5 Spoke	NC
☐	129RW90	20" Lite Spoke Wheel	\$224.00
☐	129RW91	22" Lite Spoke Wheel	\$224.00
☐	129RW92	24" Lite Spoke Wheel	NC
☐	129RW93	25" Lite Spoke Wheel	NC
☐	129RW94	26" Lite Spoke Wheel	\$224.00
☐	129RW06	22" Ultra Lightweight Spoke Wheel	\$279.00
☐	129RW10	24" Ultra Lightweight Spoke Wheel	\$279.00
☐	129RW28	22" Spinergy SPOX Black	\$944.00
☐	129RW95	24" Spinergy SPOX Black	\$944.00
☐	129RW96	25" Spinergy SPOX Black	\$944.00
☐	129RW97	26" Spinergy SPOX Black	\$944.00
☐	129RW32	22" Spinergy LX Black	\$1,185.00
☐	129RW98	24" Spinergy LX Black	\$1,185.00
☐	129RW99	25" Spinergy LX Black	\$1,185.00
☐	129RW40	26" Spinergy LX Black	\$1,185.00
☐	129RW50	24" Mountain	\$279.00

Rear Wheel Tire

See Availability Matrix on PAGE 4.

☐	129RT1	Pneumatic	NC
☐	129RT2	E2212 E2213 Pneumatic w/ Airless Insert N/A with Xtender	\$96.00
☐	129RT3	Full Poly	NC
☐	129RT38	E2218 SHOX Performance Poly	\$106.00
☐	129RT39	E2211 E2212 High-Pressure Clincher Black	\$106.00
☐	129RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$157.00
☐	129RT27	E2211 E2212 Knobby V- Trak by Kenda	\$106.00
☐	129RT12	E2211 E2212 Mountain Tire Only available in 24"	\$106.00

Handrim

Plastic is not available on all wheels in all sizes. Natural Fits available in 22" 24" and 25". Natural Fit 5 spoke MAG only in 24". Handrim Selection N/A w/ Xtender.

See matrix figure 7 in appendix page 6 for reference.

☐	129HR1	Aluminum Anodized	NC
☐	129HR4	Plastic Coated N/A w/ 25" & 26" Lite Spoke Wheels	\$123.00
☐	129HR22	Natural Fit - Standard Grip	\$526.00
☐	129HR25	Natural Fit - No Thumb Grip	\$455.00
☐	129HR23	Natural Fit - Super Grip	\$585.00
☐	129HR30	Natural Fit LT - Standard Grip	\$526.00
☐	129HR31	Natural Fit LT - Super Grip	\$585.00
☐	129HR32	Natural Fit LT - No Thumb Grip	\$455.00
☐	129HR26	The Surge - Oval w/Gription Strip	\$472.00
☐	129HR27	The Surge LT - Oval w/Gription Strip	\$508.00
☐	129HR28	Q-Grip Neoprene Coated	\$414.00
☐	129HR24	Spinergy Flex Rim Only available on Spinergy LX wheels 24" & 25"	\$782.00
☐	129HR8	8 Oblique Projections N/A w/ SHOX tires, only available with 24" Lite Spoke	\$363.00
☐	129HR7	Omit Handrim N/A w/ Ultra Lightweight Wheels	NC

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR WHEELS / SEAT HEIGHT (CONT)									


Axle			
☐	129AX1	Quick Release Steel	NC
☐	129AX13	Quad Quick Release Only available on 24" Mountain rear wheels	\$141.00

Axle Plates			
Selection not required for XTR Suspension			
☐	129AL40	Standard Short (COG 0.5" to 4") Some limitations with STFH	NC
☐	129AL43	Standard Long (COG 0.5" to 4")	NC
☐	129AL100	E0959 Axle plate amputee short	\$196.00
☐	129AL101	E0959 Axle plate amputee long	\$196.00
☐	129AL38	E1015 Suspension (COG -0.5 to +3) N/A with XTR Suspension Frame	\$441.00
☐	129AL39	Xtender (COG 1"-4")	NC

Rear Suspension Rating			
Selection not required for XTR Suspension. N/A with Xtender			
☐	129SV04	From 0 weight to 100 lbs (Yellow)	NC
☐	129SV05	From 101 lbs weight to 190 lbs (Red)	NC

Center of Gravity			
Measure from front of back post to center of rear axle. XTR Suspension Frame will be factory set at 1.5", Xtender 1"- 4"			
☐	129CC29	-1.5"	NC
☐	129CC11	-1"	NC
☐	129CC10	-.5"	NC
☐	129CC12	0"	NC
☐	129CC1	.5"	NC
☐	129CC2	1"	NC
☐	129CC3	1.5"	NC
☐	129CC4	2"	NC
☐	129CC5	2.5"	NC
☐	129CC6	3"	NC
☐	129CC7	3.5"	NC
☐	129CC8	4"	NC

Tire	Wheel				
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain
Pneumatic	20,22,24	20,22,24,26		22,24,26	
Pneumatic w/ Airless Insert	20,22,24	20,22,24		22,24	
Full Poly	20,22,24	20,22,24			
SHOX Performance Poly	20,22,24	20,22,24,25		24,25	
High-Pressure Clincher Blk		22,24,26	22,24	22,24,26	
Schwalbe Marathon		24,25,26		24,25,26	
Knobby V Trak by Kenda		24		24	
Mountain Tire					24
20" = 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590mm					



5. WHEEL LOCKS

Wheel Lock			
Not all wheel locks available in all seat depths, wheel sizes and COG settings. XTR suspension only available w/ high mount push			
☐	129WL46	High Mount Push	NC
☐	129WL50	High Mount Pull	NC
☐	129WL68	Ergo Scissor	NC
☐	129WL64	Compact	\$157.00
☐	129WL10	Omit Wheel Lock	NC
☐	129WL12	Do Not Mount	NC

Extension Handle			
Only available w/ High Mount Push or Pull Wheel Locks			
☐	129EH7	E0961 Extension handles (Pair)	\$84.00

Grade Aid			
Only avail w/ high mount push or pull. N/A w/ clincher, poly, SHOX, Schwalbe, or Xtender			
☐	129GA1	E0974 Grade Aid	\$168.00

6. FOOTREST OPTIONS

Footrest			
Angle Adjustable Flip Back N/A w/ 19-20 wide & w/ 0" inset. Some Flip back configurations won't allow 2" ground clearance and will need to raise up. Flip back not available 1" inset & 5"- 6" casters. See matrix in Appendix			
☐	129F229	Aluminum Tubular	NC
☐	129F274	Aluminum Tubular w/ Snap On Cover	\$34.00
☐	129F276	K0040 Angle Adjustable Flip Back	\$207.00
☐	129F51	K0040 Aluminum Platform Footrest	\$106.00
☐	129F278	K0040 High Mount Platform	\$106.00

Adaptive Platform Spacer			
Only available with Platform footrest & Angle Adj. Flip Back.			
☐	129AD1	Adaptive Platform 2" Spacer	\$106.00
☐	129AD2	Adaptive Platform 4" Spacer	\$106.00
☐	129AD3	Adaptive Platform 6" Spacer	\$106.00

Footrest Options			
☐	129DI03	K0038 Leg Strap	\$45.00
☐	129DI04	K0038 Extra Leg Strap	\$45.00
☐	129DI05	Gel Hanger Cover	\$118.00
☐	129DI07	Neoprene Impact Guards	\$101.00
☐	129DI02	Wear Pads	\$56.00
☐	129DI06	Position Plates (N/A for Tubular Footrest)	\$56.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST									

Back Type			
☐	129BT3	Non-Folding Angle Adjustable	NC
☐	129BT82	Folding w/ Double Lock Down Angle Adj	NC

Height Adjustable Backrest			
☐	129B110	10" - 14" Back Height	NC
		☐ 10 ☐ 11 ☐ 12 ☐ 13 ☐ 14	
☐	129B127	12" - 16" Back Height	NC
		☐ 12 ☐ 13 ☐ 14 ☐ 15 ☐ 16	
☐	129B111	14" - 18" Back Height	NC
		☐ 14 ☐ 15 ☐ 16 ☐ 17 ☐ 18	
☐	129B137	16" - 20" Back Height	NC
		☐ 16 ☐ 17 ☐ 18 ☐ 19 ☐ 20	

Push Handles			
☐	129PH6	Integral Push Handles	NC
☐	129PH5	Fold-Down Push Handles N/A w/ 10-14" Back Height or Standard Upholstery.	\$111.00
☐	129PH2	Omit Push Handles	NC

JAY J3 Back			
☐	129JB01	JAY J3 Back	See OF

Back Upholstery			
☐	129BU1	Standard Upholstery	NC
☐	129BU88	Standard Upholstery EXO	\$56.00
☐	129BU70	E2611 Tension Adj Ballistic Nylon *Reference 186BU70 for billing	\$352.00
☐	129BU71	E2611 Tension Adj 3DX Vented *Reference 186BU71 for billing	\$370.00
☐	129BU89	E2611 Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$370.00
☐	129BU14	Omit Upholstery (Use w/ JAY Backs)	-\$15.00

Military Service Patches			
☐	129VU1	US Marines Patch	\$111.00
☐	129VU2	US Army Patch	\$111.00
☐	129VU3	US Navy Patch	\$111.00
☐	129VU4	US Air Force Patch	\$111.00
☐	129VU12	US Coast Guard	\$111.00
☐	129VU5	Purple Heart	\$111.00
☐	129VU8	Purple Heart with Marines Patch	\$173.00
☐	129VU9	Purple Heart with Army Patch	\$173.00
☐	129VU10	Purple Heart with Navy Patch	\$173.00
☐	129VU11	Purple Heart with Air Force Patch	\$173.00
☐	129VU13	Purple Heart with Coast Guard	\$173.00

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away			
Includes angle adj receiver hardware. Some limitations may exist w/ integral push handles, some RSHs and wheel sizes			
☐	129AA18	E0973 Padded Swing-Away STD (9" to 14") Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"	\$252.00

Armrest - Height Adjustable			
☐	129AA11	E0973 Single Post Ht. Adjust. Std (8.5"-13.5")	\$319.00
☐	129AA29	E0973 Single Post Ht. Adjust. Low (7"-11")	\$319.00

Armpad - For Height Adjustable			
Recomend using classic style for tray tables			
☐	129AM23	Desk Length (10") Classic	NC
☐	129AM21	Desk Length (10") Contour	NC
☐	129AM24	Full Length (14") Classic	NC
☐	129AM22	Full Length (14") Contour	NC

Side Guards			
Some interference may exist w/ padded SA arms at lowest settings			
☐	129SG33	K0108 Plastic - Standard	\$162.00
☐	129SG29	K0108 Plastic - Low	\$162.00
☐	129SG35	K0108 Aluminum N/A w/ push, pull or compact wheel lock. N/A w/ XTR suspension	\$274.00
☐	129SG36	K0108 Aluminum w/ Fender N/A w/ push, pull or compact wheel locks. N/A w/ XTR suspension	\$308.00
☐	129SG12	K0108 Carbon Fiber N/A w/ push, pull or compact wheel locks	\$514.00

Color For Aluminum Side Guard			
☐	129GE01	Black	NC
☐	129GE02	Matching Frame - Standard Colors	NC
☐	129GE03	Matching Frame - Personalized Color	\$56.00

9. KEY OPTIONS AND ACCESSORIES

Anti Tips			
Anti Tip Single Side N/A 12W-14W			
Anti Tip Single Side N/A with XTR, Xtender and Suspension Axles.			
☐	129AT1	E0971 Anti Tips (Pair)	\$168.00
☐	129AT3	Omit Rear Anti Tip	NC
☐	129AT14	E0971 Anti Tip Single Side - Flip Up ☐ Left ☐ Right (Facing from rear)	\$196.00

Pack Accessories			
☐	129A1	Black Utility Backpack	\$84.00
☐	129A2	Back Pack Black/Charcoal Kids	\$67.00
☐	1291A01	Seat Pouch Black	\$56.00

Positioning Belt			
☐	129A826	E0978 Belt 1.5" Auto Buckle Universal	\$78.00
☐	129A802	E0978 Belt 2" Aircraft Buckle	\$123.00
☐	129A803	E0978 Belt 2" Aircraft Buckle Padded	\$133.00
☐	129A804	Positioning Belt Clamps Only	\$56.00

Spoke Guards			
Only available with Lite spoke. N/A w/ Xtender Power Assist			
☐	1294A01	K0065 Quickie Black Spoke Guard	\$128.00
☐	1294A02	K0065 Quickie Clear Spoke Guard	\$128.00

Miscellaneous			
☐	129A2100	E2207 Crutch Holder Only available with Integral or Fold Down push handle. Not available with XTR Suspension	\$218.00
☐	1298A22	Bolt-On Push Handles	\$128.00
☐	1298A21	Tool Kit	\$61.00
☐	1298A05	Gel-Pad swing away arm	\$168.00
☐	129A999	Quickie Water Bottle & Holder	\$61.00

Upper Extremity Support Tray			
☐	129A29	E0950 X-Small (Suggested for 10"-12")	\$285.00
☐	129A30	E0950 Small (Suggested for 13"-14")	\$285.00
☐	129A31	E0950 Medium (Suggested for 15"-17")	\$285.00
☐	129A32	E0950 Large (Suggested for 18"-20")	\$285.00

10. COLOR

Standard Frame Colors

☐	129C140	Aztec Gold	NC
☐	129C2	Black (Not avail as base w/Zebra)	NC
☐	129C141	Black Cherry	NC
☐	129C10	Black Opal	NC
☐	129C38	Blue Opal	NC
☐	129C36	Candy Blue	NC
☐	129C14	Candy Purple	NC
☐	129C35	Candy Red	NC
☐	129C142	Electric Blue	NC
☐	129C143	Evergreen	NC
☐	129C44	Glossy White	NC
☐	129C126	Glow	NC
☐	129C144	Green Apple	NC
☐	129C147	Hot Sparkle Pink	NC
☐	129C148	Mauve Pink	NC
☐	129C103	Rootbeer	NC
☐	129C114	Sparkle Silver	NC
☐	129C145	Sunrise Orange	NC
☐	129C113	Titanium Color Paint	NC
☐	129C8	Yellow	NC

Matte Frame Colors

☐	129C116	Matte Black	NC
☐	129C151	Matte Black Cherry	NC
☐	129C150	Matte Electric Blue	NC
☐	129C149	Matte Evergreen	NC
☐	129C152	Matte Purple	NC

Personalized Frame Color

☐	129C112	Desert Camouflage	\$441.00
☐	129C43	Mossy Oak Camouflage	\$441.00
☐	129C42	Stars & Stripes	\$441.00
☐	129C108	Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Blck, Matte Black, Blck Opal or Titanium)	\$441.00

Anodized Color Package

Package includes only forks, backcane bracket, axle sleeve clamp, footrest clamp, Axle plate plug.

Frog Legs, Xtender, suspension axle plates black only

☐	129EB02	Black	NC
☐	129EB06	Blue	\$179.00
☐	129EB05	Lime Green	\$179.00
☐	129EB04	Purple	\$179.00
☐	129EB03	Red	\$179.00
☐	129EB01	Ti Grey	\$179.00
☐	129EB13	Gold	\$179.00

Built For Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011



11. APPENDIX

Figure 6

		Footrest-Caster Fork Matrix			Tubular/Flat			Platform			AA Flipback		
		Caster Size	Fork	Inset	Inset	Inset	Inset	Inset	Inset	Inset	Inset	Inset	Inset
75° Frame Angle	3"	All Forks									N/A		
	4"	All Forks									N/A		
	5"	All Multi/Carb/Sngl									N/A	N/A	
	6"	Froglegs									N/A	N/A	
80° Frame Angle	3"	All Forks									N/A		
	4"	All Multi/Carb/Sngl									N/A	N/A	
	5"	Froglegs									N/A	N/A	
	6"	All Multi/Carb/Sngl									N/A	N/A	
85° Frame Angle	3"	All Forks									N/A		
	4"	All Multi/Carb/Sngl									N/A	N/A	
	5"	Froglegs									N/A	N/A	
	6"	All Multi/Carb/Sngl (tall only)									N/A	N/A	
88° Frame Angle	3"	All Forks									N/A		
	4"	All Multi/Carb/Sngl									N/A	N/A	
	5"	Froglegs									N/A	N/A	
	6"	All Multi/Carb/Sngl (tall only)									N/A	N/A	

Go

N/A = No Go

#

= "Soft" Go with limited seat angle (FSH-RSH). Min seat angle listed in cell.



Sunrise Medical (US) LLC
2842 Business Park Ave. • Fresno, CA 93727 • USA

Wheels

Handrims	5 Spoke Mag	Lite Spoke	Ultra Lightwgt	Spinergy LX	Spinergy SPOX	Mountain
Aluminum Anodized	20,22,24	20,22,24,25,26	22,24	22,24,25,26	22,24,25,26	24
Plastic Coat	20,22,24	22,24		24,26	24,26	
Natural Fit - Standard Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit - No Thumb Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit - Super Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - Std Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - Super Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - No Thumb Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Surge - Oval w/Gription Strip		22,24,25,26		22,24,25,26	22,24,25,26	
Surge LT - Oval w/Gription Strip		22,24,25,26		22,24,25,26	22,24,25,26	
Q Grip Neoprene Coated		24,25		24,25	24,25	
Spinergy Flex Rim				24,25		
8 Oblique Projection		24				
Omit Handrim	20,22,24			22,24,25,26	22,24,25,26	24