



Quickie ACCESS

February 2023

EIR14

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
1. FRAME									

Base Model			
☐	208WM105	E1161	Quickie Access 14" sizes use E1233
☐	208WM110	E1161	Quickie Access - Transit included 14" sizes use E1233

Pre-Set Tilt Angle			
Certain wheel combinations will limit the seat angle			
☐	208TA11	E1161	-5° to 20° Degree Tilt
☐	208TA12	E1161	0° to 25° Degree Tilt

Tilt Release Options			
☐	208TT8		Foot Release NA with 19"-22" Depth
☐	208TT12		Backrest trigger cable release

Transit			
☐	208TO45		Transit Option
Frame Options			
☐	208FS23		Standard Base and Seat Frame (Maximum Weight: 300 lbs)

User's Weight			
Gas spring selection depends User's weight selection.			
☐	208WG26		Up to 150 lbs
☐	208WG27		151 - 180 lbs
☐	208WG28		181 - 210 lbs
☐	208WG29		211 - 240 lbs
☐	208WG30		241 - 270 lbs
☐	208WG31		271 - 300 lbs



ACCESS

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
2. SIZE / SEATING									


Frame Width				
21" & 22" NA w/ removable stroller handle extension or back upholsteries.				
Small Frame 14" - 17"				
☐	208W5		14"	NC
☐	208W6		15"	NC
☐	208W7		16"	NC
☐	208W8		17"	NC
Large Frame 18" - 22"				
☐	208W9		18"	NC
☐	208W10		19"	NC
☐	208W11	E2201	20"	\$559.00
☐	208W12	E2201	21"	\$559.00
☐	208W13	E2201	22"	\$559.00

Seat Depth				
See Figure 1 Dimension A for measurement				
Foot release tilt option NA w/19"-22"				
Short 14" - 16"				
☐	208SD7		14"	NC
☐	208SD8		15"	NC
☐	208SD9		16"	NC
Medium 16" - 18"				
☐	208SD9		16"	NC
☐	208SD10		17"	NC
☐	208SD11		18"	NC
Long 18" - 20"				
☐	208SD11		18"	NC
☐	208SD12		19"	NC
☐	208SD13	E2203	20"	\$614.00
Extra Long 20" - 22"				
☐	208SD13	E2203	20"	\$614.00
☐	208SD14	E2203	21"	\$614.00
☐	208SD33	E2204	22"	\$1,006.00

Seat Type				
☐	208ST35	E2231	Standard Seat Pan <i>For growth ranges, see Figure 1.5</i>	\$224.00
☐	208ST54		Omit seat Pan	NC

JAY Cushions (Optional)				
To order any JAY cushions on-chair, please see the JAY price list.				
Cushion name or part number must be filled out with desired width and depth.				
☐	208JC20	JAY Bundled Cushion Cushion Name/Part #:		Price List
		Width:	Depth:	

Custom Seating				
Add potential NoGos with Custom Seating (pending review)				
☐	208QC01		JAY Custom (AES/JCF)	
☐	208QC02		JAY Surefit	
☐	208QC03		JAY Fit	



Frame Width	Seat Pan Growth
14"	14" - 16"
15"	15" - 16"
16"	16" - 18"
17"	17" - 18"
18"	18" - 20"
19"	19" - 20"
20"	20" - 22"
21"	21" - 22"
22"	NA

3. CASTERS / SEAT HEIGHTS / REAR WHEELS

Front Seat Height								
12"-12.5" with 4" casters are limited to a 1.5" floor clearance.								
Rear Wheel Size	CASTER SIZE							
	4"	5"	6"	7"	8" x 1"	8"x 1.5"	8" x 2"	
12" - 24"	☐ 12.0"	☐ 12.0"	☐ 12.0"					
	☐ 12.5"	☐ 12.5"	☐ 12.5"					
	☐ 13.0"	☐ 13.0"	☐ 13.0"	☐ 13.0"				
	☐ 13.5"	☐ 13.5"	☐ 13.5"	☐ 13.5"		☐ 13.5"	☐ 13.5"	
	☐ 14.0"	☐ 14.0"	☐ 14.0"	☐ 14.0"	☐ 14.0"			
	☐ 14.5"	☐ 14.5"	☐ 14.5"	☐ 14.5"		☐ 14.5"	☐ 14.5"	
	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	☐ 15.0"	
	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"	
	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"	
	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"	
	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	
	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"	
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"	
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"	
	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.0"			
	☐ 19.5"	☐ 19.5"	☐ 19.5"	☐ 19.5"				
	☐ 20.0"	☐ 20.0"	☐ 20.0"	☐ 20.0"		☐ 20.0"		

Axle Plate				
☐	208AL33		Standard Axle Plate	NC
Casters				
☐	208CW3		4" Low Profile Polyurethane N/A w/ HD package	NC
☐	208CW72		4" x 1" Aluminum Polyurethane	\$123.00
☐	208CW31		4" Micro Lighted	\$78.00
☐	208CW53		5" x 1" Polyurethane	NC
☐	208CW7		5" x 1.5" Semi Pneumatic	\$96.00
☐	208CW45		6" x 1" Polyurethane	NC
☐	208CW35		6" x 1.5" Semi Pneumatic	\$96.00
☐	208CW16		7" Polyurethane	NC
☐	208CW46		8" x 1" Polyurethane	NC
☐	208CW39		8"x 1.5" Semi Pneumatic	\$96.00
☐	208CW37	E2214 K0108	8"x 2" Pneumatic w/Insert	\$252.00

Rear Wheel				
☐	208RW13		12" Mag Wheel	NC
☐	208RW14		16" Mag Wheel	NC
☐	208RW92		24" Lite Spoke Wheel	NC
☐	208RW103		20" Mag Wheel 5 Spoke	NC
☐	208RW104		22" Mag Wheel 5 Spoke	NC
☐	208RW102		24" Mag Wheel 5 Spoke	NC
☐	208RW150		24" Drum brake wheel - 36 spokes	NC

Rear Wheel Tire				
☐	208RT3		Full Polyurethane	NC
☐	208RT4		Low Polyurethane	NC
☐	208RT1	E2211 E2212	Pneumatic	NC
☐	208RT2	E2213 E2211	Pneumatic Airless Insert	\$96.00

Rear Tire & Wheel Matrix				
Tire	Rear Wheel			
	5 Spoke Mag	Spoke	Mag Wheel	Drum Brake
Pneumatic	20,22,24	24	12,16	24
Pneumatic w/ Airless Insert	20,22,24	24	12,16	24
Full Poly	20,22,24	24	N/A	24
Low Poly	N/A	N/A	12	N/A
12" = 203mm 16" = 305mm 20" = 451mm 22" = 501mm 24" = 540mm				

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

3. CASTERS / SEAT HEIGHTS / REAR WHEELS (CONT)

Handrim			
<i>12" and 16" wheels not available with hand rims. See matrix for availability</i>			
☐	208HR1	Aluminum Anodized	NC
☐	208HR4	Plastic Coated	\$123.00
☐	208HR22	Natural Fit - Standard Grip	\$526.00
☐	208HR23	Natural Fit - Super Grip	\$585.00
☐	208HR25	Natural Fit - No Thumb Grip	\$455.00
☐	208HR30	Natural Fit LT - Standard Grip	\$526.00
☐	208HR31	Natural Fit LT - Super Grip	\$585.00
☐	208HR32	Natural Fit LT - No Thumb Grip	\$455.00
☐	208HR7	Omit Handrim	NC

Axle			
☐	208AX1	Quick Release Axles	NC
☐	208AX17	Threaded Axles <i>N/A w/ 24" spoke wheels & 5 Spoke Mags</i>	NC

Rear Wheel & Handrim Matrix

Handrim	Rear Wheel		
	5 Spoke Mag	Spoke	Drum Brake
Aluminum Anodized	20,22,24	24	24
Plastic Coat	20,22,24	24	24
Natural Fit - Standard Grip	24	24	
Natural Fit - No Thumb Grip	24	24	
Natural Fit - Super Grip	24	24	
Natural Fit LT - Standard Grip	24	24	
Natural Fit LT - Super Grip	24	24	
Natural Fit LT - No Thumb Grip	24	24	

4. WHEEL LOCKS

Wheel Lock			
☐	208WL8	Push to Lock	NC
☐	208WL41	E2206 Attendant Wheel Lock (Foot Lock) NA w/ ext handles. NA w/19"-22" Depth.	\$90.00
☐	208WL79	E228 Drum Brake	\$559.00
☐	208WL10	Omit Wheel Lock	NC

Extension Handle			
☐	208EH7	E0961 6" Pair	\$84.00

Hub Lock			
<i>Only available with smaller wheels</i>			
☐	208LK01	12" Hub Lock	\$111.00
☐	208LK02	16" Hub Lock	\$111.00

5. HANGER AND FOOTREST

Hangers			
<i>Footrests mount to extension tubes</i>			
☐	208H145	60° Swing In/Swing Out	NC
☐	208H139	70° Swing In/Swing Out	NC
☐	208H141	80° Swing In/Swing Out <i>NA w/ 8" casters</i>	NC
☐	208H160	E0990 Elevating Hanger - Ext Mount.	\$352.00
☐	208H159	K0053 Articulating Hanger	\$413.00
☐	208H32	Omit Hanger	-\$15.00

Extension Tubes			
☐	208E1	Short	NC
☐	208E2	Medium <i>NA with 12.5" seat height</i>	NC
☐	208E5	Long <i>NA with 12.5" seat height</i>	NC
☐	208E16	High Mount <i>N/A with ELR/ALR hangers. Only available with adult footplates</i>	\$84.00

Footrest Footplates			
<i>Heel loop or leg strap not included. See Footplate Options Section</i>			
☐	208F1	Composite	STD
☐	208F271	K0040 E0951 Composite Angle Adj-(16"-22"w)	\$202.00
☐	208F5	Aluminum <i>Only available w/ELR & ALR</i>	NC
☐	208F3	K0040 Adult Angle Adjustable - Ext Mount	\$212.00
☐	208F7	Omit Footplate and Ext. Tube	-\$15.00

Footrest Options			
☐	208HL1	E0951 Heel Loop	\$45.00
☐	208HL5	K0038 Leg Strap	\$45.00
☐	208IG2	Impact Guard Neoprene <i>Only available w/ 60° or 70° ext. mount hanger</i>	\$101.00

JAY Positioning			
☐	208K3	E1020 Left JAY Lower Extremity S/A support	\$441.00
		☐ 208K23 Pad 6" x 8" Left	NC
		☐ 208K25 Pad 8" x 10" Left	NC
☐	208K4	E1020 Right JAY Lower Extremity S/A support	\$441.00
		☐ 208K22 Pad 6" x 8" Right	NC
		☐ 208K24 Pad 8" x 10" Right	NC
☐	208K2	Omit Hanger ☐ Left ☐ Right	-\$15.00


Articulating / Elevating Hanger Cover			
☐	208DI47	Neoprene Cover	\$380.00
☐	208DI48	Gel Pad Cover	\$447.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
6. BACKREST / SEATING									

Backpost Style - (Choose ONE Option - A, B, C)

A. Stroller Handle Backpost (Fixed)			
☐	208BP9	Low (18")	NC
☐	208BP10	Medium (21")	NC
☐	208BP11	Tall (24")	NC

B. Height Adjustable Back w/ Adj. Stroller Handles			
Select desired Push Handle			
☐	208PH8	Standard Push Handle	NC
☐	208PH9	Long Push Handle	NC
AND select desired Back Cane			
☐	208B200	Standard 15"-21"	NC
☐	208B201	Tall 18"-24"	NC



C. Reclining Backrest (0-50°)			
<i>NA w/ 22" Deep Frame.</i>			
☐	208BT21	E1225 Reclining Backrest	\$989.00
	☐	208B5 Medium	NC
	☐	208B9 Tall	NC
☐	208BP64	Angle Adj stroller handles Height adj range 19"- 23"	\$168.00
	☐	208PH8 Standard push handle	NC
	☐	208PH9 Long push handle	NC

Backrest Options and Accessories			
☐	208BP13	Removable Extension NA w/ 21"-22" frame width, height adj back with adj stroller handles	\$96.00
☐	208BT101	Dynamic Backrest NA w/ Reclining Backrest or 19"-22" seat depths	\$1,671.00


JAY J3 Back			
☐	208JB01	JAY J3 Back NA w/ O2 Holder	See OF

Backrest Upholstery			
<i>ONLY Available with Low Fixed, Medium Fixed, Standard Height Adjustable (15"-21" back height) stroller handles and in 14"-20" widths. N/A with Dynamic, Reclining Backrest, O2 holder, IV Pole, JAY Anterior Thoracic Support and JAY Padded Trunk Strap Accessories. Towel bar included. EXO Tension Adj upholstery binding is <u>Standard black</u>.</i>			
☐	208BU70	E2611 Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$352.00
☐	208BU89	E2611 Backrest Tension Adj EXO *Reference 186BU70 for billing	\$370.00
Colors: ☐ Silver ☐ Orange ☐ Red ☐ Blue			\$15.00

7. ARMREST

Armrest			
☐	208AA65	E0973 Flip Back Height Adj (LITE 8.5"-13.5")	\$352.00
☐	208AA11	E0973 Single Post Height Adjustable - Std	\$319.00
☐	208AA29	E0973 Single Post Height Adjustable - Low	\$319.00
☐	208AA7	E0973 Length Adjustable Locking Flip-Up NA w/ Contour Armpad	\$308.00
☐	208AA1	E0973 Adjustable Locking Flip-Up	\$296.00
☐	208AA14	Omit Armrest	NC

Armpad			
<i>Ergo armpad only available with Single post Height adj armrest</i>			
☐	208AM23	Desk Length Classic	NC
☐	208AM21	Desk Length Contour	NC
☐	208AM24	Full Length Classic	NC
☐	208AM22	Full Length Contour	NC
☐	208AM9	Foam Short <i>Option for Length</i> Adj Locking Flip Up	NC
☐	208AM10	Foam Long <i>Option for Length</i> Adj Locking Flip Up	NC
☐	208AM27	Ergo armpad	NC

Armrest Matrix for Single Post						
Armrest Size	Transfer Tube Length D	Receiver	Classic Pad			
			H Range	Pad length L		
				Desk	Full	
Std	8"	Std	10-3/4" to 15-3/4"	10-1/4"	14"	
		Low	8-3/4" to 13-3/4"			
Low	5-1/4"	Std	9" to 13"			
		Low	7-1/4" to 11-1/4"			
			Contour Pad			
			Receiver	H Range	Pad length L	
					Desk	Full
			Std	10" to 15"	10-1/2"	14-1/4"
			Low	8-1/4" to 13-1/4"		
			Std	8-1/2" to 12-1/2"		
			Low	6-3/4" to 10-3/4"		
			Receiver	H Range	Ergo Pad	
					Pad length L	
			Std	10-1/2" to 15-1/2"	11-3/4"	
			Low	8-3/4" to 13-3/4"		
Std	9" to 13"					
Low	7-1/4" to 11-1/4"					

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

8. KEY OPTIONS AND ACCESSORIES

Anti Tips

Not ordered as separate line item.

☐	208AT1	E0971	Anti Tips (Pair)	\$168.00
--------------------------	--------	-------	------------------	----------

Positioning Belt

☐	208ZZ25	E0978	JAY 2PT Belt 2" Push Btn No Pad	\$46.00
☐	208ZZ26	E0978	JAY 2PT Belt 2" Push Btn 9.5" Pad	\$97.00
☐	208ZZ14	E0978	JAY 2PT Belt 1.5" Push Btn No Pad	\$46.00
☐	208ZZ15	E0978	JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$97.00
☐	208A117	E0978	2" Aircraft Buckle w/Clamps	\$123.00
☐	208A118	E0978	2" Aircraft Buckle Pad w/Clamp	\$133.00
☐	208A120	K0108	Positioning Belt Clamps Only	\$56.00

Wheelchair Tray Table

Only available if armrests are included. NA with JAY Arm Troughs.

☐	208A30	E0950	Small (13"-14")	\$285.00
☐	208A31	E0950	Medium (15"-17")	\$285.00
☐	208A32	E0950	Large (18"-20")	\$285.00
☐	208A33	E0950	X-Large (21"-22")	\$285.00

JAY Anterior Thoracic Support

Style

☐	208JY01	E0960	Classic Dynamic	\$170.00
☐	208JY02	E0960	Classic Structure	\$170.00
☐	208JY03	E0960	Contour Dynamic	\$170.00
☐	208JY04	E0960	Contour Structure	\$170.00
☐	208JY05	E0960	Center Opening with Buckle	\$170.00

Size

☐	208JZ01		Small	NC
☐	208JZ02		Medium	NC
☐	208JZ03		Large	NC
☐	208JZ04		X-Large	NC

JAY Padded Trunk Strap

☐	208JA01	E0960	Padded Trunk Strap - Small	\$158.00
☐	208JA02	E0960	Padded Trunk Strap - Medium	\$158.00
☐	208JA03	E0960	Padded Trunk Strap - Large	\$158.00
☐	208JA04	E0960	Padded Trunk Strap - XL	\$158.00

Miscellaneous Accessories

☐	208A1		Black Utility Backpack	\$84.00
☐	208A8		Black Seat Pouch	\$56.00
☐	208LS14	K0038	Extra Leg Strap Black	\$34.00
☐	208A601	E2208	O2 Holder NA w/ Height Adj Stroller Handles. NA w/ Reclining Backrest	\$179.00
☐	208A701	K0105	IV Pole	\$308.00
☐	208A999		Quickie Water Bottle & Holder	\$61.00

9. COLOR

Standard Frame Color

Color options only to available to hangers, seat tubes on seat. Front and back extensions on the base

☐	208C140	Aztec Gold	NC	☐	208C44	Glossy White	NC
☐	208C2	Black	NC	☐	208C126	Glow	NC
☐	208C141	Black Cherry	NC	☐	208C144	Green Apple	NC
☐	208C10	Black Opal	NC	☐	208C147	Hot Sparkle Pink	NC
☐	208C38	Blue Opal	NC	☐	208C148	Mauve Pink	NC
☐	208C36	Candy Blue	NC	☐	208C103	Rootbeer	NC
☐	208C14	Candy Purple	NC	☐	208C114	Sparkle Silver	NC
☐	208C35	Candy Red	NC	☐	208C145	Sunrise Orange	NC
☐	208C142	Electric Blue	NC	☐	208C113	Titanium	NC
☐	208C143	Evergreen	NC	☐	208C8	Yellow	NC

Matte Frame Color

☐	208C116		Matte Black	NC
☐	208C151		Matte Black Cherry	NC
☐	208C150		Matte Electric Blue	NC
☐	208C149		Matte Evergreen	NC
☐	208C152		Matte Purple	NC

Personalized Frame Color

☐	208C112		Desert Camouflage	\$441.00
☐	208C43		Mossy Oak Camouflage	\$441.00
☐	208C42		Stars & Stripes	\$441.00
☐	208C108		Zebra (Choose above base color) Std with White unless above base paint is chosen NA w/ Black	\$441.00

Anodized Color Package

Anodize package include Forks, axle plate, back plate.

☐	208EB07		Black	NC
☐	208EB08		Blue	\$179.00
☐	208EB13		Gold	\$179.00
☐	208EB09		Lime Green	\$179.00
☐	208EB10		Purple	\$179.00
☐	208EB11		Red	\$179.00
☐	208EB12		Ti Grey	\$179.00

Additional Options Through Built for Me


To inquire about any unique personalization for Quickie products, contact: Phone: (800) 333-4000, Fax: (800) 333 9011



Lower Leg Length Matrix

Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers.

Hanger	Extension	Seat Pan to Top of Footrest	
		Min	Max
60° Swing In/Out	Short	8 1/2	14
	Med	8 1/2	15 1/2
	Long	9 1/2	17 1/2
70° Swing In/Out	Short	9 1/2	14 1/2
	Med	10	16 1/2
	Long	11 1/2	18 1/2
80° Swing In/Out	Short	10	14
	Med	13	16
	Long	13	18
ALR	Short	9	13
	Med	9	14 1/2
	Long	10	16 1/2
ELR	Short	9	13
	Medium	9	14 1/2
	Long	10	16 1/2



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA



JAY® Signature Series Backrest and Accessories

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.com


When ordering mark the part number desired and it will be added to the order as a separate line item .

JAY BACKREST

JAY J3 Posterior (PA) BACK - 2.2" Contour Depth

	Shell Height		Shell Width							
Mid Thoracic (MT)	12.0" (S)	J3PAMTSxxFS	☐ 12	\$649	☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
		J3PAMTSxxTH	☐ 12	\$649	☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
Upper Thoracic (UT)	16.5" (S)	J3PAUTSxxFS	☐ 12	\$649	☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
		J3PAUTSxxTH	☐ 12	\$649	☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
	18.0" (M)	J3PAUTMxxFS			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
		J3PAUTMxxTH			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
	19.5" (T)	J3PAUTTxxFS			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
		J3PAUTTxxTH			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
Shoulder Height (SH)	21.0" (S)	J3PASHSxxFS			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
		J3PASHSxxTH			☐ 14	\$649	☐ 16	\$649	☐ 18	\$649
	22.5" (M)	J3PASHMxxFS					☐ 16	\$649	☐ 18	\$649
		J3PASHMxxTH					☐ 16	\$649	☐ 18	\$649
	24.0" (T)	J3PASHTxxFS						☐ 18	\$649	☐ 20
		J3PASHTxxTH						☐ 18	\$649	☐ 20


FS Mounting Hardware



JAY J3 Posterior Lateral (PL) BACK - 3" Contour Depth

Mid Thoracic (MT)	15.0" (T)	J3PLMTTxxFS	☐ 12	\$748	☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
		J3PLMTTxxTH	☐ 12	\$748	☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
Upper Thoracic (UT)	16.5" (S)	J3PLUTSxxFS	☐ 12	\$748	☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
		J3PLUTSxxTH	☐ 12	\$748	☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
	18.0" (M)	J3PLUTMxxFS			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
		J3PLUTMxxTH			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
	19.5" (T)	J3PLUTTxxFS			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
		J3PLUTTxxTH			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
Shoulder Height (SH)	21.0" (S)	J3PLSHSxxFS			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
		J3PLSHSxxTH			☐ 14	\$748	☐ 16	\$748	☐ 18	\$748
	22.5" (M)	J3PLSHMxxFS					☐ 16	\$748	☐ 18	\$748
		J3PLSHMxxTH					☐ 16	\$748	☐ 18	\$748
	24.0" (T)	J3PLSHTxxFS						☐ 18	\$748	☐ 20
		J3PLSHTxxTH						☐ 18	\$748	☐ 20

TH Mounting Hardware



JAY J3 Posterior Deep (PD) BACK - 6" Contour Depth

Mid Thoracic (MT)	15.0" (T)	J3PDMTTxxFS	☐ 12	\$889	☐ 14	\$889	☐ 16	\$889	☐ 18	\$889
		J3PDMTTxxTH	☐ 12	\$889	☐ 14	\$889	☐ 16	\$889	☐ 18	\$889
Upper Thoracic (UT)	18.0" (M)	J3PDUTMxxFS			☐ 14	\$889	☐ 16	\$889	☐ 18	\$889
		J3PDUTMxxTH			☐ 14	\$889	☐ 16	\$889	☐ 18	\$889
Shoulder Height (SH)	21.0" (S)	J3PDSHSxxFS			☐ 14	\$889	☐ 16	\$889	☐ 18	\$889
		J3PDSHSxxTH			☐ 14	\$889	☐ 16	\$889	☐ 18	\$889

JAY J3 Back Whitmyer Headrest

AXYS Mount and Pad

☐	J3HEDX8	E0955/E1028	AXYS Headrest with 8" Plush Pad	\$430.00
☐	J3HEDX10	E0955/E1028	AXYS Headrest with 10" Plush Pad	\$430.00
☐	J3HEDX14	E0955/E1028	AXYS Headrest with 14" Plush Pad	\$444.00
☐	J3HEDXAAPN	E0955/E1028	AXYS with Adj-a-Plush Narrow Pad	\$542.00
☐	J3HEDXAAPS	E0955/E1028	AXYS with Adj-a-Plush Std Pad	\$542.00

LINX Mount and Pad

☐	J3HEDL8	E0955/E1028	LINX Headrest with 8" Plush Pad	\$499.00
☐	J3HEDL10	E0955/E1028	LINX Headrest with 10" Plush Pad	\$499.00
☐	J3HEDL14	E0955/E1028	LINX Headrest with 14" Plush Pad	\$513.00
☐	J3HEDLAAPN	E0955/E1028	LINX with Adj-a-Plush Narrow Pad	\$611.00
☐	J3HEDLAAPS	E0955/E1028	LINX with Adj-a-Plush Std Pad	\$611.00

JAY J3 Back Laterals Swing Away Laterals With Pad - Pair

☐	JLATSA44PR	E0956/E1028	SwingAway Lat with 4"x4" Pad (PR)	\$872.00	☐	JLATSA64PR	E0956/E1028	SwingAway Lat with 6"x4" Pad (PR)	\$872.00
☐	JLATSA45PR	E0956/E1028	SwingAway Lat with 4"x5" Pad (PR)	\$872.00	☐	JLATSA65PR	E0956/E1028	SwingAway Lat with 6"x5" Pad (PR)	\$872.00
☐	JLATSA46PR	E0956/E1028	SwingAway Lat with 4"x6" Pad (PR)	\$872.00	☐	JLATSA66PR	E0956/E1028	SwingAway Lat with 6"x6" Pad (PR)	\$872.00



JAY® Signature Series Cushions

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.com.

When ordering mark the part number desired and it will be added to the order as a separate line item.

JAY CUSHIONS

JAY Fusion Cushion				
☐	JFUSION1416	E2622	JAY Fusion 14" x 16"	\$487.00
☐	JFUSION1616	E2622	JAY Fusion 16" x 16"	\$487.00
☐	JFUSION1618	E2622	JAY Fusion 16" x 18"	\$487.00
☐	JFUSION1620	E2622	JAY Fusion 16" x 20"	\$487.00
☐	JFUSION1717	E2622	JAY Fusion 17" x 17"	\$487.00
☐	JFUSION1816	E2622	JAY Fusion 18" x 16"	\$487.00
☐	JFUSION1818	E2622	JAY Fusion 18" x 18"	\$487.00
☐	JFUSION1820	E2622	JAY Fusion 18" x 20"	\$487.00
☐	JFUSION2018	E2622	JAY Fusion 20" x 18"	\$487.00
☐	JFUSION2020	E2622	JAY Fusion 20" x 20"	\$487.00
☐	JFUSION2220	E2623	JAY Fusion 22" x 20"	\$580.00

JAY Fusion Cushion w/ Single Valve Air Insert				
☐	JFUSION1416AS	E2622	JAY Fusion 14" x 16"	\$548.00
☐	JFUSION1616AS	E2622	JAY Fusion 16" x 16"	\$548.00
☐	JFUSION1618AS	E2622	JAY Fusion 16" x 18"	\$548.00
☐	JFUSION1620AS	E2622	JAY Fusion 16" x 20"	\$548.00
☐	JFUSION1717AS	E2622	JAY Fusion 17" x 17"	\$548.00
☐	JFUSION1816AS	E2622	JAY Fusion 18" x 16"	\$548.00
☐	JFUSION1818AS	E2622	JAY Fusion 18" x 18"	\$548.00
☐	JFUSION1820AS	E2622	JAY Fusion 18" x 20"	\$548.00
☐	JFUSION2018AS	E2622	JAY Fusion 20" x 18"	\$548.00
☐	JFUSION2020AS	E2622	JAY Fusion 20" x 20"	\$548.00
☐	JFUSION2220AS	E2623	JAY Fusion 22" x 20"	\$639.00

JAY Fusion Cushion w/ Dual Valve Air Insert				
☐	JFUSION1416AD	E2622	JAY Fusion 14" x 16"	\$553.00
☐	JFUSION1616AD	E2622	JAY Fusion 16" x 16"	\$553.00
☐	JFUSION1618AD	E2622	JAY Fusion 16" x 18"	\$553.00
☐	JFUSION1620AD	E2622	JAY Fusion 16" x 20"	\$553.00
☐	JFUSION1717AD	E2622	JAY Fusion 17" x 17"	\$553.00
☐	JFUSION1816AD	E2622	JAY Fusion 18" x 16"	\$553.00
☐	JFUSION1818AD	E2622	JAY Fusion 18" x 18"	\$553.00
☐	JFUSION1820AD	E2622	JAY Fusion 18" x 20"	\$553.00
☐	JFUSION2018AD	E2622	JAY Fusion 20" x 18"	\$553.00
☐	JFUSION2020AD	E2622	JAY Fusion 20" x 20"	\$553.00
☐	JFUSION2220AD	E2623	JAY Fusion 22" x 20"	\$644.00

JAY Union Cushion				
☐	JU1416	E2607	JAY UNION 14" x 16"	\$498.00
☐	JU1616	E2607	JAY UNION 16" x 16"	\$498.00
☐	JU1618	E2607	JAY UNION 16" x 18"	\$498.00
☐	JU1620	E2607	JAY UNION 16" x 20"	\$498.00
☐	JU1717	E2607	JAY UNION 17" x 17"	\$498.00
☐	JU1816	E2607	JAY UNION 18" x 16"	\$498.00
☐	JU1818	E2607	JAY UNION 18" x 18"	\$498.00
☐	JU1820	E2607	JAY UNION 18" x 20"	\$498.00
☐	JU2018	E2607	JAY UNION 20" x 18"	\$498.00
☐	JU2020	E2607	JAY UNION 20" x 20"	\$498.00
☐	JU2220	E2608	JAY UNION 22" x 20"	\$609.00

JAY Union Reduced Profile Cushion				
☐	JU1416RP	E2607	JAY Union RP 14" x 16"	\$498.00
☐	JU1616RP	E2607	JAY Union RP 16" x 16"	\$498.00
☐	JU1618RP	E2607	JAY Union RP 16" x 18"	\$498.00
☐	JU1620RP	E2607	JAY Union RP 16" x 20"	\$498.00
☐	JU1717RP	E2607	JAY Union RP 17" x 17"	\$498.00

JAY J2 Cushion				
☐	2108-MJ	E2622	JAY J2 14" x 16"	\$502.00
☐	2105-MJ	E2622	JAY J2 15.5" x 16"	\$502.00
☐	2106-MJ	E2622	JAY J2 15.5" x 18"	\$502.00
☐	2107-MJ	E2622	JAY J2 15.5" x 20"	\$530.00
☐	2117-MJ	E2622	JAY J2 17" x 17"	\$502.00
☐	2100-MJ	E2622	JAY J2 18" x 16"	\$502.00
☐	2101	E2622	JAY J2 18" x 18"	\$502.00
☐	2103-MJ	E2622	JAY J2 18" x 20"	\$531.00
☐	2118-MJ	E2622	JAY J2 20" x 18"	\$531.00
☐	2102-MJ	E2622	JAY J2 20" x 20"	\$614.00

JAY J2 Deep Contour Cushion				
☐	2446	E2622	JAY J2 Deep Contour 14" x 16"	\$514.00
☐	2466	E2622	JAY J2 Deep Contour 16" x 16"	\$514.00
☐	2468	E2622	JAY J2 Deep Contour 16" x 18"	\$514.00
☐	2460	E2622	JAY J2 Deep Contour 16" x 20"	\$526.00
☐	2477	E2622	JAY J2 Deep Contour 17" x 17"	\$514.00
☐	2486	E2622	JAY J2 Deep Contour 18" x 16"	\$514.00
☐	2488	E2622	JAY J2 Deep Contour 18" x 18"	\$514.00
☐	2480	E2622	JAY J2 Deep Contour 18" x 20"	\$537.00
☐	2408	E2622	JAY J2 Deep Contour 20" x 18"	\$537.00
☐	2400-MJ	E2622	JAY J2 Deep Contour 20" x 20"	\$649.00
☐	2428	E2623	JAY J2 Deep Contour 22" x 18"	\$649.00

JAY X2 Cushion				
☐	JX21416	E2622	JAY X2 14" x 16"	\$480.00
☐	JX21616	E2622	JAY X2 16" x 16"	\$480.00
☐	JX21618	E2622	JAY X2 16" x 18"	\$480.00
☐	JX21620	E2622	JAY X2 16" x 20"	\$480.00
☐	JX21717	E2622	JAY X2 17" x 17"	\$480.00
☐	JX21816	E2622	JAY X2 18" x 16"	\$480.00
☐	JX21818	E2622	JAY X2 18" x 18"	\$480.00
☐	JX21820	E2622	JAY X2 18" x 20"	\$480.00
☐	JX22020	E2622	JAY X2 20" x 20"	\$480.00

JAY Ion Cushion				
☐	JION1416	E2607	JAY ION 14" x 16"	\$443.00
☐	JION1616	E2607	JAY ION 16" x 16"	\$443.00
☐	JION1618	E2607	JAY ION 16" x 18"	\$443.00
☐	JION1620	E2607	JAY ION 16" x 20"	\$443.00
☐	JION1717	E2607	JAY ION 17" x 17"	\$443.00
☐	JION1816	E2607	JAY ION 18" x 16"	\$443.00
☐	JION1818	E2607	JAY ION 18" x 18"	\$443.00
☐	JION1820	E2607	JAY ION 18" x 20"	\$443.00
☐	JION2018	E2607	JAY ION 20" x 18"	\$443.00
☐	JION2020	E2607	JAY ION 20" x 20"	\$443.00
☐	JION2220	E2608	JAY ION 22" x 20"	\$554.00

☐	JU1816RP	E2607	JAY Union RP 18" x 16"	\$498.00
☐	JU1818RP	E2607	JAY Union RP 18" x 18"	\$498.00
☐	JU1820RP	E2607	JAY Union RP 18" x 20"	\$498.00
☐	JU2018RP	E2607	JAY Union RP 20" x 18"	\$498.00
☐	JU2020RP	E2607	JAY Union RP 20" x 20"	\$498.00