

**Pacific Rehab, Inc**  
 36805 N Nevermind Trail, Box 5406  
 Carefree AZ, 85377-5406

**Meywalk Mk 4**  
**Order Form 2022**

Toll Free: 888.222.9040  
 Fax: 480-575-7907

**Account Information**

Company Name: _____ Address: _____ City, State Zip: _____ Phone: _____ Fax: _____ Email: _____ Purchase Order: _____	Mark For: _____ Ship To: _____ Address: _____ State Zip: _____ Phone: _____
--	---

**Order Information**

Meywalk Mk4 Medium	Meywalk Mk4 Large
Seat Height*: 23.6 - 33" ( Inseam )	Seat Height*: 32" - 40" ( Inseam )
Trunk Support Height*: 39.8 - 52.3" ( Axilla )	Trunk Support Height*: 48" - 59" ( Axilla )
Ring Size**: 41.5" - 47" Circumference ( 105cm )	Ring Size**: 47" - 51" Circumference ( 120cm )
Weight Capacity: 285lbs (130kg)	Weight Capacity: 285lbs
Frame Color: Red	Frame Color: Black


011200  Meywalk MK4 Medium 4,961.00 Ea      011250  Meywalk MK4 Large 5,175.00 Ea

\*To insure correct sizing, please include the required measurements below.  
 \*\*Other ring sizes available: 70, 90, 105, 120, 130cm circumference.

**Mk 4 Accessories**

11147 <input type="checkbox"/> Antitip support, for Medium	295 Pr	10192 <input type="checkbox"/> Drag Brakes	295 Pr
11148 <input type="checkbox"/> Antitip support, for Large	295 Ea	10193 <input type="checkbox"/> Non Reverse Brake	375 Pr
10192 <input type="checkbox"/> Leg Divider	890 Ea	10189 <input type="checkbox"/> Castor Swivel Locks	145 Pr
10193 <input type="checkbox"/> Lateral Leg Guides	925 Pr		
13075 <input type="checkbox"/> Hip Pads	1005 Pr		

**Measuring and sizing the Meywalk MK4**



**\*Required User Information**

Axilla Height: \_\_\_\_\_

Inseam Height: \_\_\_\_\_

User Weight: \_\_\_\_\_

Notes: