



CUSTOMER CONTACT DETAILS

Customer Name			Ship to Customer Address ☐		
Advisor			Shipping Name		
Address			Shipping Address		
State	ZIP	Phone	State	ZIP	Phone
Email / Fax			Email / Fax		

ADDITIONAL INFORMATION

Current Wheelchair	Notes
Customer Weight lbs Customer Height in	
Customer Disability	

RGK CONTACT AUTHORISATION

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK ☐

Signature	Date
-----------	------

RGK INTERNAL

SERIAL NO	PROD DATE
-----------	-----------

GUIDE

☒ Base Product

☐ Required Option

☐ Optional Extra

☐ Critical Information


☐ Information Needed

PRODUCT STYLE

☒ **FXTXZ011** **OCTANE FX**
\$5,493.00 K0005

☐ XA070 Open design

Max User Weight 110kg (242lbs)




CAD APPROVAL

☐ XW010 CAD drawing for approval
\$169.00

KEY MEASUREMENTS


- ☒ XY010 Seat width - Inside of sideguards
Min - 250mm | Max - 460mm
☐
- ☒ XY020 Seat depth
Min - 250mm | Max - 460mm
☐
- ☐ XY030 Extra frame length
Provides additional leg support and hand grip
☐
- ☒ XY050 Front frame height
Min - 480mm | Max - 550mm
☐
- ☒ XY060 Rear frame height
Min - 360mm | Max - 550mm
☐
- ☒ XY070 Occupied frame length
Min = XY020 + XY030 + 140mm | Max - 625mm
☐
- ☒ XY072 Caster position
☐
- ☐ Optimize measurements
- ☐ XA020 Ergonomic seat
Standard = 120mm
☐ Specify length



FRONT FRAME AND FOOTREST

FRONT FRAME OPTIONS


- ☐ XY080 Front frame taper
- ! Specify width
- ☐ Optimize measurements
XY010 Seat width - 100mm



FOOTREST DESIGN (CHOOSE ROUND OR SQUARE)


- ☐ XY100 Round design

Caster Wheel Size	Max Height Achievable	With XB200
3"	110mm	165mm
4"	125mm	175mm
5"	135mm	185mm



- ☐ XY105 Square design
Non certified Sub4 option

Caster Wheel Size	Max Height Achievable	With XB200
3"	135mm	185mm
4"	145mm	195mm
5"	155mm	205mm



- ☐ XB200 High level footrest
\$122.00 Only with XB010, XB080, XB150

- ☐ XB200 High level footrest
\$122.00 Only with XB010, XB080, XB150


FOOTREST MEASUREMENTS

- XY115 Front footrest height

!

- XY120 Rear footrest height
STD = 20mm lower than front height

!



- XY125 Footrest length
STD = 20mm
160mm and over requires heavy duty upgrade XB150

!

FOOTREST STYLE AND COVER

- ☐ FXTXB080 Tubular footrest welded w/ aluminum plate
Supplied with spacer for +20mm adjustment
- ☐ FXTXB120 Clamp-on aluminum plate footrest
Angle / height adjustable
- ☐ FXTXB110 Tubular footrest welded w/ Carbon fiber plate
\$116.00 Supplied with spacer for +20mm adjustment
- ☐ 3k carbon fiber ☐ Spinerger CX6 Cast Carbon Fiber

FOOTREST OPTIONS

- ☐ XB165 Rubber grip tape
\$32.00 3M safety walk

CASTER WHEELS, FORKS AND WHEEL LOCKS

CASTER WHEELS

75mm (3")

- ☐ XC150 Aluminum center 3" x 1.4"
\$111.00 Silver hub | Black urethane tire

100mm (4")

- | | |
|---|--|
| ☐ XC060 Soft roll 4" x 1.3"
Black | ☐ XC110 Aluminum Center 4" x 1"
\$111.00 Silver hub Black urethane tire |
| ☐ XC130 Aluminum center 4" x 1.4"
\$111.00 Choose color below | ☐ XC140 Spinergy CX6 Cast Carbon Fiber soft roll 4" x 1.5"
\$150.00 Black only |
- ☐ Silver hub | Black tire ☐ Green hub | Black tire
☐ Black hub | Black tire ☐ Blue hub | Black tire
☐ Red hub | Black tire ☐ Purple hub | Black tire
☐ Gold hub | Black tire

125mm (5")

- | | |
|---|---|
| ☐ XC070 Soft roll 5" x 1.3"
Black | ☐ XC120 Aluminum Center 5" x 1"
\$111.00 Silver hub Black urethane tire |
| ☐ XC100 Pneumatic aluminum center 5" x 1"
\$124.00 Grey | ☐ XC140 Aluminum center 5" x 1.4"
\$150.00 Silver hub Black urethane tire |

150mm (6")


- | | |
|--|---|
| ☐ XC080 Soft roll 6" x 1.3" | ☐ XC145 Aluminum Center 6" x 1.4"
\$111.00 Silver hub Black urethane tire |
|--|---|

CASTER FORKS

- | | |
|--|--|
| ☐ XC190 Low ride caster fork
STD when no option selected
NOTE: Color follows accent | ☐ XC210 Froglegs suspension caster fork
\$403.00 |
|--|--|
- ☐ Black ☐ Silver

CASTER PIN WIDTH

- ☒ XY330 Caster pin width
XY010 (Seat width) + 70mm



WHEEL LOCK OPTIONS (CHOOSE ONE)

- | | |
|--|---|
| ☐ XD020 Underseat scissor | ☐ XD050 Push to lock |
| ☐ XD030 Compact scissor
\$87.00 | ☐ XD046 Out Front push to lock
\$150.00 |
| ☐ XD035 Compact scissor with lightweight fixing
\$103.00 Only with fixed center of gravity | ☐ XD010 Omit wheel locks |
| ☐ XD045 Out Front Compact scissor
\$150.00 | |

SIDEGUARDS (WITHOUT OR WITH FENDER)

SIDEGUARD WITHOUT FENDER

Aluminum Options

- ☐ XE080 Lift out aluminum
\$353.00
- ☐ FXTXE130 Fold flat aluminum
\$307.00

Carbon Fiber Options

- ☐ XE100 Lift out carbon fiber
\$479.00
- ☐ FXTXE150 Fold flat carbon fiber
\$496.00
- Spinergy CX6 Cast
Carbon Fiber panel
OR
Standard 3k weave
carbon fiber

ABS Plastic Options

- ☐ XE120 Lift out ABS plastic
\$229.00
- ☐ FXTXE170 Fold flat ABS plastic
\$307.00

Armrest Options

- ☐ XF010 Lift out armrest with plastic sideguard
\$395.00
Indicate fixed height height
! mm
Height (Top of seatrail to top of armrest)
- ☐ XF020 Lift out armrest with aluminum sideguard
\$395.00
Indicate fixed height
! mm
Height (Top of seatrail to top of armrest)

SIDEGUARD MEASUREMENTS WITHOUT FENDER

Optimised Measurements

- ☐ 25mm above tire
Follow tire
- ☐ Custom above tire
Follow tire

! mm



- ☐ Level with tire
Follow tire

Custom Measurements

- ☐ Sideguard length

 mm

- ☐ Sideguard height

 mm


SIDEGUARD WITH FENDER

Aluminum Options

- ☐ XE090 Lift out aluminum with fender
\$383.00
- ☐ FXTXE140 Fold flat aluminum with fender
\$337.00

Carbon Fiber Options

- ☐ XE110 Lift out carbon fiber with fender
\$558.00
- ☐ FXTXE160 Fold flat carbon fiber with fender
\$604.00

SIDEGUARD MEASUREMENTS WITH FENDER




Fender Height

- ☐ 10mm above tire
STD
- ☐ Custom above tire
Min - 10mm

! mm

Fender Position

- ☐ Center of tire
STD if no option selected
- ☐ Outside of tire



BACKREST

BACKREST HEIGHT

- ☐ XY230 Fixed height backrest
Max - 400mm
- ① mm
Specify height



- ☐ XY245 Adjustable height backrest
\$116.00
Max - 400mm | +/-25mm
- ① mm
Specify height




BACKREST DESIGN

- ☒ XY360 Backrest angle
Max -10/+10 Deg


ERGONOMIC OPTIONS

- ☒ XY360 Backrest angle
Max -7/+7 Deg
- ① ☐ ☐ ☐ ☐ ☐
-5° -3° 0° +3° +5°
- deg
Custom




- ☐ XG110 Inward taper
\$70.00
Inside of backrest posts
- ① mm
Specify width or per side

- ☐ XG100 Ergonomic backrest
\$70.00
- ① mm
Vertical height
- ① mm
Rearward position



- ☐ XG120 Outward flare
\$70.00
Inside of backrest posts
- ① mm
Specify width or per side



PUSH HANDLE OPTIONS

- ☐ XG070 Screw in push handles
\$136.00
N/A with XY245
- ☐ XG080 Flip down push handles
\$208.00
N/A with XY245

- ☐ XG085 Flip down push handles height adjustable
\$229.00
N/A with XY245
- ☐ XG097 Integrated push handle
\$228.00
Quickie® design

CAMBER BAR (FIXED OR ADJUSTABLE CENTER OF GRAVITY)

ADJUSTABLE CENTER OF GRAVITY

- ☐ XH110 Square back adjustable center of gravity
25-150mm | Max seat width 500mm

mm
Specify measurement



XH110

- ☐ XH120 Round back adjustable center of gravity
60-150mm

mm
Specify measurement




XH120

FIXED CENTER OF GRAVITY

- ☐ FXXH090 V-style fixed center of gravity
\$289.00 25-125mm

mm
Specify measurement



FXAXH090

ADJUSTABLE CENTER OF GRAVITY CAMBER ANGLE

- ☒ XY390 Camber angle

☐ 0° ☐ 1° ☐ 2° ☐ 3° ☐ 5°

FIXED CENTER OF GRAVITY CAMBER ANGLE

- ☒ XY390 Camber angle

☐ 0° ☐ 1° ☐ 2° ☐ 3° ☐ 5° Custom

CAMBER BAR OPTIONS

- ☐ FXTXK227 Standard quick release axle
120 or 130mm

- ☐ FXTXK226 100mm Titanium quick release axle
\$158.00 Fixed Axle only

- ☐ XK220 Titanium 130mm quick release axle (pair) N/A
\$99.00 with XK010, XK070, FXTXK226

- ☐ XK240 Quad release axle nuts (pair)
\$136.00 N/A with FXTXK226

ANTI-TIP OPTIONS

- ☐ XJ010 Left side swing away anti-tip
\$124.00 Only with square back frame, N/A with XH040

- ☐ XJ020 Right side swing away anti-tip
\$124.00 Only with square back frame, N/A with XH040

- ☐ XJ100 Single sport anti-tip
\$383.00 Center mounted



SERIAL NO

PROD DATE

OCTANEFX

REAR WHEELS

REAR WHEEL SIZE

☐ XY420 24"(540mm) Wheels☐ XY430 25"(590mm) Wheels

REAR WHEEL GAP

☒ XY460 Rear wheel gap specify in millimeters

 ☐ 15 ☐ 20 ☐ 25 Custom

XY460



REAR WHEEL OPTIONS

☐ XK010 Fusion 16 spoke
STD option Silver hub | Black spokes☐ XK060 Spinergy SPOX 18 spoke
\$493.00☐ XK320 Spinergy XLX 12 spoke
\$690.00 Cross/radial spoked☐ XK050 Spinergy LX 12 spoke
\$655.00☐ XK070 Spinergy BLXL 18 blade spoke
\$690.00 N/A with silver hub☐ XK330 Spinergy CX6 Cast Carbon Fiber
\$1,404.00 18 spoke Black hub

SPINERGY COLOR OPTIONS

 ☐ Black spokes ☐ White spokes ☐ Yellow spokes
☐ Red spokes ☐ Blue spokes ☐ Pink spokes
☐ Purple spokes ☐ Orange spokes ☐ Green spokes

☐ XK180 Custom spoke design
\$153.00 Specify color

Inside

Outside

SPINERGY ANODIZED OPTIONS

☐ XK190 Black hub
STD option☐ XK200 Silver hub
\$80.00 NA with BLXL / CLX☐ XK201 Red hub / rims
\$162.00 LX, XLX, BLXL only☐ XK205 Custom hub / rims
\$213.00 Specify color☐ XK202 Purple hub / rims
LX, XLX, BLXL only☐ XK203 Blue hub / rims
\$162.00 LX, XLX, BLXL only☐ XK204 Gold hub / rims
\$213.00 LX, XLX, BLXL only

SUPPLY WITHOUT WHEELS

☐ SBTXX040 Supply without wheels
-\$161.00

Wheel type

Tire type

HANDRIM POSITION AND OPTIONS

☐ XY480 Narrow mount handrim☐ XY485 Wide mount handrim☐ XL090 Remove excess mounting tabs
\$89.00 Recommended for small wheel gaps

HANDRIM STYLES

☐ XL010 Aluminum hard anodized 6 tab
\$295.00 STD if no option selected☐ XL050 Stainless steel 6 tab
\$116.00☐ XL140 Ergo gripp 6 tab wide
\$258.00☐ XL020 Tetra grip 6 tab
\$248.00☐ XL100 Para grip 6 tab
\$248.00☐ XL160 Elipse 6 tab
\$295.00 Replacement for Natural-Fit Surge☐ XL040 Titanium 6 tab
\$210.00

TIRE OPTIONS

☐ XM010 Kenda
 ☐ Grey ☐ Blue ☐ Red☐ XM040 Shox solid airless
\$97.00 Black☐ XM060 Schwalbe One lightest option
\$83.00 24", 25" only☐ XM020 Schwalbe Marathon Plus
\$60.00☐ XM050 Kenda K885 mountain bike tire
\$23.00 24", 25" only☐ XM130 Kenda black
\$16.00

INNER TUBE VALVE OPTIONS

☐ XY500 Schraeder inner tube valve
STD if no option selected☐ XY490 Presta inner tube valve

WHEEL ACCESSORIES

☐ XK300 Offroad wheel set 22" rim with handrim, MTB tire and axle pins
\$564.00 Replaces 24" wheel with XM010☐ XK310 Offroad wheel set 24" rim with handrim, MTB tire and axle pins
\$564.00 Replaces 25" wheel with XM010



SERIAL NO

PROD DATE

OCTANE^{FX}

FINISH OPTIONS

FRAME AND COMPONENT FINISH

☐ XN010 Powder coated frame

Frame	Comp.	GLOSS
☐	☐	
☐	☐	White gloss 100
☐	☐	Grey gloss 132
☐	☐	Cream gloss 135
☐	☐	Red gloss 108
☐	☐	Yellow gloss 120
☐	☐	Forest green gloss 143
☐	☐	Green gloss 107
☐	☐	Dark blue gloss 151
☐	☐	Royal blue gloss 151
☐	☐	Blue gloss 174
☐	☐	Duckegg blue gloss 141
☐	☐	Pink gloss 188
☐	☐	Orange gloss 180
☐	☐	Purple gloss 157

☐ FXTN030 \$208.00 Beadblast frame
All parts black unless XN160 chosen☐ FXTN050 \$574.00 Polished frame
All parts black unless XN170 chosen☐ XN150 Powder coated components
Frame color if not specified

Frame	Comp.	MATT
☐	☐	
☐	☐	White matt 101
☐	☐	Grey matt 131
☐	☐	Light grey matt 133
☐	☐	Navy blue matt 152

☐ XN200 \$150.00 Custom paint color
State color / RAL code
Frame Comp.☐ XN200 \$150.00 Custom paint color
State color / RAL code
Frame Comp.

Frame	Comp.	SPECIAL
☐	☐	Chrome mirror 300
☐	☐	Sparkle silver 170
☐	☐	Metallic black matt 104
☐	☐	Metallic grey matt 130
☐	☐	Black sandpaper 201
☐	☐	Grey sandpaper 202
☐	☐	White sandpaper 200

☐ XN160 \$115.00 Beadblast components☐ XN170 \$169.00 Polished components

ACCENT COLOR

Accent color includes forks and decals

☐ Black accent ☐ Red accent ☐ Purple accent
☐ Silver accent ☐ Blue accent☐ XN250 \$150.00 Custom accent color

UPHOLSTERY

SEAT UPHOLSTERY

☐ XP040 Seat sling strap adjustable☐ FXTXP045 \$39.00 Sub4 Seat sling strap adjustable☐ XR050 \$295.00 Carbon fiber seat plate

BACKREST UPHOLSTERY

☐ XP050 Airtech backrest strap adjustable☐ FXTXP035 \$81.00 Airtech Sub4 backrest strap adjustable☐ XP060 \$302.00 Backrest upgrade to full leather
Black☐ XP070 \$105.00 Backrest upgrade to leather corners
Black

BACKREST UPHOLSTERY STITCHING COLOR

☐ Silver ☐ Black ☐ Blue ☐ Orange
☐ Red ☐ Light blue ☐ Purple

BACKREST UPHOLSTERY LOGO COLOR

☐ White - STD ☐ Black ☐ Blue ☐ Custom
☐ Red ☐ Light blue ☐ Purple

BACKREST ACCESSORIES

☐ XP110 \$78.00 Backrest with zip pocket
N/A with Sub4 backrest☐ XP170 \$119.00 Custom backrest embroidery
Provide .jpeg file and positioning



SERIAL NO

PROD DATE

OCTANEFX

CUSHION

CUSHION FOAM OPTIONS

- ☐ XP080
\$73.00
Cushion foam
Select density and thickness below
- ☐ Soft - CMS040 ☐ Medium - CMS050 ☐ Hard - 9lb recon
- ☐ 50mm ☐ 75mm ☐ 100mm

CUSHION COVER UPGRADES

- ☐ XP085
Nylon cushion cover
Black
- ☐ XP090
\$56.00
Airlite cushion cover
Black
- ☐ XP100
\$185.00
Leather cushion cover
Black

JAY® PRODUCTS

- ☐ JAY® Cushion
Price dependent on size / option
- ☐ JAY® Backrest
Price dependent on size / option
- ☐ JAY® Whitmeyer® Headrest
Price dependent on size / option
- ☐ JAY® Order form attached
- ☐ JAY® Order form attached
- ☐ JAY® Order form attached

Notes

ACCESSORIES

- ☐ XR080
\$636.00
Frontwheel assembly
Only with tubular style footrest
- ☐ XP120
\$63.00
Frame protectors Neoprene zip
fasten Black
- ☐ XP420
\$42.00
Frame protectors Nylon velcro fasten
Black
- ☐ XR010
\$33.00
Seat pouch velcro attach
Black
- ☐ XR200
\$56.00
Rucksack with RGK logo
Black
- ☐ XR100
\$293.00
FX carry bag
- ☒ XR090
Tool kit
- ☒ XS060
Calf Strap

