

Ti

Grade 9

**Customer Contact Details**

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

OtherVAT Exemption Declaration (UK Only):

I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.

Signature: _____ Date: _____

RGK Contact Authorisation:

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK: ☐

Signature: _____ Date: _____


Guide:

☒ Base product☐ Required option☐ Optional Extra

Critical Information



Information Needed

PRODUCT STYLE☒ **SBTX2010** **OCTANE SUB4**☐ **XA070** Open design
OPEN**Max User Weight 110kg (242lbs)****CAD APPROVAL**☐ **XW010** CAD drawing for approval
CAD**KEY MEASUREMENTS**☒ **XY010** Seat width (*Inside of sideguards*)**SW**(Max - 440mm)
(Min - 250mm)☒ **XY020** Seat length**SL**(Max - 440mm)
(Min - 250mm)☐ **XY030** Extra seat length (*Provides additional leg support & hand grip*)**ESL** mm


(Max - 440mm total seat)

☒ **XY050** Front frame height**FH**(Max - 520mm)
(Min - 380mm)☒ **XY060** Rear frame height**RH**(Max - 500mm)
(Min - 360mm)☒ **XY070** Frame length**FL**(Min = XY020 + XY030 + 140mm)
(Max - 625mm)☒ **XY072** Castor position**CP** mm☐ Optimise measurements☐ **XA020** Ergonomic seat**ERGO SEAT**

Specify length

 mm

(Standard = 120mm)




FRONT FRAME AND FOOTREST

Front Frame Options

- ☐ XY080 Seat taper design
TAPER

i Specify width
mm


(When combined with XY091, dimension is inbetween the centre of the main bends)



- ☐ XY091 Parallel taper design
PARALLEL

i Specify width


(When combined with XY080, dimensions must be smaller by 20mm)



Footrest Design (Choose one)


- ☐ XY100 Round design
ROUND

Castor Wheel Size	Max Height Achievable
4"	110



- ☐ XY105 Square design
***Non certified Sub4 option**
SQUARE

Castor Wheel Size	Max Height Achievable
4"	130mm



Footrest Measurements

- ☒ XY115 Front footrest height
FH



- ☒ XY120 Rear footrest height
RH




Std = 20 lower than front height

- ☒ XY125 Footrest length
FL



Maximum length 140mm (Certified weight)
(160mm and over requires heavy duty upgrade XB150)



Footrest Style

- ☐ XB010 Tubular aluminium footrest
(Height adjustable)
TUBULAR ALUMINIUM

- ☐ XB020 Tubular titanium footrest
(Height adjustable)
(*Non-certified option*)
TUBULAR TITANIUM

Footrest Cover

- ☐ XB045 Lightweight carbon fibre cover
(With weight saving cut-out)
L/W CARBON



- ☐ Footplate width Sub4 certified
(Std if no option selected)
MAX PLATE WIDTH
- ☐ Footplate width Standard (optimised plate width)
***non certified Sub4 option**
CARBON

CASTOR WHEELS, FORKS & WHEEL LOCKS

Castor Wheels

- ☐ XC135 Aluminium centre 4"x1.25"
Choose colour below:
FROG LEG




- | | |
|---|---|
| <input type="checkbox"/> Silver hub & black tyre (Std)
SILVER/BLK | <input type="checkbox"/> Green hub & black tyre
GREEN/BLK |
| <input type="checkbox"/> Black hub & black tyre
BLK/BLK | <input type="checkbox"/> Red hub & black tyre
RED/BLK |
| <input type="checkbox"/> Purple hub & black tyre
PURPLE/BLK | <input type="checkbox"/> Blue hub & black tyre
BLUE/BLK |

Castor Forks

- ☒ XC190 Low ride castor fork
LOW RIDE
Note: colour follows accent

Castor Pin Width

- ☒ XY330 Castor pin width
CPW
☒ Optimise measurements
(Minimum width for clear castor rotation and maximum wheelbase length)



Wheel Lock Options

- | | |
|---|--|
| <input type="checkbox"/> SBTXD036 Compact scissor with lightweight fixing
COMPACT L/W | <input type="checkbox"/> XD010 Omit wheel locks
OMIT |
|---|--|

SIDEGUARDS

Sideguard without Fender


- ☐ XE055 Screw on lightweight carbon fibre
SCREW ON L/W CARBON

Sideguard Measurements

- ☐ Individual above tyre (Follow tyre - max. 25mm)


 mm

- ☐ Level with tyre (Follow tyre)



Sideguard with Fender

- ☐ XE065 Screw on lightweight carbon fibre w/fender
SCREW ON L/W CARBON FENDER



Sideguard Measurements

- ☐ 10mm above tyre (Standard)





- ☐ Individual above tyre (Min = 10mm)


 mm

- ☒ Centre of tyre



BACKREST			
Backrest Height - Fixed			
<input checked="" type="checkbox"/> XY230	Backrest height B/REST HT Specify height <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div>	<input type="checkbox"/> XG010	Fixed Backrest (non-folding) FIXED
Ergonomic Options			
<input checked="" type="checkbox"/> XY360	Backrest angle ANGLE <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">-5 <input type="checkbox"/></div> <div style="text-align: center;">-3 <input type="checkbox"/></div> <div style="text-align: center;">0 <input type="checkbox"/></div> <div style="text-align: center;">+3 <input type="checkbox"/></div> <div style="text-align: center;">+5 <input type="checkbox"/></div> </div> <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div> <div style="text-align: center; font-size: small;">Individual (Max -10/+20 Deg) Deg</div>	<input type="checkbox"/> XG100	Ergonomic backrest ERGO B/REST Vertical height <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div> mm Rearward position <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div> mm
<input type="checkbox"/> XG110	Inward taper (inside of backrest posts) INWARD Specify width or per side <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div> mm	<input type="checkbox"/> XG120	Outward flare (inside of backrest posts) OUTWARD Specify width or per side <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div> mm
CAMBER BAR			
Camber Angle		Centre of Gravity	
<input checked="" type="checkbox"/> XY390	Camber angle <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">1 <input type="checkbox"/></div> <div style="text-align: center;">2 <input type="checkbox"/></div> <div style="text-align: center;">3 <input type="checkbox"/></div> </div>	<input type="checkbox"/> SBTXH100	Single post fixed centre of gravity (25-125mm) SINGLE POST FIXED Specify length <div style="border: 1px solid black; height: 20px; width: 100%; margin-top: 5px;"></div>
FRAME FINISH OPTIONS			
Frame and Footrest			
<input type="checkbox"/> SBTXN030	Beadblast frame BEADBLAST	<input type="checkbox"/> SBTXN050	Polished frame POLISHED
Accent Colour (Choose One)			
Note: Accent colour includes forks & decals.			
<input type="checkbox"/> Black accent BLK <input type="checkbox"/> Red accent RED <input type="checkbox"/> Purple accent PURPLE	<input type="checkbox"/> Silver accent SILVER <input type="checkbox"/> Blue accent BLUE <input type="checkbox"/> Carbon accent CARBON (Castor caps only, black for rest)	<input type="checkbox"/> XN250 Individual Accent Colour INDIVIDUAL <div style="border-bottom: 1px solid black; width: 100%; height: 20px; margin-top: 5px;"></div>	
UPHOLSTERY			
Seat Upholstery			
<input type="checkbox"/> XP045	Sub4 Seat sling strap adjustable SUB4 STRAP ADJ		
Backrest Upholstery			
<input type="checkbox"/> XP035	Airtech Sub4 backrest strap adjustable AIRTECH SUB4 STRAP ADJ	<input type="checkbox"/> XP050	Airtech backrest strap adjustable *Non Sub4 option AIRTECH STRAP ADJ
Backrest Upholstery Stitching Options			
<input type="checkbox"/> Silver SILVER <input type="checkbox"/> Red RED	<input type="checkbox"/> Black BLACK <input type="checkbox"/> Light blue LIGHT BLUE	<input type="checkbox"/> Blue BLUE <input type="checkbox"/> Purple PURPLE	<input type="checkbox"/> Orange ORANGE
Backrest Upholstery Logo Colour			
<input type="checkbox"/> White (Std) WHITE STD <input type="checkbox"/> Red RED	<input type="checkbox"/> Black BLACK <input type="checkbox"/> Light blue LIGHT BLUE	<input type="checkbox"/> Blue BLUE <input type="checkbox"/> Purple PURPLE	<input type="checkbox"/> Individual: <div style="border-bottom: 1px solid black; width: 100%; height: 20px; margin-top: 5px;"></div>

REAR WHEELS (* non certified weight options)	
Rear Wheel Size	
<input type="checkbox"/> XY420 24" Wheels	<input type="checkbox"/> XY430 25" Wheels
Rear Wheel Gap	
<input checked="" type="checkbox"/> XY460 Rear wheel gap	
Rear Wheel Options	
<input type="checkbox"/> XK075 Spinergy BLXL 18 blade spoke (N/A with silver hub)	<input type="checkbox"/> SBTXK150 Spinergy CLX 18 carbon blade spoke £432,00 (N/A with silver hub)
<input checked="" type="checkbox"/> XK226 Titanium 100mm quick release axle pins	
Spinergy Colour Options: 	
<input type="checkbox"/> Black spokes BLK	<input type="checkbox"/> White spokes WHITE
<input type="checkbox"/> Red spokes RED	<input type="checkbox"/> Blue spokes BLUE
<input type="checkbox"/> Purple spokes PURPLE	<input type="checkbox"/> Orange spokes ORANGE
<input type="checkbox"/> Yellow spokes YELLOW	<input type="checkbox"/> Pink spokes PINK
<input type="checkbox"/> Green spokes GREEN	<input type="checkbox"/> XK180 Individual spoke design (Specify colour)
Inside: _____	
Outside: _____	
Spinergy Anodise Options	
<input type="checkbox"/> XK190 Black hub (std option) BLK	<input type="checkbox"/> XK200 Silver hub (n/a with BLXL/CLX) SILVER
<input type="checkbox"/> XK202 Purple hub/rims (LX, XLX & BLXL only) PURPLE	<input type="checkbox"/> XK203 Blue hub/rims (LX, XLX & BLXL only) BLUE
<input type="checkbox"/> XK201 Red hub/rims (LX, XLX & BLXL only) RED	<input type="checkbox"/> XK204 Gold hub/rims (LX, XLX & BLXL only) GOLD
<input type="checkbox"/> XK205 Individual hub/rims Select colour: _____	
Supply Without Wheels	
<input type="checkbox"/> SBTXK040 Supply without Wheels OMIT	 Wheel type: <input type="text"/>
Tyre type: <input type="text"/>	
Handrim Position / Options	
<input type="checkbox"/> XY480 Narrow set handrim NARROW	<input type="checkbox"/> XY485 Wide set handrim WIDE
<input type="checkbox"/> XL090 Remove excess mounting tabs (recommended for small wheel gaps) REMOVE EXCESS	
Handrim Styles	
<input type="checkbox"/> XL010 Aluminium hard anodised 6 tab HARD ANO 6 TAB	<input type="checkbox"/> XL040 Titanium 6 tab TITANIUM
<input type="checkbox"/> XL160 Elipse 6 tab (replacement for Surge) ELIPSE	<input type="checkbox"/> XL050 Stainless steel 6 tab STAINLESS STEEL
<input type="checkbox"/> XL020 Tetra grip 6 tab TETRA	<input type="checkbox"/> XL140 Ergo gripp 6 tab wide ERGO
<input type="checkbox"/> XL100 Para grip 6 tab PARA	
Tyre Options	
<input type="checkbox"/> XM060 Schwalbe One (Lightest option) ONE	
<input type="checkbox"/> XM010 Kenda KENDA	 <input type="checkbox"/> Grey <input type="checkbox"/> Blue <input type="checkbox"/> Red
<input type="checkbox"/> XM020 Schwalbe Marathon Plus £42,00 MARATHON PLUS	<input type="checkbox"/> XM130 Kenda black KENDA BLK
<input type="checkbox"/> XM040 Shox solid airless SHOX	<input type="checkbox"/> XM050 Kenda K885 mountain bike tyre (24", 25" only) K885 KENDA
Inner Tube Valve Options	
<input type="checkbox"/> XY500 Schraeder inner tube valve (Standard if no option selected)	<input type="checkbox"/> XY490 Presta inner tube valve

* Indicates items that are not part of the certified weight of the Tiga Sub4

CUSHION (non certified weight options)**

Cushion foam options

- ☐ XP080 Cushion foam **CUSHION**
(select density and size below)
- | | | |
|--|--|---|
| <input type="checkbox"/> Soft (CMS040) | <input type="checkbox"/> Medium (CMS050) | <input type="checkbox"/> Hard (9lb recon) |
| <input type="checkbox"/> 50mm | <input type="checkbox"/> 75mm | <input type="checkbox"/> 100mm |

Cushion cover upgrades


- ☐ XP085 Nylon cushion cover
(Black)
NYLON COVER
- ☐ XP090 Airlite cushion cover
(Black)
AIRLITE COVER
- ☐ XP100 Leather cushion cover
(Black)
LEATHER COVER

Jay Products

- | | |
|--|------------------------------------|
| <input type="checkbox"/> JAY Cushion
(Price dependant on size / option)
JAY Cushion | <div>Refer to Jay Order Form</div> |
| <input type="checkbox"/> JAY Backrest
(Price dependant on size / option)
JAY Backrest | <div>Refer to Jay Order Form</div> |
| <input type="checkbox"/> JAY Whitmyer Headrest
(Price dependant on size / option)
JAY Whitmyer Headrest | <div>Refer to Jay Order Form</div> |


FRAME OPTIONS (non certified weight options)**

- ☐ XR060 Welded calf strap mounting loops
(To locate calf strap on front frame)
CALF STRAP LOOPS
Specify height from floor
- ☐ Optimise measurements

**ACCESSORIES (** non certified weight options)**

- ☐ XR080 Frontwheel assembly
(Only with XB010, TGAXB080, XB150)
FRONTWHEEL
- ☐ XP120 Frame protectors Neoprene zip fasten (Black)
FRAME PROTECTORS NEOPRENE
- ☐ XP420 Frame protectors Nylon velcro fasten (Black)
FRAME PROTECTORS NYLON
- ☐ XR010 Seat pouch velcro attach (Black)
SEAT POUCH VELCRO
- ☐ XR200 Rucksack with RGK logo (Black)
RUCKSACK
- ☐ XR020 Underseat large storage bag (Black)
BLACK BOX
- ☒ XR090 Tool kit
- ☒ XS065 Sub4 Calf Strap

** indicates options that, when chosen, will not enable RGK to guarantee a Sub4kg product weight

OTHER (** non certified weight options)			
Footrest			
<input type="checkbox"/> XB200	High level footrest HIGH LEVEL F/REST	Indicate front footrest height: <div style="border: 1px solid black; width: 100px; height: 20px; display: flex; align-items: center; justify-content: center;">mm</div>	
Castor Wheel Size			
<input type="checkbox"/> XC150	Aluminium centre 3"x1.4" (Silver hub and black tyre)	FROG LEG	<input type="checkbox"/> XC145
<input type="checkbox"/> XC140	Aluminium centre 5"x1.4" (Silver hub and black tyre)	FROG LEG	Aluminium centre 6"x1.4" (Silver hub and black tyre)
Wheel Lock Options			
<input type="checkbox"/> XD020	Underseat scissor	SCISSOR	<input type="checkbox"/> XD050
			Push to lock PUSH TO LOCK
Sideguards			
<input type="checkbox"/> XE030	Screw on aluminium, indicate size on page 4 SCREW ON ALUMINIUM		
Backrest Options			
<input type="checkbox"/> SBTXG020	Folding backrest with 2 position snap lock and angle adjust FOLDING B/REST		
Push handle Options			
<input type="checkbox"/> XG070	Screw in push handles SCREW IN PH	<input type="checkbox"/> XG080	Flip down push handles FLIP DOWN PH
Centre of Gravity			
<input type="checkbox"/> SBTXH130	Square back adjustable centre of gravity (25-125mm) SQUARE BACK ADJ COG	<input type="checkbox"/> SBTXH140	Round back adjustable centre of gravity (60-125mm) ROUND BACK ADJ COG
Anti-Tip Options			
<input type="checkbox"/> XJ100	Single sport anti tip SINGLE SPORT (required to change to Adjustable axle)		
Seat Plate			
<input type="checkbox"/> XR050	Carbon fibre seat plate CARBON SEAT PLATE		
Backrest Upholstery Leather Options (only with XP050)			
<input type="checkbox"/> XP070	Backrest upgrade to leather corners (Black) (only with Airtech upholstery XP050) LEATHER CORNERS		
Backrest Accessories (only with XP050)			
<input type="checkbox"/> XP110	Backrest with zip pocket ZIP POCKET	<input type="checkbox"/> XP170	Individual backrest embroidery (Please provide .jpeg file and positioning) EMBROID
** indicates options that, when chosen, will not enable RGK to guarantee a Sub4kg product weight			

NOTES