

ORDER FORM




Q200R with Rehab Seat

EIPW35

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

HCPC	K0848 / K0856
Weights starting at: (with Batteries)	213 lbs
Power Seating	22° Tilt
Seat Type	Rehab Seat w/ Fold Down Back and 12° Angle Adjustment
Seat Size Range	Width 15" to 21" Depth 16" to 20"
Seat-to-Floor Height Range	17" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Rear-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 23", Length 42"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 15.5 miles
Drive Motors	2-pole, 6 mph / 10 kph
Drive Wheels	12" w/ Flat Free Inserts 14" w/ Pneumatic Inserts 14" w/ Flat Free Inserts
Front Casters	9" w/ Pneumatic Inserts 9" w/ Flat Free Inserts
Charger	8A Dual Mode (Compact), 10A
Battery Size	Group 22
Color	7 color accents
Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Rear Wheel Drive Base <i>End Item Number EIPW35 (Rehab Seat)</i>			\$
350PC302	K0848	Quickie Q200 R NPC No Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,345.00
350PC303	K0856	Quickie Q200 R SPC Single Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,345.00
User Weight			\$
350WG23	–	Up to 300 lbs.	NC
Motor Package			\$
350PF107	–	2-pole 6 mph / 10 kph	NC
Drive Wheels			\$
350DW77	–	Drive Wheel, 12", Solid, Knobby, Black	NC
350DW78	–	Drive Wheel, 14", Pneumatic, Knobby, Black	\$268.00
350DW79	–	Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			\$
350CW121	–	Caster Wheels 9", Pneumatic, Black	\$268.00
350CW122	–	Caster Wheels 9", Solid, Black	NC
Color Accents			\$
350EZ20	–	Apple Green	NC
350EZ22	–	Black	NC
350EZ21	–	Blue	NC
350EZ23	–	Orange	NC
350EZ24	–	Red	NC
350EZ19	–	White	NC
350EZ25	–	Silver	NC
ISO Transit Tie Downs			\$
350T050	–	ISO Transit Package	NC
LED Lights			\$
350LI04	–	LED Lights and Indicators	\$1,288.00
350LI02	–	Omit LED Lights	NC
350LI08	–	LED Fender Lights	\$521.00

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Seat Heights <i>(Seat height is measured from front of seat pan to floor with 0° frame angle.)</i>			\$
350HT55	–	Seat Height 17" (N/A with Tilt)	NC
350HT56	–	Seat Height 18"	NC
350HT57	–	Seat Height 19" (N/A with Tilt)	NC
Batteries and Chargers			\$
Group 22 Battery			
350AG31	E2361	GP22 Flat Top Pair w/Straps <i>(N/A w/ Tilt)</i>	\$564.00
350AG35	E2361	GP22/50Ah Flat Top Pair w/Straps	\$626.00
Chargers			
350BJ30	–	10A Dual Mode (Compact)	\$102.00
350BJ19	–	8A Dual Mode (Compact)	NC
SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options			\$
Seat Tilt			
350PA56	E1002	Posterior 22° Tilt	\$7,246.00
Pre Set Seat Frame Angle <i>(No Power Seating Options, Increases Seat to Floor Height)</i>			\$
350FA10	–	0° <i>(Adds 0" to STFH)</i>	NC
350FA3	–	3° <i>(Adds 1" to STFH)</i>	NC
350FA4	–	6° <i>(Adds 2" to STFH) (N/A w/ 19" STFH)</i>	NC
Backrest Frame Type			\$
350BF27	–	Tool Angle Adjustable Back, 90° to 102°, Stroller Handle	NC
350BF30	–	Tool Angle Adjustable Back, 90° to 102°, Trailing Push Handles	NC
Backrest Frame Option			\$
350BF31	–	Fold Down, Quick Release Backrest	NC
Back Cane Height			\$
350B266	–	Tall, 23"	NC
Back Cushion Type			\$
350BT148	–	Lite Sling Back, with Tension Adjustable Straps	NC
350BT135	–	JAY® J3® (Section A)	NC
350BT133	–	Omit Back Cushion	NC

SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
A - JAY® J3 Back			\$
JAY® J3 Back Contour			
350JP01	E2613	J3 PA Back for 14" to 20" Wide	\$575.00
350JP11	E2615	J3 PL Back for 14" to 20" Wide	\$671.00
350JP13	E2620	J3 PD Back for 14" to 20" Wide <i>(18" and 21" shell height only)</i>	\$805.00
JAY® J3 Shell Height			
350JR14	–	16.5" Tall Shell <i>(Upper Thoracic Short)</i>	NC
350JR15	–	18" Tall Shell <i>(Upper Thoracic Med)</i>	NC
350JR16	–	19.5" Tall Shell <i>(Upper Thoracic Tall)</i>	NC
350JR17	–	21" Tall Shell <i>(Shoulder Height)</i>	NC
350JR27	–	22.5" Tall Shell <i>(Shoulder Height)</i>	NC
350JR28	–	24" Tall Shell <i>(Shoulder Height)</i>	NC
JAY® Anterior Thoracic Support			
350K10	E0960	Medium <i>(ATS8302)</i>	\$163.00
350K11	E0960	Large <i>(ATS8303)</i>	\$163.00
350K12	E0960	Extra Large <i>(ATS8304)</i>	\$163.00
Lateral Thoracic Swing-Away Supports			
350JW19	E0956	Lateral Thoracic (LxH) 100 mm x 100 mm / 4" x 4" Pad - Right	\$135.00
350JW25	E0956	Lateral Thoracic (LxH) 100 mm x 100 mm / 4" x 4" Pad - Left	\$135.00
350JW20	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Right	\$135.00
350JW26	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Left	\$135.00
350JW21	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$135.00
350JW27	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$135.00
350JW22	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$135.00
350JW28	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$135.00
350JW23	E0956	Lateral Thoracic (LxH) 152 mm x 127 mm / 6" x 5" Pad - Right	\$135.00
350JW29	E0956	Lateral Thoracic (LxH) 152 mm x 127 mm / 6" x 5" Pad - Left	\$135.00
350JW24	E0956	Lateral Thoracic (LxH) 152 mm x 152 mm / 6" x 6" Pad - Right	\$135.00
350JW35	E0956	Lateral Thoracic (LxH) 152 mm x 152 mm / 6" x 6" Pad - Left	\$135.00
Lateral Thoracic Swing-Away Hardware			
350JW31	E1028	Lateral Thoracic Swing Away Bracket - Right	\$282.00
350JW32	E1028	Lateral Thoracic Swing Away Bracket - Left	\$282.00
Seat Type			\$
350ST80	–	JAY® <i>(Comes with Solid Seat Pan)</i> (See Section B)	NC
350ST34	E2607	SEDEO Cushion <i>(Sizes 16", 18", & 20" Wide)</i> <i>(Comes with Solid Seat Pan)</i>	\$424.00
350ST85	–	Solid Seat Pan	NC
350ST84	–	Seat Sling	NC

SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Headrests *(Requires Pad, Upholstery and Mount)*

\$

Headrest Pad

350HE36	E0955	SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" <i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>	\$173.00
350HE37	E0955	SEDEO Pad Flexible, (L x H) 400 x 110 mm / 16" x 4.5" <i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>	\$279.00
350HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$220.00
350HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$220.00
350HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$230.00
350HE45	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Small <i>(Narrow)</i> <i>(Adjustable from 6" to 14" Wide x 3" Tall)</i>	\$338.00
350HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium <i>(Adjustable from 7" to 16" Wide x 4" Tall)</i>	\$338.00
350HE47	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Large <i>(Standard)</i> <i>(Adjustable from 8" to 18" Wide x 5" Tall)</i>	\$338.00
350HE48	–	OMIT Headrest Pad	NC

Headrest Mount

350MH32	E1028	SEDEO Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO Pad)</i>	\$215.00
350MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$177.00
350MH34	E1028	Whitmyer® Linx ² Headrest Mount with Quick Release	\$215.00
350MH35	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$113.00
350MH17	–	OMIT Headrest Mount	NC

Seat Size

\$

Seat Pan Width *(Seat width is determined by inside of armpads.)*

350R78	–	15" <i>(16" seat pan width)</i>	NC
350R53	–	16" <i>(16" seat pan width)</i>	NC
350R54	–	17" <i>(16" seat pan width)</i>	NC
350R55	–	18" <i>(19" seat pan width)</i>	NC
350R56	–	19" <i>(19" seat pan width)</i>	NC
350R57	–	20" <i>(19" seat pan width)</i>	NC
350R58	–	21" <i>(19" seat pan width)</i>	NC

Seat Pan Depth *(Front of Seat Pan to Back Structure.)*

350RF51	–	16"	NC
350RF45	–	17"	NC
350RF46	–	18"	NC
350RF47	–	19"	NC
350RF48	–	20"	NC

B - JAY® Seat Cushion

\$

350CU59	E2607	JAY® Ion® <i>(Sizes 16" - 21" Wide)</i>	\$424.00
350CU61	E2607	JAY® Union® <i>(Sizes 16" - 21" Wide)</i>	\$478.00
350CU64	E2622	JAY® Fusion® <i>(Sizes 16" - 21" Wide)</i>	\$473.00
350CU104	E2622	JAY® Fusion® with Cryo® Technology	\$579.00

JAY® Bundling

\$

Note: To order any JAY® cushion not listed above please see the JAY® price list.

350JC20	–	JAY® Bundled Cushion	See Price List
		Cushion name or part number	Width
			Depth

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Armrest Structure			\$
350EG42	–	Single Post Removable <i>(Only available with Single Power, SPC)</i>	NC
350EG43	E0973	Single Post Removable	\$161.00
Side Guards (Pair)			\$
350SG49	E0953	Side Guards <i>(Depth Adjustable)</i>	\$188.00
Arm Pads			\$
350AM15	–	2" x 14" Armpad	NC
Swing-Away Legrest (Pair)			\$
350H177	–	SEDEO PRO Legrests, Swing Away, 70° <i>(14" to 21" Leg Length)</i>	NC
350H178	E0990	SEDEO PRO Legrests, Manual Elevating, Swing Away <i>(Includes Padded Calf Pads) (14.5 to 21.5" Leg Length)</i>	\$376.00
350H179	E1010	SEDEO PRO Legrests, Power Elevating, Swing Away <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$2,308.00
350H180	–	SEDEO PRO Hanger Receivers Only	NC
Swing Away Legrest Calf Supports			
350DI41	–	Calf Support, Depth Adjustable, Integral Foam <i>(7" Wide x 10" Tall, N/A 16-17" W)</i>	\$317.00
350DI42	–	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$424.00
Footplates			
350F370	K0040	Angle Adjustable Flip Up Footplates	\$182.00
Accessories			
350F364	K0108	Footplates, Coupled	\$268.00
350DI44	E0951	Heel Loops	\$54.00
350DI40	–	Calf Strap <i>(Only For 70° Swing-Away Legrests)</i>	\$107.00
350KG4	–	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$215.00
SEDEO Lower Extremity Swing Away Support			
350VD01	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$424.00
350VD02	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$424.00
Pelvic Positioning Belts			
350A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
350A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
350A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$81.00
Accessories			\$
Backpack			
350A1	–	Black Utility	\$59.00
Medial Thigh Support (N/A with Sling Seat)			
3508A80	E0957	Medial Thigh Support	\$182.00
3508A81	E2207	Pommel Swing-Away Hardware	\$268.00
Misc. Accessories			
3508A52	–	Ram Mount Phone	\$92.00
3508A53	–	Ram Mount Tablet	\$376.00
3508A54	–	Ram Mount Cup Holder	\$92.00
350A611	E2207	Crutch Holder	\$135.00
3508A58	–	Rear View Mirror	\$92.00
3509A10	–	Utility Hooks	\$56.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
VR-2 Joystick			\$
350J127	–	VR-2 Joystick	NC
R-net Joystick (Remote)			\$
350J121	–	Standard LED Joystick (w/o Screen Display)	NC
350J138	–	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	\$268.00
350J118	–	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,132.00
350J119	–	Omit Remote Control	–\$25.00
Joystick Mounting			\$
350JM34	–	Fixed on Armrest	NC
350JM35	E1028	Parallel Swing-Away	\$317.00
Joystick Mount (Location)			\$
350HM16	–	Right Hand Armrest	NC
350HM17	–	Left Hand Armrest	NC
Display Language Selection R-net			\$
350L1	–	English	NC
350L3	–	French	NC
350L7	–	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick			\$
350JH13	E2323	Ball Handle	\$135.00
350JH14	E2323	Foam Ball Handle	\$135.00
350JH15	E2323	T-Handle	\$135.00
350JH16	E2323	Goalpost	\$135.00
Electronic Programmings			\$
350PR10	–	R-net PC Programmer w/ Dongle	\$564.00
350PR11	–	VR-2 Hand Held Programmer	\$370.00
R-net Power Modules			\$
350ER12	–	Power Module, VR-2 90A, Non-Expandable	NC
350ER13	–	Power Module, R-net, 90A, Expandable	\$1,282.00
350EX15	E2377	Expandable 90A Controller - R-net	\$676.00
350EX16	E2313	Harness for 90A Expandable Controller	\$606.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Seating Control

\$

Powered Seating Control

350E040	–	Omit Thru-Drive <i>(Requires Switch Control and Standard LED Joystick Requires Ctrl+5) (N/A w/ VR-2 electronics)</i>	NC
350E037	E2310	Thru-drive Control for 1 Actuator	\$2,061.00
350E038	E2311	Thru-drive Control for 2+ Actuators	\$3,731.00

Switch Control *(Available with Omit Thru-Drive only) (Requires use of Omni2 or Joystick Jacks)*

350E033	–	Dual Toggle for 1 Actuator	NC
350E034	–	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
350E035	–	Switch-It Mini 1 button for 1 Actuator	NC
350E036	–	Switch-It Mini 2 button for 1 Actuator	NC
350E041	–	Ctrl+5 Switch Box for 2 Actuators or more	NC

Additional Power Seating Control

350QG11	–	Ctrl+5 Switch Box with Free Assignable Functions	\$1,342.00
---------	---	--	------------

Dual Toggle Mounting Location

350TN1	–	Toggle Mount, Right	NC
350TN2	–	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

350EY19	–	Right Armrest	NC
350EY20	–	Left Armrest	NC
350EY21	–	Back Cane, Right	\$161.00
350EY24	–	Back Cane, Left	\$161.00

Ctrl+5 Switch Type

350EY22	–	Push Buttons	NC
350EY23	K0108	Toggle Switches	\$322.00

Attendant Controls

\$

350A05	E2331	Dual Control Module for Attendant	\$1,396.00
350A06	E2331	R-net Color LCD (CJSM1), Attendant Control Only	\$1,396.00

Attendant Mounting Location

\$

350DN5	–	Right Back Cane	NC
350DN6	–	Left Back Cane	NC

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

350PM12	–	Include Assignable Buttons	NC
---------	---	----------------------------	----

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Specialty Input Controls

\$

R-net Omni2

350ON14	–	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$1,906.00
350ON18	–	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	\$416.00

Omni2 Mount (If used w/ Standard Joystick, it must be mounted opposite side)

350DS14	E1028	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Right	\$532.00
350DS15	E1028	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Left	\$532.00

Proportional Hand Input Devices

350EU14	E2312	Switch-It® MicroGuide	\$3,537.00
350EU15	E2312	Switch-It® MicroPilot	\$3,537.00

Hand Input Device Mounting Position (Required)

350EM05	–	Link-It to the Fixed Joystick Mount	\$317.00
350EM06	E1028	Link-It to the Parallel Swing- Away Joystick Mount	\$532.00
350EM08	–	Omit Mount	NC

Hand Input Device Mounting Platforms (Optional)

350MG16	E0950	SunDisk Hand Tray <i>(only available with MicroSeries joysticks)</i>	\$317.00
350MG15	E0950	Bullet Hand Tray	\$317.00

Hand Input Device Mounting Location (Required)

350IM09	–	Right	NC
350IM10	–	Left	NC

Proportional Head Array Input Devices

350EU25	E2328	Switch-It® DualPro Proportional Head Array - Facial Pads	\$6,135.00
350EU26	E2328	Switch-It® DualPro Proportional Head Array - Spot Pads	\$6,135.00

Head Array Mounting (Required)

350IJ13	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$177.00
350IJ14	E1028	Whitmyer® Linx ² Headrest Mount with Quick Release	\$215.00
350IJ15	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$113.00



Q200 R with Rehab Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

