

# **SUPERSTAND HLT**

Pediatric Standing System

Order Form / Price List  
Effective 9-1-2018



*Complete Positioning*



*Easy Loading*



*Portability*



# SUPERSTAND HLT

Pediatric Standing System

Order Form/Price List

Effective 9-1-2018

Client: \_\_\_\_\_ Age: \_\_\_\_\_  
 Height: \_\_\_\_\_ Weight: \_\_\_\_\_ Dx: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Facility: \_\_\_\_\_  
 Therapist: \_\_\_\_\_ Phone: \_\_\_\_\_

PO: \_\_\_\_\_ Date: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_

Standard Features include: Horizontal high-loading surface; Prone, Vertical, Supine or Multi-Position configurable; hip and chest pads; hands free standing angle adjustments; height and depth adj knee system; lock out block and 3" locking casters.

Specifications subject to change without notice.

**Superstand HLT EI**  
 Client Height: Supine 20" - 38" Prone 22" - 40"  
 Base: 19" wide x 28" long. Loading Height 27"  
 Client weight limit 65 pounds

**Superstand HLT Size 1**  
 Client Height: Supine 28" - 48" Prone 30" - 50"  
 Base: 21" wide x 37 3/4" long. Loading Height 28"  
 Client weight limit 125 pounds

**Superstand HLT Size 2**  
 Client Height: Supine 40" - 60" Prone 40" - 62"  
 Base: 25" wide x 45" long. Loading Height 29"  
 Client weight limit 150 pounds

**Folding Base:** Includes quick release folding base, mast mounting hardware, gas cylinder and 4 plastic locking casters.  
 Base Size:  **EI \$ 895.00**  **Size 1 \$ 1,095.00**  **Size 2 \$ 1,295.00**  
 Unit Configuration / Setup:  **Supine**  **Prone**  **Multi-Position**  **Custom**  
 Color:  **Candy Red**  **Candy Blue**  **Candy Purple**  **Neon Orange**  **Neon Pink**  **Shocker Green**  **Black**  **White**  **Razzberry**

**Mast**  
 **EI Mast:** EI Mast Assembly fits 20" - 38" client, includes chest and hip support pad knee supports and solid foot platform.....\$ 1,425.00  
 **Size One Mast:** Assembly fits 28" - 48" client, includes chest and hip support pad knee supports and solid foot platform.....\$ 1,625.00  
 **Size Two Mast:** Assembly fits 40" - 60" client, includes chest and hip support pad knee supports and solid foot platform.....\$ 1,950.00  
 \*\*\*Note: EI Mast fits in all bases, Size One mast fits in Size one and Size Two base, Size Two mast fits in Size Two base only.  
 Upholstery Color:  **Black Vinyl**  **Blue Vinyl**  **Black**  **Red Vinyl**  **Custom**

**Lateral Positioning Supports**

<input type="checkbox"/> <b>Lateral Chest Supports EI:</b> 3 1/2" deep by 3" high.....	\$ 219.00
<input type="checkbox"/> <b>Lateral Chest Supports Size 1:</b> 4 1/2" deep by 4" high.....	\$ 219.00
<input type="checkbox"/> <b>Lateral Chest Supports Size 2:</b> 4 1/2" deep by 6" high.....	\$ 219.00
<input type="checkbox"/> <b>Lateral Hip Supports EI:</b> 3 1/2" deep by 3" high.....	\$ 219.00
<input type="checkbox"/> <b>Lateral Hip Supports Size 1:</b> 4 1/2" deep by 5 1/2" high.....	\$ 199.00
<input type="checkbox"/> <b>Lateral Hip Supports Size 2:</b> 4 1/2" deep by 6 1/2" high.....	\$ 199.00
<input type="checkbox"/> <b>Sliding Adj. Tracks on Size 1 Non-Contoured Chest pad or on Size 2 Standard Chest pad***:</b> .....	\$ 185.00
<input type="checkbox"/> <b>Sliding Adj. Tracks on Hip pad***:</b> Sliding tracks supply maximum adj. for lateral pads for height and angle.....	\$ 185.00

\*\*\*Note: Sliding adjustment tracks not available on EI Unit.

**Chest Pad lateral mounting brackets:**  90° Standard  1" Offset  2" Offset  
**Hip Pad lateral mounting brackets:**  90° Standard  1" Offset  2" Offset

<b>EI Lateral Support Width Adjustment range</b> 90° 6" - 10" ~ 1" offset 4" - 8" ~ 2" offset 3" - 6"	<b>Size 1 Lateral Support Width Adjustment range</b> 90° 8" - 12" ~ 1" offset 6" - 10" ~ 2" offset 4" - 8"	<b>Size 2 Lateral Support Width Adjustment range</b> 90° 10" - 14" ~ 1" offset 9" - 12" ~ 2" offset 6" - 11"
--	---	---

**Straps and Supports**

<input type="checkbox"/> <b>Chest Belt:</b> <input type="checkbox"/> <b>EI</b> <input type="checkbox"/> <b>Size 1</b> <input type="checkbox"/> <b>Size 2</b> Black nylon velcro chest belt.....	\$ 75.00
<input type="checkbox"/> <b>Padded Chest Belt:</b> <input type="checkbox"/> <b>EI</b> <input type="checkbox"/> <b>Size 1</b> <input type="checkbox"/> <b>Size 2</b> Black nylon chest belt with pad.....	\$ 129.00
<input type="checkbox"/> <b>Body Support Vest:</b> <input type="checkbox"/> <b>EI</b> <input type="checkbox"/> <b>Size 1</b> <input type="checkbox"/> <b>Size 2</b> Quick release vest for upper body support.....	\$ 219.00
<input type="checkbox"/> <b>Soft Hip Support:</b> <input type="checkbox"/> <b>EI</b> <input type="checkbox"/> <b>Size 1</b> <input type="checkbox"/> <b>Size 2</b> Soft strap for hip support while in Prone.....	\$ 195.00
<input type="checkbox"/> <b>Single Patient Use Body Support Vest EI size only:</b> Unlined light duty washable vest (denim only).....	\$ 139.00
<input type="checkbox"/> <b>Single Patient Use Soft Hip Support EI size only:</b> Unlined light duty washable prone hip support strap (denim only).....	\$ 145.00

**Harness Colors:**  **Black Denim**  **Blue Denim**  **Black Vinyl**  **Blue Vinyl**  **Red Vinyl**  **Customer supplied**  **Custom**

**Head Support:**

<input type="checkbox"/> <b>Headrest:</b> Height and depth adjustable padded headrest (includes 1/2 Square Tube Mount).....	\$ 295.00
<input type="checkbox"/> <b>1/2" Square Tube Mount:</b> Height and depth adjustable mount for any headrests that use a 1/2" horizontal mounting bar.....	\$ 149.00
<input type="checkbox"/> <b>i2i Flip Back Mount:</b> Flip back mounting hardware for the i2i headrest (hardware only).....	\$ 349.00
<input type="checkbox"/> <b>1/2" Square Tube Flip Back Mount:</b> Flip back mount for any headrest using a 1/2" horizontal mounting bar.....	\$ 349.00
<input type="checkbox"/> <b>Chinrest:</b> Height adjustable chin support.....	\$ 139.00
<input type="checkbox"/> <b>Padded Head Strap:</b> Sheepskin-lined strap with Velcro closure.....	\$ 99.00

**PRIME ENGINEERING** 4202 W. SIERRA MADRE AVE., FRESNO, CA 93722 1-800-827-8263 Fax 800-800-3355  
 A Division of Axiom Industries Inc. E-mail: info@primeengineering.com Web Site: www.primeengineering.com  
 All specifications subject to change without notice Doc # 9-18\_SS\_HLT\_PLOF

**Lower Extremity Positioning**

- Independent Footplates:** Individually height adjustable footplates with heel loops.....\$ 259.00
- Multi-Position Foot System:** Individually adj. for planter-flexion, dorsi-flexion, rotation, height, with heel loops and foot straps.....\$ 475.00
- Size-Adjustable Knee Pads (pair):** Allows knee pad width adjustment from 2 1/2" - 4 1/2" in width (not available in EI).....\$ 365.00
- Contracture Knee System:** Supports contracted Lower Extremitys with a fully adjustable system (on System) .....\$ 799.00
- D-Ring Footstraps:** 1" Velcro D-Ring straps.....\$ 75.00
- Shoe Holders (pair):** Formed ABS shoe holders with ankle and toe straps.....\$ 219.00
  - Size:  Sm (6" long x 4 1/4" wide)  Med (8" long x 4 7/8" wide)  Lrg (10 1/2" long x 6" wide)
- Height Adjustable Pommel:** Assists in maintaining lower extremity position.....\$ 225.00
- Pommel W/ Rigid Hip Pad:** Height adjustable pommel with height and depth adjustable rigid hip pad.....\$ 495.00
- ABduction Leg System:** Includes rotating knees and fully adjustable shoe holders.....\$ 525.00
- ABduction Leg System with Pommel:** Includes rotating knees, fully adj shoe holders and pommel with rigid pad.....\$ 695.00

**Trays**

- Opaque Upper Extremity Positioning Tray:** Multi-Positional tray with height, depth and angle adjustment.....\$ 499.00
- Clear Upper Extremity Positioning Tray:** Multi-Positional tray with height, depth and angle adjustment.....\$ 649.00
- Tray Options:**  Opaque overlay.....\$ 179.00  Clear overlay.....\$ 175.00  Tray Easel.....\$ 249.00  Edging.....\$ 129.00
- Activity Bowl:** Tray cut-out with 7" bowl (clear tray only)..... Add \$ 195.00
- Padded Tray Covers:**  Full size removable padded cover.....\$ 225.00  Wrap around non removable elbow padding.....\$ 229.00

**Casters**

- Steel Caster Upgrade:** High performance easy rolling steel casters that lock for both swivel and roll.....\$ 199.00



**Standard Superstand HLT Support Pad Specifications:**

Specifications subject to change without notice.

**EI:** Chest Pad: 8 1/2" W x 7" H ~ Chest Laterals: 3 1/2" D x 3" H ~ Hip Pad: 8 1/2" W x 5" H ~ Hip Laterals: 3 1/2" D x 3" H ~ Knee Pads: 2" W x 3" H  
**Size 1:** Chest Pad: 10" W x 8 1/2" H ~ Chest Laterals: 4 1/2" D x 4" H ~ Hip Pad: 10" W x 8" H ~ Hip Laterals: 4 1/2" D x 5 1/2" H ~ Knee Pads: 4" W x 3 1/2" H  
**Size 2:** Chest Pad: 12 1/2" W x 15" H ~ Chest Laterals: 4 1/2" D x 6" H ~ Hip Pad: 12 1/2" W x 9" H ~ Hip Laterals: 4 1/2" D x 6 1/2" H ~ Knee Pads: 4 1/4" W x 5" H  
**Headrest Specifications:** **EI** and **Size 1** 6 1/2" W x 3 1/2" H      **Size 2** 8" W x 4 1/2" H

**Custom Configurations:** We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. We encourage you to call us when you need something special for your client. Supply the information in the following area's to us know your needs and we will assist in every way that we can.


**Client Measurements:**

A: Thorax Width \_\_\_"    B: Hip Width \_\_\_"    C: Hip-to-Knee \_\_\_"    D: Knee-to-Heel \_\_\_"    E: Shoe Size \_\_\_"    F: Shoulder Height \_\_\_"    G: Height \_\_\_"  
 H: Inseam \_\_\_"    I: Weight \_\_\_ lbs.    **Note:** All Measurements Assume Symmetrical Standing and/or Sitting Posture.


**Custom Pads:** Use this area for custom-sized pads on your Superstand HLT. All custom pads must be quoted for price and manufacturing time.

Chest Pad: \_\_\_" Wide x \_\_\_" High  Standard Shape Pad with Top Contour  Custom Shape, describe below ~ Chest Laterals \_\_\_" Wide x \_\_\_" High  
 Hip Pad: \_\_\_" Wide x \_\_\_" High  Standard Shape Pad  Custom Shape, describe below ~ Hip Laterals \_\_\_" Wide x \_\_\_" High

**Custom Tray Surface:** \_\_\_" Wide X \_\_\_" Deep Standard Tray is 22 1/2"w x 19"d for size one and 28"w x 24"d for size two.



Include any information for any tray cutouts or modifications desired in this area or submit a paper or cardboard template. Other information may be added in the area below.



Special Requests / Additional Information: \_\_\_\_\_

# SUPERSTAND HLT MPS Pre Configured Packages\*\*\*

Effective 9-1-2018

## Multi-Position All-In-One Standing System Pre-configured units.

These Fully Configured Units supply all of the components needed for a complete standing program with many clients with the ease of a single order item with minor modifications available.

These units include: High-loading platform ht; Folding base and components to support Prone, Vertical, Supine or Multi-Position configurations. Packages include: Hip pad with laterals, Chest pad with laterals, Support vest; Soft hip pad; Headrest; Chinrest; Multi-Position foot system; Opaque tray; Hands free standing angle adjustments; Height and depth adjustable knee pads; Lock out block and Upgraded steel locking casters. This packaged system can be accessorized and quoted to fit your client's needs with the items listed on this page.



Systems pictured above include optional clear tray top.

<input type="checkbox"/> <b>Superstand HLT EI</b> \$4,695.00 Client Height: Supine 20" - 38" Prone 22" - 40" Base: 19" wide x 28" long. Loading Height 27" Client weight limit 65 pounds	<input type="checkbox"/> <b>Superstand HLT Size 1</b> \$4,995.00 Client Height: Supine 28" - 48" Prone 30" - 50" Base: 21" wide x 37 3/4" long. Loading Height 28" Client weight limit 125 pounds	<input type="checkbox"/> <b>Superstand HLT Size 2</b> \$5,595.00 Client Height: Supine 40" - 60" Prone 40" - 62" Base: 25" wide x 45" long. Loading Height 29" Client weight limit 150 pounds
---	--	--

\*\*\*Note, Superstand HLT packages are configured to be a complete system package, configurable items are listed on this page only. Superstand HLT packages are not available for Custom Configuration. When more configurations or custom items are needed use the Build To Suit section of this order form.  
 Package specifications subject to change without notice.

**Frame Color:**  Candy Red  Candy Blue  Candy Purple  Neon Orange  Neon Pink  Shocker Green  Black  White  Razzberry  
**Upholstery Color:**  Red  Blue  Black **Harness Color:**  Black Denim  Blue Denim  Black Vinyl  Blue Vinyl  Red Vinyl

**Lateral Support Positioning Options for pre-configured packages.** Lateral Chest and Hip Supports are included in pre-configured packages.

<input type="checkbox"/> Sliding Adj. Tracks on Size 1 Non-Contoured Chest pad or on Size 2 Standard Chest pad:.....	\$ 185.00
<input type="checkbox"/> Sliding Adj. Tracks on Hip pad: Sliding tracks supply maximum adj. for hip lateral pads for height and angle .....	\$ 185.00

\*\*\*Note: Sliding adjustment tracks not available on EI Unit.

<b>Chest Pad lateral mounting brackets:</b>	<input type="checkbox"/> 90° Standard	<input type="checkbox"/> 1" Offset	<input type="checkbox"/> 2" Offset
<b>Hip Pad lateral mounting brackets:</b>	<input type="checkbox"/> 90° Standard	<input type="checkbox"/> 1" Offset	<input type="checkbox"/> 2" Offset

<b>EI Lateral Support Width Adjustment range</b> 90° 6" -10" - 1" offset 4" - 8" - 2" offset 3" - 6"	<b>Size 1 Lateral Support Width Adjustment range</b> 90° 8" - 12" - 1" offset 6" - 10" - 2" offset 4" - 8"	<b>Size 2 Lateral Support Width Adjustment range</b> 90° 10" - 14" - 1" offset 9" - 12" - 2" offset 6" - 11"
---	---	---

**Straps and Support Options for pre-configured units:** Support vest and a soft hip strap are already included in pre-configured packages.

<input type="checkbox"/> Chest Belt: Black nylon velcro chest belt.....	\$ 75.00
<input type="checkbox"/> Padded Chest Belt: Black nylon chest belt with pad.....	\$ 129.00
<input type="checkbox"/> Single Patient Use Body Support Vest EI size only: Unlined light duty washable vest (denim only).....	\$ -45.00
<input type="checkbox"/> Single Patient Use Soft Hip Support EI size only: Unlined light duty washable prone hip support strap (denim only).....	\$ -35.00

**Head Support Options for pre-configured units:** Pre-configured packages include a standard headrest and a chinrest, EI unit does not include a chinrest.

<input type="checkbox"/> 1/2" Square Tube Mount only: *Omits the std headrest, ht and depth adj. mount for any 1/2" tube mounting bar.....	\$ -75.00
<input type="checkbox"/> i2i Flip Back Mount: Flip back mounting hardware for the i2i headrest (hardware only) *Omits the standard headrest.....	\$ 279.00
<input type="checkbox"/> 1/2" Square Tube Flip Back Mount: Flip back mount for any 1/2" horizontal bar *Omits the std. Headrest.....	\$ 199.00

**Lower Extremity Positioning Options for pre-configured units:**

<input type="checkbox"/> Size-Adjustable Knee Pads (pair): Allows knee pad width adj. from 2 1/2" - 4 1/2" in width (not available in EI).....	\$ 365.00
<input type="checkbox"/> D-Ring Footstraps: 1" Velcro D-Ring straps.....	\$ 75.00
<input type="checkbox"/> Shoe Holders (pair): Formed ABS shoe holders with ankle and toe straps.....	\$ 219.00
Size: <input type="checkbox"/> Sm (6" long x 4 1/4" wide) <input type="checkbox"/> Med (8" long x 4 7/8" wide) <input type="checkbox"/> Lrg (10 1/2" long x 6" wide)	
<input type="checkbox"/> Height Adjustable Pommel: Assists in maintaining lower extremity position.....	\$ 245.00
<input type="checkbox"/> Pommel W/ Rigid Hip Pad: Height adjustable pommel with height and depth adjustable rigid hip pad.....	\$ 495.00
<input type="checkbox"/> ABduction Leg System: Includes rotating knees and fully adjustable shoe holders.....	\$ 525.00
<input type="checkbox"/> ABduction Leg System with Pommel: Includes rotating knees, fully adj shoe holders and pommel with rigid pad....	\$ 695.00

**Tray Options for pre-configured units:**

<input type="checkbox"/> Opaque Upper Extremity Positioning Tray: Multi-Positional tray with height, depth and angle adjustment.....	\$ Included
<input type="checkbox"/> Clear Upper Extremity Positioning Tray Upgrade: Changes tray top from opaque to clear acrylic tray top.....	\$ 149.00
Tray Options: <input type="checkbox"/> Opaque overlay.....	\$ 179.00
<input type="checkbox"/> Clear overlay.....	\$ 179.00
<input type="checkbox"/> Tray Easel.....	\$ 249.00
<input type="checkbox"/> Edging.....	\$ 129.00
<input type="checkbox"/> Activity Bowl: Tray cut-out with 7" bowl added to center of tray.....	Add \$ 195.00
Padded Covers: <input type="checkbox"/> Full size removable padded cover.....	\$ 225.00
<input type="checkbox"/> Wrap around non removable elbow padding....	\$ 229.00