


Superlight Handrim



Superlight Handrim with Micro Grip



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel. <p>
Not available with Heavy Duty Option

Item Number	Part Number	Description	Quantity	Retail
2b, 3	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy, Topolino and MAXX Spoke wheels</i>	1	\$15.45
1a	101161	Handrim Superlight Aluminum Anodized 24"	1	\$91.67
1b	505585	Handrim Superlight Aluminum Anodized 24" w/ Micro Grip	1	\$128.38
1c	101160	Handrim Superlight Aluminum Anodized 25"	1	\$91.67
1d	505586	Handrim Superlight Aluminum Anodized 25" w/ Micro Grip	1	\$129.49
2a	100536	M5 x 12 BTN SHCS ZC <i>Used with Superlight Spoke</i>	6	\$2.06
2b	100669	M5 x 10 BTN SHCS BLKZC w/Patch	6	\$1.55

Item Number	Part Number	Description	Quantity	Retail
		<i>Used with Spinergy, Topolino and MAXX Spoke wheels</i>		
3	100657	M5 Nylock Nut Blk Zn <i>Used with Spinergy, Topolino and MAXX Spoke wheels</i>	6	\$0.52
2a	502170	Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Spoke</i>	1	\$15.45

Handrim Hardware Chart

Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer	Screw	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835	100206		101106					
20" Ki Spoke	200530	6				200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349		100577	200547	100793	200202		
25" Ki Spoke	200533					200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460	
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963					
22" 5-Spoke X Core	101962						101898		101964				
24" 5-Spoke X Core	100960						100975		100976	200546			
24" 5-Spoke X Core	100960						101893				100768	100769	
20" Spinerger Spox	See Spinerger Spox Page	6			Screw: 100669 Nut: 100657	103125		103179					
22" Spinerger Spox						100827		100808		100889	100888		
24" Spinerger Spox						100766	101161	100615		100830	100828		
25" Spinerger Spox						100767	101160	100765		101464	101460		
26" Spinerger Spox						101477		101148		200200	100950		
22" Spinerger LX	See Spinerger LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888		
24" Spinerger LX						100766	101161	100615		100830	100828	200213	
25" Spinerger LX						100767	101160	100765		101464	101460	200208	
26" Spinerger LX						101477		101148		200200	100950		
18" MAXX Spoke	109285	6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835 Superlight: Screw: 100669 Nut: 100657	100206		101106					
20" MAXX Spoke	108244					200536		200542					
22" MAXX Spoke	105135					100560		100576	100569	200538	200201		
24" MAXX Spoke	105136					200349	101161	100577	200547	100793	200202		
25" MAXX Spoke	107436					200350		101870	200548	200539	200540		
26" MAXX Spoke	107437					200351		100578	200549	100907	101454		