


# Focus CR Caster Forks and Stems



Item Number	Part Number	Description	Quantity	Retail
4a, 5, 6, 13	503816	Fork Assembly 5" Std Stem Focus/Liberty <i>Used with standard weight capacity</i>	1	\$57.50
4a, 5, 7, 13	503814	Fork Assembly 6" Std Stem Focus/Liberty <i>Used with standard weight capacity</i>	1	\$57.50
4a, 5, 8, 13	503817	Fork Assembly 6" Wide Std Stem Focus/Liberty <i>Used with standard weight capacity</i>	1	\$80.53
4b, 5, 8, 13, 14	506368	Fork Assembly 6" Wide Heavy Duty Std Stem Focus <i>Used with heavy duty weight capacity</i>	1	\$208.82
4a, 5, 9, 13	503815	Fork Assembly 7" Std Stem Focus/Liberty <i>Used with standard weight capacity</i>	1	\$57.50

Item Number	Part Number	Description	Quantity	Retail
4a, 5, 10, 13	503818	Fork Assembly 7" Wide Std Stem Focus/Liberty <i>Used with standard weight capacity</i>	1	\$80.53
4b, 5, 10, 13, 14	503122	Fork Assembly 7" Wide Heavy Duty Std Stem Focus <i>Used with heavy duty weight capacity</i>	1	\$208.82
1, 3, 5, 6, 13	506369	Fork Assembly 5" Tall Stem Focus <i>Used with standard and heavy duty weight capacity</i>	1	\$121.87
1, 3, 4b, 5, 13	506370	Fork Assembly 6" Tall Stem Focus <i>Used with standard and heavy duty weight capacity</i>	1	\$121.87
1, 3, 5, 8, 13	506371	Fork Assembly 6" Wide Tall Stem Focus <i>Used with standard and heavy duty weight capacity</i>	1	\$141.64
1, 3, 5, 9, 13	506372	Fork Assembly 7" Tall Stem Focus <i>Used with standard and heavy duty weight capacity</i>	1	\$121.87
1, 3, 5, 10, 13	506373	Fork Assembly 7" Wide Tall Stem Focus <i>Used with standard and heavy duty weight capacity</i>	1	\$141.64
2, 11a, 12a	502653	Caster Attachment Hardware Kit CAT/FCS	1	\$7.73
2, 11b, 12b	502654	Caster Attachment Hardware Kit Wide Fork CAT/FCS	1	\$16.95
1	002476	Spacer, 3/4" OD x 1/2" ID x 1.23	1	\$2.10
2	100659	M8 Nylock Nut Blk Zn	1	\$1.03
1, 3, 5, 13	506374	Caster Stem Tall, Adult Tilt <i>Used with standard and heavy duty weight capacity</i>	1	\$26.40
4a, 5, 13	503819	Caster Stem Standard, Focus <i>Used with standard weight capacity</i>	1	\$17.47
4b, 5, 13	503820	Caster Stem Standard HD, Focus <i>Used with heavy duty weight capacity</i>	1	\$106.28
5	100529	M12 x 1.25mm Lock Nut Black	2	\$1.31
5, 6	503823	Caster Fork 5" Focus	1	\$44.07
5, 7	503824	Caster Fork 6" Focus	1	\$44.07
5, 8	503825	Caster Fork 6" Wide Focus	1	\$63.84
5, 9	503826	Caster Fork 7" Focus	1	\$44.07
5, 10	503827	Caster Fork 7" Wide Focus	1	\$63.84
11a	101829	M8 x 65 HHCS Blk Zn	1	\$3.09
11b	101831	M8 x 80 HHCS Blk Zn <i>Used with wide forks</i>	1	\$3.01

Item Number	Part Number	Description	Quantity	Retail
12a	100792	Caster Spacer D8 x D12 x 10mm <i>Used with all casters except 6x2 and 8x2</i>	2	\$2.10
12b	001246	Caster Spacer Black 6 x 2 Poly .472ODX.319IDX.370LG <i>Only used with 6x2 casters</i>	2	\$5.73
13	100682	Washer 1/2" .518 x .875 x .047 F/W Black Zinc	3	\$1.03
14	002593	Caster Stem Washer HD Tilt <i>Used with standard stem in combination with heavy duty option; Use qty 3 of 002593 and qty 1 of 100682</i>	3	\$18.89