


# **(Discontinued 1) Focus CR Stroller Handle Backrest**



*This Focus CR Stroller Handle Backrest was effective prior to October 2015.*

Item Number	Part Number	Description	Quantity	Retail
1, 4, 8, 14 - 16	501094	Backrest Hardware Kit Std Focus	1	\$165.85
1, 8, 14 - 16, 18, 19	501406	Backrest Hardware Kit, Std Non Folding, Pair Focus	1	\$321.72
1	001071	Back Plate Std TILT	2	\$83.11
4	001075	Stop Pin Backrest Tilt	1	\$6.61

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>	<b>Retail</b>
8	000055	Saddle 1"	4	\$1.55
9	100804	Tube End Plug 1"	1	\$1.55
14	101781	M6 x 50mm Partial Thread BHCS Blk ZC	4	\$1.55
15	100658	M6 Nylock Nut Blk Zn	4	\$1.03
16	100746	M6 Flat Washer Blk Zn	4	\$1.55
18	101815	M8 x 50 BHCS Blk Zn	1	\$3.61
19	100659	M8 Nylock Nut Blk Zn	1	\$1.03