


Catalyst Seat Upholstery



SCREW (ITEM 4) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	10
16" - 17"	12
18" - 19"	14
20"	16

The seat upholstery width measures approximately 2.5" less than the chair size

Item Number	Part Number	Description	Quantity	Retail
1a	100093	Seat Upholstery Cat 1414	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1b	100094	Seat Upholstery Cat 1415	1	\$92.70
1c	100095	Seat Upholstery Cat 1416	1	\$92.70
1d	100096	Seat Upholstery Cat 1417	1	\$92.70
1e	100097	Seat Upholstery Cat 1418	1	\$92.70
1f	100896	Seat Upholstery Cat 1419	1	\$92.70
1g	100136	Seat Upholstery Cat 1420	1	\$92.70
1h	100098	Seat Upholstery Cat 1514	1	\$92.70
1i	100099	Seat Upholstery Cat 1515	1	\$92.70
1j	100100	Seat Upholstery Cat 1516	1	\$92.70
1k	100101	Seat Upholstery Cat 1517	1	\$92.70
1l	100102	Seat Upholstery Cat 1518	1	\$92.70
1m	100103	Seat Upholstery Cat 1519	1	\$92.70
1n	100104	Seat Upholstery Cat 1520	1	\$92.70
1o	100105	Seat Upholstery Cat 1614	1	\$92.70
1p	100106	Seat Upholstery Cat 1615	1	\$92.70
1q	100107	Seat Upholstery Cat 1616	1	\$92.70
1r	100108	Seat Upholstery Cat 1617	1	\$92.70
1s	100109	Seat Upholstery Cat 1618	1	\$92.70
1t	100110	Seat Upholstery Cat 1619	1	\$92.70
1u	100111	Seat Upholstery Cat 1620	1	\$92.70
1v	100112	Seat Upholstery Cat 1714	1	\$92.70
1w	100113	Seat Upholstery Cat 1715	1	\$92.70
1x	100114	Seat Upholstery Cat 1716	1	\$92.70
1y	100115	Seat Upholstery Cat 1717	1	\$92.70
1z	100116	Seat Upholstery Cat 1718	1	\$92.70
1aa	100117	Seat Upholstery Cat 1719	1	\$92.70
1ab	100118	Seat Upholstery Cat 1720	1	\$92.70
1ac	100119	Seat Upholstery Cat 1814	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1ad	100120	Seat Upholstery Cat 1815	1	\$92.70
1ae	100121	Seat Upholstery Cat 1816	1	\$92.70
1af	100122	Seat Upholstery Cat 1817	1	\$92.70
1ag	100123	Seat Upholstery Cat 1818	1	\$92.70
1ah	100124	Seat Upholstery Cat 1819	1	\$92.70
1ai	100125	Seat Upholstery Cat 1820	1	\$92.70
1aj	101558	Seat Upholstery Cat 1914	1	\$92.70
1ak	100906	Seat Upholstery Cat 1915	1	\$92.70
1al	100126	Seat Upholstery Cat 1916	1	\$92.70
1am	100127	Seat Upholstery Cat 1917	1	\$92.70
1an	100128	Seat Upholstery Cat 1918	1	\$92.70
1ao	100129	Seat Upholstery Cat 1919	1	\$92.70
1ap	100130	Seat Upholstery Cat 1920	1	\$92.70
1aq	101559	Seat Upholstery Cat 2014	1	\$92.70
1ar	101557	Seat Upholstery Cat 2015	1	\$92.70
1as	100131	Seat Upholstery Cat 2016	1	\$92.70
1at	100132	Seat Upholstery Cat 2017	1	\$92.70
1au	100133	Seat Upholstery Cat 2018	1	\$92.70
1av	100134	Seat Upholstery Cat 2019	1	\$92.70
1aw	100135	Seat Upholstery Cat 2020	1	\$92.70
1ax	102526	Seat Upholstery Cat 2114	1	\$92.70
1ay	102527	Seat Upholstery Cat 2115	1	\$92.70
1az	101450	Seat Upholstery Cat 2116	1	\$92.70
1ba	101554	Seat Upholstery Cat 2117	1	\$92.70
1bb	100855	Seat Upholstery Cat 2118	1	\$103.00
1bc	100856	Seat Upholstery Cat 2119	1	\$103.00
1bd	100857	Seat Upholstery Cat 2120	1	\$103.00
1be	200352	Seat Upholstery Cat 2214	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1bf	200353	Seat Upholstery Cat 2215	1	\$92.70
1bg	101556	Seat Upholstery Cat 2216	1	\$92.70
1bh	101555	Seat Upholstery Cat 2217	1	\$92.70
1bi	100247	Seat Upholstery Cat 2218	1	\$103.00
1bj	100858	Seat Upholstery Cat 2219	1	\$103.00
1bk	100261	Seat Upholstery Cat 2220	1	\$103.00
2a	100230	Seat Rail 14" Deep CAT	2	\$16.73
2b	100231	Seat Rail 15" Deep CAT	2	\$16.73
2c	100232	Seat Rail 16" Deep CAT	2	\$16.73
2d	100233	Seat Rail 17" Deep CAT	2	\$16.73
2e	100234	Seat Rail 18" Deep CAT	2	\$16.73
2f	100235	Seat Rail 19" Deep CAT	2	\$16.73
2g	100236	Seat Rail 20" Deep CAT	2	\$16.73
3a	000059	Lifting Strap 14" Wide CAT	1	\$6.26
3b	000060	Lifting Strap 15" Wide CAT	1	\$6.26
3c	000061	Lifting Strap 16" Wide CAT	1	\$6.26
3d	000062	Lifting Strap 17" Wide CAT	1	\$6.26
3e	000063	Lifting Strap 18" Wide CAT	1	\$6.26
3f	000064	Lifting Strap 19" Wide CAT	1	\$6.26
3g	000058	Lifting Strap 20" Wide CAT	1	\$6.26
3h	100248	Lifting Strap 21" Wide CAT	1	\$6.27
3i	100875	Lifting Strap 22" Wide CAT	1	\$6.27
4	100734	M5 x 16 PHMS ZC	16	\$1.55