


(Discontinued) Catalyst 5Vx Vertical Axle Plate



Item Number	Part Number	Description	Quantity	Retail
2	100756	Axle Receiver 1/2"	1	\$17.33
3	100583	Axle Plate Washer 18mm	1	\$1.03
4	100514	Axle Receiver Nut M18 x D12 x 10mm	2	\$1.51
5	100674	M6 x 45 CLS10.9 DIN931 HHCS ZC	2	\$1.55
6	100746	M6 Flat Washer Blk Zn	2	\$1.55
7	100658	M6 Nylock Nut Blk Zn	2	\$1.03
8	100860	M6 x 35 HHCS P/THD ZC	2	\$1.55