


# Focus CR Wheel Lock



Item Number	Part Number	Description	Quantity	Retail
1a, 2b, 3a, 4 - 10	501486	Wheel Lock Assembly Push Left Focus	1	\$47.71
1b, 2a, 3a, 4 - 10	501487	Wheel Lock Assembly Push Right Focus	1	\$47.71
1a, 1b, 2a, 2b, 3a, 4 - 10	501113	Wheel Lock Assembly Push PAIR Focus	1	\$95.41
2b, 3a, 4 - 10, 11a, 12-15	501488	Wheel Lock Assembly Push w/ Ext Left Focus	1	\$79.87
2a, 3a, 4 - 10, 11b, 12-15	501489	Wheel Lock Assembly Push w/ Ext Right Focus	1	\$79.87

Item Number	Part Number	Description	Quantity	Retail
2a, 2b, 3a, 4 - 10, 11a, 11b, 12-15	501114	Wheel Lock Assembly Push w/ Ext PAIR Focus	1	\$159.73
1c, 2b, 3b, 4 - 10	501490	Wheel Lock Assembly Pull Left Focus	1	\$47.71
1d, 2a, 3b, 4 - 10	501491	Wheel Lock Assembly Pull Right Focus	1	\$47.71
1c, 1d, 2a, 2b, 3b, 4 - 10	501115	Wheel Lock Assembly Pull PAIR Focus	1	\$95.41
2b, 3b, 4 - 10, 11c, 12-15	501492	Wheel Lock Assembly Pull w/ Ext Left Focus	1	\$79.87
2a, 3b, 4 - 10, 11d, 12-15	501493	Wheel Lock Assembly Pull w/ Ext Right Focus	1	\$79.87
2a, 2b, 3b, 4 - 10, 11c, 11d, 12-15	501116	Wheel Lock Assembly Pull w/ Ext PAIR Focus	1	\$159.73
2a	001138	Plate Wheel Lock Right Tilt	1	\$16.00
2b	001139	Plate Wheel Lock Left Tilt	1	\$16.00
3a	003769	Wheel Lock Mount Push Tilt <i>Replaced 001140</i>	1	\$16.48
3b	003770	Wheel Lock Mount Pull Tilt <i>Replaced 001238</i>	1	\$31.89
4	001141	Weld Nut NDM 06028 Zinc	1	\$2.58
5	101824	M6 x 30 HHCS Blk Zn	1	\$2.06
6	101823	M6 x 20 HHCS Blk Zn	1	\$1.03
7	100658	M6 Nylock Nut Blk Zn	2	\$1.03
8	101853	M6 Washer Blk Zn	1	\$0.52
9	100746	M6 Flat Washer Blk Zn	1	\$1.55
10	101835	M6 Lock Washer Blk Zn	1	\$0.52
12-15	100854	Brake Extension Replacement	1	\$11.95
13	101117	M5 x 22 BTN SHCS BLK ZC	1	\$1.55
14	101856	Flat Washer, M5, Blk Zn	1	\$0.52
15	100657	M5 Nylock Nut Blk Zn	1	\$0.52