


Catalyst Reclining Backrest



Item Number	Part Number	Description	Quantity	Retail
1a, 5, 7, 8b, 9, 10a, 11, 24, 25, 39, 40	503327	Recline Back Tube Subassembly, Short 22", Right, Catalyst	1	\$139.44
1a, 5, 7, 8a, 9, 10b, 11, 24, 25, 39, 40	503328	Recline Back Tube Subassembly, Short 22", Left, Catalyst	1	\$139.44
1b, 5, 7, 8b, 9, 10a, 11, 24, 25, 39, 40	503329	Recline Back Tube Subassembly, Tall 26", Right, Catalyst	1	\$141.96
1b, 5, 7, 8a, 9, 10b, 11, 24, 25, 39, 40	503330	Recline Back Tube Subassembly, Tall 26", Left, Catalyst	1	\$141.96

Item Number	Part Number	Description	Quantity	Retail
9, 10a, 10b, 11, 12a, 30, 39, 40	503331	Recline Removable Stroller Handle Assembly 14 - 17W Catalyst	1	\$119.11
9, 10a, 10b, 11, 12b, 30, 39, 40	503332	Recline Removable Stroller Handle Assembly 18 - 22W Catalyst	1	\$119.11
1a, 2, 3a, 3b, 4a, 4b, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12a, 13-15, 16a, 17a, 17b, 18- 40	503333	Recline Back Kit, Short Height, 14-17" Width, Standard Gas Spring PAIR, Catalyst	1	\$1,127.85
1a, 2, 3a, 3b, 4a, 4b, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12b, 13-15, 16a, 17a, 17b, 18- 40	503334	Recline Back Kit, Short Height, 18-22" Width, Standard Gas Spring PAIR, Catalyst	1	\$1,127.85
1b, 2, 3a, 3b, 4a, 4b, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12a, 13-15, 16b, 17a, 17b, 18- 40	503335	Recline Back Kit, Tall Height, 14-17" Width, Standard Gas Spring PAIR, Catalyst	1	\$1,127.85
1b, 2, 3a, 3b, 4a, 4b, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12b, 13-15, 16b, 17a, 17b, 18- 40	503336	Recline Back Kit, Tall Height, 18-22" Width, Standard Gas Spring PAIR, Catalyst	1	\$1,127.85
1a, 2, 3a, 3b, 4c, 4d, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12a, 13-15, 16a, 17a, 17b, 18- 40	503337	Recline Back Kit, Short Height, 14-17" Width, Heavy Duty Gas Spring, Catalyst	1	\$1,127.85
1a, 2, 3a, 3b, 4c, 4d, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12b, 13-15, 16a, 17a, 17b, 18- 40	503338	Recline Back Kit, Short Height, 18-22" Width, Heavy Duty Gas Spring, Catalyst	1	\$1,127.85
1b, 2, 3a, 3b, 4c, 4d, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12a, 13-15, 16b, 17a, 17b, 18- 40	503339	Recline Back Kit, Tall Height, 14-17" Width, Heavy Duty Gas Spring, Catalyst	1	\$1,127.85

Item Number	Part Number	Description	Quantity	Retail
1b, 2, 3a, 3b, 4c, 4d, 5-7, 8a, 8b, 9, 10a, 10b, 11, 12b, 13-15, 16b, 17a, 17b, 18-40	503340	Recline Back Kit, Tall Height, 18-22" Width, Heavy Duty Gas Spring, Catalyst	1	\$1,127.85
2	003217	Backplate, Catalyst Recline, Outside	1	\$27.40
3a	003218	Backplate, Catalyst Recline, Inside Left	1	\$33.17
3b	003219	Backplate, Catalyst Recline, Inside Right	1	\$33.17
4a	001264	Gas Spring Lock 92mm 20kg, Left (Standard Weight Capacity)	1	\$149.21
4b	001265	Gas Spring Lock 92mm 20kg, Right (Standard Weight Capacity)	1	\$149.21
4c	001375	Gas Spring Lock 92mm 40kg, Left (Heavy Duty)	1	\$149.21
4d	001376	Gas Spring Lock 92mm 40kg, Right (Heavy Duty)	1	\$149.21
6	001267	Plug, Back Tube, Lower, Recline	1	\$2.16
9	003222	Recline Stroller Handle Adaptor	1	\$15.57
10a, 11	501090	Stroller Handle w/ Foam, Right	1	\$50.17
10b, 11	501091	Stroller Handle w/ Foam, Left	1	\$50.17
12a	002578	Growth Strut, Stroller, 10-20 Widths <i>Used on chair widths 14 - 17</i>	1	\$25.22
12b	002579	Growth Strut, Stroller, 18-24 Widths <i>Used on chair widths 18 - 22</i>	1	\$27.59
13	003638	Trigger Mount, No Lockout	1	\$3.09
14	001123	Trigger Horizontal	1	\$7.71
15	001124	Pin with Lanyard	1	\$18.53
16a	001268	Cable, Recline, Short	1	\$91.13
16b	001269	Cable, Recline, Long	1	\$92.59
17a	001298	Cover, Gas Spring, Left	1	\$3.89
17b	001299	Cover, Gas Spring, Right	1	\$3.89
18	003223	Recline Backplate Stop Pin	3	\$2.58
19	001271	1" Saddle, .323 Hole	2	\$1.59
20	000612	1" Half Saddle Washer	4	\$1.88

Item Number	Part Number	Description	Quantity	Retail
21	102081	Bushing Flange .319 IDx.394 OD	2	\$0.55
22	101861	Nylon Washer 5/16 x 1/2 x .03	2	\$0.52
23	101076	1/4" E Ring	1	\$1.55
26	101814	M8 x 25 BHCS Blk Zn	1	\$1.03
27	101592	M6 x 45 CLS10.9 DIN7991 FHCS BLKZC	2	\$2.58
28	100670	M6 x 45 BTN SHCS BLK ZC	3	\$1.55
29	101820	M5 x 22 HHCS Blk Zn	1	\$5.67
30	101843	M5 x 20 SHCS Blk Zn	1	\$0.52
31	100734	M5 x 16 PHMS ZC	1	\$1.55
32	102328	M4 x 25 Pan Head Phil Blk Zn	1	\$1.03
33	100659	M8 Nylock Nut Blk Zn	1	\$1.03
34	100658	M6 Nylock Nut Blk Zn	5	\$1.03
35	100657	M5 Nylock Nut Blk Zn	1	\$0.52
36	102329	M4 Nylock Nut Blk Zn	1	\$1.03
37	100746	M6 Flat Washer Blk Zn	3	\$1.55
38	001117	Cable Clip 1"	1	\$3.75
39	100669	M5 x 10 BTN SHCS BLKZC w/Patch	2	\$1.55
40	100533	Detent Button	1	\$3.52