




Catalyst Tension Adjustable Back Upholstery - Standard Backposts



Wings
003094
003095

Wings
003096
003097



Backrest Range and Height	Wings	Total Strap Qty	Total D-Ring Qty
Short 8.5"	003094 003095	3	0
Short 9.5" - 10.5"	003094 003095	4	2
Short 11.5"	003094 003095	5	4
Medium 12.5"	003096 003097	5	0
Medium 13.5" - 14.5"	003096 003097	6	2
Medium 15.5"	003096	7	4
Tall 16.5" - 17.5"	003096 003097	7	4
Tall 18.5" - 19.5 "	003096 003097	8	6

Assemblies contain the cover, amount of straps and d-rings used for the tallest height listed in the assembly. <p>Tension Straps (6a-6d) are based on seat width and used with the combination of Fold Down Push Handle and Tension Adjustable Back Upholstery

Item Number	Part Number	Description	Quantity	Retail
1a	002497	Backrest Cover 13" Rigid / 14" Cat x 9-12"	1	\$130.00
1b	002498	Backrest Cover 13" Rigid / 14" Cat x 11-14"	1	\$130.00
1c	002499	Backrest Cover 13" Rigid / 14" Cat x 13-16"	1	\$130.00
1d	002500	Backrest Cover 13" Rigid / 14" Cat x 15-18"	1	\$130.00
1e	003043	Backrest Cover 13" Rigid / 14" Cat x 17-20"	1	\$150.00
1f	002501	Backrest Cover 14" Rigid / 15" Cat x 9-12"	1	\$140.00
1g	002502	Backrest Cover 14" Rigid / 15" Cat x 11-14"	1	\$140.00
1h	002503	Backrest Cover 14" Rigid / 15" Cat x 13-16"	1	\$140.00
1i	002504	Backrest Cover 14" Rigid / 15" Cat x 15-18"	1	\$140.00
1j	003044	Backrest Cover 14" Rigid / 15" Cat x 17-20"	1	\$150.00
1k	002505	Backrest Cover 15" Rigid / 16" Cat x 9-12"	1	\$140.00
1l	002506	Backrest Cover 15" Rigid / 16" Cat x 11-14"	1	\$140.00
1m	002507	Backrest Cover 15" Rigid / 16" Cat x 13-16"	1	\$140.00
1n	002508	Backrest Cover 15" Rigid / 16" Cat x 15-18"	1	\$140.00
1o	003045	Backrest Cover 15" Rigid / 16" Cat x 17-20"	1	\$150.00
1p	002509	Backrest Cover 16" Rigid / 17" Cat x 9-12"	1	\$140.00
1q	002510	Backrest Cover 16" Rigid / 17" Cat x 11-14"	1	\$140.00
1r	002511	Backrest Cover 16" Rigid / 17" Cat x 13-16"	1	\$140.00
1s	002512	Backrest Cover 16" Rigid / 17" Cat x 15-18"	1	\$140.00
1t	003046	Backrest Cover 16" Rigid / 17" Cat x 17-20"	1	\$150.00
1u	003048	Backrest Cover 17" Rigid / 18" Cat x 9-12"	1	\$150.00
1v	003049	Backrest Cover 17" Rigid / 18" Cat x 11-14"	1	\$150.00
1w	003050	Backrest Cover 17" Rigid / 18" Cat x 13-16"	1	\$150.00
1x	003051	Backrest Cover 17" Rigid / 18" Cat x 15-18"	1	\$150.00
1y	003052	Backrest Cover 17" Rigid / 18" Cat x 17-20"	1	\$150.00
1z	003054	Backrest Cover 18" Rigid / 19" Cat x 9-12"	1	\$150.00
1aa	003055	Backrest Cover 18" Rigid / 19" Cat x 11-14"	1	\$150.00
1ab	003056	Backrest Cover 18" Rigid / 19" Cat x 13-16"	1	\$150.00

Item Number	Part Number	Description	Quantity	Retail
1ac	003057	Backrest Cover 18" Rigid / 19" Cat x 15-18"	1	\$150.00
1ad	003058	Backrest Cover 18" Rigid / 19" Cat x 17-20"	1	\$150.00
1ae	003060	Backrest Cover 19" Rigid / 20" Cat x 9-12"	1	\$160.00
1af	003061	Backrest Cover 19" Rigid / 20" Cat x 11-14"	1	\$160.00
1ag	003062	Backrest Cover 19" Rigid / 20" Cat x 13-16"	1	\$160.00
1ah	003063	Backrest Cover 19" Rigid / 20" Cat x 15-18"	1	\$160.00
1ai	003064	Backrest Cover 19" Rigid / 20" Cat x 17-20"	1	\$160.00
1aj	003066	Backrest Cover 20" Rigid / 21" Cat x 9-12"	1	\$160.00
1ak	003067	Backrest Cover 20" Rigid / 21" Cat x 11-14"	1	\$160.00
1al	003068	Backrest Cover 20" Rigid / 21" Cat x 13-16"	1	\$160.00
1am	003069	Backrest Cover 20" Rigid / 21" Cat x 15-18"	1	\$160.00
1an	003070	Backrest Cover 20" Rigid / 21" Cat x 17-20"	1	\$160.00
1ao	003072	Backrest Cover 21" Rigid / 22" Cat x 13-16"	1	\$175.00
1ap	003073	Backrest Cover 21" Rigid / 22" Cat x 15-18"	1	\$175.00
1aq	003074	Backrest Cover 21" Rigid / 22" Cat x 17-20"	1	\$175.00
2a	003095	Tension Adjustable Upholstery Wing 5.5" Height Left	1	\$50.00
2b	003097	Tension Adjustable Upholstery Wing 9.5" Height Left	1	\$60.00
3a	003094	Tension Adjustable Upholstery Wing 5.5" Height Right	1	\$50.00
3b	003096	Tension Adjustable Upholstery Wing 9.5" Height Right	1	\$60.00
4	100973	Back Upholstery D-Ring	2	\$10.00
5a	003487	Back Upholstery Strap 15.3" <i>Used on 14" chair widths</i>	3	\$12.00
5b	002532	Back Upholstery Strap 17.3" <i>Used on 15" chair widths</i>	3	\$12.00
5c	002533	Back Upholstery Strap 19.3" <i>Used on 16" chair widths</i>	3	\$12.00
5d	002534	Back Upholstery Strap 21.3" <i>Used on 17" chair widths</i>	3	\$12.00
5e	002535	Back Upholstery Strap 23.3" <i>Used on 18" chair widths</i>	3	\$12.00

Item Number	Part Number	Description	Quantity	Retail
5f	002536	Back Upholstery Strap 25.3" <i>Used on 19" chair widths</i>	3	\$15.00
5g	002537	Back Upholstery Strap 27.3" <i>Used on 20" chair widths</i>	3	\$15.00
5h	002538	Back Upholstery Strap 29.3" <i>Used on 21" chair widths</i>	3	\$15.00
5i	002539	Back Upholstery Strap 31.3" <i>Used on 22" chair widths</i>	3	\$15.00
6a	001750	Tension Strap Folding Handle 14-15	1	\$19.81
6b	001751	Tension Strap Folding Handle 16-17	1	\$20.92
6c	001752	Tension Strap Folding Handle 18-19	1	\$22.16
6d	001753	Tension Strap Folding Handle 20-22	1	\$23.28
7	100734	M5 x 16 DIN7985 PHMS ZC	2	\$1.50
1a, 2a, 3a, 4, 5a	505518	Back Upholstery Assembly Tension Adjustable 14W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1c, 2b, 3b, 4, 5a	505519	Back Upholstery Assembly Tension Adjustable 14W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1e, 2b, 3b, 4, 5a	505520	Back Upholstery Assembly Tension Adjustable 14W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1f, 2a, 3a, 4, 5b	505521	Back Upholstery Assembly Tension Adjustable 15W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1h, 2b, 3b, 4, 5b	505522	Back Upholstery Assembly Tension Adjustable 15W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1j, 2b, 3b, 4, 5b	505523	Back Upholstery Assembly Tension Adjustable 15W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1k, 2a, 3a,4, 5c	505524	Back Upholstery Assembly Tension Adjustable 16W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1m, 2b, 3b, 4, 5c	505525	Back Upholstery Assembly Tension Adjustable 16W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1o, 2b, 3b, 4, 5c	505526	Back Upholstery Assembly Tension Adjustable 16W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1p, 2a, 3a, 4, 5d	505527	Back Upholstery Assembly Tension Adjustable 17W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1r, 2b, 3b, 4, 5d	505528	Back Upholstery Assembly Tension Adjustable 17W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1t, 2b, 3b,4, 5d	505529	Back Upholstery Assembly Tension Adjustable 17W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1u, 2a, 3a, 4, 5e	505530	Back Upholstery Assembly Tension Adjustable 18W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1w, 2b, 3b, 4, 5e	505531	Back Upholstery Assembly Tension Adjustable 18W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1y, 2b, 3b, 4, 5e	505532	Back Upholstery Assembly Tension Adjustable 18W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1z, 2a, 3a, 4, 5f	505533	Back Upholstery Assembly Tension Adjustable 19W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53

Item Number	Part Number	Description	Quantity	Retail
1ab, 2b, 3b, 4, 5f	505534	Back Upholstery Assembly Tension Adjustable 19W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1ad, 2b, 3b, 4, 5f	505535	Back Upholstery Assembly Tension Adjustable 19W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1ae, 2a, 3a, 4, 5g	505536	Back Upholstery Assembly Tension Adjustable 20W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1ag, 2b, 3b, 4, 5g	505537	Back Upholstery Assembly Tension Adjustable 20W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1ai, 2b, 3b, 4, 5g	505538	Back Upholstery Assembly Tension Adjustable 20W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1aj, 2a, 3a, 4, 5h	505539	Back Upholstery Assembly Tension Adjustable 21W, 8.5" - 11.5" Height, Standard Backposts	1	\$379.53
1al, 2b, 3b, 4, 5h	505540	Back Upholstery Assembly Tension Adjustable 21W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1an, 2b, 3b, 4, 5h	505541	Back Upholstery Assembly Tension Adjustable 21W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86
1ao, 2b, 3b, 4, 5i	505543	Back Upholstery Assembly Tension Adjustable 22W, 12.5" - 15.5" Height, Standard Backposts	1	\$495.61
1aq, 2b, 3b, 4, 5i	505544	Back Upholstery Assembly Tension Adjustable 22W, 16.5" - 19.5" Height, Standard Backposts	1	\$518.86