


# Focus Pneumatic Tire

---



Item Number	Part Number	Description	Quantity	Retail
1a	100911	Tire 20" x 1 3/8" Pneumatic	1	\$24.00
1b	100634	Tire 22" x 1 3/8" Pneumatic	1	\$24.00
1c	100635	Tire 24" x 1 3/8" Pneumatic (K103)	1	\$25.00

## Tire Compatibility Chart

Wheel		Tire										
Wheel	Wheel Part Number	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke	200529	101587	101992	100918	101587 / 101092							
20" Ki Spoke	200530	100911	100917	102085	100911 / 100913	102974	102957					
22" Ki Spoke	200531	100634	100621	100632	100634 / 100503	103165	102961					101475
24" Ki Spoke	200532	100635	100622	100633	100635 / 100504	102976	102958	102982	100864	103176	102920	101476
25" Ki Spoke	200533					102977	102959		101357	103177		100166
26" Ki Spoke	200534	100636	100806		100636 / 100505	102978	102960		100865			
24" Superlight	101159	100635					102958	102982	100864	103176	102920	
25" Superlight	101158						102959		101357	103177		
12" Mag <sup>1</sup>	N/A	001424	001426		001425							
16" Mag <sup>1</sup>	N/A	001159	001161	003453	001160							
20" 5-Spoke X Core	101961	100911	100917	102085	100911 / 100913							
22" 5-Spoke X Core	101962	100634	100621	100632	100634 / 100503							
24" 5-Spoke X Core	100960	100635	100622	100633	100635 / 100504			102982			102920	
20" Spinergy Spox	See Spinergy Spox Page	100911	100917	102085	100911 / 100913	102974	102957					
22" Spinergy Spox		100634	100621	100632	100634 / 100503	103165	102961					101475
24" Spinergy Spox		100635	100622	100633	100635 / 100504	102976	102958	102982	100864	103176	102920	101476
25" Spinergy Spox						102977	102959		101357	103177		100166
26" Spinergy Spox		100636	100806			100636 / 100505	102978	102960		100865		
22" Spinergy LX	See Spinergy LX Page	100634	100621	100632	100634 / 100503		102961					101475
24" Spinergy LX		100635	100622	100633	100635 / 100504	102976	102958	102982	100864	103176	102920	101476
25" Spinergy LX						102977	102959		101357	103177		100166
26" Spinergy LX		100636	100806			100636 / 100505	102978	102960		100865		
24" Red Topolino	105045	100635					102958	102982	100864	103176	102920	
24" Yellow Topolino	Discontinued	100635					102958	102982	100864	103176	102920	
25" Red Topolino	105047						102959		101357	103177		
25" Yellow Topolino	105048						102959		101357	103177		

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly

## Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part No.	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Marathon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke	200529	101587	101588	101137																	
20" Ki Spoke	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke	200532	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
24" Superlight	101159	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	100635	100589	100596				102982	100589	100596					102920	100589	100596				
20" Spinerger Spox	See Spinerger Spox Page	100911	100912	101137	102957	102844	101137														
22" Spinerger Spox		100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinerger Spox		100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinerger Spox					102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Spinerger Spox		100636	100590	100608	102960	100582	100597				100865		100582	100597							
22" Spinerger LX	See Spinerger LX Page	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinerger LX		100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinerger LX					102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Spinerger LX		100636	100590	100608	102960	100582	100597				100865		100582	100597							
24" Red Topolino	105045	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
24" Yellow Topolino	DXD	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Red Topolino	DXD				102959	101457	101452				101357	103177	101457	101452							
25" Yellow Topolino	105048				102959	101457	101452				101357	103177	101457	101452							