


Catalyst 5Vx Side Frame Assembly



Upper Pivot Assemblies are based on seat width

Item Number	Part Number	Description	Quantity
2a, 16 - 18	500024	Upper Pivot Assembly 14" Width <i>This kit also contains the 14" slider link (not pictured)</i>	1
2b, 16 - 18	500025	Upper Pivot Assembly 15" Width <i>This kit also contains the 15" slider link (not pictured)</i>	1
2b, 16 - 18	500026	Upper Pivot Assembly 16" Width <i>This kit also contains the 16" slider link (not pictured)</i>	1

Item Number	Part Number	Description	Quantity
2b, 16 - 18	500027	Upper Pivot Assembly 17" Width <i>This kit also contains the 17" slider link (not pictured)</i>	1
2b, 16 - 18	500028	Upper Pivot Assembly 18" Width <i>This kit also contains the 18" slider link (not pictured)</i>	1
2b, 16 - 18	500029	Upper Pivot Assembly 19" Width <i>This kit also contains the 19" slider link (not pictured)</i>	1
2b, 16 - 18	500030	Upper Pivot Assembly 20" Width <i>This kit also contains the 20" slider link (not pictured)</i>	1
2b, 16 - 18	501479	Upper Pivot Assembly 21/22" Width <i>This kit also contains the 21/22" slider link (not pictured)</i>	1
5, 6, 7a, 8	501005	Seat Hook and Support Hardware Kit, Standard, Catalyst <i>This kit contains quantity of 1 each</i>	1
5, 6, 7b, 8	114528	Seat Hook and Support Hardware Kit, Standard, for Heavy Duty Catalyst <i>Used with Heavy Duty / 350 lbs Weight Capacity. This kit contains quantity of 1 each</i>	1
1	000011	Frame Connector Upper CAT	1
3	000009	Frame Connector Lower CAT	1
4	000001	Lower Frame Pivot CAT	1
5	000003	Seat Hook, Catalyst	2
6	000004	Seat Hook Support, Catalyst	2
7a	100668	M6 x 20 CLS10.9 BTN SHCS ZC w/Patch <i>Used with Standard Weight Capacity</i>	2
7b	101180	M6 x 20 SHCS BLK ZC <i>Used with Heavy Duty / 350 lbs Weight Capacity</i>	1
8	100831	1/2" Washer OD .040" Thick	3
9	000065	M6 x 23.68 Threaded Barrel <i>Additional 2 are needed with chairs 18"-20" deep</i>	2
10	000066	M6 x 27 Threaded Barrel	2
11	100667	M6 x 12 BTN SHCS BLKZC w/Patch	8
12	000336	Latch Block	1
13	101819	M6 x 40 FHCS Blk Zn	1
14	100658	M6 Nylock Nut Blk Zn	2
15	100704	M6 x 90 CLS8.8 DIN931 HHCS ZC	1

Item Number	Part Number	Description	Quantity
19	100746	M6 Flat Washer Blk Zn	1