

## Little Wave Handrims - Natural Fit



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel. Tab mounted handrims must be installed prior to the tire being installed.

| Item Number | Part Number | Description                                                                                                                 | Quantity |
|-------------|-------------|-----------------------------------------------------------------------------------------------------------------------------|----------|
| 1a          | 107168      | Handrim Natural Fit Standard 6 Rivet Nut 22" w/ Micro Grip<br><i>Used on 22" Ki Spoke, 22" MAXX Spoke, and 22" MAXX Mag</i> | 1        |
| 1c          | 100768      | Handrim Natural Fit Standard 5 Rivet Nut 24"<br><i>Used on 24" 5-Spoke Mag</i>                                              | 1        |
| 1d          | 100889      | Handrim Natural Fit Standard 6 Tab 22" (3R22T1)<br><i>Used on 22" Spinergy</i>                                              | 1        |
| 2a          | 107177      | Handrim Natural Fit LT 6 Rivet Nut 22" w/ Micro Grip<br><i>Used on 22" Ki Spoke, 22" MAXX Spoke, and 22" MAXX Mag</i>       | 1        |
| 2d          | 100888      | Handrim Natural Fit LT 6 Tab 22"<br><i>Used on 22" Spinergy</i>                                                             | 1        |
| 2d          | 107182      | Handrim Natural Fit LT 6 Tab 22" w/ Micro Grip<br><i>Used on 22" Spinergy</i>                                               | 1        |
| 4b          | 100829      | Natural Fit Thumb Grip Super 24"                                                                                            | 1        |
| 5-7         | 502169      | Handrim Hardware Kit, 6 Rivet Nut, Natural Fit<br><i>Used with Ki Spoke and MAXX Spoke Wheels</i>                           | 1        |
| 10-11       | 502174      | Handrim Hardware Kit, 6 Tab, EACH<br><i>Used with Spinergy</i>                                                              | 1        |
| 8b-9        | 117379      | Handrim Hardware Kit, MAXX Mag, EACH<br><i>Used with MAXX Mag wheels</i>                                                    | 1        |

| Item Number | Part Number | Description                                                                                                           | Quantity |
|-------------|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 5           | 100698      | Handrim Tab<br><i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i> | 6        |
| 6           | 100835      | 10/32 x 1 3/4" Pan Head Screw<br><i>Used with Ki Spoke and MAXX Spoke Wheels</i>                                      | 6        |
| 7           | 100653      | Handrim Spacer Spoke<br><i>Used with Ki Spoke and MAXX Spoke Wheels</i>                                               | 6        |
| 8a          | 100629      | Handrim Spacer Alum Ano Riv<br><i>Used with 5-Spoke X Core Wheels</i>                                                 | 5        |
| 8b          | 101756      | Handrim Spacer Plastic Coated Rivet<br><i>Used with MAXX Mag Wheels</i>                                               | 6        |
| 9           | 101893      | 10-32 x 2 PHTH ZN<br><i>Used with 5-Spoke X Core (Qty 5 needed) and MAXX Mag wheels (Qty 6 needed)</i>                | 1        |
| 10          | 100669      | M5 x 10 BTN SHCS BLKZC w/Patch<br><i>Used with Spinergy</i>                                                           | 6        |
| 11          | 100657      | M5 Nylock Nut Blk Zn<br><i>Used with Spinergy</i>                                                                     | 6        |

| Handrim Hardware Chart |                |                    |                                                              |                                              |                                                                  |                   |                |                  |                  |                    |                  |                      |                  |          |  |
|------------------------|----------------|--------------------|--------------------------------------------------------------|----------------------------------------------|------------------------------------------------------------------|-------------------|----------------|------------------|------------------|--------------------|------------------|----------------------|------------------|----------|--|
| Wheel                  |                | Handrim Connection |                                                              |                                              |                                                                  | Handrim           |                |                  |                  |                    |                  |                      |                  |          |  |
| Wheel                  | Wheel Part No. | Connection Points  | Rivet/Tab                                                    | Spacer*                                      | Screw*                                                           | Aluminum Anodized | Superlight     | Plastic Coated   | Projection       | Ergonomic Standard | Ergonomic LT     | Natural Fit Standard | Natural Fit LT   | Flex Rim |  |
| 22" Spoke              | 116730         | 5 or 6             | 116732 (Used with 5 Rivet)<br><br>100698 (Used with 6 Rivet) | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection<br>Ergonomic: 122063           | 101898 (5 Rivet)  |                | 101964 (5 Rivet) | 112824 (6 Rivet) | 113082 (6 Rivet)   | 113085 (6 Rivet) | 200538 (6 Rivet)     | 200201 (6 Rivet) |          |  |
| 24" Spoke              | 116728         | 5 or 6             | Not used w/ Superlight Handrims                              |                                              | Natural Fit: 100835<br><br>Superlight: 100669<br><br>Nut: 100657 | 100975 (5 Rivet)  | 101161 (6 Tab) | 100976 (5 Rivet) | 112825 (6 Rivet) | 113083 (6 Rivet)   | 113086 (6 Rivet) | 100793 (6 Rivet)     | 200202 (6 Rivet) |          |  |
| 18" MAXX Spoke         | 109285         | 3                  | 100698 (Not used w/ Superlight Handrim)                      | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection, Ergonomic: 122063             | 100206            |                |                  |                  |                    |                  |                      |                  |          |  |
| 20" MAXX Spoke         | 108244         | 6                  |                                                              |                                              | Natural Fit: 100835                                              | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Spoke         | 105135         |                    |                                                              |                                              | Superlight: 100669                                               | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Spoke         | 105136         |                    |                                                              |                                              | Nut: 100657                                                      | 200349            | 101161         | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 25" MAXX Spoke         | 107436         |                    |                                                              |                                              | 200350                                                           |                   | 112822         | 200548           | 113084           | 113087             | 200539           | 200540               |                  |          |  |
| 26" MAXX Spoke         | 107437         |                    |                                                              |                                              | 200351                                                           |                   | 112823         | 200549           |                  |                    | 100907           | 101454               |                  |          |  |
| 20" MAXX Mag           | 111853         | 6                  |                                                              | Aluminum, Plastic Coated, Projection: 100629 | Aluminum, Plastic Coated, Projection, Ergonomic: 122064          | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Mag           | 111854         |                    |                                                              | Natural Fit & Ergonomic: 101756              | Natural Fit: 101893                                              | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Mag           | 111855         |                    |                                                              |                                              |                                                                  | 105137            |                | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 24" Superlight         | 101159         | 6                  |                                                              |                                              | 100536                                                           | 100754            | 101161         | 100836           |                  | 113080             | 113081           | 100830               | 100828           |          |  |
| 25" Superlight         | 101160         |                    |                                                              |                                              |                                                                  |                   |                | 101160           | 101091           |                    |                  |                      | 101464           | 101460   |  |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 122065 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 122063 for the screw on MAXX Mag wheels.

Effective 4/19/2024

| Handrim Hardware Chart (Continued) |                         |                    |           |         |                    |                   |            |                |            |                    |              |                      |                |                       |
|------------------------------------|-------------------------|--------------------|-----------|---------|--------------------|-------------------|------------|----------------|------------|--------------------|--------------|----------------------|----------------|-----------------------|
| Wheel                              |                         | Handrim Connection |           |         |                    | Handrim           |            |                |            |                    |              |                      |                |                       |
| Wheel                              | Wheel Part No.          | Connection Points  | Rivet/Tab | Spacer* | Screw*             | Aluminum Anodized | Superlight | Plastic Coated | Projection | Ergonomic Standard | Ergonomic LT | Natural Fit Standard | Natural Fit LT | Flex Rim              |
| 20" Spinerger Spox                 | See Spinerger Spox Page | 6                  |           |         | Ergonomic: 100536  | 103125            |            | 103179         |            |                    |              |                      |                |                       |
| 22" Spinerger Spox                 |                         |                    |           |         | All Others: 100669 | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger Spox                 |                         |                    |           |         |                    | 100766            | 101161     | 100615         |            |                    |              | 100830               | 100828         |                       |
| 25" Spinerger Spox                 |                         |                    |           |         | Nut: 100657        | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger Spox                 |                         |                    |           |         |                    | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |
| 22" Spinerger LX                   | See Spinerger LX Page   | 6                  |           |         | Ergonomic: 100536  | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger LX                   |                         |                    |           |         | All Others: 100669 | 100766            | 101161     | 100615         |            | 113080             | 113081       | 100830               | 100828         | See Spinerger LX Page |
| 25" Spinerger LX                   |                         |                    |           |         |                    | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger LX                   |                         |                    |           |         | Nut: 100657        | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |

| Handrim Hardware Chart (Discontinued Wheels) |                |                    |           |                                           |                                           |                   |            |                |            |                      |                |          |
|----------------------------------------------|----------------|--------------------|-----------|-------------------------------------------|-------------------------------------------|-------------------|------------|----------------|------------|----------------------|----------------|----------|
| Wheel                                        |                | Handrim Connection |           |                                           |                                           | Handrim           |            |                |            |                      |                |          |
| Wheel                                        | Wheel Part No. | Connection Points  | Rivet/Tab | Spacer*                                   | Screw*                                    | Aluminum Anodized | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim |
| 18" Ki Spoke                                 | 200529         | 3                  | 100698    | 100653                                    | Aluminum, Plastic Ctd, Projection: 122063 | 100206            |            | 101106         |            |                      |                |          |
| 20" Ki Spoke                                 | 200530         | 6                  |           |                                           |                                           | 200536            |            | 200542         |            |                      |                |          |
| 22" Ki Spoke                                 | 200531         |                    |           |                                           |                                           | 100560            |            | 100576         | 100569     | 200538               | 200201         |          |
| 24" Ki Spoke                                 | 200532         |                    |           |                                           |                                           | 200349            |            | 100577         | 200547     | 100793               | 200202         |          |
| 25" Ki Spoke                                 | 200533         |                    |           |                                           | Natural Fit: 100835                       | 200350            |            | 101870         | 200548     | 200539               | 200540         |          |
| 26" Ki Spoke                                 | 200534         |                    |           |                                           |                                           | 200351            |            | 100578         | 200549     | 100907               | 101454         |          |
| 20" 5-Spoke X Core                           | 101961         | 5                  |           | Aluminum, Projection, Natural Fit: 100629 | 100724                                    | 101897            |            | 101963         |            |                      |                |          |
| 22" 5-Spoke X Core                           | 101962         |                    |           |                                           |                                           | 101898            |            | 101964         |            |                      |                |          |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           |                                           | 100975            |            | 100976         | 200546     |                      |                |          |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           | Plastic Ctd: 101756                       | 101893            |            |                |            |                      | 100768         | 100769   |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 122065 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 122063 for the screw on MAXX Mag wheels.

Effective 4/19/2024