

<b>Bill To</b> (Ultimate Medical Supplier Only)	
Purchase Order #:	_____
Quote #:	_____
Account #:	_____
Ordered By:	_____
Supplier Company:	_____
Address:	_____
City, State, Zip:	_____
Phone: _____	Email: _____
Payer Source:	_____
<b>Name of ATP (Required):</b> _____	

<b>Ship To</b>	
Name:	_____
Address:	_____
City, State, Zip:	_____
Phone:	_____
Tag For:	_____

<b>Evolv XT User Specifications</b>	
Approximate Height Range	6'0"-6'10"
Weight Limit	350 lbs.
Seat Depth Range (seat pivot to flat back)	20"-25"
Seat to Foot Plate Range (from seat pivot)	17"-23"

<b>Base</b>	<i>Note: The DEFAULT items will be added if no selection is made.</i>	
<b>EasyStand® Evolv® XT</b>		<b>\$3,873</b>
<b>Price Includes:</b>		
PNG50209-1 EasyStand Evolv XT Base		\$2,409
PNG50377 Extra Tall Configuration		\$1,464

<b>Lift Mechanisms-must select one</b>		
<input type="checkbox"/> PNG50320	Manual Hydraulic Actuator <b>DEFAULT</b>	N/C
<b>Manual Hydraulic Actuator Handle Location</b>		
<input type="checkbox"/> PNG50059	Right Side <b>DEFAULT</b>	N/C
<input type="checkbox"/> PNG50126	Left Side	N/C
<b>Manual Hydraulic Actuator Handle Accessories</b>		
<input type="checkbox"/> PNG50328	Adjustable/Removable Actuator Handle	\$257
<input type="checkbox"/> PNG50044	Quad Grip Handle Extension	\$278
<input type="checkbox"/> PNG50070	T-Style Handle Extension	\$174
<input type="checkbox"/> PNG50321	Pow'r Up Lift (battery powered)	\$1,419

**Front Frames and Trays**  
**MUST SELECT ONE FRAME AND ONE TRAY IN THE SAME FRAME SECTION**

<b>Shadow Tray Frame-back required*</b>		
<input type="checkbox"/> PNG51054	Column Front for Shadow Tray with Wheels	\$571
<input type="checkbox"/> PNG51056	Column Front for Shadow Tray with Swivel Casters	\$618
<b>Shadow Tray*</b>		
<input type="checkbox"/> PNG50515	Black Molded Adjustable Shadow Tray-Right <small>(7"-12"H from seat, 4"-20"D from back)</small>	\$1,593
<input type="checkbox"/> PNG50537	Black Molded Adjustable Shadow Tray-Left <small>(7"-12"H from seat, 4"-20"D from back)</small>	\$1,593
<input type="checkbox"/> PNG50517	Clear Adjustable Shadow Tray-Right <small>(6"-11"H from seat, 4"-20"D from back)</small>	\$1,778
<input type="checkbox"/> PNG50538	Clear Adjustable Shadow Tray-Left <small>(6"-11"H from seat, 4"-20"D from back)</small>	\$1,778

<b>Column Frame*</b>		
<input type="checkbox"/> PNG51055	Column with Wheels <b>DEFAULT</b>	\$647
<input type="checkbox"/> PNG51057	Column with Swivel Casters	\$708
<b>Column Tray*</b>		
<input type="checkbox"/> PNG50247-1	Black Molded Tray <b>DEFAULT</b>	\$660
<input type="checkbox"/> PNG50248-1	Clear Tray	\$849
<input type="checkbox"/> PNG50194-1	Oversized Clear Tray	\$1,005
<input type="checkbox"/> PNG50368-1	Black Molded Angle Adjustable Tray	\$1,052
<input type="checkbox"/> PNG50366-1	Clear Angle Adjustable Tray	\$1,234

<b>Mobile Frame</b>		
<input type="checkbox"/> PNG51058	Column with Mobile Option	\$2,718
<b>Mobile Tray</b>		
<input type="checkbox"/> PNG50372-1	Narrow Clear Tray	\$849
<input type="checkbox"/> PNG50327-1	Narrow Clear Angle Adjustable Tray	\$1,234

<b>Tray Accessories</b>		
<input type="checkbox"/> PNG50042	Hand Grips (for clear tray only, pair)*	\$213
<input type="checkbox"/> PNG30295	Padded Tray Cover (not available with oversized clear trays or mobile trays)*	\$255

<b>Chest Pad-must select one</b>		
<input type="checkbox"/> PNG30001	Medium Chest Pad <b>DEFAULT</b>	N/C
<input type="checkbox"/> PNG30000	Large Contoured Chest Pad	\$161
<input type="checkbox"/> PNG50331	Large Contoured Adjustable Chest Pad	\$231

<b>Knee Pads-must select one</b>		
<input type="checkbox"/> PNG50460-1	1-Piece Knee Pad (pair) <b>DEFAULT</b>	\$371
<input type="checkbox"/> PNG50384-1	Independent Knee Pads (pair)	\$917
<input type="checkbox"/> PNG50385-1	Independent Roho® Knee Pads (pair)	\$1,054

\* not available with mobile frame

**Foot Plates-must select one**

- PNG50278-1 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$373

**Foot Plate Accessories**

- PNG30030 Foot Straps-15"L (length over top of foot, pair) \$88
- PNG50505 Foot Straps-21"L (length over top of foot, pair) \$99
- PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair) \$142
- PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair) \$153

**Seats-must select one**

- PNG30006-1 Planar Seat **DEFAULT** \$398
- PNG30270 Transfer Seat (replaces standard planar seat) \$454
- PNG50318 Rotating Seat (replaces standard planar seat) \$664
- PNG30023 Seat with Roho® Insert (replaces standard planar seat) \$772

**Seat Accessories**

- PNG50500 Easy-Adjust Seat Depth-**CANNOT be ordered off-unit** \$123  
(roller bearing allows for smoother seat depth adjustment)

**Pelvic Support**

- PNG50161 Hip Supports (10"-15" range, pair) \$409
- PNG50198 Hip Supports (12"-17" range, pair, limited adjustment with shadow tray) \$409
- PNG50518 Hip Supports (17"-22" range, pair, limited adjustment with shadow tray) \$409
- PNG30029 Velcro® Positioning Belt-26"-55" (back required) \$88
- PNG30306 Velcro® Positioning Belt-38"-67" (back required) \$99
- PNG30028 Positioning Belt with Airline Style Buckle-26"-55" (back required) \$191

**Backs**

- PNG50066 Flat Back-15"H \$750
- PNG50067 Flat Back-19"H \$791
- PNG50062 Contoured Back-15"H \$825
- PNG50068 Contoured Back-19"H \$866
- PNG50171 Removable Flat Back-15"H \$860
- PNG50173 Removable Flat Back-19"H \$901
- PNG50172 Removable Contoured Back-15"H \$935
- PNG50174 Removable Contoured Back-19"H \$978

**Back Accessories**

- P600002 Back Angle Adjustment Wheel (back required) \$169
- P80246 Angle Locator \$44

**Notes:**

---



---



---



---



---

**Chest Straps-back required**

- PNG50192 Chest Strap-30"-46" \$196
- PNG50193 Chest Strap-34"-54" \$206
- PNG50504 Chest Strap-46"-66" \$217

**Upper Body Support-back required**

**ATTENTION - ALL OPTIONS BELOW REQUIRE THE ACCESSORIES MOUNTING BRACKET**

- PNG50034 Accessories Mounting Bracket \$104

**Chest Vests-back and accessories mounting bracket required**

- PNG11212 High Mount Chest Vest Bracket \$89
- PNG30170 X-Style Chest Vest-11"Lx9.5"W \$334
- PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$341
- PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$334
- PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$341
- PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W \$386
- PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W \$393
- PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W \$386
- PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W \$393

**Lateral Supports-back and accessories mounting bracket required**

- PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair) \$409
- PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair) \$409
- PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair) \$409
- PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair) \$409

**Head Supports-back and accessories mounting bracket required**

- PNG50240 Head Support-Short \$445
- PNG50037 Head Support-Tall \$455

**Other Accessories**

- PNG50040 Push Handles (back required, pair) \$174
- P80809 Tool Pouch \$31

\* not available with mobile

**(A) Seat Depth Range: 20"-25"**  
**(seat pivot to flat back)**

Measure the distance between the back of the buttocks and the back of the knee.

**(B) Seat to Foot Plate Range: 17"-23"**  
**(from seat pivot)**

Measure the distance between the foot/shoe and the knee bend.



**AltimateMedical**  
creating products · changing lives

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

webquotes@altimatemedical.com easystand.com

**Email form to webquotes@altimatemedical.com**