



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

☐ Quote ☐ PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: ☐ Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Effective 07/18/2022

Rogue ALX K0005, E1235 **\$2070**
250 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

See reference A on page 6

- ☐ 14" N/C
☐ 15" N/C
☐ 16" N/C
☐ 17" N/C
☐ 18" N/C
☐ 19" \$125
☐ 20" E2201 \$500
☐ 21" E2201 \$500
☐ 22" E2201 \$500

Seat Depth

See reference B on page 6

- ☐ 14" N/C
☐ 15" N/C
☐ 16" N/C
☐ 17" N/C
☐ 18" N/C
☐ 19" \$75
☐ 20" E2203 \$550

Weight Capacity

- ☐ 250 lb Std
☐ 350 lb (Heavy Duty) \$750
Requires 6" casters

Transit

- ☐ Transit System \$320
125 lb maximum weight capacity
☐ Omit Transit N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height
Height Adjustable	N/C	☐ Short	☐ 11 ² ☐ 12 ☐ 13 ☐ 14
		☐ Medium	☐ 13 ☐ 14 ☐ 15 ☐ 16 ☐ 17
		☐ Tall	☐ 17 ☐ 18 ☐ 19 ☐ 20
Fixed Height ¹	\$85	☐ Fixed	☐ 9 ² ☐ 10 ² ☐ 11 ² ☐ 12 ³ ☐ 13 ☐ 14 ☐ 15 ☐ 16 ☐ 17 ☐ 18

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

³ Fold Down Push Handle requires Height Adjustable Rigidizer Bar

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 6.

- ☐ 0" (Performance) N/C
☐ 1" N/C
☐ 2" (Classic) N/C

Adjustable Depth

Adds the capability to grow seat in depth

- ☐ No Rear Adjustment N/C
☐ 1/2" N/C
☐ 1" \$85
☐ 2" \$165
☐ 3" \$165

Long Sling

Available only with 1", 2" and 3" Adjustable Depth

- ☐ Long Sling N/C
Sling will be maximum 2" longer than seat depth

Seat Option

- ☐ Standard Seat Sling Std
☐ Solid Seat Pan E2231 \$265
Seat Pan is 1" longer when Long Sling option is selected

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 6.

- ☐ 90° \$125
☐ 85° N/C
☐ 80° N/C
☐ 75° N/C
☐ 70° \$125

Push Handle

- ☐ Omit Push Handle N/C
☐ Standard Push Handle \$85
☐ Low Profile Push Handle \$150
2" shorter than standard push handles
☐ Bolt-On Push Handle \$125
Handles will be 1" - 2" above back height
☐ Fold Down Push Handle \$100
Not available with 350 lb Weight Capacity

DCN0929.1 Pricing and specifications subject to change without notice

Backrest continued

Rigidizer Bar

- ☐ Adjustable Height Rigidizer Bar _____ \$75
☐ Non-Adjustable Height Rigidizer Bar _____ N/C

Locking Backrest

- ☐ Locking Fold Down Back _____ \$95
Rogue style locking back. Not available with 350 lb Weight Capacity

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

- ☐ Standard ☐ +3° ☐ +6° ☐ +9°
☐ -3° ☐ -6° ☐ -9° ☐ -12°

Back Upholstery

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable

Not available with 350 lb Weight Capacity

- ☐ Standard Cover _____ Std
☐ Air Mesh Breathable Cover _____ \$75

Tension Adjustable E2611

- ☐ Standard Cover _____ \$320
☐ Air Mesh Breathable Cover _____ \$390

- ☐ Omit Back Upholstery _____ - \$10

Military Patch

US Military Service Patch _____ \$90

- ☐ Army ☐ Navy ☐ Marines ☐ Air Force ☐ Coast Guard

☐ US Purple Heart Patch _____ \$90

Purple Heart with Military Service Patch Combo _____ \$145

- ☐ Army ☐ Navy ☐ Marines ☐ Air Force ☐ Coast Guard

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend.

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube.

2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 6.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 6.

Footrest Taper	Caster Size	Front Seat Height														
☐ Straight	3"	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5
	4"	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5	
	5"	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5		
	6"	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5			
☐ 1" Taper	3"	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5						
	4"	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5							
	5"	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5								
	6"	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5									
☐ 2" Taper	3"	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5						
	4"	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5							
	5"	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5								
	6"	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5									

Caster

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	☐	☐		
	.75" Lighted Rollerblade	\$70	☐	☐		
	1" Poly	N/C		☐	☐	☐
	1" Poly Aluminum	\$110		☐	☐	
	1" Pneumatic	\$80				☐
	1.5" Poly	\$55		☐	☐	☐
	1.5" Soft Roll Aluminum	\$185	☐	☐	☐	☐

Fork Options

- ☐ Standard Fork _____ Std
☐ Slotted Fork _____ \$75
☐ Single Sided Fork _____ \$170

Not available with 6" Caster Size

☐ Frog Legs Suspension Fork E1015 _____ \$450

Not available with 350 lb Weight Capacity

- ☐ User Weight 80 - 120 lbs
☐ User Weight 120 - 160 lbs
☐ User Weight 160 - 220 lbs
☐ User Weight 220 - 250 lbs

☐ Frog Legs II Suspension System _____ \$450

Lightweight shock absorbing fork. Not available with 350 lb Weight Capacity.

☐ Caster Pin Locks K0073 _____ \$185

Available only with Slotted Fork

Footrest

Footplate

- ☐ Tubular "Open" Footplate _____ Std
Not available w/footrest widths greater than 17.5". See Reference G on page 6.
- ☐ Tubular Footplate with ABS Cover _____ \$30
Not available w/footrest widths greater than 15.5". See Reference G on page 6.
- ☐ Aluminum Angle Adjustable K0040 _____ \$95
- ☐ Aluminum Angle Adjustable with Riser K0040, K0037 _____ \$145
Not available with 350 lb Weight Capacity
- ☐ High Mount Aluminum Angle Adjustable K0040, K0037 _____ \$95
Not available with 350 lb Weight Capacity
- ☐ Angle Adjustable Flip Under K0040, K0037 _____ \$185
Not available with footrest widths less than 9.5". Not available with 350 lb Weight Capacity
- ☐ High Mount Angle Adjustable Flip Under K0040, K0037 _____ \$270
Not available with footrest widths greater than 15.5" or 350 lb Weight Capacity.
See Reference G on page 6.

Footrest Protector

- ☐ Footrest Protector _____ \$25
Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest.
Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 6.

☐ 8"	☐ 8.5"	☐ 9"	☐ 9.5"	☐ 10"	☐ 10.5"
☐ 11"	☐ 11.5"	☐ 12"	☐ 12.5"	☐ 13"	☐ 13.5"
☐ 14"	☐ 14.5"	☐ 15"	☐ 15.5"	☐ 16"	☐ 16.5"
☐ 17"	☐ 17.5"	☐ 18"	☐ 18.5"	☐ 19"	☐ 19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels.
Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 6.

Wheel Size	Rear Seat Height													
22"	☐ 13.5 ¹	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20
24"	☐ 14.5 ¹	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20		
25"	☐ 15 ¹	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20			
26"	☐ 15.5 ¹	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20				

¹ Available only with Low Profile tires

Camber Bar Tower

- ☐ Standard Tower _____ N/C
- ☐ Mx Tower _____ \$100

Standard Tower allows rear seat height adjustment up or down in 1/4" increments. 1" adjustment available; 1/2" in lowest end of range.
Mx Tower required with select low seat heights or when large rear seat height adjustment range is required. Provides 3" adjustment range in 1/4" increments.

Rear Wheel

		Rear Wheel Size				
		Price	22"	24"	25"	26"
Rear Wheel Type	MAXX Performance Spoke	\$125	☐	☐	☐	☐
	MAXX Performance Spoke	N/C	☐	☐	☐	☐
	Mag	N/C	☐	☐	☐	☐
	Spinergy Spox	\$845	☐	☐	☐	☐
	Spinergy LX	\$1050	☐	☐	☐	☐

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color:

- ☐ Black ☐ Red ☐ White ☐ Blue ☐ Yellow
☐ Pink ¹ ☐ Green ¹ ☐ Purple ¹ ☐ Orange ¹

¹ Available only with Spinergy LX

Handrim Tab:

Available only with Aluminum Anodized and Plastic Coated Handrims

- ☐ Short Tab _____ \$25
- ☐ Long Tab _____ N/C

Rear Tire Type

Full Profile

- ☐ Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- ☐ Full Poly 1 3/8" E2218 _____ \$55
Not available on 25" or 26" wheel sizes
- ☐ Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$85
Not available on 25" wheel size
- ☐ Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
- ☐ Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size
- ☐ Mountain 2" Black E2211, E2212 _____ \$125
Available only on 22", 24" and 25" MAXX Performance and Spinergy wheels

Low Profile

- ☐ Low Poly 1" E2218 _____ N/C
Not available on 25" or 26" wheel sizes or 350 lb Weight Capacity
- ☐ Sentinel Solid Tire 1" E2218 _____ \$95
Not available on 22" Mag wheel
- ☐ High Pressure 1" E2211, E2212 _____ \$80
Not available on Mag wheel
- ☐ Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$140
Not available on Mag wheel or 22" wheel size

Rear Seat Height and Wheels continued

Handrim

Standard Handrims

- ☐ Aluminum Anodized _____ Std
- ☐ Plastic Coated _____ \$112
- ☐ Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with MAXX Performance or 24" Mag wheel. Not Available with 25" wheels.
- ☐ Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- ☐ Natural Fit LT Regular Grip _____ \$465
- ☐ Natural Fit Standard Regular Grip _____ \$445
- ☐ Natural Fit LT Super Grip _____ \$515
- ☐ Natural Fit Standard Super Grip _____ \$495
- ☐ Natural Fit LT Omit Grip _____ \$395
- ☐ Natural Fit Standard Omit Grip _____ \$375

Handrim Option

- ☐ 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axle

- ☐ Quick Release _____ Std
- ☐ Quad Quick Release _____ \$185
Not available with Mag or Spinergy wheels

Wheel Locks

Wheel Lock

- ☐ Push to Lock _____ Std
- ☐ Push to Lock Flush Mount _____ N/C
- ☐ Pull to Lock _____ N/C
- ☐ Under Seat Scissor Lock _____ N/C
- ☐ Low Profile Scissor Lock _____ \$65
Not available with Low Poly, Full Poly or Sentinel Solid tires
- ☐ Short Thro Scissor Lock _____ \$110
- ☐ Omit Wheel Locks _____ N/C

Armrests

Armrest

Height Adjustable T-Arm Tall (11.5" - 15") E0973 **New!**

1.5" minimum Rear Wheel Spacing

- ☐ Desk _____ \$325
- ☐ Full _____ \$325

Height Adjustable T-Arm (8.5" - 13.5") E0973

1.5" minimum Rear Wheel Spacing

- ☐ Desk _____ \$325
- ☐ Full _____ \$325

Height Adjustable T-Arm Low (6.5" - 10.5") E0973

1.5" minimum Rear Wheel Spacing. If requested configuration is unavailable, Standard Height T-Arms will be used.

- ☐ Desk _____ \$325
- ☐ Full _____ \$325

Center of Gravity Preset

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 6.

Amputee Positions

☐ -1.5"	☐ -1.25"	☐ -1"	☐ -.75"
☐ -.5"	☐ -.25"	☐ 0"	☐ .25"

Standard Positions

☐ .5"	☐ .75"	☐ 1"	☐ 1.25"
☐ 1.5"	☐ 1.75"	☐ 2"	☐ 2.25"
☐ 2.5"	☐ 2.75"	☐ 3"	☐ 3.25"
☐ 3.5"	☐ 3.75"	☐ 4"	☐ 4.25"

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 6.

- ☐ .75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75"

Rear Wheel Camber

- ☐ 0° ☐ 2° ☐ 4° ☐ 6° ☐ 8°

Wheel Lock Options

- ☐ Extension Handles E0961 _____ \$77
Available only with Push to Lock or Pull to Lock Wheel Locks
- ☐ Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Direct Mount

☐ Composite Side Guards - Direct Mount Receiver K0108 \$150

Aluminum Straight K0108

☐ Black Matte \$195

☐ Match Frame Color \$245

Aluminum Adjustable Fender with Foam Pad K0108

Not available with Adjustable Depth greater than 1/2"

☐ Black Matte \$275

☐ Match Frame Color \$325

☐ Carbon Fiber Adjustable Fender K0108 \$485

Not available with Adjustable Depth greater than 1/2"

Clamp-On Mount

☐ Composite Side Guards - Clamp-on Receiver K0108 \$175

Aluminum Straight K0108

☐ Black Matte \$225

☐ Match Frame Color \$275

Aluminum Adjustable Fender with Foam Pad K0108

☐ Black Matte \$305

☐ Match Frame Color \$355

☐ Carbon Fiber Adjustable Fender K0108 \$515

Fabric Side Guard

Not available when omitting Back Upholstery

☐ Fabric Side Guards K0108 \$95

Accessories

Anti-Tip

☐ Rear Anti-Tips E0971 \$150

☐ Omit Rear Anti-Tips N/C

Calf Strap

☐ Velcro Adjustable K0038 \$42

☐ Padded Velcro Adjustable K0038 \$68

☐ Bodypoint Padded AeroMesh K0038 \$118

Pelvic Positioning Belt

☐ 1 1/2" Auto Buckle E0978 \$80

☐ 1 1/2" Padded Auto Buckle E0978 \$170

☐ 2" Airline Buckle E0978 \$110

Bodypoint Evoflex - Pelvic Stabilizer E0978 \$188

Size: ☐ Small (16-24") ☐ Medium (19-29") ☐ Large (25-37")

Belt Mounting Kit

☐ Ki Mobility Belt Mounting Clamp \$35

☐ Bodypoint Belt Mounting Bracket Kit \$108

Other Accessories

☐ Tool Kit \$45

☐ Removable Underseat Pouch \$50

☐ Ki Mobility Backpack \$65

☐ Luggage Carrier \$85

☐ Black Neoprene Impact Guards \$90

☐ Under Chair Carrier and Removable Bag \$85

Not available with scissor locks

☐ Clear Spoke Guards K0065 \$115

Not available with Mag wheel

☐ Black Spoke Guards K0065 \$115

Not available with 20" wheel size or Mag wheel

☐ Cane and Crutch Holder \$195

Frame Color

Gloss Finish

☐ Black N/C

☐ Blue Skies N/C

☐ Bubblegum Pink N/C

☐ Candy Blue N/C

☐ Candy Purple N/C

☐ Candy Red N/C

☐ Copper N/C

☐ Granny Smith N/C

☐ Metallic Orange N/C

☐ Silver N/C

☐ Yellow N/C

Matte Finish

☐ Charcoal Gray N/C

☐ Emerald Green N/C

☐ Iridescent Burgundy N/C

☐ Iridescent Navy N/C

☐ Military Green N/C

☐ Retro Red N/C

☐ Royal Blue N/C

☐ Shadow Black N/C

☐ Smooth White N/C

References

A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"
21"	18.5"	16.5"	14.5"
22"	19.5"	17.5"	15.5"

H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2





Diagram 3



Axiom Cushions

Weight capacity: 350 lbs

Axiom G Cushion and Extra Outer Cover

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$123	\$77	\$123	\$77	\$123	\$77	\$123	\$77	\$123	\$77	\$123	\$77	\$123	\$77	\$204	\$77
	HCPCS	E2601		E2601		E2601		E2601		E2601		E2601		E2601		E2602	
Depth	14"	☐	☐														
	15"	☐	☐	☐	☐	☐	☐										
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐
	22"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐

Axiom S Cushion and Extra Outer Cover

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$216	\$88	\$216	\$88	\$216	\$88	\$216	\$88	\$216	\$88	\$216	\$88	\$216	\$88	\$323	\$88
	HCPCS	E2603		E2603		E2603		E2603		E2603		E2603		E2603		E2604	
Depth	14"	☐	☐														
	15"	☐	☐	☐	☐	☐	☐										
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐
	22"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐

Axiom P Cushion and Extra Outer Cover

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$364	\$98	\$364	\$98	\$364	\$98	\$364	\$98	\$364	\$98	\$364	\$98	\$364	\$98	\$472	\$98
	HCPCS	E2605		E2605		E2605		E2605		E2605		E2605		E2605		E2606	
Depth	14"	☐	☐														
	15"	☐	☐	☐	☐	☐	☐										
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐		
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐
	22"					☐	☐			☐	☐	☐	☐	☐	☐	☐	☐

Axiom Cushions continued

Weight capacity: 350 lbs

Axiom SP Visco Cushion and Extra Inner / Outer Cover

		Width																	
		14"		15"		16"		17"		18"		19"		20"		21"		22"	
		Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer
Retail		\$441	\$108	\$441	\$108	\$441	\$108	\$441	\$108	\$441	\$108	\$441	\$108	\$441	\$108	\$539	\$108	\$539	\$108
HCPCS	E2607			E2607		E2607		E2607		E2607		E2607		E2607		E2607		E2608	
Depth	14"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	15"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	22"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Axiom SP Fluid Cushion and Extra Inner / Outer Cover

		Width																	
		14"		15"		16"		17"		18"		19"		20"		21"		22"	
		Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer	Cushion	Cover Inner/Outer
Retail		\$482	\$108	\$482	\$108	\$482	\$108	\$482	\$108	\$482	\$108	\$482	\$108	\$482	\$108	\$590	\$108	\$590	\$108
HCPCS	E2607			E2607		E2607		E2607		E2607		E2607		E2607		E2607		E2608	
Depth	14"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	15"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	22"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Solid Seat Insert

Select size in chart below. Weight capacity: 350 lbs.

Solid Seat Insert \$62

		Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Depth	14"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	15"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	16"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	17"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	21"	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	22"	☐	☐	☐	☐	☐	☐	☐	☐	☐	

Star Stabil-Air Cushions

Star Stabil-Air Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$549		\$69	\$80	\$549		\$69	\$80	\$549		\$69	\$80	\$549		\$69	\$80	\$675	\$91 \$103
HCPCS	E2622					E2622				E2622				E2622				E2623	
Depth	13.5"	☐	☐	☐															
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	17.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"							☐	☐	☐	☐	☐	☐	☐	☐	☐			

Star Stabil-Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$549		\$69	\$80	\$549		\$69	\$80	\$549		\$69	\$80	\$549		\$69	\$80	\$675	\$91 \$103
HCPCS	E2622					E2622				E2622				E2622				E2623	
Depth	13.5"	☐	☐	☐															
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	17.25"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Star Standard Air Cushions

Star Standard Air Cushion (2" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$475		\$69	\$80	\$475		\$69	\$80	\$475		\$69	\$80	\$475		\$69	\$80	\$580	\$91 \$103
HCPCS	E2622					E2622				E2622				E2622				E2623	
Depth	13.5"	☐	☐	☐															
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	17.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"							☐	☐	☐	☐	☐	☐	☐	☐	☐			

Star Standard Air Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"					
		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion		Cover		Cushion	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$475		\$69	\$80	\$475		\$69	\$80	\$475		\$69	\$80	\$475		\$69	\$80	\$475	\$69 \$80
HCPCS	E2622					E2622				E2622				E2622				E2622	
Depth	13.5"																		
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	17.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"							☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Star Standard Air Cushions continued

Star Standard Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
	Retail	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$580	\$91	\$103
	HCPCS	E2622			E2622			E2622			E2622			E2622			E2623		
Depth	13.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐									
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐			
	17.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"							☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Star Standard Air Cushions Dual Valve

Star Standard Air Cushion Dual Valve (4" Tall Cells) and Extra Standard / Incontinent Cover


		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
	Retail	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80
	HCPCS	E2624			E2624			E2624			E2624			E2624		
Depth	13.5"	☐	☐	☐												
	15.5"	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	17.25"				☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19.25"							☐	☐	☐	☐	☐	☐			
	21"										☐	☐	☐	☐	☐	☐

Axiom Back


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width					
		12"	14"	16"	18"	20"	22"
Retail		\$574	\$574	\$574	\$574	\$574	\$815
HCPC		E2613	E2613	E2613	E2613	E2613	E2614
Height	Low	8"	☐ AXPB1208	☐ AXPB1408	☐ AXPB1608	☐ AXPB1808	☐ AXPB2008
		10"		☐ AXPB1410	☐ AXPB1610	☐ AXPB1810	☐ AXPB2010
		12"	☐ AXPB1212	☐ AXPB1412	☐ AXPB1612	☐ AXPB1812	☐ AXPB2012
	Mid	14"		☐ AXPB1414	☐ AXPB1614	☐ AXPB1814	☐ AXPB2014
		16"	☐ AXPB1216	☐ AXPB1416	☐ AXPB1616	☐ AXPB1816	☐ AXPB2016
		18"		☐ AXPB1418	☐ AXPB1618	☐ AXPB1818	☐ AXPB2018



Axiom Posterior Back - Low




Axiom Posterior Back - Mid


		Axiom Posterior Back Measurement Reference			
		Back Height	A	B	Shell Contour Depth
Height	Low	8"	10"	6.5"	2.25"
		10"	12"	8.5"	
		12"	14"	10.5"	
	Mid	14"	16"	8.5"	
		16"	18"	9.5"	
		18"	20"	10.5"	

Axiom Posterior Lateral Back (PL)


		Width					
		12"	14"	16"	18"	20"	22"
Retail		\$656	\$656	\$656	\$656	\$656	\$892
HCPC		E2615	E2615	E2615	E2615	E2615	E2616
Height	Low	8"	☐ AXPL1208	☐ AXPL1408	☐ AXPL1608	☐ AXPL1808	☐ AXPL2008
		10"		☐ AXPL1410	☐ AXPL1610	☐ AXPL1810	☐ AXPL2010
		12"	☐ AXPL1212	☐ AXPL1412	☐ AXPL1612	☐ AXPL1812	☐ AXPL2012
	Mid	14"		☐ AXPL1414	☐ AXPL1614	☐ AXPL1814	☐ AXPL2014
		16"	☐ AXPL1216	☐ AXPL1416	☐ AXPL1616	☐ AXPL1816	☐ AXPL2016
		18"		☐ AXPL1418	☐ AXPL1618	☐ AXPL1818	☐ AXPL2018
	Tall	20"	☐ AXPL1220	☐ AXPL1420	☐ AXPL1620	☐ AXPL1820	☐ AXPL2020
		22"		☐ AXPL1422	☐ AXPL1622	☐ AXPL1822	☐ AXPL2022



Axiom Posterior Lateral Back - Low



Axiom Posterior Lateral Back - Mid



Axiom Posterior Lateral Back - Tall


Axiom Posterior Lateral Back Measurement Reference					
		Back Height	A	B	Shell Contour Depth
Height	Low	8"	10"	6.5"	3.65"
		10"	12"	8.5"	
		12"	14"	10.5"	
	Mid	14"	16"	8.5"	
		16"	18"	9.5"	
		18"	20"	10.5"	
	Tall	20"	22"	11.5"	
		22"	24"	12.5"	

Axiom Back continued


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Lateral Back (LB)


		Width						
		12"	14"	16"	18"	20"	22"	
		Retail HCPC	\$759 E2620	\$759 E2620	\$759 E2620	\$759 E2620	\$759 E2620	\$841 E2621
Height	Low	8"	☐ AXLB1208	☐ AXLB1408	☐ AXLB1608	☐ AXLB1808	☐ AXLB2008	
		10"		☐ AXLB1410	☐ AXLB1610	☐ AXLB1810	☐ AXLB2010	☐ AXLB2210
		12"	☐ AXLB1212	☐ AXLB1412	☐ AXLB1612	☐ AXLB1812	☐ AXLB2012	
	Mid	14"		☐ AXLB1414	☐ AXLB1614	☐ AXLB1814	☐ AXLB2014	☐ AXLB2214
		16"	☐ AXLB1216	☐ AXLB1416	☐ AXLB1616	☐ AXLB1816	☐ AXLB2016	
		18"		☐ AXLB1418	☐ AXLB1618	☐ AXLB1818	☐ AXLB2018	☐ AXLB2218
	Tall	20"	☐ AXLB1220	☐ AXLB1420	☐ AXLB1620	☐ AXLB1820	☐ AXLB2020	
		22"		☐ AXLB1422	☐ AXLB1622	☐ AXLB1822	☐ AXLB2022	☐ AXLB2222



Axiom Lateral Back - Low



Axiom Lateral Back - Mid




Axiom Lateral Back - Tall


		Axiom Lateral Back Measurement Reference			
		Back Height	A	B	Shell Contour Depth
Height	Low	8"	10"	6.5"	3.65"
		10"	12"	8.5"	
		12"	14"	10.5"	
	Mid	14"	16"	8.5"	
		16"	18"	9.5"	
		18"	20"	10.5"	
	Tall	20"	22"	11.5"	
		22"	24"	12.5"	

Axiom Deep Lateral Back (DL)

			Width					
			12"	14"	16"	18"	20"	22"
			Retail	\$759	\$759	\$759	\$759	\$759
HCPC			E2620	E2620	E2620	E2620	E2620	E2621
Height	Mid	14"	☐ AXDL1214	☐ AXDL1414	☐ AXDL1614	☐ AXDL1814	☐ AXDL2014	☐ AXDL2214
		16"		☐ AXDL1416	☐ AXDL1616	☐ AXDL1816	☐ AXDL2016	
		18"	☐ AXDL1218	☐ AXDL1418	☐ AXDL1618	☐ AXDL1818	☐ AXDL2018	☐ AXDL2218
	Tall	20"	☐ AXDL1220	☐ AXDL1420	☐ AXDL1620	☐ AXDL1820	☐ AXDL2020	
		22"		☐ AXDL1422	☐ AXDL1622	☐ AXDL1822	☐ AXDL2022	☐ AXDL2222



Axiom Deep Lateral Back - Mid



Axiom Deep Lateral Back - Tall


Axiom Deep Lateral Back Measurement Reference									
		Back Height	A	B ¹					Shell Contour Depth
Height	Mid	14"	16"	8.25"	8.25"	8.25"	8.25"	8.25"	6"
		16"	18"	9.25"	9.25"	9.25"	9.25"	9.25"	
		18"	20"	10.25"	10.25"	10.25"	10.25"	10.25"	
	Tall	20"	22"	9.5"	9.75"	10"	10.25"	10.25"	
		22"	24"	10.75"	11"	11.25"	11.5"	11.5"	

¹ Measurement based on back width

Axiom Back Accessories

Headrest and Hardware

Choose headrest and hardware combination in chart below. Not available with Low Backs.



				Hardware					
				Standard Straight Back Mount Hardware with Ball		Multi Axis Removable Hardware with Horizontal Bar		Headrest Hardware with Ball	
				HEDHDW101 / HEDHDW117		HEDHDW201 / HEDHDW217		HEDHDW401 / HEDHDW417	
	Description	Item Code	HCPC	Kio Code	Retail	Kio Code	Retail	Kio Code	Retail
Headrest Pad	Comfort Plus Headrest 8"	HEDPAD408	E0955 E1028	☐ HDR001	\$284	☐ HDR009	\$321	☐ HDR017	\$398
	Comfort Plus Headrest 10"	HEDPAD410	E0955 E1028	☐ HDR002	\$284	☐ HDR010	\$321	☐ HDR018	\$398
	Comfort Plus Headrest 14"	HEDPAD414	E0955 E1028	☐ HDR003	\$284	☐ HDR011	\$321	☐ HDR019	\$398
	Small All Positioning Headrest Pad and Ring Only	HEDPAD201	E0955 E1028	☐ HDR004	\$284	☐ HDR012	\$321	☐ HDR020	\$398
	Medium All Positioning Headrest Pad and Ring Only	HEDPAD202	E0955 E1028	☐ HDR005	\$284	☐ HDR013	\$321	☐ HDR021	\$398
	Large All Positioning Headrest Pad and Ring Only	HEDPAD203	E0955 E1028	☐ HDR006	\$284	☐ HDR014	\$321	☐ HDR022	\$398
	Tri-Comfort Plus Headrest 14"	HEDPAD502	E0955 E1028	☐ HDR007	\$296	☐ HDR015	\$333	☐ HDR023	\$409
	Tri-Comfort Plus Headrest 19"	HEDPAD503	E0955 E1028	☐ HDR008	\$296	☐ HDR016	\$333	☐ HDR024	\$409

Headrest Hardware Options

	Description	Item Code	HCPC	Retail (Each)
☐	Headrest Flip Down Mount <i>Only Available with HEDHDW401/HEDHDW417</i>	HDR100	E1028	\$175
☐	Universal Mounting Plate <i>For use with other manufacturer headrests</i>	UMP	K0108	\$41

Axiom Back Accessories continued

Lateral Support

Laterals for Posterior Back					Item Code	HCPC	Pair	Each
Pair	Left	Right	Description					
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back		LF3x4P	E0956	\$313	\$157
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back		LF4x5P	E0956	\$313	\$157
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back		LS3x4P	E0956, E1028	\$667	\$334
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back		LS4x5P	E0956, E1028	\$667	\$334
Laterals for Posterior Lateral Back					Item Code	HCPC	Pair	Each
Pair	Left	Right	Description					
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back		LF3x4L	E0956	\$313	\$157
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back		LF4x5L	E0956	\$313	\$157
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back		LS3x4L	E0956, E1028	\$667	\$334
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back		LS4x5L	E0956, E1028	\$667	\$334

Axiom Lateral Support Modifications

Extended lead times will apply. Modified Lateral Supports are not eligible for refund.

	Description	Item Code	HCPC	Retail
☐	Modification for Chest Strap attached to Lateral Supports	AXMLATMT	K0108	N/C

Lumbar Pads

Available only with Axiom Lateral Back and Axiom Deep Lateral back

	Description	HCPC	Retail
☐	Lumbar Pad	K0108	\$50

Privacy Flap

	Description	HCPC	Retail
☐	Privacy Flap	K0108	\$47

Shoulder Guides

	Description	Item Code	HCPC	Retail
☐	Shoulder Strap Guides and Mount <i>Selection of a modification from the Axiom Back Modifications required in order to mount Shoulder Guides. Not available with Low Backs.</i>	AXBSG	K0108	\$77

Axiom Back Accessories continued

Posture Support

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	HCPC	Retail
☐	Dynaform Chest Support: 2-Piece with Hook & Loop, 3" High Pad, Small	ATRP2P14941	E0960	\$121
☐	Dynaform Chest Support: 2-Piece with Hook & Loop, 4" High Pad, Medium	ATRP2P14942	E0960	\$121
☐	Dynaform Chest Support: 2-Piece with Hook & Loop, 5" High Pad, Large	ATRP2P14943	E0960	\$121
☐	Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	ATRD2P14924	E0960	\$100
☐	Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	ATRD2P14927	E0960	\$100
☐	Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	ATRD2P14928	E0960	\$100
☐	AirLogic Anterior Trunk Support: Standard, Structure, Small	ATRP4P14282	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Standard, Structure, Medium	ATRP4P14283	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Standard, Structure, Large	ATRP4P14284	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Standard, Structure, X-Large	ATRP4P14285	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	ATRD4P14257	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	ATRD4P14258	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	ATRD4P14259	E0960	\$164.17
☐	AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	ATRD4P14260	E0960	\$164.17

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Axiom Back Modifications

Choose modification to add holes to the back shell in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	HCPC	Retail
☐	Modification for 4 Point Harness	AXM4PT	K0108	N/C
☐	Modification for 2 Point Chest Strap	AXM2PT	K0108	N/C
☐	Modification for Shoulder Guides	AXMSG	K0108	N/C

Painted Back Shell Modification

If no selection is made, standard Shadow Black will be used with no additional lead time

	Gloss Finish	Item Code	Retail
☐	Black	AXMBLK	N/C
☐	Blue Skies	AXMSKY	N/C
☐	Bubblegum Pink	AXMPNK	N/C
☐	Candy Blue	AXMCBL	N/C
☐	Candy Purple	AXMPUR	N/C
☐	Candy Red	AXMRED	N/C
☐	Copper	AXMCOP	N/C
☐	Granny Smith	AXMGSM	N/C
☐	Metallic Orange	AXMMOR	N/C
☐	Silver	AXMSIL	N/C
☐	Teal	AXMTEA	N/C
☐	Yellow	AXMYLW	N/C

	Matte Finish	Item Code	Retail
☐	Charcoal Gray	AXMCGM	N/C
☐	Emerald Green	AXMEGM	N/C
☐	Iridescent Burgundy	AXMIBM	N/C
☐	Iridescent Navy	AXMINM	N/C
☐	Military Green	AXMAGM	N/C
☐	Retro Red	AXMRRM	N/C
☐	Royal Blue	AXMRBM	N/C
☐	Smooth White	AXMSWM	N/C