



ORDER FORM

Print pages 1-25 to exclude pictures.


Q500H with SEDEO® PRO™
and JAY® Syncra™ Seating

EIPW32

Customer	Ship To
Account #	Name
Date	Address
PO#	City
Buyer	State Zip
Mark For	Phone

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

HCPC	K0848 / K0850 / K0856 / K0858 / K0861 K0862 / K0868 / K0877 / K0884
Weights Starting At: (w/out Batteries)	NPC: 343 lbs. SPC: 367 lbs. MPC: 381 lbs.
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction / 12" seat elevator
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Drive Base Specifications

Drive Wheel Position	Hybrid (mid/rear) drive
User Weight Capacity	300 lbs. / 350 lbs. HD
Base Size	Width 24.5", Length 36"
Minimum Turning Radius	25.5" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A / VR-2 90A
Driving Range	Up to 17.7 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 5 mph HD, 8 mph
Drive Wheels	14" or 13" w/flat-free inserts (black) 14" or 13" w/pneumatic inserts (black)
Casters	9" front / 4" rear (black)
Charger	8A or 10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Color	7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Hybrid-Wheel Drive Base End Item Number EIPW32 (PRO™ Seat)			\$
320PC248	K0848	Quickie Q500 H NPC No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,262.00
320PC249	K0856	Quickie Q500 H SPC Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,262.00
320PC250	K0861	Quickie Q500 H MPC Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q500 H Base, R-net 90amp Expandable Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,262.00
Group 3 Heavy Duty Hybrid-Wheel Drive Base End Item Number EIPW32 (PRO™ Seat)			\$
320PC251	K0850	Quickie Q500 H NPC HD No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$10,630.00
320PC252	K0858	Quickie Q500 H SPC HD Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$10,630.00
320PC253	K0862	Quickie Q500 H MPC HD Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q500 H Base, R-net 90amp Expandable Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$10,630.00
Group 4 Standard Duty Hybrid-Wheel Drive Base End Item Number EIPW32 (PRO™ Seat)			\$
320PC254	K0868	Quickie Q500 H NPD No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, R-net 120amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$13,156.00
320PC255	K0877	Quickie Q500 H SPD Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 H Base, R-net 120amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$13,156.00
320PC256	K0884	Quickie Q500 H MPD Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q500 H Base, R-net 120amp Electronics, Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$13,156.00

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
User Weight			\$
320WG22	—	Up to 250 lbs.	NC
320WG03	—	251 to 300 lbs.	NC
320WG15	—	301 to 350 lbs.	NC
Motor Package			\$
320PF98	—	4-pole 8 kph / 5 mph <i>(Standard for Group 3 HD K0850, K0858, K0862)</i>	NC
320PF99	—	4-pole 10 kph / 6 mph <i>(Standard for Group 3 K0848, K0856, K0861)</i>	NC
320PF100	—	4-pole 12.5 kph / 8 mph <i>(Standard for Group 4 K0868, K0877, K0884)</i>	NC
320PF101	—	4-pole 12.5 kph / 8 mph <i>(Upgrade for Group 3 K0848, K0856, K0861)</i>	\$1,235.00
Tracking			\$
320MI01	K0108	SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,498.00
Drive Wheels  P. 27			\$
320DW58	—	Drive Wheel, 14", Dual Density, Knobby, Black	NC
320DW67	—	Drive Wheel, 14", Pneumatic, Knobby, Black	\$268.00
320DW69	—	Drive Wheel, 13", Dual Density, Knobby, Black	NC
320DW76	—	Drive Wheel, 13", Pneumatic, Knobby, Black	\$268.00
Caster Wheels  P. 27			\$
320CW57	—	Caster Wheels 9", Solid, Black	NC
320CW56	—	Caster Wheels 9", Pneumatic, Black	\$209.00
Wheel Locks			\$
320WL42	—	Wheel Locks	\$113.00
Color Accents			\$
320EZ20	—	Apple Green	NC
320EZ22	—	Black	NC
320EZ21	—	Blue	NC
320EZ23	—	Orange	NC
320EZ24	—	Red	NC
320EZ19	—	White	NC
320EZ25	—	Silver	NC

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Service Decals			\$
320D10	—	Army	NC
320D11	—	Air Force	NC
320D13	—	Marines	NC
320D14	—	Navy	NC
320D15	—	Purple Heart	NC
ISO Transit Tie Downs			\$
320T050	—	ISO Transit Package	NC
LED Lights  P. 27			\$
320LI04	—	LED Lights and Indicators	\$1,288.00
320LI02	—	Omit LED Lights	NC
Seat Heights			\$
320HT103	—	16.5"	NC
320HT55	—	17"	NC
320HT104	—	17.5"	NC
320HT56	—	18"	NC
320HT105	—	18.5"	NC
320HT57	—	19"	NC
Batteries and Chargers  P. 27			\$
Group 24 Battery Options (Recommended for max daily range)			
320AG25	E2363	GP24 Flat Top Pair w/Strap MK Gel	\$1,047.00
320AG26	E2363	GP24 Flat Top Pair w/Strap Brightway AGM	\$789.00
Group 34 Battery Options (Required for lower seat-to-floor heights)			
320AG27	E2359	GP34 Flat Top Pair w/Straps MK Gel	\$1,009.00
320AG29	E2359	GP34 Flat Top Pair w/Straps Brightway AGM	\$762.00
Omit Batteries			
320AG24	—	Battery Harnesses Only <i>(Omit Batteries)</i>	NC
Chargers			
320BJ30	—	10A Dual Mode (Compact)	\$102.00
320BJ19	—	8A Dual Mode	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options  P. 27			\$
Seat Elevator			
320PI07	E2300	Seat Elevator, Powered, 300 mm / 12" (Not available with HD or 8 mph motors)	\$3,489.00
Seat Tilt			
320PA35	E1002	Posterior 50° Tilt	\$7,783.00
320PA36	E1002	Posterior 50° HD Tilt	\$8,857.00
Seat Recline			
320PA37	E1004	150° Powered Recline w/ 4" Shear Reduction	\$7,783.00
320PA39	E1004	150° HD Powered Recline w/ 4" Shear Reduction	\$8,857.00
Seat Tilt/Recline			
320PA41	E1007	Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	\$12,614.00
320PA43	E1007	Posterior 50° HD Tilt and 150° HD Recline Powered, w/ 4" Shear Reduction	\$13,688.00
Power Recline Backrest Preset Home Angle			
320QK1	—	90°	NC
320QK2	—	95° (Default Position)	NC
320QK3	—	100°	NC
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			\$
320FA10	—	0° (Adds 0" to STFH)	NC
320FA3	—	3° (Adds 1" to STFH)	NC
320FA4	—	6° (Adds 2" to STFH)	NC
320FA5	—	9° (Adds 3" to STFH)	NC
Backrest Frame Type			\$
320BF10	—	Tool Angle Adjustable Back, 90°-120°	NC
320BF11	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Right (Adjustment lever mounted on opposite side of joystick)	\$628.00
320BF12	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Left (Adjustment lever mounted on opposite side of joystick)	\$628.00
Backrest Frame Option			\$
320BF13	—	Flip Down, Quick Release	\$135.00
Back Cushion Type			\$
320BT134	—	SEDEO® Back Cushion (Choose Seating in Section A)	NC
320BT138	—	JAY® Syncra™ (Choose Seating in Section B)	NC
320BT135	—	JAY® J3® (Choose Seating in Section C)	NC
320BT136	—	JAY® J3® Interface Only (Allows JAY Your Way modifications)	NC
320BT137	—	SEDEO® PRO™ Aftermarket Back Adapter (Height adjustable round tubes for traditional mounting)	\$483.00
320BT133	—	OMIT (For personalized mounting. No headrest mount)	NC

POWERED SEATING					
OPTION	HCPC	DESCRIPTION	RETAIL		
A - SEDEO® Backrest Dimensions					
Back Width					
320BN28	—	400mm / 16"	NC		
320BN29	—	460mm / 18"	NC		
320BN30	—	500mm / 20"	NC		
320BN31	—	560mm / 22"	NC		
Back Shell Height			Tool Angle Adjust.	150° Recline	
320B255	—	460mm / 18"	22"-23"	23"-25"	Installed Back Height
320B256	—	560mm / 22"	24"-27"	26"-29"	
320B257	—	660mm / 26"	28"-30"	30"	
Installed Back Height - (From Pan to Top of Back)					
320JR32	—	22"	NC		
320JR33	—	23"	NC		
320JR34	—	24"	NC		
320JR35	—	25"	NC		
320JR36	—	26"	NC		
320JR37	—	27"	NC		
320JR38	—	28"	NC		
320JR39	—	29"	NC		
320JR40	—	30"	NC		
Lateral Thoracic Support Pads					
320KD31	E0956	Lateral, Swing-Away, Adjustable, Long Pad 150 x 130 mm / 6" x 5" , Right	\$140.00		
320KD32	E0956	Lateral, Swing-Away, Adjustable, Long Pad 150 x 130 mm / 6" x 5" , Left	\$140.00		
320KD33	E0956	Lateral, Swing-Away, Adjustable, Short Pad 90 x 120 mm / 3.5" x 4.5" , Right	\$140.00		
320KD34	E0956	Lateral, Swing-Away, Adjustable, Short Pad 90 x 120 mm / 3.5" x 4.5" , Left	\$140.00		
Lateral Thoracic Swing-Away Hardware					
320KE45	E1028	Lateral Swing-Away Hardware, Right	\$247.00		
320KE46	E1028	Lateral Swing-Away Hardware, Left	\$247.00		

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
B - JAY® Syncra™ Seating			
B - JAY® Syncra™ Back  P. 26			
320JP14	E2620	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	\$1,047.00
320JP15	E2621	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	\$1,155.00
B - JAY® Syncra™ Back Width			
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>			
320JQ07	—	16" Wide (Adjustment range 16 to 17" Wide)*	NC
320JQ08	—	18" Wide (Adjustment range 18 to 19" Wide)*	NC
320JQ09	—	20" Wide (Adjustment range 20 to 21" Wide)*	NC
320JQ10	—	22" Wide (Adjustment range 22 to 23" Wide)	NC
<i>*Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.</i>			
B - JAY® Syncra™ Back Shell Options and Height Ranges			
<i>JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.</i>			
JAY® Syncra™ Shell Height			
320JR42	—	16" Tall Shell (22" to 24" Back Height)	NC
320JR43	—	18" Tall Shell (24" to 26" Back Height)	NC
320JR44	—	20" Tall Shell (26" to 28" Back Height)	NC
320JR45	—	22" Tall Shell (28" to 30" Back Height)	NC
JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)			
320JR50	—	22" (16" Shell Height)	NC
320JR51	—	23" (16" Shell Height)	NC
320JR52	—	24" (16" Shell Height)	NC
320JR53	—	24" (18" Shell Height)	NC
320JR54	—	25" (18" Shell Height)	NC
320JR55	—	26" (18" Shell Height)	NC
320JR56	—	26" (20" Shell Height)	NC
320JR57	—	27" (20" Shell Height)	NC
320JR58	—	28" (20" Shell Height)	NC
320JR67	—	28" (22" Shell Height)	NC
320JR68	—	29" (22" Shell Height)	NC
320JR69	—	30" (22" Shell Height)	NC
B - JAY® Syncra™ Back Upholstery			
320BU96	—	Syncra™ Premium Leatherette	NC
320BU97	—	Syncra™ Stretch	NC
320BU98	—	Syncra™ 3DX Microclimatic	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

B - JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

\$

JAY® Lateral Thoracic Support Pads

320JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$135.00
320JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$135.00
320JW36	E0956	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$135.00
320JW39	E0956	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	\$135.00
320JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$135.00
320JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$135.00

JAY® Lateral Thoracic Support Pad Upholstery

320NU01	—	Syncra™ Premium Leatherette	NC
320NU02	—	Syncra™ Stretch	NC
320NU03	—	Syncra™ 3DX Microclimatic	NC

JAY® Lateral Thoracic Swing-Away Mount

320JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$268.00
320JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$268.00

B - JAY® Syncra™ Fixed Lateral Support Accessory

\$

JAY® Fixed Lateral Support Accessory Hardware (JAY® Fixed Lateral Support is height adjustable)

320JW33	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	\$135.00
320JW34	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	\$135.00

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

C - JAY® J3® Back**C - JAY® J3® Back - Contour** **P. 27** **\$**

320JP01	E2613	J3® PA Back for 14"-20" Wide	\$575.00
320JP02	E2614	J3® PA Back for 22" Wide	\$854.00
320JP11	E2615	J3® PL Back for 14"-20" Wide	\$671.00
320JP13	E2620	J3® PD Back for 14"-20" Wide (18" and 21" shell height only)	\$805.00

C - JAY® J3® Back Width **\$**

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

320JQ02	—	16"	NC
320JQ03	—	18"	NC
320JQ04	—	20"	NC
320JQ05	—	22"	NC

C - JAY® J3® Back Shell Options and Height Ranges **\$**

JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown, based on the back frame type and the amount of usable hardware adjustability.

JAY® J3® Shell Height			Tool Angle Adjust.	150° Recline	Installed Back Height	
320JR41	—	15" tall shell (Upper Thoracic Med) for 22" Wide	22"-23"	22"-23"		NC
320JR14	—	16.5" tall shell (Upper Thoracic Short)	21"-23"	22"-24"		NC
320JR15	—	18" tall shell (Upper Thoracic Med)	22"-24"	23"-25"		NC
320JR16	—	19.5" tall shell (Upper Thoracic Tall)	24"-26"	24"-26"		NC
320JR18	—	18" tall shell (Shoulder Height) for 22" wide	24"	24"		NC
320JR17	—	21" tall shell (Shoulder Height)	24"-26"	25"-28"		NC
320JR27	—	22.5" tall shell (Shoulder Height)	26"-28"	27"-30"		NC
320JR28	—	24" tall shell (Shoulder Height)	27"-29"	28"-30"		NC

JAY® J3® Installed Back Height (From Pan to Top of Back)

320JR29	—	21"	NC
320JR19	—	22"	NC
320JR20	—	23"	NC
320JR21	—	24"	NC
320JR22	—	25"	NC
320JR23	—	26"	NC
320JR24	—	27"	NC
320JR25	—	28"	NC
320JR30	—	29"	NC
320JR31	—	30"	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

C - JAY® J3® Back Spine-Align Fit Kit 2 P. 27 \$

320JS02	K0108	410 mm / 16"	\$149.00
320JS03	K0108	460 mm / 18"	\$149.00
320JS04	K0108	510 mm / 20"	\$149.00

C - JAY® Anterior Thoracic Support P. 27 \$

Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)

320K10	E0960	Medium (ATS8302)	\$163.00
320K11	E0960	Large (ATS8303)	\$163.00
320K12	E0960	Extra Large (ATS8304)	\$163.00

C - JAY® J3® Lateral Thoracic Swing-Away Supports \$
Lateral Thoracic Support Pads

320JW19	E0956	Lateral Thoracic (LxH) 100 mm x 100 mm / 4" x 4" Pad - Right	\$135.00
320JW25	E0956	Lateral Thoracic (LxH) 100 mm x 100 mm / 4" x 4" Pad - Left	\$135.00
320JW20	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Right	\$135.00
320JW26	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Left	\$135.00
320JW21	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$135.00
320JW27	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$135.00
320JW22	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$135.00
320JW28	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$135.00
320JW23	E0956	Lateral Thoracic (LxH) 152 mm x 127 mm / 6" x 5" Pad - Right	\$135.00
320JW29	E0956	Lateral Thoracic (LxH) 152 mm x 127 mm / 6" x 5" Pad - Left	\$135.00
320JW24	E0956	Lateral Thoracic (LxH) 152 mm x 152 mm / 6" x 6" Pad - Right	\$135.00
320JW30	E0956	Lateral Thoracic (LxH) 152 mm x 152 mm / 6" x 6" Pad - Left	\$135.00

Lateral Thoracic Swing-Away Hardware

320JW31	E1028	Lateral Thoracic Swing-Away Bracket - Right	\$282.00
320JW32	E1028	Lateral Thoracic Swing-Away Bracket - Left	\$282.00

Headrests (Requires Pad, Upholstery, & Mount) P. 29 \$
Headrest Pad

320HE36	E0955	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$173.00
320HE37	E0955	SEDEO® Pad Flexible, (LxH) 400 x 110 mm / 16" x 4.5" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$279.00
320HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$220.00
320HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$220.00
320HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$230.00
320HE45	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Small (Narrow) (Adjustable from 6 to 14" Wide x 3" Tall)	\$338.00
320HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$338.00
320HE47	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Large (Standard) (Adjustable from 8 to 18" Wide x 5" Tall)	\$338.00
320HE48	—	OMIT Headrest Pad	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Headrest Pad Upholstery			
320NU10	—	SEDEO® Upholstery	NC
320NU04	—	Syncra™ Premium Leatherette	NC
320NU05	—	Syncra™ Stretch	NC
320NU06	—	Syncra™ 3DX Microclimatic	NC
320NU11	—	Whitmyer® Lycra	NC
Headrest Mount			
320MH32	E1028	SEDEO® Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO® Pad)</i>	\$215.00
320MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$177.00
320MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$215.00
320MH35	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$113.00
320MH36	E1028	Whitmyer® AXYS Offset Headrest Mount with Quick Release	\$135.00
320MH37	E1028	Whitmyer® Cobra FB Headrest Mount with Quick Release	\$322.00
320MH17	—	OMIT Headrest Mount	NC
Push Handles			\$
320PH27	—	Height Adjustable Push Handles	\$322.00
Seat Frame Rails			\$
320SN14	—	Seat Rail Short <i>(For effective seat depth 14" to 18")</i>	NC
320SN15	—	Seat Rail Long <i>(For effective seat depth 17" to 20")</i>	NC
Seat Cushion Type			\$
320ST79	—	SEDEO® (See Section D)	NC
320ST81	—	JAY® Syncra™ Cushion (See Section E)	NC
320ST80	—	JAY® (See Section F)	NC
320ST75	—	Omit	NC
C - Seat Size			\$
Seat Pan Width <i>(Odd sizes not available on SEDEO® cushions)</i>			
320R53	—	16"	NC
320R54	—	17"	NC
320R55	—	18"	NC
320R56	—	19"	NC
320R57	—	20"	NC
320R58	—	21" <i>(Not available with JAY® Syncra™ Cushions)</i>	NC
320R59	—	22"	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
Seat Pan Depth (Front of seat pan to front of back cushion – effective depth)			
320RF52	–	14"	NC
320RF53	–	15"	NC
320RF51	–	16"	NC
320RF45	–	17"	NC
320RF46	–	18"	NC
320RF47	–	19"	NC
320RF48	–	20"	NC
320RF49	–	21" (Not available with JAY® Syncra™ Cushions)	NC
320RF50	–	22"	NC

D - SEDEO® Cushions P. 27 \$

320CU82	E2607	SEDEO® Cushion (Sizes 16", 18", & 20" Wide) (Not available with 14" Depth)	\$424.00
320CU83	E2608	SEDEO® Cushion (Size 22" Wide)	\$532.00

E - JAY® Syncra™ Cushion P. 26 (JAY Your Way modifications available through Customer Service (800) 333-4000) \$

320CU87	E2607	JAY® Syncra™ Cushion with Visco Foam (16" to 20" wide)	\$445.00
320CU88	E2608	JAY® Syncra™ Cushion with Visco Foam (22" wide)	\$526.00
320CU89	E2622	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (16" to 20" wide)	\$580.00
320CU90	E2623	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (22" wide)	\$671.00
320CU91	E2622	JAY® Syncra™ Cushion with Air Insert - Single Valve (16" to 20" wide)	\$532.00
320CU92	E2623	JAY® Syncra™ Cushion with Air Insert - Single Valve (22" wide)	\$623.00
320CU93	E2622	JAY® Syncra™ Cushion with Air Insert - Dual Valve (16" to 20" wide)	\$537.00
320CU94	E2623	JAY® Syncra™ Cushion with Air Insert - Dual Valve (22" wide)	\$628.00

JAY® Syncra™ Cushion Cover Upholstery

320YR73	–	Syncra™ Premium Leatherette	NC
320YR74	–	Syncra™ Stretch	NC
320YR75	–	Syncra™ 3DX Microclimatic	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
F - JAY® Seat Cushion		 P. 28	\$
320CU55	E2622	JAY® X2® Cushion <i>(Sizes 16" - 20" Wide)</i>	\$457.00
320CU56	E2622	JAY® J2® Cushion <i>(Sizes 16" - 20" Wide)</i>	\$467.00
320CU57	E2622	JAY® J2® Deep Contour Cushion <i>(Sizes 16" - 20" Wide)</i>	\$483.00
320CU58	E2623	JAY® J2® Deep Contour Cushion <i>(Size 22" Wide)</i>	\$596.00
320CU59	E2607	JAY® Ion® <i>(Sizes 16" - 21" Wide)</i>	\$424.00
320CU60	E2608	JAY® Ion® <i>(Size 22" Wide)</i>	\$532.00
320CU61	E2607	JAY® Union® <i>(Sizes 16" - 21" Wide)</i>	\$478.00
320CU62	E2608	JAY® Union® <i>(Size 22" Wide)</i>	\$585.00
320CU63	E2607	JAY® Union® Reduced Profile <i>(Sizes 16" - 20" Wide)</i>	\$478.00
320CU64	E2622	JAY® Fusion® <i>(Sizes 16" - 21" Wide)</i>	\$473.00
320CU65	E2623	JAY® Fusion® <i>(Size 22" Wide)</i>	\$564.00
320CU66	E2622	JAY® Fusion® Air Insert Single Valve <i>(Sizes 16" - 21" Wide)</i>	\$532.00
320CU67	E2623	JAY® Fusion® Air Insert Single Valve <i>(Size 22" Wide)</i>	\$623.00
320CU68	E2622	JAY® Fusion® Air Insert Dual Valve <i>(Sizes 16" - 21" Wide)</i>	\$537.00
320CU69	E2623	JAY® Fusion® Air Insert Dual Valve <i>(Size 22" Wide)</i>	\$628.00
320CU70	E2622	JAY® Fusion® Reduced Profile <i>(Sizes 16" - 21" Wide)</i>	\$473.00
320CU71	E2623	JAY® Fusion® Reduced Profile <i>(Size 22" Wide)</i>	\$564.00
320CU72	E2622	JAY® Fusion® Reduced Profile Air Insert Single Valve <i>(Sizes 16" - 21" Wide)</i>	\$532.00
320CU73	E2623	JAY® Fusion® Reduced Profile Air Insert Single Valve <i>(Size 22" Wide)</i>	\$623.00
320CU74	E2622	JAY® Fusion® Reduced Profile Air Insert Dual Valve <i>(Sizes 16" - 21" Wide)</i>	\$537.00
320CU75	E2623	JAY® Fusion® Reduced Profile Air Insert Dual Valve <i>(Size 22" Wide)</i>	\$628.00
320CU76	E2624	JAY® Fusion® Positioning <i>(Sizes 16" - 21" Wide)</i>	\$494.00
320CU77	E2625	JAY® Fusion® Positioning <i>(Size 22" Wide)</i>	\$585.00
320CU78	E2624	JAY® Fusion® Positioning Air Insert Single Valve <i>(Sizes 16" - 21" Wide)</i>	\$553.00
320CU79	E2625	JAY® Fusion® Positioning Air Insert Single Valve <i>(Size 22" Wide)</i>	\$644.00
320CU80	E2624	JAY® Fusion® Positioning Air Insert Dual Valve <i>(Sizes 16" - 21" Wide)</i>	\$558.00
320CU81	E2625	JAY® Fusion® Positioning Air Insert Dual Valve <i>(Size 22" Wide)</i>	\$649.00

F - JAY® Solid Seat	\$
----------------------------	----

320CU84	—	JAY® Solid Seat	\$204.00
---------	---	-----------------	----------

F - JAY® Bundling	\$
--------------------------	----

Note: To order any JAY® cushion not listed above please see the JAY® price list.

320JC20	—	JAY® Bundled Cushion	See Price List
		Cushion name or part number	Width Depth

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Knee and Thigh Swing-Away Support Pads

\$

Lateral Knee Support Pads and Hardware *(must pick pad and mount)*

320XY05	E0953	Lateral Knee Support Pad, Right (4" round)	\$188.00
320XY06	E0953	Lateral Knee Support Pad, Left (4" round)	\$188.00
320XX08	E1028	Lateral Swing-Away Bracket Knee - Right	\$247.00
320XX09	E1028	Lateral Swing-Away Bracket Knee - Left	\$247.00

Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*

300KF22	E0953	Lateral Thigh Support Pad Short, Right (4.8" H x 4.8" L / 120mm x 120mm)	\$135.00
300KF23	E0953	Lateral Thigh Support Pad Short, Left (4.8" H x 4.8" L / 120mm x 120mm)	\$135.00
300KF24	E0953	Lateral Thigh Support Pad Long, Right (4.8" H x 7" L / 120mm x 180mm)	\$135.00
300KF25	E0953	Lateral Thigh Support Pad, Long, Left (4.8" H x 7" L / 120mm x 180mm)	\$135.00
300KF26	E1028	Lateral Thigh Removable Mount - Right	\$268.00
300KF27	E1028	Lateral Thigh Removable Mount - Left	\$268.00

JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount)

\$

JAY® Lateral Thigh Support Pads

320KF16	E0953	Lateral Thigh Support Pad Standard, Right (8.25" x 3.25")	\$135.00
320KF17	E0953	Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")	\$135.00
320KF18	E0953	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$135.00
320KF19	E0953	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$135.00

JAY® Lateral Thigh Support Pad Upholstery

320NU07	—	Syncra™ Premium Leatherette	NC
320NU08	—	Syncra™ Stretch	NC
320NU09	—	Syncra™ 3DX Microclimatic	NC

JAY® Lateral Thigh Removable Mount

320KF20	E1028	Lateral Thigh Removable Mount - Right	\$268.00
320KF21	E1028	Lateral Thigh Removable Mount - Left	\$268.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Armrest Structure  P. 28			\$
320EG35	—	Reclining, Flip Back, Long (<i>Height Adjustable 9" to 15"</i>)	NC
320EG36	—	Reclining, Flip Back with 2-point post locking arm (<i>Height Adjustable 11" to 15" x .625"</i>)	\$268.00
320EG32	—	Single Post, Flip Back (<i>N/A with thigh support</i>), (<i>Height Adjustable 10" to 15"</i>)	NC
320EG33	—	Single Post, Removable (<i>N/A with thigh support</i>), (<i>Height Adjustable 10" to 15"</i>)	NC
Arm Pads  P. 28			\$
SEDEO® Armrest Pad			
320AB43	—	300 mm / 12", Integral Foam, Right	NC
320AE43	—	300 mm / 12", Integral Foam, Left	NC
320AB44	—	400 mm / 16", Integral Foam, Right	NC
320AE44	—	400 mm / 16", Integral Foam, Left	NC
320AB45	—	Wide, 450 mm / 18", Integral Foam, Right (<i>Increases armpad width by 1.5". May prevent flip-back.</i>)	NC
320AE45	—	Wide, 450 mm / 18", Integral Foam, Left (<i>Increases armpad width by 1.5". May prevent flip-back.</i>)	NC
Comfort Arm Pad (Foam)			
320AB46	—	2" Wide Right	\$81.00
320AE46	—	2" Wide Left	\$81.00
320AB47	—	4" Wide Right	\$81.00
320AE47	—	4" Wide Left	\$81.00
Comfort Arm Pad (Gel)			
320AB48	—	2" Wide Right	\$161.00
320AE48	—	2" Wide Left	\$161.00
320AB49	—	4" Wide Right	\$161.00
320AE49	—	4" Wide Left	\$161.00
AES Trough/Elbow Block Arm Pads (Gel)			
320AB50	E2209	3.5" x 13" Right	\$268.00
320AE50	E2209	3.5" x 13" Left	\$268.00
320AB51	E2209	4" x 14.5" Right	\$268.00
320AE51	E2209	4" x 14.5" Left	\$268.00
Side Guards (Pair)			\$
320SG43	E0953	Side Guards (9"L x 5.5"H) (<i>N/A with Reclining Armrests</i>)	\$188.00
320SG44	E0953	Side Guards Integral Foam (11"L x 5.5"H) (<i>N/A with Reclining Armrests</i>)	\$220.00
320SG45	E0953	Side Guards Wedge Shaped, Padded (18"L x 5"H, 3" thick in front, 1.5" thick in rear) (<i>Only with Single Post Removable Armrests</i>)	\$220.00
320SG46	E0953	Side Guards Mounted On Seat Rail Bracket (<i>Only with Reclining Armrests</i>) (<i>N/A with Thigh Supports</i>)	\$242.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
SEDEO® PRO™ Swing-Away Legrest (Pair)			\$
320H177	–	Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length)	NC
320H178	E0990	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$376.00
320H179	E1010	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$2,308.00
320H180	–	Hanger Receivers Only, SEDEO® PRO™	NC
320DI76	–	Omit Legrests	NC
Swing-Away Legrest Footplates			
320F369	K0040	Angle Adjustable Flip-Up Footplates	\$182.00
Swing-Away Legrest Calf Supports			
320DI41	–	Calf Support, Depth Adjustable, Integral Foam (7" Wide x 10" Tall, NA 16-17" W)	\$317.00
320DI42	–	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$424.00
Swing-Away Legrest Options			
320F364	K0108	Footplates, Coupled	\$268.00
320DI44	E0951	Heel Loops (Required)	\$54.00
320DI40	–	Calf Strap (Only For 70° Swing-Away Legrests)	\$107.00
320KG4	–	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$215.00
SEDEO® Lower Extremity Swing-Away Support			
320VD01	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$424.00
320VD02	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$424.00
Pelvic Positioning Belts  P. 29			
320A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
320A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
320A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$81.00
Accessories			\$
Ventilator Option			
320VT12	E1029	Vent Tray Fixed (Pulmonetic LTV Models)	\$499.00
O2 Tank Option			
320A718	E2208	O2 Tank Holder Mounted to Backrest (D and E Styles)	\$349.00
Wheelchair Tray Table			
3207A09	E0950	Tray, Swing-Away	\$376.00
Backpack			
320A1	–	Black Utility	\$59.00
Medial Thigh Support			
3208A80	E0957	Medial Thigh Support	\$182.00
3208A81	E1028	Pommel Swing Away Hardware	\$268.00
Misc. Accessories			
3208A59	–	Transfer Bar (Pair)	\$376.00
3208A52	–	Ram Mount Phone	\$92.00
3208A53	–	Ram Mount Tablet	\$376.00
3208A54	–	Ram Mount Cup Holder	\$92.00
3208A58	–	Rear View Mirror	\$92.00
320A611	–	Crutch Holder	\$135.00
3208A73	–	Bag For Armrest	\$92.00
3208A75	–	Reflectors (Backrest shroud mounted)	\$89.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
Joystick (Remote)  P. 29			\$
320J127	—	VR-2 Joystick	NC
320J121	—	Standard LED Joystick (w/o Screen Display)	NC
320J116	—	Color Joystick (2.3" Color LCD Screen Display, CJS1M1)	NC
320J117	—	Heavy Duty Color Joystick (2.3" Color LCD Screen Display, Fixed Joystick Mount Only)	\$1,235.00
320J118	—	High Resolution Color Joystick, (CJS1M2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,132.00
320J119	—	Omit Remote Control (Required with arm mounted specialty controls)	—\$25.00
Joystick Mounting			\$
320JM34	—	Fixed on Armrest	NC
320JM35	E1028	Parallel Swing-Away on Armrest	\$317.00
320JM36	E1028	Link-It Mid Line Mount on Armrest (Swing-Away)	\$913.00
Joystick Mount (Location)			\$
320HM16	—	Right Hand Armrest	NC
320HM17	—	Left Hand Armrest	NC
Display Language Selection R-net			\$
320L1	—	English	NC
320L3	—	French	NC
320L7	—	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick  P. 29			\$
320JH13	E2323	Ball Handle	\$135.00
320JH14	E2323	Foam Ball Handle	\$135.00
320JH15	E2323	T-Handle	\$135.00
320JH16	E2323	Goalpost	\$135.00
Electronic Programmers  P. 30			\$
320PR10	—	R-net PC Programmer w/ Dongle	\$564.00
320PR11	—	VR-2 Hand Held Programmer	\$370.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Power Modules	P. 30	\$
---------------	---	----

Non-Expandable

320ER12	—	Power Module, VR-2, 90A, Non-Expandable <i>(Only available with 6mph motors)</i>	NC
320ER18	—	Power Module, R-net, 90A, Non-Expandable <i>(Only available with 5mph and 8mph)</i>	NC
320ER11	—	Power Module, R-net, 120A, Non-Expandable <i>(Only available with Group 4)</i>	NC

Expandable

320ER13	—	Power Module, R-net, 90A, Expandable <i>(Required for Specialty Controls or 3+ Actuators)</i>	\$1,229.00
320EX17	E2377	Expandable 90A Controller - R-net	\$649.00
320EX18	E2313	Harness for 90A Expandable Controller	\$580.00
320ER10	—	Power Module, R-net, 120A, Expandable <i>(Only available with Group 4)</i> <i>(Required for Specialty Controls or 3+ Actuators)</i>	\$1,337.00
320EX09	E2377	Expandable 120A Controller - R-net	\$703.00
320EX10	E2313	Harness for 120A Expandable Controller	\$318.00

Seating Control	P. 30	\$
-----------------	---	----

Powered Seating Control

320E040	—	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
320E037	E2310	Thru-drive Control for 1 Actuator	\$2,061.00
320E038	E2311	Thru-drive Control for 2+ Actuators	\$3,731.00

Switch Control *(Available with Omit Thru-Drive only) (1 Actuator for Tilt or Elevate only)*

320E033	—	Dual Toggle for 1 Actuator	NC
320E034	—	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
320E035	—	Switch-It Mini 1 button for 1 Actuator	NC
320E036	—	Switch-It Mini 2 button for 1 Actuator	NC
320E041	—	Ctrl+5 Switch Box for 2 Actuators or more	NC

Additional Power Seating Control *(Available with Thru-Drive only)*

320QG11	—	Ctrl+5 Switch Box With Free Assignable Functions	\$1,342.00
320QG12	—	Ctrl+5 Switch Box With Free Assignable Functions <i>(Multi-power only)</i>	NC

Dual Toggle Mounting Location

320TN1	—	Toggle Mount, Right	NC
320TN2	—	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

320EY19	—	Right Armrest	NC
320EY20	—	Left Armrest	NC
320EY21	—	Back Cane, Right	\$161.00
320EY24	—	Back Cane, Left	\$161.00

Ctrl+5 Switch Type

320EY22	—	Push Buttons	NC
320EY23	K0108	Toggle Switches	\$322.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
Attendant Controls  P. 30			\$
3200A05	E2331	Dual Control Module for Attendant	\$1,396.00
3200A06	E2331	R-net Color LCD (CJSM1) Control Attendant Only	\$1,396.00
Attendant Mounting Location			\$
320DN5	—	Right Back Cane	NC
320DN6	—	Left Back Cane	NC
R-net Programming			\$
Assignable Buttons			
<i>Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC</i>			
320PM12	—	Include Assignable Buttons	NC
R-net Omni2 Options  P. 30 <i>For Expandable Controls, choose all that apply.</i>			\$
3200N14	—	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$1,906.00
3200N15	—	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	NC
Omni2 Mount <i>If used w/ Standard Joystick, it must be mounted opposite side.</i>			\$
320DS14	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$532.00
320DS15	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$532.00

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Hand Control Packages  P. 31 *Requires Standard Joystick Omit / Only one hand control device can be ordered.*

\$

Proportional Input Devices

320EU14	E2312	Switch-It MicroGuide	\$3,537.00
320EU15	E2312	Switch-It MicroPilot	\$3,537.00
320EU16	K0108	Switch-It TouchDrive 2	\$4,128.00
320EU17	E2328	Switch-It VersaGuide	\$1,948.00
320EU18	E2328	Switch-It VersaGuide EZ	\$2,303.00

Mounting Position (Required)

320EM04	E1028	Link-It to the Backrest <i>(Not Available w/ Omit Back Cushion Type)</i>	\$1,020.00
320EM05	E1028	Link-It to the Fixed Joystick Mount	\$317.00
320EM06	E1028	Link-It to the Parallel Swing-Away Joystick Mount	\$532.00
320EM07	E1028	Link-It to the Midline Swing-Away	\$1,127.00
320EM08	E1028	Omit Mount	NC

Mounting Platforms (Optional) *Not available w/ Link-It to Backrest*

320MG16	E0950	SunDisk Hand Tray <i>(Only available w/ MicroSeries Joysticks)</i>	\$317.00
320MG15	E0950	Bullet Hand Tray	\$317.00

Mounting Location (Required)

320IM09	—	Right	NC
320IM10	—	Left	NC

Switch Options

—	—	See Switch Selections <i>(Maximum 2)</i>	—
---	---	--	---

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Head Control Packages P. 30 *(Select headrest mount on next page)*

\$

Only one head control device can be ordered.

Proportional Input Devices

320EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads	\$6,135.00
320EU44	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads - Dedicated Mode Outer Left	\$6,436.00
320EU45	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads - Dedicated Mode Outer Right	\$6,436.00
320EU46	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads - Dedicated Mode & Reverse	\$6,651.00
320EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads	\$6,135.00
320EU27	E2328	Switch-It RIM Proportional Head Control	\$6,329.00

Switched Input Devices

320EU28	E2330	Switch-It 3-Switch Head Array - Short Curved Wings	\$4,718.00
320EU29	E2330	Switch-It 3-Switch Head Array - Short Straight Wings	\$4,718.00
320EU30	E2330	Switch-It 3-Switch Head Array - Long Curved Wings	\$4,718.00
320EU31	E2330	Switch-It 3-Switch Head Array - Long Straight Wings	\$4,718.00
320EU42	E2330	Switch-It 4-Switch Head Array - Long Straight Wings Only - Mode Outer Left	\$4,933.00
320EU43	E2330	Switch-It 4-Switch Head Array - Long Straight Wings Only - Mode Outer Right	\$4,933.00
320EU33	E2330	Switch-It 5-Switch Head Array - Long Straight Wings Only - Mode Outer Left	\$5,073.00
320EU47	E2330	Switch-It 5-Switch Head Array - Long Straight Wings Only - Mode Outer Right	\$5,073.00

Switch Options

—	—	See Switch Selections <i>(Maximum 2)</i>	—
---	---	--	---

Sip 'N Puff

\$

320ID73	E2330	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i>	\$5,362.00
320ID74	K0108	Switch-It OptiStop with Sip 'N Puff	\$2,088.00

Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)*

\$

320IJ13	E1028	Whitmyer® Linx Headrest Mount w/ Quick Release	\$177.00
320IJ14	E1028	Whitmyer® Linx-2 Headrest Mount w/ Quick Release	\$215.00
320IJ15	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$113.00
320IJ16	E1028	Whitmyer® AXYS Offset Headrest Mount with Quick Release	\$135.00
320IJ17	E1028	Whitmyer® Cobra FB Headrest Mount with Quick Release	\$322.00

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Chin Control Packages P. 31

\$

Proportional Input Devices

320EU34	E2312	Switch-It MicroGuide <i>(Switches mounted to the headrest only.)</i>	\$3,537.00
320EU35	E2312	Switch-It MicroPilot <i>(Switches mounted to the headrest only.)</i>	\$3,537.00
320EU36	E2328	Switch-It VersaGuide	\$1,948.00
320EU37	E2328	Switch-It VersaGuide EZ	\$2,303.00

Mounting Position

320EM09	E1028	Link-It To Chin Harness, Small (150mm / 6")	\$1,342.00
320EM10	E1028	Link-It To Chin Harness, Large (220mm / 8")	\$1,342.00
320EM11	E1028	Link-It To Headrest (Swing-Away)	\$1,342.00
320EM53	E1028	Link-It To Manual Chin Boom - Right Side	\$4,074.00
320EM54	E1028	Link-It To Manual Chin Boom - Left Side	\$4,074.00
320EM16	E1028	Omit Mount	NC

Input Device Mounting Location

320IM13	—	Right	NC
320IM14	—	Left	NC

Switch Options

—	—	See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i>	—
---	---	---	---

Safety Devices

\$

320SA06	—	Switch-It Remote Stop Switch - XLR Connector	\$1,004.00
---------	---	--	------------

Switched Driving Packages

\$

320SW44	E2328	Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. <i>(Pro Spot proportional switches are required if proportional driving is desired.)</i> <i>(For non-proportional driving, single switches can be used.)</i>	\$3,725.00
---------	-------	---	------------

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Switches P. 32

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections

Switch / Quantity

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	\$
2.5" Single Switch - Black							\$209.00
2.5" Single Switch - Blue							\$209.00
2.5" Single Switch - Green							\$209.00
2.5" Single Switch - Red							\$209.00
2.5" Single Switch - Yellow							\$209.00
Egg Switch - Black							\$118.00
Egg Switch - Blue							\$118.00
Egg Switch - Green							\$118.00
Egg Switch - Red							\$118.00
Egg Switch - Yellow							\$118.00
MicroLite <i>(Includes color labels)</i>							\$177.00
Mini Button 1 - Blue <i>(Velcro)</i>							\$113.00
Mini Button 1 - Green <i>(Velcro)</i>							\$113.00
Mini Button 1 - Red <i>(Velcro)</i>							\$113.00
1" Single Switch - Black							\$296.00
1" Single Switch - Blue							\$296.00
1" Single Switch - Green							\$296.00
1" Single Switch - Orange							\$296.00
1" Single Switch - Red							\$296.00
1" Single Switch - Yellow							\$296.00
Electronic Switches							
Flex Switch (Headrest mount only)							\$230.00
Zero Touch Proximity Sensor							\$209.00
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red							\$510.00
ProSpot Proportional Switch - Green							\$510.00
ProSpot Proportional Switch - Blue							\$510.00
ProSpot Proportional Switch - Yellow							\$510.00
ProSpot Proportional Switch - Orange							\$510.00

Section B – Switch Mounting Position

Link-It to the Backrest <i>(Up to 1 per side)</i>							\$698.00
Link-It to the Headrest - Gooseneck & Rigid <i>(Up to 2 per side)</i>							\$698.00
Link-It to the Fixed Joystick Mount <i>(Up to 2 per side)</i>							\$532.00
Link-It to the Parallel Swing-Away Joystick Mount <i>(Up to 2 per side)</i>							\$532.00
Link-It to the Seat Rail <i>(Up to 2 per side)</i> <i>(Not avail w/knee & thigh supports)</i>							\$424.00
Omit Mount <i>(includes Velcro)</i>							NC

Section C – Switch Mounting Location

Right side of chair							\$263.00
Left side of chair							\$263.00

Joystick Integrated Switches (Pair)

320J003	–	Proximity Switches on Flex Mount, Pair <i>(Used with VersaGuide/VersaGuide EZ)</i>	\$1,068.00
---------	---	--	------------

Additional Switches (Velcro)

320SA08	K0108	Mini 2 - Button <i>(Stereo Plug)</i>	\$258.00
320SA09	K0108	Mini 2 - Button <i>(Dual Mono Plug)</i>	\$258.00
320SA10	K0108	Mini 4 - Button In-Line <i>(DB9 connector)</i>	\$585.00

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

POWER BASE

320DW58



14" Drive Wheel, Knobby

320CW55



9" Caster Wheels

320LI04



LED Lights

320BJ30



Charger

320PA41 – 320PA43



Power Seating: Tilt, Recline

320PI07



Power Seating: Elevate

SEATING AND POSITIONING

320JP01 & 320JP02



JAY® J3® Back: 2.2"

320JP11



JAY® J3® Back: 3"

320JP13



JAY® J3® Back: 6"

320K10 – 320K12



JAY® Anterior Thoracic Support

320JS02 – 320JS04



JAY® J3® Back, Spine-Align Kit 2

320ST79



SEDEO® Cushion


SEATING AND POSITIONING




JAY® Ion® Cushion



JAY® Union® Cushion



JAY® Fusion® Cushion w/Air




JAY® Fusion® Cushion w/Fluid



JAY® X2® Cushion



JAY® J2® Cushion



JAY® J2® Deep Contour Cushion



SEDEO® PRO™ Reclining Flip Back Armrest



SEDEO® PRO™ Single Post Removable Armrest



Comfort Reclining Armpad, Foam




SEDEO® Reclining Armrest Pad



Comfort Reclining Armpad, Gel

SEATING AND POSITIONING

320A837 – 320A840



Positioning Belt, Aircraft Style

320HE36



SEDEO® Headrest, Small

320HE37



SEDEO® Headrest, Flexible

320MH32



SEDEO® Multiposition Headrest Mount

320HE38 – 320HE44



Whitmyer® Plush Headrest

320HE45 – 320HE47



Whitmyer® Adjust-a-Plush Headrest

320MH33



Whitmyer® Linx Headrest Mount

320MH34



Whitmyer® Linx-2 Headrest Mount

ELECTRONICS

320J116



Color Joystick (CJSM1)

320J117



Heavy Duty Color Joystick

320J118



High Resolution Color Joystick (CJSM2+BT)

320JH16




Goalpost Joystick Handle

320JH15



T-Handle Joystick Handle

320JH14



Foam Ball Joystick Handle

320JH13



Ball Joystick Handle

ELECTRONICS



R-net PC Programmer



R-net Power Module



Attendant Control



Ctrl+5®, Push Button



Ctrl+5®, Toggle



R-net Omni2



Switch-It™ Head Array



Switch-It™ Dual Pro® Proportional Head Array



Switch-It™ RIM Proportional Head Control

ELECTRONICS

320EU14 & 320EU34



Switch-It™ MicroGuide

320EU15 & 320EU35



Switch-It™ MicroPilot

320EU17 & 320EU36



Switch-It™ VersaGuide

320EU16



Switch-It™ TouchDrive 2 (TD2)

320EU20



All Round Joystick, Light

320EU19



All Round Joystick, Standard

320EM04




Link-It to Backrest

320EM05




Link-It to Fixed Joystick Mount

320EM06




Link-It to Parallel Swing-Away Joystick Mount

320EM07



Link-It to Midline Swing-Away Joystick Mount

320EM09



Link-It to Chin Harness

320EM11



Link-It to Headrest, Swing-Away

320EM53




Link-It to Manual Chin Boom

320EM55




Link-It to Power Chin Boom

ELECTRONICS



2.5" Single Switch
black, blue, green, red, yellow




Egg Switch
black, blue, green, red, yellow




MicroLite



1" Single Switch
black, blue, green, orange, red, yellow



Twister Switch
black, green, red



Mini 2



Mini 4



Mini 1 Button (Mono Plug)
blue, green, red



Q500 H with SEDEO® PRO™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

