



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

Little Wave XP

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025

Little Wave XP (K0005, E1235) \$3465

Little Wave Width Growable Frame 165 lb weight capacity Maximum Seat Width of 19"

## Frame

### Seat Width

|     |     |
|-----|-----|
| 10" | N/C |
| 11" | N/C |
| 12" | N/C |
| 13" | N/C |
| 14" | N/C |
| 15" | N/C |
| 16" | N/C |

### Seat Depth

|     |     |
|-----|-----|
| 10" | N/C |
| 11" | N/C |
| 12" | N/C |
| 13" | N/C |
| 14" | N/C |
| 15" | N/C |
| 16" | N/C |

### Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.

|                  |     |
|------------------|-----|
| 0" (Performance) | N/C |
| 1"               | N/C |
| 2" (Classic)     | N/C |

### Transit Option

|                                |       |
|--------------------------------|-------|
| Transit System                 | \$246 |
| 125 lb maximum weight capacity |       |
| Omit                           | N/C   |

### Adjustable Depth

Adds the capability to grow seat in depth

|                    |     |
|--------------------|-----|
| No Rear Adjustment | N/C |
| 1/2"               | N/C |
| 1"                 | N/C |
| 2"                 | N/C |
| 3"                 | N/C |

### Long Sling

Available only with 1", 2" and 3" Adjustable Depth

|                                                 |     |
|-------------------------------------------------|-----|
| Long Sling                                      | N/C |
| Sling will be maximum 2" longer than seat depth |     |

### Seat Option

|                                                          |       |
|----------------------------------------------------------|-------|
| Standard Seat Sling                                      | N/C   |
| Aluminum Solid Seat Pan (E2231)                          | \$272 |
| Seat Pan is 1" longer when Long Sling option is selected |       |

### Front Frame Bend

Angle between front seat tube and ground. See reference D on page 7.

|     |     |
|-----|-----|
| 90° | N/C |
| 85° | N/C |
| 80° | N/C |
| 75° | N/C |
| 70° | N/C |

## Backrest

| Back Post / Height Range / Back Height |       |                |                  |                  |                  |                  |
|----------------------------------------|-------|----------------|------------------|------------------|------------------|------------------|
| Back                                   | Price | Height         | Back Height      |                  |                  |                  |
| Height Adjustable                      | N/C   | Short          | 11" <sup>2</sup> | 12"              | 13"              | 14"              |
|                                        |       | Medium         | 13"              | 14"              | 15"              | 16"              |
|                                        |       | Tall           | 17"              | 18"              |                  |                  |
| Fixed Height <sup>1</sup>              | \$67  | Fixed Straight | 9" <sup>2</sup>  | 10" <sup>2</sup> | 11" <sup>2</sup> | 12" <sup>2</sup> |
|                                        |       |                | 14"              | 15"              | 16"              | 17"              |
| Offset <sup>1</sup>                    | N/C   | Offset         | 11" <sup>2</sup> | 12" <sup>2</sup> | 13" <sup>2</sup> | 14" <sup>2</sup> |
|                                        |       |                | 16"              | 17"              |                  |                  |

<sup>1</sup> Fixed Height. Adjustable Height Rigidizer Bar Included.

<sup>2</sup> Not available with Fold Down Push Handle

| Rigidizer Bar                       |      |
|-------------------------------------|------|
| Adjustable Height Rigidizer Bar     | \$15 |
| Included with Offset Back Posts     |      |
| Non Adjustable Height Rigidizer Bar | N/C  |

| Seat Back Angle |      |
|-----------------|------|
| Standard        | +3°  |
| -9°             | -12° |
|                 | +6°  |
|                 | +9°  |
|                 | -3°  |
|                 | -6°  |
|                 | N/C  |

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

### Offset Back Post Option

Selection required with Offset Back Posts. Not available with Height Adjustable or Fixed Height Back Posts. Each back post is offset 1" creating 2" overall offset.

|                                        |     |
|----------------------------------------|-----|
| 2" Inward Offset                       | N/C |
| 2" Outward Offset                      | N/C |
| Not available with Swing Away Armrests |     |

### Push Handle

|                                                                |       |
|----------------------------------------------------------------|-------|
| Omit Push Handle                                               | N/C   |
| Standard Push Handle                                           | N/C   |
| Not available on Offset Back Posts with 11" - 14" Back Heights |       |
| Bolt-On Push Handle                                            | \$215 |
| Handles will be 1" - 2" above back height                      |       |
| Fold Down Push Handle                                          | \$123 |

### Stroller Handles

Stroller Handles are removable and height adjustable

|                                                   |       |
|---------------------------------------------------|-------|
| Ergo Stroller Handle                              | \$190 |
| Single Stroller Handle                            | \$169 |
| Double Stroller Handle                            | \$338 |
| Not available with Push Handles or 10" Seat Width |       |

### Back Upholstery

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

|                                                           |        |
|-----------------------------------------------------------|--------|
| Padded Velcro Adjustable with Standard Cover              | N/C    |
| Padded Velcro Adjustable with Air Mesh Breathable Cover   | \$79   |
| Not available on chair widths less than 12"               |        |
| Tension Adjustable with Standard Cover (E2611)            | \$359  |
| Tension Adjustable with Air Mesh Breathable Cover (E2611) | \$408  |
| Not available on chair widths less than 12"               |        |
| Omit                                                      | (\$10) |

## Front Seat Height and Casters

| Caster Size / Footrest Taper / Front Seat Height |  |
|--------------------------------------------------|--|
|--------------------------------------------------|--|

Select Front Seat Height and Footrest Taper in row corresponding to desired Caster Size. 6" Casters not available with Seat Depths less than 13".

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. Taper is not available with 70° Front Frame Bend. See reference E on page 7.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

| Footrest Taper      | Caster Size | Front Seat Height |       |       |       |       |       |       |       |       |       |       |       |       |       |       |
|---------------------|-------------|-------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| No Taper (Straight) | 3"          | 11.5"             | 12"   | 12.5" | 13"   | 13.5" | 14"   | 14.5" | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" |
|                     | 4"          | 12.5"             | 13"   | 13.5" | 14"   | 14.5" | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" |
|                     | 5"          | 13.5"             | 14"   | 14.5" | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |
|                     | 6"          | 14"               | 14.5" | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |       |
| 1"                  | 3"          | 14.5"             | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |       |       |
|                     | 4"          | 15.5"             | 16"   | 16.5" | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |       |       |       |       |
|                     | 5"          | 16.5"             | 17"   | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |       |       |       |       |       |       |
|                     | 6"          | 17"               | 17.5" | 18"   | 18.5" | 19"   | 19.5" | 20"   | 20.5" |       |       |       |       |       |       |       |



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# Little Wave XP

## Casters

|             |                          | Caster Size |    |    |    |    |
|-------------|--------------------------|-------------|----|----|----|----|
|             |                          | Price       | 3" | 4" | 5" | 6" |
| Caster Type | .75" Rollerblade         | N/C         |    |    |    |    |
|             | .75" Lighted Rollerblade | \$77        |    |    |    |    |
|             | 1.25" Soft Roll          | N/C         |    |    |    |    |
|             | 1" Poly Aluminum         | \$128       |    |    |    |    |
|             | 1" Pneumatic             | \$87        |    |    |    |    |
|             | 1.5" Soft Roll           | \$67        |    |    |    |    |
|             | 1.5" Soft Roll Aluminum  | \$210       |    |    |    |    |

## Fork Options

|                                      |       |
|--------------------------------------|-------|
| Standard Fork                        | N/C   |
| Single Sided Fork                    | \$178 |
| <i>Not available with 6" Casters</i> |       |
| Frog Legs II Suspension System       | \$498 |

## Footrest

### Footplates

|                                                                                      |       |
|--------------------------------------------------------------------------------------|-------|
| Tubular "Open" Footplate                                                             | N/C   |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i> |       |
| Tubular Footplate with ABS Cover                                                     | \$36  |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i> |       |
| Aluminum Angle Adjustable (K0040)                                                    | \$97  |
| Aluminum Angle Adjustable with Riser (K0040)                                         | \$190 |
| High Mount Aluminum Angle Adjustable (K0040, K0037)                                  | \$164 |
| Angle Adjustable Flip Under (K0040, K0037)                                           | \$210 |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i> |       |
| High Mount Angle Adjustable Flip Under (K0040, K0037)                                | \$318 |

### Seat to Footrest Length

*See reference H on p. 7 should be at least 2.5" shorter than FSH*

|       |     |       |     |       |     |
|-------|-----|-------|-----|-------|-----|
| 5.5"  | 6"  | 6.5"  | 7"  | 7.5"  | 8"  |
| 8.5"  | 9"  | 9.5"  | 10" | 10.5" | 11" |
| 11.5" | 12" | 12.5" | 13" | 13.5" | 14" |
| 14.5" | 15" | 15.5" | 16" | 16.5" | 17" |

## Rear Seat Height and Wheels

### Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.

| Wheel Size | Rear Seat Height   |     |       |     |       |     |       |     |       |     |       |     |       |     |       |     |
|------------|--------------------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|
| 18"        | 11.5" <sup>1</sup> | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" |
| 20"        | 12.5" <sup>1</sup> | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" |
| 22"        | 13.5" <sup>1</sup> | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" |       |     |
| 24"        | 14.5" <sup>1</sup> | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" |       |     |       |     |

<sup>1</sup> Available only with Low Profile tires

### Wheels

|            |               | Rear Wheel Size |     |     |     |     |
|------------|---------------|-----------------|-----|-----|-----|-----|
|            |               | Price           | 18" | 20" | 22" | 24" |
| Rear Wheel | MAXX Spoke    | N/C             |     |     |     |     |
|            | Mag           | \$56            |     |     |     |     |
|            | Spinergy Spox | \$975           |     |     |     |     |

### Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes   Red Spokes   Blue Spokes   White Spokes \_\_\_\_\_ N/C

### Handrims

#### Standard Handrims

|                                                                                                                                  |       |
|----------------------------------------------------------------------------------------------------------------------------------|-------|
| Aluminum Anodized                                                                                                                | N/C   |
| Plastic Coated                                                                                                                   | \$127 |
| Projection (8 Vertical Projections)                                                                                              | \$367 |
| <i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i> |       |
| Omit Handrim                                                                                                                     | N/C   |

#### Ergonomic Handrims

*Not available with 20" or 26" wheels*

|                       |       |
|-----------------------|-------|
| Lite Regular Grip     | \$477 |
| Standard Regular Grip | \$456 |
| Lite Extra Grip       | \$528 |
| Standard Extra Grip   | \$507 |
| Lite Omit Grip        | \$405 |
| Standard Omit Grip    | \$354 |



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

Little Wave XP

## Rear Tire Type

|                                                             |       |
|-------------------------------------------------------------|-------|
| Full Poly 1 3/8" (E2218)                                    | \$77  |
| Low Poly 1" (E2218)                                         | \$79  |
| Pneumatic with Airless Insert 1 3/8" (E2211, E2213)         | \$99  |
| Street Pneumatic 1 3/8" (E2211, E2212)                      | N/C   |
| High Pressure 1" (E2211, E2212)                             | \$100 |
| <i>Not available on Mag wheel or 18" wheel size</i>         |       |
| Schwalbe Marathon Plus Evolution 1" (E2211, E2212)          | \$164 |
| <i>Available only on 24" MAXX Spoke or Spinerger wheels</i> |       |

## Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest, 1.25" minimum with Height Adjustable Pediatric Desk Armrest, 1" minimum with 0° Rear Wheel Camber. See reference K on page 7. No selection needed with Wheelchair Drive Ready option.

0.75" 1" 1.25" 1.5" 1.75" N/C

## Rear Wheel Camber

0° 2° 4° 6° 8° N/C

## Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 3). Requires Rear Anti-Tips. Not available with 5th Wheel or Spinerger wheels.

|                 |     |
|-----------------|-----|
| e-motion Ready  | N/C |
| twion Ready     | N/C |
| 0 Degree Camber | N/C |
| 3 Degree Camber | N/C |

## Wheel Locks

### Wheel Lock

|                                                                                                                                                                                                      |       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| Push to Lock                                                                                                                                                                                         | N/C   |
| Pull to Lock                                                                                                                                                                                         | N/C   |
| Low Profile Scissor Lock                                                                                                                                                                             | \$67  |
| Drum Brake (E2228)                                                                                                                                                                                   | \$574 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> |       |
| Dual - Drum Brake with Push to Lock (E2228)                                                                                                                                                          | \$677 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> |       |
| Dual - Drum Brake with Pull to Lock (E2228)                                                                                                                                                          | \$677 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> |       |
| Omit                                                                                                                                                                                                 | N/C   |

## Handrim Options

### Handrim Mount

If no selection is made, Standard Mount will be used. Reduces overall width by .75.

Close Mount N/C

3M Micro Grip \$67

*Not available on Plastic Coated or Projection Handrims*

## Center of Gravity presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 7.

| Amputee Positions  |        |     |       |
|--------------------|--------|-----|-------|
| -1.5"              | -1.25" | -1" | -.75" |
| -.5"               | -.25"  | 0"  | .25"  |
| Standard Positions |        |     |       |
| .5"                | .75"   | 1"  | 1.25" |
| 1.5"               | 1.75"  | 2"  | 2.25" |
| 2.5"               | 2.75"  | 3"  | 3.25" |
| 3.5"               | 3.75"  | 4"  | 4.25" |

## Axles

No selection needed with Drum Brake

Quick Release N/C

Quad Quick Release \$194

*Not available with Mag or Spinerger wheels*

## Wheel Lock Options

Extension Handles (E0961) \$84

*Available only with Push to Lock or Pull to Lock Wheel Locks*

Grade Aids (E0974) \$154

*Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert tires*



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

Little Wave XP

## Armrests

| Armrest Type                                                    | Length | Price |
|-----------------------------------------------------------------|--------|-------|
| Height Adjustable T-Arm Tall (11.5"-15") E0973 <sup>1,3</sup>   | Desk   | \$318 |
|                                                                 | Full   |       |
| Height Adjustable T-Arm (8.5"-13.5") E0973 <sup>1,3</sup>       | Desk   | \$325 |
|                                                                 | Full   |       |
| Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>2,3,4</sup> | Desk   | \$318 |
|                                                                 | Full   |       |

1 Not available with Seat Depths less than 13"

2 Not available with Seat Depths less than 12"

3 1.5" minimum Rear Wheel Spacing

4 If requested configuration is unavailable, Standard Height T-Arm will be used

5 Requires Fixed Height Backrest and Adjustable Height Rigidizer Bar

6 Not available with back heights less than 12"

7 Best fit with Seat Depths less than 14"

8 1.25" minimum Rear Wheel Spacing

9 Not available with Arm Pads

10 Not available with Back Heights less than 13"

11 Not available with Offset

12 Not available with Transit System AND more than 2" Adjustable Depth

| Armrest Type                                                       | Price |
|--------------------------------------------------------------------|-------|
| Height Adjustable Pediatric Desk (6.5" - 10") E0973 <sup>7,8</sup> | \$231 |
| Tubular Flip Up <sup>5,9,10,11</sup>                               | \$51  |
| Swing Away (8.5"-12.5") E0973 <sup>9,11,12</sup>                   | \$244 |
| Omit Armrests                                                      | N/C   |

### Arm Pads

Not available with Tubular Flip Up or Swing Away Armrests. Full Arm Pads not available with Height Adjustable Pediatric Desk Armrests.

|                                                                         |      |
|-------------------------------------------------------------------------|------|
| Standard Arm Pad                                                        | N/C  |
| Waterfall Arm Pad                                                       | \$26 |
| Not available with Angle Adjustable Locking Extendable Flip Up Armrests |      |
| Foam Grip                                                               | N/C  |
| Available only with Angle Adjustable Locking Extendable Flip Up Armrest |      |

## Side Guards

Not available with T-Arms or Height Adjustable Pediatric Desk Armrests

### Direct Mount Side Guards

|                                                                         |       |
|-------------------------------------------------------------------------|-------|
| Composite Side Guards - Adult - Direct Mount Receiver (K0108) _____     | \$147 |
| 12" minimum Seat Depth                                                  |       |
| Composite Side Guards - Pediatric - Direct Mount Receiver (K0108) _____ | \$159 |

### Clamp-On Mount Side Guards

|                                                                                             |       |
|---------------------------------------------------------------------------------------------|-------|
| Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black Matte (K0108) _____       | \$313 |
| Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Match Frame Color (K0108) _____ | \$364 |
| Not available with Hydrographics                                                            |       |

## Accessories

### 5th Wheel

Rotating removable rear caster. See reference L on page 7.

|                                                                        |       |
|------------------------------------------------------------------------|-------|
| Standard 5th Wheel (E0971) _____                                       | \$256 |
| Dynamic 5th Wheel (K0108) _____                                        | \$302 |
| Soft Elastomer _____                                                   | N/C   |
| Both elastomers are included; the selected elastomer will be installed |       |
| Firm Elastomer _____                                                   | N/C   |
| Both elastomers are included; the selected elastomer will be installed |       |

### Anti-Tip

|                              |       |
|------------------------------|-------|
| Rear Anti-Tips (E0971) _____ | \$172 |
| Omit _____                   | N/C   |

### Calf Strap

|                                         |       |
|-----------------------------------------|-------|
| Velcro Adjustable (K0038) _____         | \$43  |
| Padded Velcro Adjustable (K0038) _____  | \$70  |
| Padded Airmesh Calf Strap (K0038) _____ | \$121 |

### Pelvic Positioning Belt

|                                                                 |       |
|-----------------------------------------------------------------|-------|
| 1" Auto Buckle Pelvic Positioning Belt (E0978) _____            | \$88  |
| 1" Padded Auto Buckle Pelvic Positioning Belt (E0978) _____     | \$179 |
| 1 1/2" Auto Buckle Pelvic Positioning Belt (E0978) _____        | \$89  |
| 1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978) _____ | \$179 |
| Small Flex-Mount Positioning Belt (16-24") (E0978) _____        | \$204 |
| Medium Flex-Mount Positioning Belt (19-29") (E0978) _____       | \$204 |
| Large Flex-Mount Positioning Belt (25-37") (E0978) _____        | \$204 |

### Belt Mounting Kit

|                                       |       |
|---------------------------------------|-------|
| Ki Mobility Belt Mounting Clamp _____ | \$37  |
| Belt Mounting Bracket _____           | \$115 |

### Other Accessories

|                                                        |       |
|--------------------------------------------------------|-------|
| Tool Kit _____                                         | \$62  |
| Removable Underseat Pouch _____                        | \$63  |
| Pediatric Ki Mobility Backpack _____                   | \$73  |
| Ki Mobility Backpack _____                             | \$81  |
| Best fit on chairs more than 10" wide                  |       |
| Luggage Carrier _____                                  | \$88  |
| Black Neoprene Impact Guards _____                     | \$97  |
| Clear Spoke Guards (K0065) _____                       | \$123 |
| Not available with Mag wheel                           |       |
| Black Spoke Guards (K0065) _____                       | \$125 |
| Not available with 18" or 20" wheel sizes or Mag wheel |       |
| Cane and Crutch Holder _____                           | \$210 |



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



Frame Color

Gloss Finish

Touch-up paint is available for purchase through the parts manual

|                    |     |
|--------------------|-----|
| Black              | N/C |
| Ruby Red           | N/C |
| Radiant Orange     | N/C |
| Yellow             | N/C |
| Electric Green     | N/C |
| Rich Teal          | N/C |
| Glacier Blue       | N/C |
| Blue Steel         | N/C |
| Deep Sapphire Blue | N/C |
| Purple Plum        | N/C |
| Opalescent Violet  | N/C |
| Dark Pink          | N/C |
| Light Pink         | N/C |
| Champagne          | N/C |

Matte Finish

Touch-up paint is available for purchase through the parts manual

|                 |     |
|-----------------|-----|
| Shadow Black    | N/C |
| Iridescent Navy | N/C |
| Bright Red      | N/C |
| Tactical Green  | N/C |
| Dry Earth Tan   | N/C |
| Platinum        | N/C |
| White           | N/C |

## References

### **A. Seat Width** See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

### **B. Seat Depth** See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

### **C. Sling Position** See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

### **D. Front Frame Bend** See Diagram 3

Angle between front seat tube and ground.

### **E. Footrest Taper** See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

### **F. Front Seat Height** See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

### **G. Footrest Width** See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

| Seat Width          | Straight                    | 1" Taper                    |
|---------------------|-----------------------------|-----------------------------|
| Outside Measurement | 2" Narrower than Seat Width | 4" Narrower than Seat Width |
| 10"                 | 7.5"                        |                             |
| 11"                 | 8.5"                        |                             |
| 12"                 | 9.5"                        |                             |
| 13"                 | 10.5"                       | 7.5"                        |
| 14"                 | 11.5"                       | 8.5"                        |
| 15"                 | 12.5"                       | 10.5"                       |
| 16"                 | 13.5"                       | 11.5"                       |

### **H. Seat to Footrest Length** See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

### **I. Rear Seat Height** See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

### **J. Center of Gravity Preset** See Diagram 3

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

### **K. Rear Wheel Spacing** See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

### **L. 5th Wheel** See Diagram 4

Rotating removable 5th wheel allows for aggressive center of gravity positioning of the rear wheels.

Diagram 1




Diagram 2




Diagram 3




Diagram 4





Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Axiom for Kids P - Positioning Cushion

**Weight Capacity 165 lbs / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.**

| Width  |     | 8"         | 9"         | 10"        | 11"        | 12"        | 13"        | 14"        | 15"        | 16"        |
|--------|-----|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| Retail |     | \$397      |            |            |            |            |            |            |            |            |
| HCPCS  |     | E2605      | E2605      | E2605      | E2605      | E2605      | E2605      | E2605      | E2605      | E2605      |
| Depth  | 8"  | PXP0808BLK | PXP0908BLK | PXP1008BLK |            |            |            |            |            |            |
|        |     | PXP0808BLU | PXP0908BLU | PXP1008BLU |            |            |            |            |            |            |
|        |     | PXP0808PUR | PXP0908PUR | PXP1008PUR |            |            |            |            |            |            |
|        | 9"  | PXP0809BLK | PXP0909BLK | PXP1009BLK | PXP1109BLK |            |            |            |            |            |
|        |     | PXP0809BLU | PXP0909BLU | PXP1009BLU | PXP1109BLU |            |            |            |            |            |
|        |     | PXP0809PUR | PXP0909PUR | PXP1009PUR | PXP1109PUR |            |            |            |            |            |
|        | 10" | PXP0810BLK | PXP0910BLK | PXP1010BLK | PXP1110BLK | PXP1210BLK | PXP1310BLK |            |            |            |
|        |     | PXP0810BLU | PXP0910BLU | PXP1010BLU | PXP1110BLU | PXP1210BLU | PXP1310BLU |            |            |            |
|        |     | PXP0810PUR | PXP0910PUR | PXP1010PUR | PXP1110PUR | PXP1210PUR | PXP1310PUR |            |            |            |
|        | 11" |            |            | PXP1011BLK | PXP1111BLK | PXP1211BLK | PXP1311BLK | PXP1411BLK |            |            |
|        |     |            |            | PXP1011BLU | PXP1111BLU | PXP1211BLU | PXP1311BLU | PXP1411BLU |            |            |
|        |     |            |            | PXP1011PUR | PXP1111PUR | PXP1211PUR | PXP1311PUR | PXP1411PUR |            |            |
|        | 12" |            |            | PXP1012BLK | PXP1112BLK | PXP1212BLK | PXP1312BLK | PXP1412BLK |            |            |
|        |     |            |            | PXP1012BLU | PXP1112BLU | PXP1212BLU | PXP1312BLU | PXP1412BLU |            |            |
|        |     |            |            | PXP1012PUR | PXP1112PUR | PXP1212PUR | PXP1312PUR | PXP1412PUR |            |            |
|        | 13" |            |            | PXP1013BLK | PXP1113BLK | PXP1213BLK | PXP1313BLK | PXP1413BLK | PXP1513BLK |            |
|        |     |            |            | PXP1013BLU | PXP1113BLU | PXP1213BLU | PXP1313BLU | PXP1413BLU | PXP1513BLU |            |
|        |     |            |            | PXP1013PUR | PXP1113PUR | PXP1213PUR | PXP1313PUR | PXP1413PUR | PXP1513PUR |            |
|        | 14" |            |            | PXP1014BLK | PXP1114BLK | PXP1214BLK | PXP1314BLK | PXP1414BLK | PXP1514BLK | PXP1614BLK |
|        |     |            |            | PXP1014BLU | PXP1114BLU | PXP1214BLU | PXP1314BLU | PXP1414BLU | PXP1514BLU | PXP1614BLU |
|        |     |            |            | PXP1014PUR | PXP1114PUR | PXP1214PUR | PXP1314PUR | PXP1414PUR | PXP1514PUR | PXP1614PUR |
|        | 15" |            |            |            |            | PXP1215BLK | PXP1315BLK | PXP1415BLK | PXP1515BLK | PXP1615BLK |
|        |     |            |            |            |            | PXP1215BLU | PXP1315BLU | PXP1415BLU | PXP1515BLU | PXP1615BLU |
|        |     |            |            |            |            | PXP1215PUR | PXP1315PUR | PXP1415PUR | PXP1515PUR | PXP1615PUR |
|        | 16" |            |            |            |            | PXP1216BLK | PXP1316BLK | PXP1416BLK | PXP1516BLK | PXP1616BLK |
|        |     |            |            |            |            | PXP1216BLU | PXP1316BLU | PXP1416BLU | PXP1516BLU | PXP1616BLU |
|        |     |            |            |            |            | PXP1216PUR | PXP1316PUR | PXP1416PUR | PXP1516PUR | PXP1616PUR |
|        | 17" |            |            |            |            |            | PXP1317BLK | PXP1417BLK | PXP1517BLK | PXP1617BLK |
|        |     |            |            |            |            |            | PXP1317BLU | PXP1417BLU | PXP1517BLU | PXP1617BLU |
|        |     |            |            |            |            |            | PXP1317PUR | PXP1417PUR | PXP1517PUR | PXP1617PUR |
|        | 18" |            |            |            |            |            | PXP1318BLK | PXP1418BLK | PXP1518BLK | PXP1618BLK |
|        |     |            |            |            |            |            | PXP1318BLU | PXP1418BLU | PXP1518BLU | PXP1618BLU |
|        |     |            |            |            |            |            | PXP1318PUR | PXP1418PUR | PXP1518PUR | PXP1618PUR |

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1075.1 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Axiom for Kids SP - Skin Protection and Positioning Cushion

**Weight Capacity 165 lbs. / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.**

| Width  |     | 8"          | 9"          | 10"         | 11"         | 12"         | 13"         | 14"         | 15"         | 16"         |
|--------|-----|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Retail |     | \$481       |             |             |             |             |             |             |             |             |
| HCPCS  |     | E2607       | E2607       | E2607       | E2607       | E2607       | E2607       | E2607       | E2607       | E2607       |
| Depth  | 8"  | PXSP0808BLK | PXSP0908BLK | PXSP1008BLK |             |             |             |             |             |             |
|        |     | PXSP0808BLU | PXSP0908BLU | PXSP1008BLU |             |             |             |             |             |             |
|        |     | PXSP0808PUR | PXSP0908PUR | PXSP1008PUR |             |             |             |             |             |             |
|        | 9"  | PXSP0809BLK | PXSP0909BLK | PXSP1009BLK | PXSP1109BLK |             |             |             |             |             |
|        |     | PXSP0809BLU | PXSP0909BLU | PXSP1009BLU | PXSP1109BLU |             |             |             |             |             |
|        |     | PXSP0809PUR | PXSP0909PUR | PXSP1009PUR | PXSP1109PUR |             |             |             |             |             |
|        | 10" | PXSP0810BLK | PXSP0910BLK | PXSP1010BLK | PXSP1110BLK | PXSP1210BLK | PXSP1310BLK |             |             |             |
|        |     | PXSP0810BLU | PXSP0910BLU | PXSP1010BLU | PXSP1110BLU | PXSP1210BLU | PXSP1310BLU |             |             |             |
|        |     | PXSP0810PUR | PXSP0910PUR | PXSP1010PUR | PXSP1110PUR | PXSP1210PUR | PXSP1310PUR |             |             |             |
|        | 11" |             |             | PXSP1011BLK | PXSP1111BLK | PXSP1211BLK | PXSP1311BLK | PXSP1411BLK |             |             |
|        |     |             |             | PXSP1011BLU | PXSP1111BLU | PXSP1211BLU | PXSP1311BLU | PXSP1411BLU |             |             |
|        |     |             |             | PXSP1011PUR | PXSP1111PUR | PXSP1211PUR | PXSP1311PUR | PXSP1411PUR |             |             |
|        | 12" |             |             | PXSP1012BLK | PXSP1112BLK | PXSP1212BLK | PXSP1312BLK | PXSP1412BLK |             |             |
|        |     |             |             | PXSP1012BLU | PXSP1112BLU | PXSP1212BLU | PXSP1312BLU | PXSP1412BLU |             |             |
|        |     |             |             | PXSP1012PUR | PXSP1112PUR | PXSP1212PUR | PXSP1312PUR | PXSP1412PUR |             |             |
|        | 13" |             |             | PXSP1013BLK | PXSP1113BLK | PXSP1213BLK | PXSP1313BLK | PXSP1413BLK | PXSP1513BLK |             |
|        |     |             |             | PXSP1013BLU | PXSP1113BLU | PXSP1213BLU | PXSP1313BLU | PXSP1413BLU | PXSP1513BLU |             |
|        |     |             |             | PXSP1013PUR | PXSP1113PUR | PXSP1213PUR | PXSP1313PUR | PXSP1413PUR | PXSP1513PUR |             |
|        | 14" |             |             | PXSP1014BLK | PXSP1114BLK | PXSP1214BLK | PXSP1314BLK | PXSP1414BLK | PXSP1514BLK | PXSP1614BLK |
|        |     |             |             | PXSP1014BLU | PXSP1114BLU | PXSP1214BLU | PXSP1314BLU | PXSP1414BLU | PXSP1514BLU | PXSP1614BLU |
|        |     |             |             | PXSP1014PUR | PXSP1114PUR | PXSP1214PUR | PXSP1314PUR | PXSP1414PUR | PXSP1514PUR | PXSP1614PUR |
|        | 15" |             |             |             |             | PXSP1215BLK | PXSP1315BLK | PXSP1415BLK | PXSP1515BLK | PXSP1615BLK |
|        |     |             |             |             |             | PXSP1215BLU | PXSP1315BLU | PXSP1415BLU | PXSP1515BLU | PXSP1615BLU |
|        |     |             |             |             |             | PXSP1215PUR | PXSP1315PUR | PXSP1415PUR | PXSP1515PUR | PXSP1615PUR |
|        | 16" |             |             |             |             | PXSP1216BLK | PXSP1316BLK | PXSP1416BLK | PXSP1516BLK | PXSP1616BLK |
|        |     |             |             |             |             | PXSP1216BLU | PXSP1316BLU | PXSP1416BLU | PXSP1516BLU | PXSP1616BLU |
|        |     |             |             |             |             | PXSP1216PUR | PXSP1316PUR | PXSP1416PUR | PXSP1516PUR | PXSP1616PUR |
|        | 17" |             |             |             |             |             | PXSP1317BLK | PXSP1417BLK | PXSP1517BLK | PXSP1617BLK |
|        |     |             |             |             |             |             | PXSP1317BLU | PXSP1417BLU | PXSP1517BLU | PXSP1617BLU |
|        |     |             |             |             |             |             | PXSP1317PUR | PXSP1417PUR | PXSP1517PUR | PXSP1617PUR |
|        | 18" |             |             |             |             |             | PXSP1318BLK | PXSP1418BLK | PXSP1518BLK | PXSP1618BLK |
|        |     |             |             |             |             |             | PXSP1318BLU | PXSP1418BLU | PXSP1518BLU | PXSP1618BLU |
|        |     |             |             |             |             |             | PXSP1318PUR | PXSP1418PUR | PXSP1518PUR | PXSP1618PUR |

### Accessories


#### Extra Cover

|                    |       |
|--------------------|-------|
| Extra Cover Black  | \$118 |
| Extra Cover Blue   | \$118 |
| Extra Cover Purple | \$118 |

#### Solid Seat Insert

|                           |      |
|---------------------------|------|
| Solid Seat Insert (E0992) | \$72 |
|---------------------------|------|


| Axiom for Kids Back Width (all models) | Fits back Frame Widths    |
|----------------------------------------|---------------------------|
| 10" (25.5 cm)                          | 9" - 11" (23 cm - 28 cm)  |
| 12" (30.5 cm)                          | 11" - 13" (28 cm - 33 cm) |
| 14" (35.5 cm)                          | 13" - 15" (33 cm - 38 cm) |
| 16" (40.5 cm)                          | 15" - 17" (38 cm - 43 cm) |



## Axiom for Kids PL - Posterior Lateral Back

Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back options and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.

|                                |                        |                         |                         |                         |
|--------------------------------|------------------------|-------------------------|-------------------------|-------------------------|
| <b>Back Width</b>              | 10" (25.5 cm)          | 12" (30.5 cm)           | 14" (35.5 cm)           | 16" (40.5 cm)           |
| <b>Fits Back Frame Widths:</b> | 9"-11" (23 cm - 28 cm) | 11"-13" (28 cm - 33 cm) | 13"-15" (33 cm - 38 cm) | 15"-17" (38 cm - 43 cm) |
| <b>Actual Back Width (F):</b>  | 10" (25.5 cm)          | 12" (30.5 cm)           | 14" (35.5 cm)           | 16" (40.5 cm)           |
| <b>Fits Trunk Widths (F):</b>  | 10" (25.5 cm)          | 12" (30.5 cm)           | 14" (35.5 cm)           | 16" (40.5 cm)           |
| <b>Contour Depth (H):</b>      | 3" (7.5 cm)            | 3" (7.5 cm)             | 3" (7.5 cm)             | 3" (7.5 cm)             |



### Back

Includes Lumbar Pad. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

|                   |        | Wheelchair Frame Width |          |          |          |
|-------------------|--------|------------------------|----------|----------|----------|
|                   |        | 10                     | 12       | 14       | 16       |
| Axiom Back Height | Retail | \$740                  |          |          |          |
|                   | HCPCS  | E2615                  | E2615    | E2615    | E2615    |
|                   | 10     | PXPL1010               |          |          |          |
|                   | 12     | PXPL1012               | PXPL1212 | PXPL1412 |          |
|                   | 14     | PXPL1014               | PXPL1214 | PXPL1414 | PXPL1614 |
|                   | 16     |                        | PXPL1216 | PXPL1416 | PXPL1616 |
|                   | 18     |                        |          |          | PXPL1618 |

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.


Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1076.3 Pricing and specifications subject to change without notice.

## Axiom for Kids AL - Posterior Back for Adjustable Laterals

*Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back options and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.*

|                                                                           |                          |                             |                             |                             |
|---------------------------------------------------------------------------|--------------------------|-----------------------------|-----------------------------|-----------------------------|
| <b>Back Width</b>                                                         | 10" (25.5 cm)            | 12" (30.5 cm)               | 14" (35.5 cm)               | 16" (40.5 cm)               |
| <b>Fits Back Frame Widths:</b>                                            | 9"-11" (23 cm - 28 cm)   | 11"-13" (28 cm - 33 cm)     | 13"-15" (33 cm - 38 cm)     | 15"-17" (38 cm - 43 cm)     |
| <b>Actual Back Width (F):</b>                                             | 8" (20cm)                | 10" (25.5 cm)               | 12" (30.5 cm)               | 14" (35.5 cm)               |
| <b>Component Lateral Width Range (J)*:</b><br>(*1" with included spacers) | 8" - 11" (20 cm - 28 cm) | 10" - 13" (25.5 cm - 33 cm) | 12" - 15" (30.5 cm - 38 cm) | 14" - 17" (35.5 cm - 43 cm) |
| <b>Contour Depth (H):</b>                                                 | 1.25" (3 cm)             | 1.25" (3 cm)                | 1.25" (3 cm)                | 1.25" (3 cm)                |



### Back

*Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.*

|                   |        | Wheelchair Frame Width |          |          |          |
|-------------------|--------|------------------------|----------|----------|----------|
|                   |        | 10"                    | 12"      | 14"      | 16"      |
| Axiom Back Height | Retail | \$870                  |          |          |          |
|                   | HCPCS  | E2291                  | E2291    | E2291    | E2291    |
|                   | 10"    | PXAL1010               |          |          |          |
|                   | 12"    | PXAL1012               | PXAL1212 | PXAL1412 |          |
|                   | 14"    | PXAL1014               | PXAL1214 | PXAL1414 | PXAL1614 |
|                   | 16"    |                        | PXAL1216 | PXAL1416 | PXAL1616 |
|                   | 18"    |                        |          |          | PXAL1618 |

### Axiom for Kids - Adjustable Laterals Trunk Supports (for PXAL Back)

*Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)*

| Description                                                                             | Item Code   | HCPCS        | Retail |
|-----------------------------------------------------------------------------------------|-------------|--------------|--------|
| 3 x 4 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)              | PLF3x4PP    | E0956        | \$305  |
| 4 x 5 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)              | PLF4x5PP    | E0956        | \$305  |
| 3 x 4 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1028)   | PLTAS3x4PP  | E0956, E1028 | \$790  |
| 4 x 5 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1028)   | PLTAS4x5PP  | E0956, E1028 | \$790  |
| 3 x 4 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1028) | PLTASC3x4PP | E0956, E1028 | \$790  |
| 4 x 5 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1028) | PLTASC4x5PP | E0956, E1028 | \$790  |

### Axiom for Kids - Adjustable Lateral Pelvic Supports (for PXAL Back)


*Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)*

| Description                                                                      | Item Code  | HCPCS        | Retail |
|----------------------------------------------------------------------------------|------------|--------------|--------|
| 3 x 4 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)            | PLPF3x4PP  | E0956        | \$305  |
| 4 x 5 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)            | PLPF4x5PP  | E0956        | \$305  |
| 3 x 4 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1028) | PLASP3x4PP | E0956, E1028 | \$790  |
| 4 x 5 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1028) | PLASP4x5PP | E0956, E1028 | \$790  |

## Axiom for Kids ALCL - Posterior Back for Contoured Adjustable Laterals

*Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back options and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.*

|                                           |                             |                             |                             |                             |
|-------------------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| <b>Back Width</b>                         | 10" (25.5 cm)               | 12" (30.5 cm)               | 14" (35.5 cm)               | 16" (40.5 cm)               |
| <b>Fits Back Frame Widths:</b>            | 9"-11" (23 cm - 28 cm)      | 11"-13" (28 cm - 33 cm)     | 13"-15" (33 cm - 38 cm)     | 15"-17" (38 cm - 43 cm)     |
| <b>Actual Back Width (F):</b>             | 8" - 11" (20 cm - 28 cm)    | 10" - 13" (25.5 cm - 33 cm) | 12" - 15" (30.5 cm - 38 cm) | 14" - 17" (35.5 cm - 43 cm) |
| <b>Contoured Lateral Width Range (J):</b> | 8" - 11" (20 cm - 28 cm)    | 10" - 13" (25.5 cm - 33 cm) | 12" - 15" (30.5 cm - 38 cm) | 14" - 17" (35.5 cm - 43 cm) |
| <b>Contour Depth (K):</b>                 | 5" - 5.5" (12.5 cm - 14 cm) | 5" - 5.5" (12.5 cm - 14 cm) | 5" - 5.5" (12.5 cm - 14 cm) | 5" - 5.5" (12.5 cm - 14 cm) |



### Back

*Requires Contoured Adjustable Lateral Trunk Support selection below. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.*

| Axiom Back Height | Wheelchair Frame Width |            |            |            |
|-------------------|------------------------|------------|------------|------------|
|                   | 10"                    | 12"        | 14"        | 16"        |
|                   | Retail                 | \$970      |            |            |
|                   | HCPCS                  | E2293      | E2293      | E2293      |
|                   | 10"                    | PXAL1010CL |            |            |
|                   | 12"                    | PXAL1012CL | PXAL1212CL | PXAL1412CL |
|                   | 14"                    | PXAL1014CL | PXAL1214CL | PXAL1414CL |
|                   | 16"                    |            | PXAL1216CL | PXAL1416CL |
|                   | 18"                    |            |            | PXAL1618CL |

### Axiom for Kids - Contoured Adjustable Lateral Trunk Supports (for PXALCL Back)

| Description                                                            | Item Code | HCPCS | Retail |
|------------------------------------------------------------------------|-----------|-------|--------|
| 3" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956) | PCL3      | E0956 | \$325  |
| 4" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956) | PCL4      | E0956 | \$325  |

*4" Contoured Lateral Wings are not available on 10" or 12" Back Heights*

## Accessories

### Axiom for Kids Back Accessories & Modifications

|                                                                              |      |
|------------------------------------------------------------------------------|------|
| Lumbar Pad (K0108)                                                           | \$53 |
| Privacy Flap (K0108)                                                         | \$53 |
| <i>Privacy Flap not available for 10" Wide backs</i>                         |      |
| Axiom DEX for Kids Hardware Extension (K0108)                                | \$95 |
| <i>Hardware Extension adds an additional 2" of anterior depth adjustment</i> |      |
| Axiom for Kids Adjustable Shoulder Guides (K0108)                            | \$97 |

### Headrests

| Description                         | Standard Straight Removable Adjustable |        | Multi-Axis Removable Adjustable Hardware |        | Posalinc Removable Adjustable Hardware |        | Posalinc Flip-Back Adjustable Hardware |        | HCPCS       |
|-------------------------------------|----------------------------------------|--------|------------------------------------------|--------|----------------------------------------|--------|----------------------------------------|--------|-------------|
|                                     | Item Code                              | Retail | Item Code                                | Retail | Item Code                              | Retail | Item Code                              | Retail |             |
| Comfort Soft Headrest 8"            | CS08SS                                 | \$347  | CS08MA                                   | \$369  | CS08PL                                 | \$519  | CS08PLF                                | \$585  | E0955/E1028 |
| Comfort Soft Headrest 10"           | CS10SS                                 | \$347  | CS10MA                                   | \$369  | CS10PL                                 | \$519  | CS10PLF                                | \$585  | E0955/E1028 |
| Comfort Soft Headrest 14"           | CS14SS                                 | \$347  | CS14MA                                   | \$369  | CS14PL                                 | \$519  | CS14PLF                                | \$585  | E0955/E1028 |
| Contour Positioning Pad, 8" X 6.25" | CP08SS                                 | \$347  | CP08MA                                   | \$369  | CP08PL                                 | \$519  | CP08PLF                                | \$585  | E0955/E1028 |
| Contour Positioning Pad, 9" X 6.5"  | CP09SS                                 | \$347  | CP09MA                                   | \$369  | CP09PL                                 | \$519  | CP09PLF                                | \$585  | E0955/E1028 |
| Contour Positioning Pad, 10" X 7.5" | CP10SS                                 | \$347  | CP10MA                                   | \$369  | CP10PL                                 | \$519  | CP10PLF                                | \$585  | E0955/E1028 |
| Tri-soft Headrest 10"               | TS10SS                                 | \$347  | TS10MA                                   | \$369  | TS10PL                                 | \$519  | TS10PLF                                | \$585  | E0955/E1028 |
| Tri-soft Headrest 14"               | TS14SS                                 | \$347  | TS14MA                                   | \$369  | TS14PL                                 | \$519  | TS14PLF                                | \$585  | E0955/E1028 |

### Headrest Hardware

*Not required with Axiom Headrests*

| Description                       | Item Code | HCPCS | Retail |
|-----------------------------------|-----------|-------|--------|
| Universal Headrest Mounting Plate | UMP       | K0108 | \$53   |



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Anterior Trunk Supports

| Description                                                                | Item Code | HCPCS | Retail |
|----------------------------------------------------------------------------|-----------|-------|--------|
| Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small         | 14941     | E0960 | \$143  |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium        | 14942     | E0960 | \$143  |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large         | 14943     | E0960 | \$143  |
| Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small  | 14924     | E0960 | \$143  |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium | 14927     | E0960 | \$143  |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large  | 14928     | E0960 | \$143  |
| Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small                | 14851     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Small                  | 14852     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Medium                 | 14853     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Large                  | 14854     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small                | 14856     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small                  | 14857     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium                 | 14858     | E0960 | \$175  |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large                  | 14859     | E0960 | \$175  |
| AirLogic Anterior Trunk Support: Standard, Structure, Small                | 14282     | E0960 | \$199  |
| AirLogic Anterior Trunk Support: Standard, Structure, Medium               | 14283     | E0960 | \$199  |
| AirLogic Anterior Trunk Support: Standard, Structure, Large                | 14284     | E0960 | \$199  |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small                  | 14257     | E0960 | \$199  |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium                 | 14258     | E0960 | \$199  |

## Gloss Finish Back Shell

*If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.*

|                    |     |
|--------------------|-----|
| Black              | N/C |
| Ruby Red           | N/C |
| Radiant Orange     | N/C |
| Yellow             | N/C |
| Electric Green     | N/C |
| Rich Teal          | N/C |
| Glacier Blue       | N/C |
| Blue Steel         | N/C |
| Deep Sapphire Blue | N/C |
| Purple Plum        | N/C |
| Opalescent Violet  | N/C |
| Dark Pink          | N/C |
| Light Pink         | N/C |
| Champagne          | N/C |

## Matte Finish Back Shell

*If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.*

|                 |     |
|-----------------|-----|
| Iridescent Navy | N/C |
| Bright Red      | N/C |
| Tactical Green  | N/C |
| Dry Earth Tan   | N/C |
| Platinum        | N/C |
| White           | N/C |