


(Discontinued) Rogue / Rogue XP One Arm Drive Blade Style



One Arm Drive Kits do not include rear tires <p>

If retrofitting to add the One Arm Drive, new rear towers are required in addition to the complete One Arm Drive system. If retrofitting from Hex One Arm Drive, you will need a new drive shaft, axles (#12 and #14), One Arm Drive spacer (#6), Axle Mount Assemblies (#16) and installation instruction DCN0504. Reference Rear Wheel Spacing order form No-Go's to determine new rear wheel spacing.

Item Number	Part Number	Description	Quantity	Retail
1-6, 7a, 8a, 9a, 10 - 12, 13a, 14a, 15, 16, 17a, 18a, 19, 20	109496	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 8.72" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13a, 14b, 15, 16, 17a, 18a, 19, 20	109497	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 8.72" Drive Shaft, Short Axle Assembly	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13b, 14a, 15, 16, 17a, 18a, 19, 20	109498	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 9.72" Drive Shaft	1	\$2,095.95

Item Number	Part Number	Description	Quantity	Retail
1-6, 7a, 8a, 9a, 10 - 12, 13c, 14a, 15, 16, 17a, 18a, 19, 20	109499	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 10.74" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13d, 14a, 15, 16, 17a, 18a, 19, 20	109500	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 11.72" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13e, 14a, 15, 16, 17a, 18a, 19, 20	109501	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 12.72" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13f, 14a, 15, 16, 17a, 18a, 19, 20	109502	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 13.82" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13g, 14a, 15, 16, 17a, 18a, 19, 20	109503	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 14.67" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13i, 14a, 15, 16, 17a, 18a, 19, 20	109505	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 16.67" Drive Shaft	1	\$2,095.95
1-6, 7a, 8a, 9a, 10 - 12, 13j, 14a, 15, 16, 17a, 18a, 19, 20	109506	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Aluminum Drive Handrim, 17.67" Drive Shaft	1	\$2,095.95
1-6, 7d, 8a, 9a, 10 - 12, 13a, 14a, 15, 16, 17a, 18a, 19, 20	109507	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 8.72" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13a, 14b, 15, 16, 17a, 18a, 19, 20	109508	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 8.72" Drive Shaft, Short Axle Assembly	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13b, 14a, 15, 16, 17a, 18a, 19, 20	109509	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 9.72" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13c, 14a, 15, 16, 17a, 18a, 19, 20	109510	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 10.74" Drive Shaft	1	\$2,306.07

Item Number	Part Number	Description	Quantity	Retail
1-6, 7d, 8a, 9a, 10 - 12, 13d, 14a, 15, 16, 17a, 18a, 19, 20	109511	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 11.72" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13e, 14a, 15, 16, 17a, 18a, 19, 20	109512	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 12.72" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13f, 14a, 15, 16, 17a, 18a, 19, 20	109513	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 13.82" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13g, 14a, 15, 16, 17a, 18a, 19, 20	109514	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 14.67" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13i, 14a, 15, 16, 17a, 18a, 19, 20	109516	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 16.67" Drive Shaft	1	\$2,306.07
1-6, 7d, 8a, 9a, 10 - 12, 13j, 14a, 15, 16, 17a, 18a, 19, 20	109517	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 20" Mag with Poly Coated Drive Handrim, 17.67" Drive Shaft	1	\$2,306.07
1-6, 7b, 8b, 9b, 10 - 12, 13a, 14a, 15, 16, 17a, 18b, 19, 20	109518	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 8.72" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13a, 14b, 15, 16, 17a, 18b, 19, 20	109519	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 8.72" Drive Shaft, Short Axle Assembly	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13b, 14a, 15, 16, 17a, 18b, 19, 20	109520	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 9.72" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13c, 14a, 15, 16, 17a, 18b, 19, 20	109521	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 10.74" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13d, 14a, 15, 16, 17a, 18b, 19, 20	109522	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 11.72" Drive Shaft	1	\$2,095.95

Item Number	Part Number	Description	Quantity	Retail
1-6, 7b, 8b, 9b, 10 - 12, 13e, 14a, 15, 16, 17a, 18b, 19, 20	109523	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 12.72" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13f, 14a, 15, 16, 17a, 18b, 19, 20	109524	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 13.82" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13g, 14a, 15, 16, 17a, 18b, 19, 20	109525	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 14.67" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13i, 14a, 15, 16, 17a, 18b, 19, 20	109527	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 16.67" Drive Shaft	1	\$2,095.95
1-6, 7b, 8b, 9b, 10 - 12, 13j, 14a, 15, 16, 17a, 18b, 19, 20	109528	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Aluminum Drive Handrim, 17.67" Drive Shaft	1	\$2,095.95
1-6, 7e, 8b, 9b, 10 - 12, 13a, 14a, 15, 16, 17a, 18b, 19, 20	109529	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 8.72" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13a, 14b, 15, 16, 17a, 18b, 19, 20	109530	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 8.72" Drive Shaft, Short Axle Assembly	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13b, 14a, 15, 16, 17a, 18b, 19, 20	109531	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 9.72" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13c, 14a, 15, 16, 17a, 18b, 19, 20	109532	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 10.74" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13d, 14a, 15, 16, 17a, 18b, 19, 20	109533	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 11.72" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13e, 14a, 15, 16, 17a, 18b, 19, 20	109534	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 12.72" Drive Shaft	1	\$2,306.07

Item Number	Part Number	Description	Quantity	Retail
1-6, 7e, 8b, 9b, 10 - 12, 13f, 14a, 15, 16, 17a, 18b, 19, 20	109535	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 13.82" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13g, 14a, 15, 16, 17a, 18b, 19, 20	109536	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 14.67" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13i, 14a, 15, 16, 17a, 18b, 19, 20	109538	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 16.67" Drive Shaft	1	\$2,306.07
1-6, 7e, 8b, 9b, 10 - 12, 13j, 14a, 15, 16, 17a, 18b, 19, 20	109539	Rogue Blade One Arm Drive Retrofit Kit, 0° Camber, 22" Mag with Poly Coated Drive Handrim, 17.67" Drive Shaft	1	\$2,306.07
	DCN0504	Instruction - Clik and Rogue Blade One Arm Drive Retrofit Installation	1	
1-6, 12, 14a, 15	108075	Blade One Arm Drive Axle Assembly	1	\$637.50
1-6, 12, 14b, 15	108079	Blade One Arm Drive Axle Assembly Short <i>Used on 11" seat width with rear wheel spacings less than 1.75" or 12" seat width with rear wheel spacing of 1.25"</i>	1	\$504.90
1	101459	1/2-20 NTE JAM NYL INS L/N BLK ZC	2	\$1.58
2	108082	1/2" Washer .030 Thick	6	\$1.02
3	104056	10-32 X 3/4" Socket Head Cap Screw	4	\$1.02
5	108083	Hub Handrim, One Arm Drive	1	\$81.46
6	107564	One Arm Drive Spacer <i>Qty 1 always used on drive side. Qty 2 used with standard 5.03" axle (part number 107574)</i>	2	\$21.21
7a	003200	One Arm Drive Handrim Aluminum 20"	1	\$262.65
7b	003199	One Arm Drive Handrim Aluminum 22"	1	\$262.65
7c	003198	One Arm Drive Handrim Aluminum 24"	1	\$262.65
7d	003203	One Arm Drive Handrim Poly Coated 20"	1	\$414.99
7e	003202	One Arm Drive Handrim Poly Coated 22"	1	\$414.99
7f	003201	One Arm Drive Handrim Poly Coated 24"	1	\$414.99
8a	101897	Handrim Aluminum Anodized 5 Rivet Nut 20"	1	\$68.29
8b	101898	Handrim Aluminum Anodized 5 Rivet Nut 22"	1	\$68.29

Item Number	Part Number	Description	Quantity	Retail
8c	100975	Handrim Aluminum Anodized 5 Rivet Nut 24"	1	\$57.78
9a	101961	Wheel X Core Mag 5-Spoke 20" Black	1	\$147.08
9b	101962	Wheel X Core Mag 5-Spoke 22" Black	1	\$141.83
9c	100960	Wheel X Core Mag 5-Spoke 24" Black	1	\$136.58
10	101717	M5 x 25 SHCS BLK ZC	4	\$0.53
11	100657	M5 Nylock Nut Blk Zn	8	\$0.53
13a	107562	One Arm Drive - Blade Drive Shaft 8.72"	1	\$562.82
13b	107563	One Arm Drive - Blade Drive Shaft 9.72"	1	\$562.82
13c	107032	One Arm Drive - Blade Drive Shaft 10.74"	1	\$562.82
13d	107565	One Arm Drive - Blade Drive Shaft 11.72"	1	\$562.82
13e	107566	One Arm Drive - Blade Drive Shaft 12.72"	1	\$562.82
13f	107567	One Arm Drive - Blade Drive Shaft 13.82"	1	\$562.82
13g	107568	One Arm Drive - Blade Drive Shaft 14.67"	1	\$562.82
13i	107570	One Arm Drive - Blade Drive Shaft 16.67"	1	\$562.82
13j	107571	One Arm Drive - Blade Drive Shaft 17.67"	1	\$562.82
14a	107574	One Arm Drive - Blade Style Axle 5.03"	1	\$92.09
14b	108073	One Arm Drive - Blade Style Axle 4.33" <i>Used on 11" seat width with rear wheel spacings less than 1.75" or 12" seat width with rear wheel spacing of 1.25"</i>	1	\$87.36
15	108085	Key, Blade One Arm Drive, 3/16" x 2"	1	\$16.29
16	108074	One Arm Drive Rogue/Clik Axle Mount 0 Degree Camber PAIR	1	\$476.38
2, 10, 11, 16, 20	110116	One Arm Drive Blade Axle Mount Assembly 0 Degree Camber with Hardware, Rogue <i>This kit is for one axle mount with attaching hardware</i>	1	\$483.70
17a, 18a, 19	505175	20" X Core Mag with One Arm Drive Hub	1	\$254.68
17a, 18b, 19	505176	22" X Core Mag with One Arm Drive Hub	1	\$254.68
17b, 18c, 19	003238	24" X Core Mag with One Arm Drive Hub	1	\$254.68
20	101734	M5 x 35 SHCS BLK ZC	4	\$1.05

Rogue/Rogue XP Blade OAD Drive Shaft Configurations Guide

Rogue Seat Width	Wheel Spacing	Part No.	Compressed Drive Shaft Length	Extended Drive Shaft Length
11"	1, 1.25	*~#107562	7.12"	8.72"
11"	1.5	*#107562	7.12"	8.72"
11"	1.75	~107562	7.12"	8.72"
12"	1, 1.25	*107562	7.12"	8.72"
12"	1.5	~107562	7.12"	8.72"
12"	1.75	~107563	8.12"	9.72"
13"	1, 1.25	107562	7.12"	8.72"
13"	1.5	107563	8.12"	9.72"
13"	1.75	~107032	9.14"	10.74"
14"	1, 1.25	107563	8.12"	9.72"
14"	1.5	107032	9.14"	10.74"
14"	1.75	~107565	10.12"	11.72"
15"	1, 1.25	107032	9.14"	10.74"
15"	1.5	107565	10.12"	11.72"
15"	1.75	~107566	11.12"	12.72"
16"	1, 1.25	107565	10.12"	11.72"
16"	1.5	107566	11.12"	12.72"
16"	1.75	~107567	12.22"	13.82"
17"	1, 1.25	107566	11.12"	12.72"
17"	1.5	107567	12.22"	13.82"
17"	1.75	~107568	13.07"	14.67"
18"	1, 1.25	107567	12.22"	13.82"
18"	1.5	107568	13.07"	14.67"
18"	1.75	~107569	14.07"	15.67"
19"	1, 1.25	107568	13.07"	14.67"
19"	1.5	107569	14.07"	15.67"
19"	1.75	~107570	15.07"	16.67"
20"	1, 1.25	107569	14.07"	15.67"
20"	1.5	107570	15.07"	16.67"
20"	1.75	~107571	16.07"	17.67"

Rogue XP Seat Width	Wheel Spacing	Part No.	Compressed Drive Shaft Length	Extended Drive Shaft Length
11"	1, 1.25	*107562	7.12"	8.72"
11"	1.5	*107562	7.12"	8.72"
11"	1.75	107562	7.12"	8.72"
12"	1, 1.25	107562	7.12"	8.72"
12"	1.5	107562	7.12"	8.72"
12"	1.75	107563	8.12"	9.72"
13"	1, 1.25	107563	8.12"	9.72"
13"	1.5	107563	8.12"	9.72"
13"	1.75	107563	9.14"	10.74"
14"	1, 1.25	107032	9.14"	10.74"
14"	1.5	107032	9.14"	10.74"
14"	1.75	107565	10.12"	11.72"
15"	1, 1.25	107565	10.12"	11.72"
15"	1.5	107565	10.12"	11.72"
15"	1.75	107566	11.12"	12.72"
16"	1, 1.25	107566	11.12"	12.72"
16"	1.5	107566	11.12"	12.72"
16"	1.75	107567	12.22"	13.82"
17"	1, 1.25	107567	12.22"	13.82"
17"	1.5	107568	12.22"	13.82"
17"	1.75	107568	13.07"	14.67"
18"	1, 1.25	107568	13.07"	14.67"
18"	1.5	107568	13.07"	14.67"
18"	1.75	107569	14.07"	15.67"
19"	1, 1.25	107569	14.07"	15.67"
19"	1.5	107569	14.07"	15.67"
19"	1.75	107570	15.07"	16.67"
20"	1, 1.25	107570	15.07"	16.67"
20"	1.5	107570	15.07"	16.67"
20"	1.75	107571	16.07"	17.67"

*For this configuration, do not use spacer 107564 and replace axle 107574 with axle 108073.

~For this configuration, tower assemblies are to be mounted outboard of the frame.

#For this configuration, camber plugs require adjustment during drive shaft installation.

NOTE: 1" and 1.25" wheel spacings are not available with Height Adjustable T-Arms.

NOTE: Tower assemblies are to be mounted outboard of the frame when T-Arms are used.