


# BACK PLATE REPLACEMENT BOLT KITS (QUICKIE IRIS)

EXAMPLE OF A LEFT SIDE  
OFFSET OUTER PLATE

\*NOTE: BOLT QUANTITIES  
ARE 2 PER SIDE



IRISMG5

## BACK PLATE BOLT KIT REFERENCE GUIDE

ACCESSORY TYPE	BOLT LENGTH AND PART NUMBER				
	2"	2 1/4"	2 1/2"	2 3/4"	3"
	221066	221067	221061	221062	221063
DUAL POST ARMREST				AB	
DUAL POST WITH TRANSIT					AB
DUAL POST FLIP-BACK ARMREST		A		B	
DUAL POST FLIP-BACK WITH TRANSIT			A		B
OTHER	AB				
OTHER WITH TRANSIT		AB			