


Rogue XP Camber Tube - Growth



Black Camber Tube kits based on seat width. Color Anodized Camber Tube kits used on all chair widths. <p>

New camber tube required when changing the degree of camber by more than 2 degrees. <p>

If Wheelchair Drive Ready option is present, replacement parts need to be processed through Alber USA.

Item Number	Part Number	Description	Quantity	Retail
1a	500720	Camber Tube Rigid 10" Kit Black	1	\$173.01
1b	500721	Camber Tube Rigid 11" Kit Black	1	\$173.01
1c	500722	Camber Tube Rigid 12" Kit Black	1	\$173.01
1d	500723	Camber Tube Rigid 13" Kit Black	1	\$173.01
1e	500724	Camber Tube Rigid 14" Kit Black	1	\$173.01
1f	500725	Camber Tube Rigid 15" Kit Black	1	\$173.01
1g	500726	Camber Tube Rigid 16" Kit Black	1	\$173.01
1h	500727	Camber Tube Rigid 17" Kit Black	1	\$173.01
1i	500728	Camber Tube Rigid 18" Kit Black	1	\$173.01

Item Number	Part Number	Description	Quantity	Retail
1j	501749	Camber Tube Rigid Blue Anodized	1	\$187.19
1k	501750	Camber Tube Rigid Bright Green Anodized	1	\$187.19
1l	501751	Camber Tube Rigid Orange Anodized	1	\$187.19
1m	501752	Camber Tube Rigid Pink Anodized	1	\$187.19
1n	501753	Camber Tube Rigid Pewter Anodized	1	\$187.19
1o	501754	Camber Tube Rigid Purple Anodized	1	\$187.19
1p	501755	Camber Tube Rigid Red Anodized	1	\$187.19
1q	501756	Camber Tube Rigid Gold Anodized	1	\$187.19
2	100667	M6 x 12 BTN SHCS BLKZC w/Patch	4	\$1.55
3	000066	M6 x 27 Threaded Barrel	2	\$4.37

Color Anodizing Samples

