

# Lightning Mobility Base

The LIGHTNING is a lightweight, highly functional Pediatric Mobility Base designed with the recreational look that compliments today's lifestyle. The LIGHTNING is an innovative vehicle constructed to be super 'tuff' and is especially convenient for transporting and storing due to it's easy to operate "Stand & Fold" feature. For function, performance and durability in a beautiful package... Nothing Beats The LIGHTNING!



## Lightning features.

Heavy Duty Upholstery includes moldable 'stays'

7000 Series Aluminum is used for aircraft structural members because of its superb strength to weight properties

Push Handle (for pushing & control of chair)

Locking Strap

Seat Depth Adjustable

Offset Arms provide additional width

Folding Actuator

Carrying Handle

Front Rigging: Swing Away, Flip-Up Adjusts in the following areas- extension, fore & aft, angle & lateral

Storage Basket

Wheel Locks ( operates both sides with single actuation )



## Accessories Available for the Lightning

**SE** Standard Accessories of the SE Model. **Upholstery Color Choices:** NAVY BLUE  RED GLOW 

Photo of the Standard SE Model.



### Accessories.



HEADREST

**SE** OPTIONAL



HEADREST EXTENSION

**SE** OPTIONAL




CANOPY

**SE** OPTIONAL



SINGLE FLAP

**SE** OPTIONAL



HARNESSES COVERS

**SE** OPTIONAL



DOUBLE FLAP

**SE** OPTIONAL



HEEL LOOPS/ANKLE CUFFS

**SE** OPTIONAL



SEAT EXTENSION KIT

OPTIONAL



SEAT DEPTH REDUCTION KIT

OPTIONAL



WORK-N-PLAY TRAY

OPTIONAL



TIE DOWN

OPTIONAL

## Foot Positioning Options



The Positioning Footrests will easily adjust between 1" & 17" extension in 1" increments. The riggings remove easily & the Positioning Pin is tethered to prevent loss.



Foot plate incorporates 50 degrees of adjustable angle for Dorsa & Plantar Flexion, & adjusts laterally over a 2" range.



Foot plate incorporates Fore & Aft adjustment over a 1.5" range.



### LIGHTNING CRASH TESTING:

**Transportation-** The Lightning meets the requirements specified in 5.3 of section 19 of Ansi/Resna WC/Vol.1 prior to May 2002.

**Point Of Impact-** The Lightning is so tough that it exhibits virtually no distortion at point of impact.



## Lightning Specifications

A RANGE OF OPTIONS & ACCESSORIES AVAILABLE.

AVAILABLE MODELS LT11, LT11SE, LT14, LT14SE, LT16 & LT16SE.

### Specification Chart:

SPECIFICATIONS	LIGHTNING 11	LIGHTNING 14	LIGHTNING 16	LIGHTNING SE
Overall Seat Width	12.5"	15.5"	16.5"	✓
Seat Depth*	10"-12"	12"-14"	14"-16"	✓
Seat Angle	30 deg	30 deg	30 deg	✓
OA Ht./Push Handle Ht.	36"	36"	36"	✓
Folded Dimensions	19.25"Wx19.25Dx25.5H	24"Wx19.25"Dx27.5H	25"Wx19.25Dx28H	✓
Weight Capacity	75 lbs	100 lbs	150 lbs	✓
Weight of Chair without Riggings	22 lbs	25 lbs	26 lbs	✓
Foldable	Yes	Yes	Yes	✓
Transit Ready	Yes	Yes	Yes	✓
Color Choices Frame	Silver	Silver	Silver	✓
SPECIFICATIONS	LIGHTNING 11	LIGHTNING 14	LIGHTNING 16	LIGHTNING SE
Removable Adjustable SwingAway Footrests	Yes	Yes	Yes	✓
Depth Adjustable Crotch Strap	Yes	Yes	Yes	✓
Heavy Duty Reinforced Upholstery	Yes	Yes	Yes	✓
SPECIFICATIONS	LIGHTNING 11	LIGHTNING 14	LIGHTNING 16	LIGHTNING SE
Headrest - Winged	Yes	Yes	Yes	STD
Headrest Extension	Yes	Yes	Yes	STD
Weather Protector	Yes	Yes	Yes	✓
Harness Covers	Yes	Yes	Yes	STD
Heel Loops & Ankle Cuffs	Yes	Yes	Yes	STD
Work-N-Play Tray	Yes	Yes	Yes	✓
Canopy	Yes	Yes	Yes	STD
Tie Down Kit	Yes	Yes	Yes	✓
Basket	STD	STD	STD	STD
2" Deeper Seat Kits	Yes	Yes	Yes	✓
Single Flap Lat. Support	Yes	Yes	Yes	STD
Double Flap Lat. Support	Yes	Yes	Yes	STD