


CR45 Transit Forks and Stems



| Item Number | Part Number | Description | Quantity |
|-------------------|-------------|--|----------|
| 4, 5, 6, 13, 14a | 123024 | Transit Fork Assembly, Black, 4" Std Stem, Tilt | 1 |
| 4, 5, 7a, 13, 14b | 123109 | Transit Fork Assembly 5" Std Stem, Tilt | 1 |
| 4, 5, 7b, 13, 14b | 123107 | Transit Fork Assembly 6" Std Stem, Tilt | 1 |
| 4, 5, 8, 13, 15 | 123110 | Transit Fork Assembly 6" Wide Std Stem, Tilt <i>Only used with 2" wide casters.</i> | 1 |
| 4, 5, 9, 13, 16 | 123108 | Transit Fork Assembly 7" Std Stem, Tilt | 1 |
| 4, 5, 10, 13, 17 | 123111 | Transit Fork Assembly 7" Wide Std Stem, Tilt | 1 |

| Item Number | Part Number | Description | Quantity |
|----------------------|-------------|--|----------|
| | | <i>Only used with 2" wide casters.</i> | |
| 1, 3, 5, 6, 13, 14a | 123120 | Transit Fork Assembly 4" Tall Stem, Tilt | 1 |
| 1, 3, 5, 7a, 13, 14b | 123127 | Transit Fork Assembly 5" Tall Stem Tilt | 1 |
| 1, 3, 5, 7b, 13, 14b | 123128 | Transit Fork Assembly 6" Tall Stem Tilt | 1 |
| 1, 3, 5, 8, 13, 15 | 123129 | Transit Fork Assembly 6" Wide Tall Stem Tilt <i>Only used with 2" wide casters.</i> | 1 |
| 1, 3, 5, 9, 13, 16 | 123130 | Transit Fork Assembly 7" Tall Stem Tilt | 1 |
| 2, 11a, 12a | 502653 | Caster Attachment Hardware Kit CAT/FCS | 1 |
| 2, 11b, 12b | 502654 | Caster Attachment Hardware Kit for 6" and 7" Wide Fork, CAT/FCS <i>Only used with 2" wide casters</i> | 1 |
| 1 | 002476 | Spacer, 3/4" OD x 1/2" ID x 1.23 | 1 |
| 2 | 100659 | M8 Nylock Nut Blk Zn | 1 |
| 1, 3, 5, 13 | 506374 | Caster Stem Tall, Adult Tilt | 1 |
| 4, 5, 13 | 503819 | Caster Stem Standard, Tilt | 1 |
| 5 | 100529 | M12 x 1.25mm Lock Nut Black | 2 |
| 5, 6, 14a | 123132 | Transit Fork 4", Tilt | 1 |
| 5, 7a, 14b | 123112 | Transit Fork 5" Tilt | 1 |
| 5, 7b, 14b | 123113 | Transit Fork 6" Tilt | 1 |
| 5, 8, 15 | 123114 | Transit Fork 6" Wide Tilt <i>Only used with 2" wide casters</i> | 1 |
| 5, 9, 16 | 123115 | Transit Fork 7" Tilt | 1 |
| 5, 10, 17 | 123116 | Transit Fork 7" Wide Tilt <i>Only used with 2" wide casters</i> | 1 |
| 11a | 101829 | M8 x 65 HHCS Blk Zn | 1 |
| 11b | 101831 | M8 x 80 HHCS Blk Zn <i>Only used with 2" wide casters</i> | 1 |
| 12a | 100792 | Caster Spacer D8 x D12 x 10mm / Close Mount Handrim Spacer <i>Used with all casters except 6x2</i> | 2 |
| 12b | 001246 | Caster Spacer Black 6 x 2 Poly .472ODX.319IDX.370LG <i>Only used with 6" x 2" casters</i> | 2 |
| 13.0 | 100682 | Washer 1/2" .518 x .875 x .047 F/W Black Zinc | 3 |