


# Ethos Height Adjustable T-Arm - Growth



If Height Adjustable T-Arms are present, depending on change in depth, seat rail saddles and attaching hardware may be required

Item Number	Part Number	Description	Quantity
1	005259	Seat Rail Saddle, Offset 1.125" <i>Use quantity 4 per side for 14" - 15" seat depths, quantity 5 per side for 16" - 18" seat depths and quantity 6 per side for 19" - 20" seat depths. If Seat Option is Omit Seat Sling, Seat Rail Saddles are not used.</i>	6
2	101117	M5 x 22 BTN SHCS BLK ZC <i>Use quantity 4 per side for 14" - 15" seat depths, quantity 5 per side for 16" - 18" seat depths and quantity 6 per side for 19" - 20" seat depths</i>	6