


Little Wave Plastic Coated Handrim



Tab mounted handrims must be installed prior to the tire being installed.

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|--|----------|---------|
| 8-10 | 505263 | Handrim Hardware Kit, 3 Rivet Nut, EACH <i>Used with 18" Ki Spoke and 18" MAXX Spoke wheels</i> | 1 | \$19.31 |
| 11 | 101106 | Handrim Plastic Coated 3 Rivet Nut 18" | 1 | \$74.28 |

| Item Number | Part Number | Description | Quantity | Retail |
|-------------|-------------|---|----------|----------|
| | | <i>Used on 18" Ki Spoke and 18" MAXX Spoke wheels</i> | | |
| 8-10 | 502168 | Handrim Hardware Kit, 6 Rivet Nut, EACH <i>Used with Ki Spoke and MAXX Spoke wheels</i> | 1 | \$38.61 |
| 7a | 200542 | Handrim Plastic Coated 6 Rivet Nut 20" <i>Used on 20" Ki Spoke and 20" MAXX Spoke wheels</i> | 1 | \$160.34 |
| 7b | 100576 | Handrim Plastic Coated 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke wheels</i> | 1 | \$143.41 |
| 7c | 112821 | Handrim Plastic Coated 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke wheels. Replaced part 100577.</i> | 1 | \$143.41 |
| 5, 6 | 502173 | Handrim Hardware Kit, 5 Rivet Nut, Plastic Coated, EACH <i>Used with X Core Mag wheels</i> | 1 | \$33.10 |
| 4a | 101963 | Handrim Plastic Coated 5 Rivet Nut 20" <i>Used on 20" 5-Spoke X Core Mag wheels</i> | 1 | \$140.71 |
| 4b | 101964 | Handrim Plastic Coated 5 Rivet Nut 22" <i>Used on 22" 5-Spoke X Core Mag wheels</i> | 1 | \$140.71 |
| 4c | 100976 | Handrim Plastic Coated 5 Rivet Nut 24" <i>Used on 24" 5-Spoke X Core Mag wheels</i> | 1 | \$132.37 |
| 2, 3 | 502174 | Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy and Topolino wheels</i> | 1 | \$16.55 |
| 1a | 100836 | Handrim Plastic Coated 6 Tab 24" <i>Used on 24" Superlight/Topolino wheels</i> | 1 | \$163.98 |
| 1b | 103179 | Handrim Plastic Coated Spinergy 6 Tab 20" <i>Used on 20" Spinergy Spox wheels</i> | 1 | \$204.08 |
| 1c | 100808 | Handrim Plastic Coated Spinergy 6 Tab 22" <i>Used on 22" Spinergy wheels</i> | 1 | \$204.08 |
| 1d | 112815 | Handrim Plastic Coated Spinergy 6 Tab 24" <i>Used on 24" Spinergy wheels. If chair was invoiced prior to 8/14/19, quantity 2 must be ordered due to a change in finish. Replaced 100615.</i> | 1 | \$194.36 |
| 2 | 100669 | M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinergy and Topolino wheels</i> | 6 | \$1.66 |
| 3 | 100657 | M5 Nylock Nut Blk Zn <i>Used with Spinergy and Topolino wheels</i> | 6 | \$0.56 |
| 5 | 101756 | Handrim Spacer Plastic Coated Riv <i>Used with X Core Mag wheels</i> | 5 | \$0.56 |
| 6 | 100724 | M5 x 55 PH TRS HD M/S ZC <i>Used with X Core Mag wheels</i> | 5 | \$1.10 |

| Item Number | Part Number | Description | Quantity | Retail |
|--------------------|--------------------|---|-----------------|---------------|
| 8 | 100653 | Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke wheels</i> | 6 | \$1.10 |
| 9 | 100698 | Handrim Tab <i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i> | 6 | \$4.17 |
| 10 | 100654 | M5 x 40 DIN7985 PH Pan HD M/S ZC <i>Used with Ki Spoke and MAXX Spoke wheels</i> | 6 | \$1.66 |

Handrim Hardware Chart

| Wheel | | Handrim Connection | | | | Handrim | | | | | | | |
|--------------------|------------------------|--------------------|--|--|--|-------------------|------------|----------------|------------|----------------------|----------------|----------|--|
| Wheel | Wheel Part No. | Connection Points | Tab | Spacer* | Screw* | Aluminum Anodized | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim | |
| 18" Ki Spoke | 200529 | 3 | 100698 | 100653 | Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835 | 100206 | | 101106 | | | | | |
| 20" Ki Spoke | 200530 | | | | | 200536 | | 200542 | | | | | |
| 22" Ki Spoke | 200531 | | | | | 100560 | | 100576 | 100569 | 200538 | 200201 | | |
| 24" Ki Spoke | 200532 | | | | | 200349 | 101161 | 100577 | 200547 | 100793 | 200202 | | |
| 25" Ki Spoke | 200533 | | | | | 200350 | | 101870 | 200548 | 200539 | 200540 | | |
| 26" Ki Spoke | 200534 | | | | | 200351 | | 100578 | 200549 | 100907 | 101454 | | |
| 24" Superlight | 101159 | 6 | | | 100536 | 100754 | 101161 | 100836 | | 100830 | 100828 | | |
| 25" Superlight | 101160 | | | | | | 101160 | 101091 | | 101464 | 101460 | | |
| 20" 5-Spoke X Core | 101961 | 5 | | Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756 | 100724 | 101897 | | 101963 | | | | | |
| 22" 5-Spoke X Core | 101962 | | | | | | 101898 | | 101964 | | | | |
| 24" 5-Spoke X Core | 100960 | | | | | | 100975 | | 100976 | 200546 | | | |
| 24" 5-Spoke X Core | 100960 | | | | | | 101893 | | | | 100768 | 100769 | |
| 20" Spinergy Spox | See Spinergy Spox Page | 6 | | | Screw: 100669 Nut: 100657 | 103125 | | 103179 | | | | | |
| 22" Spinergy Spox | | | | | | 100827 | | 100808 | | 100889 | 100888 | | |
| 24" Spinergy Spox | | | | | | 100766 | 101161 | 100615 | | 100830 | 100828 | | |
| 25" Spinergy Spox | | | | | | 100767 | 101160 | 100765 | | 101464 | 101460 | | |
| 26" Spinergy Spox | | | | | | 101477 | | 101148 | | 200200 | 100950 | | |
| 22" Spinergy LX | See Spinergy LX Page | 6 | | | Screw: 100669 Nut: 100657 | 100827 | | 100808 | | 100889 | 100888 | | |
| 24" Spinergy LX | | | | | | 100766 | 101161 | 100615 | | 100830 | 100828 | 200213 | |
| 25" Spinergy LX | | | | | | 100767 | 101160 | 100765 | | 101464 | 101460 | 200208 | |
| 26" Spinergy LX | | | | | | 101477 | | 101148 | | 200200 | 100950 | | |
| 18" MAXX Spoke | 109285 | 6 | 100698 (Not used w/ Superlight Handrim) | 100653 (Not used w/ Superlight Handrim) | Aluminum, Plastic Ctd, Projection: 100654 Natural Fit: 100835 Superlight: Screw: 100669 Nut: 100657 | 100206 | | 101106 | | | | | |
| 20" MAXX Spoke | 108244 | | | | | 200536 | | 200542 | | | | | |
| 22" MAXX Spoke | 105135 | | | | | 100560 | | 100576 | 100569 | 200538 | 200201 | | |
| 24" MAXX Spoke | 105136 | | | | | 200349 | 101161 | 100577 | 200547 | 100793 | 200202 | | |
| 25" MAXX Spoke | 107436 | | | | | 200350 | | 101870 | 200548 | 200539 | 200540 | | |
| 26" MAXX Spoke | 107437 | | | | | 200351 | | 100578 | 200549 | 100907 | 101454 | | |

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.