


# Focus CR Dynamic Back - Growth



Has your width changed? If yes, verify the current growth strut will work or select appropriately below.

Item Number	Part Number	Description	Quantity
1, 2a	503107	Growth Strut Assembly Short, Focus/Flip <i>Used on chairs 14" - 17" wide</i>	1
1, 2b	503108	Growth Strut Assembly Long, Flip/Focus <i>Used on chairs 18" - 22" wide</i>	1
1, 2c	114289	Growth Strut Assembly Extra Long, Focus <i>Used on chairs 23"-24" wide</i>	1