



## (Discontinued 1) Rogue Style Fixed Height Backrest with Non-Adjustable Height Rigidizer Bar on Rogue ALX - Growth (Formerly Tsunami)



Fixed Height with Omit Push Handles



Fixed Height with Push Handles

*This Rogue Style Backrest is not available with the Tsunami Heavy Duty option*

Item Number	Part Number	Description	Quantity
1a	108760	Fixed Height Rigid Backrest 5", 14" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1b	108762	Fixed Height Rigid Backrest 6", 14" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1c	108764	Fixed Height Rigid Backrest 7", 14" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1d	108766	Fixed Height Rigid Backrest 8", 14" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1

Item Number	Part Number	Description	Quantity
1e	108768	Fixed Height Rigid Backrest 9", 14" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1f	108770	Fixed Height Rigid Backrest 10", 14" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1g	108772	Fixed Height Rigid Backrest 11", 14" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1h	108774	Fixed Height Rigid Backrest 12", 14" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1i	108775	Fixed Height Rigid Backrest 13", 14" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1j	108776	Fixed Height Rigid Backrest 14", 14" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1k	108777	Fixed Height Rigid Backrest 15", 14" Wide <i>Used with 17" back height with omit push handles</i>	1
1l	108778	Fixed Height Rigid Backrest 16", 14" Wide <i>Used with 18" back height with omit push handles</i>	1
1m	108779	Fixed Height Rigid Backrest 5", 15" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1n	108781	Fixed Height Rigid Backrest 6", 15" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1o	108783	Fixed Height Rigid Backrest 7", 15" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1p	108785	Fixed Height Rigid Backrest 8", 15" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1q	108787	Fixed Height Rigid Backrest 9", 15" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1r	108789	Fixed Height Rigid Backrest 10", 15" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	
1s	108791	Fixed Height Rigid Backrest 11", 15" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1t	108793	Fixed Height Rigid Backrest 12", 15" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1u	108794	Fixed Height Rigid Backrest 13", 15" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1v	108795	Fixed Height Rigid Backrest 14", 15" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1w	108796	Fixed Height Rigid Backrest 15", 15" Wide <i>Used with 17" back height with omit push handles</i>	1
1x	108797	Fixed Height Rigid Backrest 16", 15" Wide <i>Used with 18" back height with omit push handles</i>	1
1y	108798	Fixed Height Rigid Backrest 5", 16" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1z	108800	Fixed Height Rigid Backrest 6", 16" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1aa	108802	Fixed Height Rigid Backrest 7", 16" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1ab	108804	Fixed Height Rigid Backrest 8", 16" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ac	108806	Fixed Height Rigid Backrest 9", 16" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1ad	108808	Fixed Height Rigid Backrest 10", 16" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1ae	108810	Fixed Height Rigid Backrest 11", 16" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1

Item Number	Part Number	Description	Quantity
1af	108812	Fixed Height Rigid Backrest 12", 16" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1ag	108813	Fixed Height Rigid Backrest 13", 16" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ah	108814	Fixed Height Rigid Backrest 14", 16" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ai	108815	Fixed Height Rigid Backrest 15", 16" Wide <i>Used with 17" back height with omit push handles</i>	1
1aj	108816	Fixed Height Rigid Backrest 16", 16" Wide <i>Used with 18" back height with omit push handles</i>	1
1ak	108817	Fixed Height Rigid Backrest 5", 17" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1al	108819	Fixed Height Rigid Backrest 6", 17" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1am	108821	Fixed Height Rigid Backrest 7", 17" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1an	108823	Fixed Height Rigid Backrest 8", 17" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ao	108825	Fixed Height Rigid Backrest 9", 17" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1ap	108827	Fixed Height Rigid Backrest 10", 17" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1aq	108829	Fixed Height Rigid Backrest 11", 17" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1ar	108831	Fixed Height Rigid Backrest 12", 17" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1as	108832	Fixed Height Rigid Backrest 13", 17" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	
1at	108833	Fixed Height Rigid Backrest 14", 17" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1au	108834	Fixed Height Rigid Backrest 15", 17" Wide <i>Used with 17" back height with omit push handles</i>	1
1av	108835	Fixed Height Rigid Backrest 16", 17" Wide <i>Used with 18" back height with omit push handles</i>	1
1aw	108836	Fixed Height Rigid Backrest 5", 18" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1ax	108838	Fixed Height Rigid Backrest 6", 18" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1ay	108840	Fixed Height Rigid Backrest 7", 18" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1az	108842	Fixed Height Rigid Backrest 8", 18" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1ba	108844	Fixed Height Rigid Backrest 9", 18" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bb	108846	Fixed Height Rigid Backrest 10", 18" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1bc	108848	Fixed Height Rigid Backrest 11", 18" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1bd	108850	Fixed Height Rigid Backrest 12", 18" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1be	108851	Fixed Height Rigid Backrest 13", 18" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1bf	108852	Fixed Height Rigid Backrest 14", 18" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bg	108853	Fixed Height Rigid Backrest 15", 18" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 17" back height with omit push handles</i>	
1bh	108854	Fixed Height Rigid Backrest 16", 18" Wide <i>Used with 18" back height with omit push handles</i>	1
1bi	108855	Fixed Height Rigid Backrest 5", 19" Wide <i>Used with 9" back height with standard or low profile push handles</i>	1
1bj	108857	Fixed Height Rigid Backrest 6", 19" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1bk	108859	Fixed Height Rigid Backrest 7", 19" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1bl	108861	Fixed Height Rigid Backrest 8", 19" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1bm	108863	Fixed Height Rigid Backrest 9", 19" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bn	108865	Fixed Height Rigid Backrest 10", 19" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1bo	108867	Fixed Height Rigid Backrest 11", 19" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1bp	108869	Fixed Height Rigid Backrest 12", 19" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1bq	108870	Fixed Height Rigid Backrest 13", 19" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1br	108871	Fixed Height Rigid Backrest 14", 19" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bs	108872	Fixed Height Rigid Backrest 15", 19" Wide <i>Used with 17" back height with omit push handles</i>	1
1bt	108873	Fixed Height Rigid Backrest 16", 19" Wide <i>Used with 18" back height with omit push handles</i>	1
1bu	108874	Fixed Height Rigid Backrest 5", 20" Wide	1

Item Number	Part Number	Description	Quantity
		<i>Used with 9" back height with standard or low profile push handles</i>	
1bv	108876	Fixed Height Rigid Backrest 6", 20" Wide <i>Used with 10" back height with standard or low profile push handles or 12" back height with fold down push handles</i>	1
1bw	108878	Fixed Height Rigid Backrest 7", 20" Wide <i>Used with 9" back height with omit push handles, 11" back height with standard or low profile push handles or 13" back height with fold down push handles</i>	1
1bx	108880	Fixed Height Rigid Backrest 8", 20" Wide <i>Used with 10" back height with omit push handles, 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1
1by	108882	Fixed Height Rigid Backrest 9", 20" Wide <i>Used with 11" back height with omit push handles, 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1
1bz	108884	Fixed Height Rigid Backrest 10", 20" Wide <i>Used with 12" back height with omit push handles, 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1
1ca	108886	Fixed Height Rigid Backrest 11", 20" Wide <i>Used with 13" back height with omit push handles, 15" back height with standard or low profile push handles or 17" back height with fold down push handles</i>	1
1cb	108888	Fixed Height Rigid Backrest 12", 20" Wide <i>Used with 14" back height with omit push handles, 16" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1
1cc	108889	Fixed Height Rigid Backrest 13", 20" Wide <i>Used with 15" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1cd	108890	Fixed Height Rigid Backrest 14", 20" Wide <i>Used with 16" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ce	108891	Fixed Height Rigid Backrest 15", 20" Wide <i>Used with 17" back height with omit push handles</i>	1
1cf	108892	Fixed Height Rigid Backrest 16", 20" Wide <i>Used with 18" back height with omit push handles</i>	1
2	100804	Tube End Plug 1"	2
3	000034	Sleeve 3/4" x 1 1/2"	2
4	000073	M5 x 24 Threaded Barrel	2
5	100669	M5 x 10 BTN SHCS BLKZC w/Patch	4