


Note: for part number information, please see attachment in informations panel.

Note: Cushion widths 16" - 20" use B size bladders and 21"-22" use C size bladders.



Obliquity Pad Ordering Guide

Cushion Width	1st Layer	2nd Layer
16"-20"	JBALCOP-BLH JBALCOP-BRH	JBALCOP-CLH JBALCOP-CRH
21"-22"	JBALCOP-CLH JBALCOP-CRH	JBALCOP-DLH JBALCOP-DRH

Note: Obliquity Pad is 1/2":  
 One layer of pads = 1/2" height  
 Two layers of pads = 1" height

Pos.	Item Number	Description	Remarks	USD
1	CJFUSION1616CFHRP	COV FUSION RP CRYO HYB 16X16		
1	CJFUSION1617CFHRP	COV FUSION RP CRYO HYB 16X17		
1	CJFUSION1618CFHRP	COV FUSION RP CRYO HYB 16X18		
1	CJFUSION1619CFHRP	COV FUSION RP CRYO HYB 16X19		
1	CJFUSION1620CFHRP	COV FUSION RP CRYO HYB 16X20		
1	CJFUSION1716CFHRP	COV FUSION RP CRYO HYB 17X16		
1	CJFUSION1717CFHRP	COV FUSION RP CRYO HYB 17X17		
1	CJFUSION1718CFHRP	COV FUSION RP CRYO HYB 17X18		
1	CJFUSION1719CFHRP	COV FUSION RP CRYO HYB 17X19		
1	CJFUSION1720CFHRP	COV FUSION RP CRYO HYB 17X20		
1	CJFUSION1721CFHRP	COV FUSION RP CRYO HYB 17X21		
1	CJFUSION1816CFHRP	COV FUSION RP CRYO HYB 18X16		
1	CJFUSION1817CFHRP	COV FUSION RP CRYO HYB 18X17		
1	CJFUSION1818CFHRP	COV FUSION RP CRYO HYB 18X18		
1	CJFUSION1819CFHRP	COV FUSION RP CRYO HYB 18X19		

Pos.	Item Number	Description	Remarks	USD
1	CJFUSION1820CFHRP	COV FUSION RP CRYO HYB 18X20		
1	CJFUSION1821CFHRP	COV FUSION RP CRYO HYB 18X21		
1	CJFUSION1822CFHRP	COV FUSION RP CRYO HYB 18X22		
1	CJFUSION1916CFHRP	COV FUSION RP CRYO HYB 19X16		
1	CJFUSION1917CFHRP	COV FUSION RP CRYO HYB 19X17		
1	CJFUSION1918CFHRP	COV FUSION RP CRYO HYB 19X18		
1	CJFUSION1919CFHRP	COV FUSION RP CRYO HYB 19X19		
1	CJFUSION1920CFHRP	COV FUSION RP CRYO HYB 19X20		
1	CJFUSION1921CFHRP	COV FUSION RP CRYO HYB 19X21		
1	CJFUSION1922CFHRP	COV FUSION RP CRYO HYB 19X22		
1	CJFUSION1924CFHRP	COV FUSION RP CRYO HYB 19X24		
1	CJFUSION2016CFHRP	COV FUSION RP CRYO HYB 20X16		
1	CJFUSION2017CFHRP	COV FUSION RP CRYO HYB 20X17		
1	CJFUSION2018CFHRP	COV FUSION RP CRYO HYB 20X18		
1	CJFUSION2019CFHRP	COV FUSION RP CRYO HYB 20X19		
1	CJFUSION2020CFHRP	COV FUSION RP CRYO HYB 20X20		
1	CJFUSION2021CFHRP	COV FUSION RP CRYO HYB 20X21		
1	CJFUSION2022CFHRP	COV FUSION RP CRYO HYB 20X22		
1	CJFUSION2118CFHRP	COV FUSION RP CRYO HYB 21X18		
1	CJFUSION2119CFHRP	COV FUSION RP CRYO HYB 21X19		
1	CJFUSION2120CFHRP	COV FUSION RP CRYO HYB 21X20		
1	CJFUSION2121CFHRP	COV FUSION RP CRYO HYB 21X21		
1	CJFUSION2122CFHRP	COV FUSION RP CRYO HYB 21X22		
1	CJFUSION2218CFHRP	COV FUSION RP CRYO HYB 22X18		
1	CJFUSION2219CFHRP	COV FUSION RP CRYO HYB 22X19		
1	CJFUSION2220CFHRP	COV FUSION RP CRYO HYB 22X20		
1	CJFUSION2221CFHRP	COV FUSION RP CRYO HYB 22X21		
1	CJFUSION2222CFHRP	COV FUSION RP CRYO HYB 22X22		
2	CFHB	CRYO FLUID HYBRID BLADDER B	Used in 16"-20" wide	500.00
2	CFHC	CRYO FLUID HYBRID BLADDER C	Used in 21"-22" wide	500.00
3	ICJFUSION1616RP	INN COV FUSION RP FLD 16X16		101.00
3	ICJFUSION1617RP	INN COV FUSION RP FLD 16X17		101.00
3	ICJFUSION1618RP	INN COV FUSION RP FLD 16X18		101.00
3	ICJFUSION1619RP	INN COV FUSION RP FLD 16X19		101.00
3	ICJFUSION1620RP	INN COV FUSION RP FLD 16X20		101.00
3	ICJFUSION1716RP	INN COV FUSION RP FLD 17X16		101.00
3	ICJFUSION1717RP	INN COV FUSION RP FLD 17X17		101.00
3	ICJFUSION1718RP	INN COV FUSION RP FLD 17X18		101.00
3	ICJFUSION1719RP	INN COV FUSION RP FLD 17X19		101.00
3	ICJFUSION1720RP	INN COV FUSION RP FLD 17X20		101.00
3	ICJFUSION1721RP	INN COV FUSION RP FLD 17X21		101.00
3	ICJFUSION1816RP	INN COV FUSION RP FLD 18X16		101.00
3	ICJFUSION1817RP	INN COV FUSION RP FLD 18X17		101.00
3	ICJFUSION1818RP	INN COV FUSION RP FLD 18X18		101.00
3	ICJFUSION1819RP	INN COV FUSION RP FLD 18X19		101.00
3	ICJFUSION1820RP	INN COV FUSION RP FLD 18X20		101.00
3	ICJFUSION1821RP	INN COV FUSION RP FLD 18X21		101.00
3	ICJFUSION1822RP	INN COV FUSION RP FLD 18X22		101.00
3	ICJFUSION1916RP	INN COV FUSION RP FLD 19X16		101.00
3	ICJFUSION1917RP	INN COV FUSION RP FLD 19X17		101.00
3	ICJFUSION1918RP	INN COV FUSION RP FLD 19X18		101.00
3	ICJFUSION1919RP	INN COV FUSION RP FLD 19X19		101.00
3	ICJFUSION1920RP	INN COV FUSION RP FLD 19X20		101.00
3	ICJFUSION1921RP	INN COV FUSION RP FLD 19X21		101.00
3	ICJFUSION1922RP	INN COV FUSION RP FLD 19X22		101.00
3	ICJFUSION2016RP	INN COV FUSION RP FLD 20X16		101.00
3	ICJFUSION2017RP	INN COV FUSION RP FLD 20X17		101.00

Pos.	Item Number	Description	Remarks	USD
3	ICJFUSION2018RP	INN COV FUSION RP FLD 20X18		101.00
3	ICJFUSION2019RP	INN COV FUSION RP FLD 20X19		101.00
3	ICJFUSION2020RP	INN COV FUSION RP FLD 20X20		101.00
3	ICJFUSION2021RP	INN COV FUSION RP FLD 20X21		101.00
3	ICJFUSION2022RP	INN COV FUSION RP FLD 20X22		101.00
3	ICJFUSION2024RP	INN COV FUSION RP FLD 20X24		101.00
3	ICJFUSION2118RP	INN COV FUSION RP FLD 21X18		101.00
3	ICJFUSION2119RP	INN COV FUSION RP FLD 21X19		101.00
3	ICJFUSION2120RP	INN COV FUSION RP FLD 21X20		101.00
3	ICJFUSION2121RP	INN COV FUSION RP FLD 21X21		101.00
3	ICJFUSION2122RP	INN COV FUSION RP FLD 21X22		101.00
3	ICJFUSION2218RP	INN COV FUSION RP FLD 22X18		152.00
3	ICJFUSION2219RP	INN COV FUSION RP FLD 22X19		152.00
3	ICJFUSION2220RP	INN COV FUSION RP FLD 22X20		152.00
3	ICJFUSION2221RP	INN COV FUSION RP FLD 22X21		152.00
3	ICJFUSION2222RP	INN COV FUSION RP FLD 22X22		152.00
3	ICJFUSION2224RP	INN COV FUSION RP FLD 22X24		152.00
4	BJFUSION1616RP	BASE JAY FUSION RPROF 16X16		200.00
4	BJFUSION1617RP	BASE JAY FUSION RPROF 16X17		200.00
4	BJFUSION1618RP	BASE JAY FUSION RPROF 16X18		200.00
4	BJFUSION1619RP	BASE JAY FUSION RPROF 16X19		200.00
4	BJFUSION1620RP	BASE JAY FUSION RPROF 16X20		200.00
4	BJFUSION1716RP	BASE JAY FUSION RPROF 17X16		200.00
4	BJFUSION1717RP	BASE JAY FUSION RPROF 17X17		200.00
4	BJFUSION1718RP	BASE JAY FUSION RPROF 17X18		200.00
4	BJFUSION1719RP	BASE JAY FUSION RPROF 17X19		200.00
4	BJFUSION1720RP	BASE JAY FUSION RPROF 17X20		200.00
4	BJFUSION1721RP	BASE JAY FUSION RPROF 17X21		200.00
4	BJFUSION1816RP	BASE JAY FUSION RPROF 18X16		200.00
4	BJFUSION1817RP	BASE JAY FUSION RPROF 18X17		200.00
4	BJFUSION1818RP	BASE JAY FUSION RPROF 18X18		200.00
4	BJFUSION1819RP	BASE JAY FUSION RPROF 18X19		200.00
4	BJFUSION1820RP	BASE JAY FUSION RPROF 18X20		200.00
4	BJFUSION1821RP	BASE JAY FUSION RPROF 18X21		200.00
4	BJFUSION1822RP	BASE JAY FUSION RPROF 18X22		200.00
4	BJFUSION1916RP	BASE JAY FUSION RPROF 19X16		200.00
4	BJFUSION1917RP	BASE JAY FUSION RPROF 19X17		200.00
4	BJFUSION1918RP	BASE JAY FUSION RPROF 19X18		200.00
4	BJFUSION1919RP	BASE JAY FUSION RPROF 19X19		200.00
4	BJFUSION1920RP	BASE JAY FUSION RPROF 19X20		200.00
4	BJFUSION1921RP	BASE JAY FUSION RPROF 19X21		200.00
4	BJFUSION1922RP	BASE JAY FUSION RPROF 19X22		200.00
4	BJFUSION2016RP	BASE JAY FUSION RPROF 20X16		200.00
4	BJFUSION2017RP	BASE JAY FUSION RPROF 20X17		200.00
4	BJFUSION2018RP	BASE JAY FUSION RPROF 20X18		200.00
4	BJFUSION2019RP	BASE JAY FUSION RPROF 20X19		200.00
4	BJFUSION2020RP	BASE JAY FUSION RPROF 20X20		200.00
4	BJFUSION2021RP	BASE JAY FUSION RPROF 20X21		200.00
4	BJFUSION2022RP	BASE JAY FUSION RPROF 20X22		200.00
4	BJFUSION2024RP	BASE JAY FUSION RPROF 20X24		200.00
4	BJFUSION2118RP	BASE JAY FUSION RPROF 21X18		200.00
4	BJFUSION2119RP	BASE JAY FUSION RPROF 21X19		200.00
4	BJFUSION2120RP	BASE JAY FUSION RPROF 21X20		200.00
4	BJFUSION2121RP	BASE JAY FUSION RPROF 21X21		200.00
4	BJFUSION2122RP	BASE JAY FUSION RPROF 21X22		200.00
4	BJFUSION2218RP	BASE JAY FUSION RPROF 22X18		250.00
4	BJFUSION2219RP	BASE JAY FUSION RPROF 22X19		250.00

Pos.	Item Number	Description	Remarks	USD
4	BJFUSION2220RP	BASE JAY FUSION RPROF 22X20		250.00
4	BJFUSION2221RP	BASE JAY FUSION RPROF 22X21		250.00
4	BJFUSION2222RP	BASE JAY FUSION RPROF 22X22		250.00
4	BJFUSION2424	BASE JFUSION2424		220.00
5A	JBALCOP-BLH	OBLIQUITY PAD LH PLA B	Left- 16 thru 20 wide	16.00
5A	JBALCOP-BRH	OBLIQUITY PAD RH PLA B	Right- 16 thru 20 wide	16.00
5B	JBALCOP-CLH	OBLIQUITY PAD LH PLA C	Left- 21 thru 22 wide	16.00
5B	JBALCOP-CRH	OBLIQUITY PAD RH PLA C	Right- 21 thru 22 wide	16.00
5C	JBALCOP-DLH	OBLIQUITY PAD LH PLA D	Left	16.00
5C	JBALCOP-DRH	OBLIQUITY PAD RH PLA D	Right	16.00