

Bill To (Altimate Medical Supplier Only)	Ship To
Purchase Order #: _____ Quote #: _____	Name: _____
Account #: _____ Ordered By: _____	Address: _____
Supplier Company: _____	Address: _____
Address: _____ City, State, Zip: _____	City, State: _____
Phone: _____ Fax: _____	Zip: _____
Email: _____ Payer Source: _____	Phone: _____
Name of ATP (Required): _____	Tag For: _____

Evolv Medium User Specifications	
Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat pivot to flat back)	14"-19"
Seat to Foot Plate Range (from seat pivot)	11.5"-18"

Base

■ **PNG50162-1 EasyStand[®] Evolv[®] Medium** \$2,182

The Base unit includes: Two 5" Front Wheels, Two 5" Rear Locking Casters, Manual Hydraulic Actuator with Handle, Planar Seat, Height-Only Adjustable Foot Plates, No-Tray Chest Pad, Transfer Handles and Black Upholstery.

Note: The default items will be added if no selection is made.

Default items include: Actuator Handle Locations and 1-Piece Knee Pad.

Choose actuator handle location: Right Side **default** Left Side

Seat Adjustment-CANNOT be ordered off-unit

PNG50499 Easy-Adjust Seat Depth (roller bearing allows for smoother seat depth adjustment) \$108

Frame Styles

PNG50025 Mobile (not available with shadow tray or swing-away front) \$1,469

PNG50314 Front Swivel Casters* \$179

PNG50483 Front Swivel Casters for Swing-Away Front* \$179

PNG50417 Swing-Away Front (tray required)* \$270

PNG50428 Swing-Away Front for Shadow Tray (tray required)* \$211

Lift Mechanisms

PNG50010 Pow'r Up Lift (replaces hydraulic actuator, battery powered) \$1,222

PNG50328 Adjustable/Removable Actuator Handle (not available with pow'r up) \$228

PNG50044 Quad Grip Handle Extension (not available with pow'r up) \$153

PNG50070 T-Style Handle Extension (not available with pow'r up) \$153

Shadow Tray-back support required*

PNG50515 Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back) \$1,415

PNG50517 Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) \$1,579

Front Mounted Trays*

PNG50247-1 Black Molded Tray \$585

PNG50248-1 Clear Tray \$753

PNG50194-1 Oversized Clear Tray \$892

PNG50368-1 Black Molded Angle Adjustable Tray \$934

PNG50366-1 Clear Angle Adjustable Tray \$1,095

Front Mounted Trays "Swing-Away"*

PNG50418-1 Black Molded Tray for Swing-Away Front \$585

PNG50456-1 Clear Tray for Swing-Away Front \$753

PNG50458-1 Oversized Clear Tray for Swing-Away Front \$892

PNG50490-1 Black Molded Angle Adjustable Tray for Swing-Away Front \$934

PNG50489-1 Clear Angle Adjustable Tray for Swing-Away Front \$1,095

PNG50491-1 Oversized Angle Adjustable Tray for Swing-Away Front \$1,232

Mobile Trays

PNG50372-1 Clear Tray for Mobile \$753

PNG50327-1 Clear Angle Adjustable Tray for Mobile \$1,095

Tray Accessories-tray required

PNG50042 Hand Grips (for clear tray only, pair)* \$153

PNG30000 Large Contoured Chest Pad \$142

PNG50331 Large Contoured Adjustable Chest Pad \$204

PNG30295 Padded Tray Cover (not available with oversized clear tray or mobile tray)* \$226

Knee Pads-must select one

PNG50460-1 1-Piece Knee Pad (pair) **default** \$328

PNG50384-1 Independent Knee Pads (pair) \$814

PNG50385-1 Independent Roho[®] Knee Pads (pair) \$936

PNG50424-1 1-Piece Knee Pad for Swing-Away Front (pair) * \$328

PNG50425-1 Independent Knee Pads for Swing-Away Front (pair) * \$814

PNG50426-1 Independent Roho[®] Knee Pads for Swing-Away Front (pair) * \$936

* not available with mobile

Foot Plates

- PNG50414 Multi-Adjustable Foot Plates (replaces standard height adjustable foot plates, pair) \$216

Foot Plate Accessories

- PNG30030 Foot Straps-15"L (length over top of foot, pair) \$46
- PNG50505 Foot Straps-21"L (length over top of foot, pair) \$46
- PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair) \$117
- PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair) \$117

Seats

- PNG30269 Transfer Seat (replaces standard planar seat) \$139
- PNG50324 Rotating Seat (replaces standard planar seat) \$241
- PNG30183 Seat with Roho® Insert (replaces standard planar seat) \$609

Seat Accessories

- PNG50161 Hip Supports (9"-15" range, pair) \$363
- PNG50198 Hip Supports (11"-17" range, pair, limited adjustment with shadow tray) \$363
- PNG50518 Hip Supports (16"-22" range, pair, limited adjustment with shadow tray) \$363

Backs

- PNG50066 Flat Back-15"H \$665
- PNG50067 Flat Back-19"H \$702
- PNG50062 Contoured Back-15"H \$732
- PNG50068 Contoured Back-19"H \$769
- PNG50171 Removable Flat Back-15"H \$763
- PNG50173 Removable Flat Back-19"H \$800
- PNG50172 Removable Contoured Back-15"H \$830
- PNG50174 Removable Contoured Back-19"H \$868

Positioning and Upper Body Support

ATTENTION - All options below require a back

- PNG30029 Velcro® Positioning Belt-26"-55" \$57
- PNG30306 Velcro® Positioning Belt-38"-67" \$57
- PNG30028 Positioning Belt with Airline Style Buckle-26"-55" \$153
- PNG50192 Chest Strap-30"-46" \$142
- PNG50193 Chest Strap-34"-54" \$142
- PNG50504 Chest Strap-46"-66" \$142
- P600002 Back Angle Adjustment Wheel \$104

ATTENTION - All options below require the accessories mounting bracket

- PNG50034 Accessories Mounting Bracket \$68
- PNG30170 X-Style Chest Vest-11"Lx9.5"W \$277
- PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$277
- PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$277
- PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$277
- PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W \$323
- PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W \$323
- PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W \$323
- PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W \$323
- PNG11212 High Mount Chest Vest Bracket \$78
- PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair) \$363
- PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair) \$363
- PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair) \$363
- PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair) \$363

Head Support

- PNG50240 Head Support-Short (back and accessories mounting bracket required) \$363
- PNG50037 Head Support-Tall (back accessories mounting bracket required) \$363

Other Accessories

- P80246 Angle Locator \$36
- PNG50040 Push Handles (back and accessories mounting bracket required, pair) \$153
- P80809 Tool Pouch \$25

* not available with mobile

Notes:

AltimateMedical

Ph: 800.342.8968 or 507.697.6393
 Fax: 877.342.8968 or 507.697.6900
 info@easystand.com easystand.com

Email form to webquotes@easystand.com